

Laboratory Equipment, Instruments and Supplies

2016 International Catalog



WORLD PRECISION INSTRUMENTS, INC.



Sarasota International Trade Center 175 Sarasota Center Blvd. Sarasota, FL 34240 P (941) 371-1003 · F (941) 377-5428 www.wpiinc.com

Dear Valued Customer,

First and foremost, thank you for your continued patronage. Our core values are driven by a desire to empower you to actualize your scientific ideas with cost-effective, quality instruments. We always strive to provide you more for less.

This year, we have once again renewed our commitment to quality and value, which is reflected in our new look. We have implemented measures (like our ISO-9001:2008 certification) throughout our home office and production facility to improve our efficiency and ensure quality for the benefit of our customers. We strive to operate with the highest standards of integrity and efficiency.

Since our earliest days, WPI has designed laboratory equipment for researchers. Our founder Harry Fein began developing electrophysiology equipment for Yale University researchers in the seventies. Nearly 50 years ago (1967) WPI was born out of Harry's passion for the advancement of science. Since then, we have worked hand-in-hand with researchers to create the equipment needed to facilitate modern science and advance biomedical research. We continue to offer the instruments that empower your scientific ideas.

Now, we are a multi-national company with offices in four countries and distributors all over the world. Originally, we specialized in electrophysiology, and now we are in several areas of study, the core being in tissue and cell biology, animal physiology and electrophysiology. Additionally, we offer a full range of quality surgical instruments, general laboratory equipment and supplies, such as microscopes and pipettes.

We look forward to serving you for many years to come. If you have questions or need assistance ordering, our customer service team is ready to help.

Best regards,

Cliff Bredenberg

WPI President and CEO

Muscle Physiology
Cell Tester ● Vertical Tissue Bath ● Muscle Tester ● Isolated Heart Syste
Surgical Instruments
Scissors ● Tweezers ● Ear Punches ● Vessel Clips ● Mouse Surgery Kit ● Electrosurgical Unit ● Cautery Instrument
Animal Surgery & Support
Anesthesia ● Analgesia ● Monitoring ● Blood Pressure Measurement ● Stereotaxic Syste
Biosensing
NO, CO, H₂O₂, O₂ and H₂S Detection ● Ion Selective Electrodes ● pH Mete
Amplifiers, Electrometers
Intracellular ● Extracellular ● Transducer ● Metal Microelectrodes ● Voltage/Current Clam
Epithelial Physiology
Epithelial Voltohmmeter ● TEER Measurement ● Ussing System ● Perfusion System
Stimulators
Digital Stimulators ● Multi-Channel ● Linear Stimulus Isolator ● High Current Isola
Laboratory Supplies 9
Adhesives ● Cables & Connectors ● Metal Wire ● Luers & Fittings ● Laboratory Bases, Rods and Clam
Glass, Holders, Electrodes 10
Capillary Glass ● Electrode Holders ● Metal Microelectrodes ● MicroFil Non-metallic Syringe Need
Microforges, Pullers 11
Pullers ● Bevelers ● Microforg
Micromanipulators, Stands 12
Manual & Motorized Micropositioners ● Magnetic Stands ● Tilting Ba
Microscopes, Cameras 14
Surgical, Stereo, Inverted Microscopes ● Cameras ● Illuminators ● Heater Controller ● Culture Dish
Pumps, Microinjection 15
Peristaltic, Syringe and Pneumatic Pumps ● Microprocessor Controlled Injector ● Microfluid
Spectroscopy 18
Spectrometers ● Sample Cells ● Light Sources ● Fiber Optic Cables ● Fiber Optic Dipping Pro
Index 20

Index \cdot Ordering Information \cdot Terms and Conditions





Cell Tester System

The **SI-CTS200** system is a revolutionary new research tool for studying single living cells. This system, which is the result of blending the latest technologies in electronics, mechanics and optics, permits researchers to investigate living systems at a new level using an innovative experimental paradigm. The Cell Tester provides researchers with the comprehensive ability to investigate and characterize the physiological, biomechanical and biophysical properties of single isolated living cells, and to extend these findings to both the sub-cellular and multi-cellular levels. The SI-CTS200 systems can be used on single living cells, small multi-cellular preparations or skinned muscle fibers and strips without modification. It is designed to sit on the stage of any standard, research-level, inverted

microscope while maintaining the optical path of the microscope for simultaneous fluorescence or confocal imaging.

The Cell Tester is an integrated system of components needed to maintain and handle single living cells, stimulate and perform pertubations of the cells, and detect, amplify, and record signals, like contractile force, from the cells.

- Optical transducer for measuring force with nN sensitivity and integrated vacuum attachment system for cell attachment. Lifetime warranty on optical heads.
- Equipped with a nanomotor for stretching and relaxing cell with nm resolution.
- A rotating cuvette system for easy alignment of cells increases productivity. Interchangeable

bath inserts provide a range of options for the handling of live cells.

- Signal conditioning amplifier system includes a force transducer amplifier with multiple gains, a position controller for moving the nanomotor used to stretch and relax cells and to open and close the microtweezers, a 2-channel temperature controller, and an antioscillation unit for eliminating the resonance frequency of the transducer and its mounting support.
- Lab-Trax-8/16 data acquisition system with MDAC software records force signals, motor positions and other data through four analog inputs. It is designed for controlling the stimulator and the position of the linear motor through two analog outputs and other devices through two digital outputs.



The manipulators tilt up and fold back to facilitate system setup. There are two programmable memory positions (home and target) for easy exchange of the 35mm dishes, providing access to the bath so you can add the live cells. The combination of these features enables high experimental throughput.



The SI-CTS200 system utilizes a unique rotating bath to dramatically improve experimental throughput. It is designed to orient cells in the XY plane so that no physical manipulation of the position of the cell itself is required prior to capture by the grabbing devices attached to the force sensor and linear actuator.

This bath has two interchangeable inserts. The first holds any 35mm glass bottom dish (WPI #FD35-100). The second is a native cuvette insert containing the live cells.

Signal Conditioning Amplifier System with CTS200 Electronics





SI-CTS200 Complete Cell Tester System

System Includes: Rotating Cuvette System; Micromanipulator System; Signal Conditioning Amplifier with four modules: Optical Transducer Amplifier; Temperature Controller; Anti-Oscillation Unit; Position Controller; data acquisition system with MDAC Software for recording, controlling stimulation and nanomotor position; Force Transducer of choice; Nanomotor; Glass Fiber Cell Mounts (1 set); MyoTak Bio-adhesive kit (5-week supply)

SI-CTS200B Cell Tester, Non-Rotating, No Micromanipulators SI-CTS200A Cell Tester, Manual Platform, No Micromanipulators

New Stretch Dependent Signaling Pathway Discovered in Cardiac Ventricular Myocytes termed "X-ROS"





A confocal microscope in linescan mode is used to record calcium spikes. The healthy cell in the top image shows a marked increase in number and intensity of calcium sparks over the duration of the stretch. The diseased muscle in the bottom image shows fewer sparks when at rest and a profound increase in sparks of short duration during the stretch.

Introduction

- Control of cardiac Ca²⁺ release is critical for the regulation of contraction and maintenance of electrical activity.
- A new signaling pathway was discovered: X-ROS signaling.
- It regulates normal Ca²⁺ release in healthy heart cells and may drive pathologic Ca2+ release in diseased cells.


Results

- Stretch triggers the generation of ROS by NOX2
- Stretch-dependent process activates NOX2 production of ROS which reversibly oxidizes nearby RyR2 ryanodine receptors.
- X-ROS oxidation tunes the sensitivity of the RyR2s, increasing the Ca2+ spark rate and enhancement of Ca2+ signaling.
- During a sustained stretch of a cardiomyocyte in a Cell Tester, a rapid elevation of ROS production subsides over the duration of stretch.
- During a repetitive cyclical stretch, a new level of steady-state ROS production is maintained.

Conclusion

The level of steady state ROS generation in the cell may be graded by diastolic length or pre-load.

REFERENCE: X-ROS signaling: rapid mechano-chemo transduction in heart. Prosser BL, Ward CW, Lederer WJ. Science. 2011 Sep 9;333(6048):1440-5. PMID:



A single skeletal muscle cell is held with microtweezers. The top image shows the unstretched cell, and the bottom shows the stretched cell.



Muscle Research System

The **SI-MKB** is the next research platform in the progression of studying living systems from the cellular level to whole organisms. It is the standard muscle research system, which can be configured to study skinned and intact muscle fibers, muscle strips and small whole muscles. The modular design of the SI-MKB Muscle Research Platform allows it to be configured for turnkey solutions for specific applications. The system is built on a solid platform making precise mechanical and optical measurements easy. Like the **Cell Tester** and the Horizontal Tissue Bath systems, the SI-MKB uses SI-KG optical force transducers and is constructed with corrosionfree materials (stainless steel, anodized aluminum and plastic).



Can be configured with specific components, like a linear motor or a stimulator, for measuring:


- Responses to electrical stimulation, including tetany
- Mechanical properties of contracting and relaxing muscle strips
- Isotonic force and the effects of constant
- Twitch amplitude and kinetics, like contraction and relaxation times and velocities
- Slack, quick stretch-release, constant velocity and eccentric contraction tests
- Vibrational studies that simulate unloaded muscle shortening
- Quick temperature-change effects on skinned and intact muscle fibers

Can be equipped with components, like a laser diode or a photometer, for measuring:

- Sarcomere length/spacing by laser diffraction at the same time muscle force is measured
- Intracellular calcium concentration/ distribution in intact muscle fibers as muscle force is measured
- ATPase activity in skinned muscle fibers as muscle force is measured

Supplied with a Labtrax 8/16 data acquisition system and MDAC software for:

- Recording signals from the force transducer and the motor position monitor, through four analog inputs
- Controlling the stimulator, or an external stimulus isolator, through an analog output
- Controlling the position of the linear motor through a second analog output
- Controlling other devices through two digita



The Sarcomere Spacing Assembly has the optical cuvette, laser diode array and camera in one unit. It also includes the electronic control module. See page 10.


SOME SYSTEMS AVAILABLE FOR INTACT MUSCLE FIBERS INCLUDE:

SI-MKBM System with a 5cm heated cuvette and a linear motor SI-MKB System with a 5cm heated cuvette and a digital micrometer

Systems include: Base plate; force transducer and stand, digital micrometer or linear motor (optional), micrometer or motor stand, cuvette and table, oxygenation system, evacuation system, signal conditioning system with transducer amplifier and temperature control module, anti-oscillation module and linear motor amplifier (if applicable), Lab-Trax 8/16 data acquisition system with MDAC software for recording force and position signals and for controlling stimulation and the movement of the linear motor.

Horizontal Tissue Bath

Research system for higher throughput of complex pharmacological/physiological assays



The new SI-H Horizontal Tissue/Organ Bath system (SI-HTB) combines the ease of use and productivity of a traditional vertical organ bath with the more advantageous features of single tissue physiology platforms.

- Two channel system for increased productivity, easily expanded to add channels
- Fully independent heating and fluid control for each channel
- Low profile/small footprint
- Variable volume, chemically inert teflon bath with shape configurations from variable to fixed 500µL-10mL
- Modular, space-saving, blade-style electronics (Control up to 4 channels with the electronics in one chassis)
- Large variety of force transducers covering mN-N forces, all with lifetime warranty
- Can be combined with automated fluid control systems
- Add a linear motor and controller (SI-MOT) to perform mechanical, electrophysiological and optical techniques, inculding isometric (standard), isotonic, eccentric and auxotonic
- Add an electrometer like the WPI Duo773 or Electro705
- Add the multi-purpose Biofluorimeter (SI-BF-100) for tissue fluorescence (calcium, NOx, ROS). See page 6.

Multiple Bath Options

The **SI-HTB** system breaks through the (large) volume limitations of the traditional organ bath, allowing volumes as low as 500µL in an inert, teflon-based bath.

The bath design allows multiple shape options for thick, long, flat and thin tissue. When pharmaceuticals are available in precious, small amounts, you will appreciate this standard feature. A wide range of transducers and tissue mounting supports complement this freedom of tissue shape, volume and size.

Low Profile

The low profile and small footprint of the bath system, combined with the modular, space-saving, chassis-mounted design of

the electronics, reduces the bench space requirement up to 4-fold when compared with standard 4-channel organ baths.

Versatile System

The SI-HTB combines advanced physiological techniques with the throughput needed in pharmacological assays in one flexible platform. Upgrades to four or more channels are easy and economical.

The motor option (SI-MOT) turns your system into a tissue work-out station with isotonic, auxotonic and eccentric force measurement capabilities. Nearly all established myomechnical tests from stretch-release to workloops and muscle fatigue are now possible in a single organ bath system. Some of these procedures require a length change solution (software/hardware) like WPI's MDAC package.

Future is Now

The solid horizontal tissue bath design is ideal for combination with electrophysiology on the same platform. Intracellular measurements can share the stable solid base of the bath system.

WPI's new fiber-optic based, multi-channel Biofluorimeter allows for tissue fluorescence measurements (calcium, NO, ROS) on the SI-HTB platform.

Now, you can design a system to meet your needs and budget. And, it is fully upgradeable in the future.

Options

WPI's 16-bit, high speed, Labview-based Muscle Data Acquisition system SI-LABTRAX-**MDAC** is perfect for this platform.

SI-HTB2 Horizontal Tissue Bath, 2-Channel System

2-Channel **SI-HTB** platform for isometric force (1) SI-KGX Force Trandsucers (2) SI-BAM21LCB Optical Force Transducer Amplifiers (2) SI-TCM2B 2-Channel Temperature Controller (1) Signal Conditioning Amplifier System Chassis (1)

SI-HTB4 Horizontal Tissue Bath, 4-Channel System

2-Channel **SI-HTB** platform for isometric force (2) SI-KGX Force Trandsucers (4) SI-BAM21LCB Optical Force Transducer Amplifiers (4) SI-TCM2B 2-Channel Temperature Controller (2) Signal Conditioning Amplifier System Chassis (1)

Horizontal Tissue Bath, 2-Channel Motorized System SI-HTB2M

2-Channel **SI-HTB** platform for isometric force (1) SI-KGX Force Trandsucers (2) SI-BAM21LCB Optical Force Transducer Amplifiers (2) SI-TCM2B 2-Channel Temperature Controller (1) SI-MOT Linear motor with controllers (SI-MOTDB) (2) SI-AOSUB Anti-oscillation Unit (2) Signal Conditioning Amplifier System Chassis (1) MDAC Data Acquistion software (1)

SI-HTB4M Horizontal Tissue Bath, 4-Channel Motorized System

SI-KGX Force Trandsucers (4) 2-Channel **SI-HTB** platform for isometric force (2) SI-BAM21LCB Optical Force Transducer Amplifiers (4) SI-TCM2B 2-Channel Temperature Controller (2) SI-MOT Linear motor with controllers (SI-MOTDB) (4) SI-AOSUB Anti-oscillation Unit (4) Signal Conditioning Amplifier System Chassis (2) MDAC Data Acquistion software (1)



Biofluorometer

Now even more reliable. simplified and affordable

The new SI-BF-100 is an LED-based fluorometer for life science applications. It is ideally suited for ratiometric calcium detection (FURA-2) and ATPase detection (via NADH fluorescence). With up to seven LED modules (wavelengths), the SI-BF-100 covers many fluorometric applications in neuroscience and cell biology.

Recent advancements in optics and LED technology simplify ratiometric calcium imaging, making this equipment more affordable. A breakthrough in WPI patented technology allows the SI-BF-100 to use wavelengths below 380nm and produce more light in those spectra. This technology significantly cuts the cost of photometric calcium imaging without sacrificing resolution

- LED light sources require less power, give off less heat and are more compact and affordable
- Sampling rates up to 1kHz (1,000 ratios/ second maximum). At lower speeds, signal averaging is used for noise reduction.
- Two auto ranging photomultiplier inputs allow you to monitor multiple wavelengths from a single emission output that can be comprised of any wavelength of light for which an LED module is available
- Using a separate reference channel, ultra-stable, continuous ratio calculations automatically compensate for LED intensity drift. This ensures less noise and produces more accurate measurements.
- Application-specific probes are available for existing tissue baths and cuvette systems.
- Ratio noise is <0.05 peak to peak, drift is less than 0.1 unit/hour
- The warm up time of less than one minute is a dramatic improvement over the common 20–60 minutes required by xenon or mercury light sources
- Replace the emission filter easily or change the LED modules to transform the SI-BF-100 into a general purpose fluorometer for many other applications



How it Works

Up till now, calcium imaging systems have been required to compensate for errors and noise introduced by the complexity of their design. The systems require mechanical filters and use expensive xenon or mercury light sources. The beauty of the SI-BF-100 is its simplicity. The elegance of its design reduces the noise introduced into the system and the errors inherent in traditional designs.

Monochromatic LED light sources using WPI patented technology eliminate the need for complex and expensive white light sources

and filter wheels. Because the LED modules can be pulsed, sampling frequencies up to 1,000

cycles per second are possible.

The LED light source emits specific excitation frequencies which travel through the probe. The excitation light can be comprised of any wavelength of light for which an LED module is available. The probe returns a single emission

output to one or two photomultiplier inputs on the front of the SI-BF-100, which are independently filtered for specific wavelengths. This design allows you to monitor multiple wavelengths from a single emission output.

The LED light source in the WPI design makes this ratiometric fluorometer more compact, energy efficient and affordable. As added benefits, the low-power light source produces much less heat, and it warms up in less than one minute!

This incredible design is not limited to calcium imaging. By simply replacing the

emission filters in front of the photomultipliers with the desired wavelength filters, your SI-BF-100

becomes a general

purpose fluorometer for any application you can imagine. Changing a filter involves removing the two screws that hold the filter carriage on the face of the SI-BF-100, swapping the filter and reinstalling the integral SMA/filter carriage.

SI-BF-100	Biofluorometer
OPTIONAL CO	DMPONENTS
M3301R	Manual Manipulator, right-handed
M3301L	Manual Manipulator, left-handed
M10	Magnetic Base

Unique fiber optic coupling probe

allows for highly efficient transfer of

light and ease of placement. Custom

probes fit your existing systems.

SI-BF-100 SPECIFICATIONS

FIBER OPTIC LIGHT INPUT/OUTPUT	.SMA terminated
BANDWIDTH	.1000 ratios/second
RATIO NOISE	.< 0.1 peak to peak

ANALOG OUTPUT RANGE0–10V (continuous, equivalent to a ratio 0–10)

CALL FOR APPLICATION

In-Vivo Applications of SI-BF-100 Using Voltage Sensitive Dyes

- Monitor electrical activity over large areas of intact brain without the limitations of microelectrode arrays
- Simultaneously monitor electrical activity and calcium concentration
- Genetically-encoded fluorescent probes make it possible to monitor neurotransmitter behavior over entire regions of intact tissue
- Integration with miniature fiber optic arrays or micro-scale endoscopes allow for measurements on animals while conscious and mobile

The new SI-BF-100 is an LED-based fluorometer for life science applications. With up to seven LED modules (wavelengths), the SI-BF-100 covers many fluorometric applications in both Neuroscience and Cell Biology. This technology significantly cuts the cost of fluorescent imaging without sacrificing resolution or quality.

Using the SI-BF-100 Biofluorometer equipped with high intensity LED modules and an appropriate fiber-optic probe, a researcher can perform many different types of analysis on intact tissue in vivo. Some potential applications include simultaneous measurement of membrane potential and calcium concentrations. With the use of genetically-encoded fluorescence probes,



WPI's SI-BF-100 Biofluorometer probes include both the excitation and emission light fibers. See page 199.

the applications are limited only by the imagination of the researcher.

The use of a probe allows for the measurement of entire areas of the intact brain without the limitations presented by the placement of arrays of microelectrodes. By removing these limitations, the researcher is able to collect data over larger groups of neurons without placing a single electrode.

Miniature fiber-optic arrays and microscale endoscopes have already been developed for use in Neuroscience. Coupling this technology with the SI-BF-100 allows the research to be conducted on conscious and mobile animals. By not requiring the animal to be attached to a traditional light source, more accurate results can be obtained.

The SI-BF-100 eliminates the need for expensive and complicated arrays of microelectrodes to conduct studies in vivo. Larger areas can now be completely analyzed at a fraction of the cost with a significant reduction in experimental complexity.


Microscopic Imaging Applications Using Voltage Sensitive Dyes

- Monitor membrane potential over small networks in cell culture
- Monitor surface membrane potential in a single cell
- Simultaneously monitor electrical activity and calcium concentration
- Integration of fluorescence microscopy onto existing patch-clamp setups
- Genetically-encoded fluorescence probes make it possible to monitor neurotransmitter behavior between cells in cell culture samples
- With the addition of a CCD, traditional spatial imaging can be conducted

The new SI-BF-100 Biofluorometer is an LED-based fluorometer for life science applications. With up to seven LED wavelength modules (three when using the high power Biofluorometer), the SI-BF-100 covers many fluorometric applications in both Neuroscience and Cell Biology. This technology significantly cuts the cost of fluorescence imaging without sacrificing resolution or quality.

Using the SI-BF-100 Biofluorometer equipped with high intensity LED modules and coupling to an inverted microscope, a researcher can perform many different types of analysis on both single cells and small cell cultures. Some potential applications include simultaneous measurement of membrane potential and calcium concentration. With the use of genetically-encoded fluorescence probes, the applications are limited only by the imagination of the researcher.

The SI-BF-100 uses LED technology, eliminating the need for expensive secondary external light sources and project-specific filter cubes. With the potential for seven different LED modules inside every Biofluorometer, more fluorophores can be analyzed simultaneously, without the need for expensive filter wheels and external timing devices.



The Biofluorometer can connect to the epifluoresence port of a microscope, and its high intensity LED light source is used for the illumination.

Signal Conditioning Amplifier System

Choose the amplifier modules you need to measure nearly anything!

- Ergonomic design
 - 8-Channel
- Small footprint
- Backplane design includes provision for configurable communication between modules

ignal conditioning in the last decade has increasingly moved in the direction of the computer, because the software preserves the raw data and is highly reliable. Operations such as integration, differentiation, filtration and even waveform generation are now efficiently handled in software. On the other hand, however, the transduction of physical signals such as bio-potentials, force, temperature, pressure or ionic concentrations must be measured with an electronic amplifier.

SIH/WPI's physiology amplifier system focuses on this idea and provides a flexible electronic platform intended to process the transduction of physical signals, displacement transducer outputs and the outputs from electrochemical free radical sensors. This platform simply focuses on the reliable transduction of the electronic signal and provides a convenient passage for the translation of real world signals to a computer for analysis.

The system consists of an 8-channel frame that includes an ultra quiet, shielded power supply. All of the module outputs are routed to rear panel connectors. If you prefer, outputs may be routed internally to the inputs of other modules. The system has a small footprint and may be stacked to provide as many channels as you need.

The SI-BMFA Power Frame is the foundation of the SI modular physiology suite. It incorporates a robust power supply that can accommodate up to eight physiology modules, which can be mixed or matched in any combination. Modules are quick and easy to install, thanks to an innovative and mechanically solid card rail system.

When the system is ordered with SI-MKB (Muscle Tester) system, the Signal Conditioning Amplifier System (chassis) is configured with an SI-BAM21-LCB. Optional modules include an SI-TCM2B Temperature Control Module, an SI- MOTDB Linear Motor Controller, an SI-PF100 Programmable Filter Mondule, the SI-SARCAM Sarcomere Spacing Module and the SI-

COLUB Constant Load Unit. The Temperature Control Module, Linear Motor Controller and Sarcomere Spacing Module require two slots each on the chassis backplane.

The system is flexible and configurable. A variety of modules are available for the Signal Conditioning Amplifier System, and you can mix and match the modules to suit your requirements.

Modules currently available include:

- SI-BAM21-LCB Optical Transducer **Amplifier**
- SI-PF100 Programmable Filter Module
- SI-SARCAM Sarcomere Spacing Module
- SI-MOTDB Linear Motor Control Module
- SI-TCM2B Temperature Control Module
- SI-COLUB Constant Load Module



SI-BMFA

Power Frame Enclosure



Optical Transducer Amp

The SI-BAM21-LCB KG Optical Force Transducer Amplifier is used in conjunction with the SI-H muscle physiology systems. The SI-BAM21-LCB powers the force transducer and converts the output of the transducer to an amplified analog voltage that is proportional to the force applied to the force transducer. The output signal can be multiplied by a factor of 1, 2, 5 or 10 to provide better resolution for a minimal change in applied force.

NOTE: An optional factory setting increases the multiplier by a factor of 10, allowing the signal to be multiplied by 10, 20, 50 and 100.

The SI-BAM21-LCB amplifier works with KG optical force transducers to:

- Generate an analog output (-10VDC to +10VDC) that is proportional to the force applied to the tissue sample.
- Supply a DC voltage that powers the KG force transducer to which it is connected.

Features

- Rapid Auto Zeroing function with fine offset adjustment
- Offset indicator LED's
- Multiple gain ranges with adjustable fine tuning for precise calibration

• Digital interface for optional Anti-Oscillation module for use in isotonic studies

Also available in a single standalone enclosure, either version provides an incredibly quiet, linear and stable transducer signal to your data recording system.

How the SI-BAM21-LCB Amplifier Works

In a typical setup, a muscle is held by a force transducer. The force transducer is connected to the SI-BAM21-LCB. As the muscle contracts or releases, the transducer converts the force into an electrical current signal which is proportional to the force applied to the transducer. The SI-BAM21-LCB converts the current signal into a voltage signal that can be displayed on the screen of the recording device.

Before initiating an experiment, the SI-BAM21-LCB must first be zeroed. This sets the baseline for measurements to follow.

The output signal is buffered and multiplied by 1, 2, 5 or 10, depending on the Gain switch setting on the front panel of the amplifier module. The X10 setting is useful when output signals are extremely small. Finally, the force proportional signal is sent through the output amplifier

The analog output has a range of –10V to +10V that drives the **LABTRAX-MDAC** data acquisition system, multimeter or oscilloscope.

NOTE: When the Signal Conditioning Amplifier System is configured at the factory for an **SI-HTB** or **SI-MKB** Muscle Tester system, the signal is routed internally from the SI-BAM21-LCB module to the SI-AOSUB module.



The stand alone SI-BAM21-LC Force Transducer Amplifier has all the capabilities of the Signal Conditioning Amplifier module (SI-BAM21-LCB).

ACCESSORIES

2851	BNC Cable
SI-KG2	0-2N Force Transducer
SI-KG2B	0-0.2N Force Transducer
SI-KG4	0-50mN Force Transducer
SI-KGxx	contact us for specialty tranducers with different ranges

See Optical Force Transducers on page 13.

Optical Transducer Amplifier SI-BAM21-LCB SI-BAM21-LC Standalone Optical Force Transducer Amplifier

SI-BAM21-LCB SPECIFICATIONS

Input Configuration Current to voltage converter

Gain 1x, 2x, 5x, 10x

Optional factory setting: 10x, 20x, 50x, 100x Input Offset Adjustment + 2 0 VDC 470 Q Output Impedance

12 VDC provided by chassis Power

± 10 VDC Output Range

Signal Conditioning Amplifier System

Sarcomere Spacing Module


The SI-SARCAM Sarcomere Spacing Module was designed for the SI-MKB and SI-HTB systems. It consists of a linear diode array detector coupled with a calibrated controller (shown below). Together with the laser diode system, it facilitates the calibrated measurement of sarcomere length in striated muscle samples. The SI-SARCAM can be combined with mechanical intact muscle parameter measurements with the SI-MKB and SI-HTB systems. It can be combined with calcium measurements and myomechanics in intact and skinned muscle fibers when you are using the SI-MKB system. The SI-SARCAM



module requires two of the eight module positions in the chassis.

How It Works

A red laser diode (λ =650nm) is used. A onedimensional (linear) CCD camera with a time resolution of 1ms scans the light intensity of the diffraction pattern beginning 6.8mm from the center beam. The camera scans only half the diffraction pattern, assuming that the pattern is symmetrical with respect to the center beam. When the scan arrives at the part of the CCD chip which is hit by the first order diffraction of the pattern, the output voltage increases to a peak level, and the first order of diffraction is captured. Then, the intensity signal decreases again. The sarcomere length is calculated from the first order diffraction distance. See the full article at



www.wpiinc. com/sarcam. A cross hair mark used for aligning the laser is etched on the front of the camera, which is mounted just behind the optical cuvette.

SI-SARCAM

Sarcomere Length Control Module

Temperature Control

The SI-H Temperature Control Unit is designed for use with the SI-H line of muscle physiology research platforms. It maintains the temperature of an SI-H cuvette up to 45°C. It is accurate to 0.1°C. The circuit is appropriate to RTD (resistive temperature device) applications. It linearly converts a temperature reading to a voltage that is displayed as a temperature on the SI-TCM2 and can



be recorded. This unit is available in a standalone model and as a module for the Signal Conditioning Amplifier System backplane.

- Controls two cuvettes simultaneously
- Uses PID control to maintain a constant temperature with ±0.1°C tolerance
- Has both high and low alarm warnings which can be user defined
- Easy to control with simple interface
- Also available as a stand-alone device call for details and price.

SI-TCM2B

2-Channel Temperature Control Module

2-Channel Temperature Control Standalone

SI-SARCAM SPECIFICATIONS

INPUT CONFIGURATION Current to voltage converter

LASER Red laser diode

LASER WAVELENGTH 650 nm

CAMERA Linear CCD Camera

POWER REQUIREMENTS 12V DC provided by the chassis

MAXIMUM POWER CONSUMPTION 1.3A AT 115V 50/60HZ,

1.8A AT 230V 50/60HZ

See the Sarcomere Spacing Assembly on page 4.

SI-TCM2 SPECIFICATIONS

INPUT CONFIGURATION Current to voltage converter

POWER REQUIREMENTS 12V DC at 2.5A 50/60Hz wall adaptor.

2.5mm ID/5.5mm OD with positive center DC barrel (included-WPI

#801513)

OPERATING TEMPERATURE RANGE Room temperature

DISPLAY PRECISION 0.1°C CONTROLLER RESOLUTION 0.1%

CUVETTE TEMPERATURE SENSOR 1000Ω RTD (1000Ω at 0°C)



Constant Load

The SI-COLUB Constant Load **Module** for performing constant load experiments, has augmented flexibility. In its primary mode (Constant Load) the unit takes an external trigger command from the force transducer to perform a constant load cycle. In addition, the module allows for a different external trigger or you can completely bypass the module without having to switch cabling.

The SI-COLUB Constant Load Module lets you maintain a constant force, muscle length or sarcomere length rather than keep the total length of the preparation constant during an isometric contraction. This is accomplished using a feedback loop.

The SI-COLUB monitors a designated parameter to determine how much force is necessary. It also monitors a feedback signal. The motor position command signal driving the motor is constantly adjusted to drive the feedback signal to the commanded setpoint.

- Offers three modes including Constant Load, Eternal Loop and Bypass
- Can be configured using MDAC software and Lab-Trax 8/16 data acquisition system

Constant Load Mode SI-COLUB External Loop Mode Force Step Bypass Mode Command Measured Force Linear Motor (from SI-AOSUB or SI-BAM21-LCB) **Bypass** External Trigger Feedback

SI-COLUB Constant Load Module

SI-COLUB SPECIFICATIONS

COMMAND REQUEST +10V **FEEDBACK** ±10V MOTOR OUTPUT +10V

POWER REQUIREMENTS 12V DC provided by the chassis

SI-PF100 Programmable Filter Module

When you use a motor, an SI-PF100 Programmable Filter is necessary to minimize the natural vibration. It is designed so you can eliminate the resonance frequency without affecting the signal of interest. It is a low pass filter set to pass signals of interest below the specified frequency. It can be calibrated from 5 to 1,000 Hz.

You may select either a Bessel or a Butterworth filter. Then, you must carefully select the cutoff frequency based on the typical resonance frequency of your force transducer and your own experimental setup.

When the Signal Conditioning Amplifier System electronics are configured at the factory with an SI-PF100 Programmable Filter, the signal is routed internally from the amplifier (SI-BAM21-LCB module) to the SI-PF100. If you prefer, the signal may be routed from the amplifier through the ports on the front panel of the Programmable Filter using a standard BNC cable.

- Bessel and Butterworth
- Programmable cutoff frequency
- Signal may be routed through the backplane or through the front panel **BNCs**





SI-PF100

Programmable Filter Module

SI-PF100 SPECIFICATIONS

POWFR 12 VDC provided by chassis

INPUT ± 10 VDC **CUTOFF FREQUENCY RANGE** 5-1000Hz

FILTER TYPES Bessel, Butterworth

Signal Conditioning Amplifier System

Linear Motor Control

The SI-H Linear Motor Controller is designed for use with the SI-H line of muscle physiology research platforms. For systems that require a linear motor, this unit provides the precision control of the motor. A linear motor is required for measuring mechanical muscle properties such as slack-test, isotonic release, constant velocity release, stretch release, vibration studies, after-loaded contractions and eccentric contractions (intact muscle). The position of the linear motor is determined by a combination of the data from the controller indicating the current position and the DC value applied to the front panel at the Position In port. The applied Position In signal can be provided by a data acquisition system (Lab-Trax 8/16). The data acquisition analog output signal is set to define the waveform and timing pattern of force to be applied to the

The Linear Motor Controller has been designed with an automatic shutoff feature. That means that the voltage driving the motor automatically shuts off if the motor draws too much current. After less than a second, the motor cycles back on again. If the setup still draws too much current, it repeats the power down cycle. This could happen if too much force is being applied to the sample. The system continues to cycle the motor off and on until the force on the motor is reduced. If this happens, the motor hums as if it is trying to work, but the motor produces no force output. If an auto-shutdown occurs, adjust the experiment and force tension on the motor.

The motor is setup and calibrated to the Linear Motor Controller that is shipped with it. The motor and controller have a maximum range of motion of ±3.5mm from the center of travel.

The Position In BNC interface provides for an external position command input. The input range is ±10VDC. The signal presented at the input will affect the length of motor travel. The motor/controller combination are calibrated for travel of 1.0mm/2.0VDC. For example, +7VDC translates to +3.5mm of travel from the center, and -7VDC equals -3.5mm of travel from the center. This input could come from the Lab-Trax-8/16 Data Acquisition system or an SI-COLUB Constant Load Module.

NOTE: Use of a linear motor with a muscle testing platform creates vibration which excites the sensor's resonant frequency and requires an Programmable Filter module (SI-PF100) for best results.

- Powers the motor and provides an output indicating the actual motor position
- Connects to Analog to Digital Converter output of the computer or data acquisition system (like Lab-Trax-8/16) to allow the programming to control the waveform and timing of the motor control



The SI-MOT motor can be used with any SI-HTB Horizontal Tissue Bath or SI-MKB Muscle Research System.



- Input range of ±10VDC
- Over current protection that automatically shuts down when the supply voltage dips below the reference value
- Linear motor position is determined by a DC value applied from the Position In port

Linear Motor for SI-HTBM system SI-MOT-MT SI-MOT-MKB Linear Motor for SI-MKB SI-MOTDB Linear Motor Controller

SI-MOTDB SPECIFICATIONS

POWER REQUIREMENTS 12V DC provided by the chassis

INPUT ±10V DC

1 mm/2 VDC current TRAVFI

MAXIMLIM TRAVEL 7 mm (±3.5 mm from center of travel)

KG Optical Force Transducers



Simple calibration

- Different models to accommodate a wide range of forces and sensitivities
- Nearly insensitive to changes in temperature and ambient light
- Extremely high level of linearity
- Virtually indestructible with normal use
- Lifetime warranty

T	R/	١N	SE	U	CE	R	SP	E	CI	FI	C	41	ΓIO	N:	S
- 1	,		1			.,1		, ,						, .	

Unloaded transducer without tissue mounting support

	Force Range	Range	Noise	Compliance (nm/mN)	Resonance Frequency	
SI-KG2	0-2 N	0-200 g	250 μΝ	150	1.3 kHz	
SI-KG2B	0-0.2 N	0-20 g	80 μΝ		590 Hz	
SI-KG4	0-50 mN	0-5 g	15 μΝ	0.5	1.2 kHz	
SI-KGxx	contact us j	for specialt	y transduc	ers with differen	t ranges.	

KG transducers are required for use with SI-BAM21-type


Tissue Mounting Hooks

Mounting hooks can be used in a variety of combinations, depending on the type of tissue to be examined.

Mounting hooks are sold in kits. Currently, there are 11 kit configurations, each available in four different sizes. The mounting hook size that is required depends on the force transducer used.



*The **97909** tissue mount adapter tube (OD:0.096", ID:0.035") allows you to use SI-KG4 size mounting hooks with SI-KG2 size force transducers.



Motor/micrometer mount pictured on the right.

Vascular hooks are available for mounting blood vessels (rings). They are normally used with a pair of blunt hooks (SI-TM8).



The micrometer and motor receive a large (SI-KG4 size) tissue mount. If a smaller tissue mount is used, the 97909 adapter is required. This adapter is included with every SI-MKB or SI-HTB system.

Ordering

When ordering tissue mounts, specify the tissue mount configuration and force transducer to be used.

	Force Transducer Mount	Micrometer/Motor Mount	Force Transducers
SI-TM1 — Papillary Muscle	Basket	Pointed Hook	Available for all force transducers
SI-TM2 — General Purpose	Pointed Hook	Pointed Hook	Available for all force transducers
SI-TM3 — Small Skeletal Muscle	Tweezer	Tendon Hook	Available for all force transducers
SI-TM4 — Trabeculae	Basket	Basket	Available for all force transducers
SI-TM5 — Papillary Muscle	Basket	Tweezer	Available for all force transducers
SI-TM6 — Muscle Rings	Blunt Hook	Blunt Hook	Available for all force transducers
SI-TM7 — General Purpose	Tweezer	Tweezer	Available for all force transducers
SI-TM8 — Muscle Rings	Blunt Hook/Vascular Hook	Blunt Hook/Vascular Hook	Available for all force transducers
SI-TM9 — Strong Skeletal Muscle	Tweezer	Spring Clip	SI-KG2, SI-KG2A, SI-KG2B Only
SI-TM10 — Strong Skeletal Muscle	Pointed Hook	Spring Clip	SI-KG2, SI-KG2A, SI-KG2B Only
SI-TM11 — Very Strong Skeletal Muscle	Pointed Hook	Screw Clamp	SI-KG2, SI-KG2A, SI-KG2B Only

BioTester

Biaxial Test System For Biomaterials

The BioTester offers you a complete system for examining any planar biological or replacement tissue-skin, ligaments, blood vessels, heart valves, sclera, membranes and scaffolds! Powerful image tracking and analysis software delivers synchronized data, image and video management.

Quick and easy sample mounting


High Resolution (integrated) CCD camera provides synchronized video tracking for live images (up to 15 frames/sec) and real-time analysis


- Image tracking and analysis software offers real-time data graphing and imaging to confirm the quality of the data collected
- Precision measurement of small samples (3mm-15mm square)
- Uniaxial or biaxial tension tests of planar tissues
- Multi-modal cyclic, simple and relaxation testing

Export data to spreadsheets or other scientific modeling software

Synchronized data, image and video management

The unique tungsten BioRakes easily pierce the toughest and most delicate soft tissue samples and provide distributed attachment sites across the geometry of the sample for uniform attachment and deformation across the edge of the sample. The sharp rakes will not damage fragile samples.





The software features live data graphing and imaging and offers a virtually limitless number of test stages and duration combinations. Previously used tests and templates can be easily edited.

CS-BIOTESTER

Biaxial Biomaterial Tester

For information on these tissue testers, see www.wpiinc.com/stretch



Uniaxial Tissue Tester

- Quality testing
- Cost-effective
- Multiple attachment options
- Available imaging and temperature controlled fluid bath

The **CS-USTRETCH** makes it possible to carry out dependable, mechanical testing of soft materials for a reasonable price.

Applications

Typical applications include hydrogels, skin, muscle, blood vessels, heart valves, ligaments, sclera and scaffolds.

Specifications

Test specimen as small as 3x3mm and as large as 20x100mm. Elongation rates can be as high as

50mm/s and applied forces can be monitored at 100Hz. Load control is also possible, and simple or complex protocols can be specified. A variety of load cells are available, with force resolutions as low as 1mN and capacities as high as 100N.

Includes imaging software similar to the BioTester.



Mount a specimen quickly and precisely using the patented BioRake attachment system.

CS-USTRETCH

Uniaxial Biomaterial Tester

Check the WPI website for complete details and specifications.



Micro-Scale Material Characterization


The MicroSquisher can be used to determine the compressive stress-strain properties of a variety of materials (hydrogel microspheres, small tissue samples, scaffolds and cell aggregates) with peak forces ranging from 1µN to 1mN.

The device can perform displacementcontrolled compression and stress relaxation testing.

Example Application

The interface tensions that exist play an important role in the organization of cells within aggregates. These properties can be determined by analyzing the force-time curve and test images from a parallel plate compression test.

The MicroSquisher image analysis module quantifies the aggregate profile, allowing cellcell and cell-medium interface tensions to be calculated



glass plate

cell aggregate (250µm diameter)

glass plate






CS-MICROSQUISHER

Micro-Scale Compression System

MechanoCulture

Mechanotransduction cell culture system



MechanoCulture allows researchers to culture cells on a deformable substrate or as part of a 3D matrix to understand how mechanical environment affects cell differentiation, mitosis and signalling.

Easy to use software

on specimens between 50 µm and 2

mm with force resolutions as small as

Displacement control is achieved by

a motorized piezo stage.

50 nN. Forces are determined from the

deflection of a flexible cantilever beam to

which one compression plate is attached.

manipulating the base of that beam using

The specimen can be tested in ambient air or in a temperature-controlled fluid

bath. An integrated camera system allows synchronized imaging at up to 5Hz.

The **MechanoCulture** software is used to specify test parameters. The test sequence can then be downloaded to the instrument. A run/pause

button is used to initiate, pause and stop the test. An LED display indicates the state of the unit, including the number of cycles remaining in the original protocol. The base unit can be stopped and disconnected from its power source without losing track of its position in the protocol.

- Uniaxial, equibiaxial, or nonequibiaxial strains
- 0 15% strain
- User downloaded test protocol
- Monolayer cells on substrate or 3D constructs
- Autoclavable
- Affordable and expandable system
- VHS videocassette-sized units fit easily in standard incubator

Other Systems Available

The MechanoCulture FX can uniaxially stretch 24 wells while capturing images on an inverted microscope. The sterile single-use silicone plates have a thin transparent bottom that has similar optical properties to a glass coverslip.

The MechanoCulture T6 can uniaxially stretch up to 6 clamp-mounted specimens from 5-80mm in length. For stimulation in tendon, ligament, and bone tissue engineering work, it can deliver up to 250N of thrust. For cardiovascular research, the system can stimulate at up to 5Hz.

CS-MECHANO-SYS Strainable Substrate for Culturing Cells **CS-MECHANO-FX** 24 Well Mechanical Stimulation for Culturing Cells CS-MECHANO-T6 6 Channel Mechanical Stimulation for Tissues

World Precision Instruments

www.wpiinc.com

Surgical Instruments

The instruments shown in this section are a sampling of the most popular tools from WPI. For a complete look at the hundreds of other items available, call today for the complete catalog.

Made in Switzerland

When we listen to our customers, we build a bridge of understanding and rapport. In that understanding, we discover the true needs of the customer. We have had numerous requests from you for affordable and good quality instruments and we are listening. World Precision Instruments is proud to introduce our own line of Swiss-made

forceps and scissors. We worked diligently with our manufacturer to make sure our products meet your expectations. With over 45 years experience, World Precision Instruments provides innovative instruments to the biomedical research community. We look forward to continuing to expand our product line to better serve your evolving

WPI Swiss Tweezers

stainless steel, 11.5cm

504502

acid-resistant, antimagnetic, 11.5cm

504503



Tips: 0.18×0.2 mm

WPI Swiss Tweezers

stainless steel, 11.5cm

504504

acid-resistant, antimagnetic, 11.5cm

504505



Tips: 0.18×0.2 mm

WPI Swiss Tweezers

stainless steel,11cm, 0.1 x 0.06mm tips

504506

acid-resistant, antimagnetic, 11cm

504507.....



WPI Swiss Tweezers

stainless steel, 11cm, 45° tips

504512

acid-resistant, antimagnetic, 11cm, 45° tips

504513



Tips: $0.06 \times 0.07 \text{ mm}$

WPI Swiss Tweezers

stainless steel, 11cm

504510

acid-resistant, antimagnetic, 11cm

504511



Tips: $0.06 \times 0.07 \text{ mm}$

WPI Swiss Tweezers

stainless steel,10.5cm, extra fine 90° tips

acid-resistant, antimagnetic, 10.5cm, extra fine 90° tips

504509.....



World Precision Instruments

Forceps by Dumont

For material information and instrument purchase guidance visit www.wpiinc.com



Material: Stainless steel Length: 12 cm (4.75 in.) 0.20 x 0..12 mm tips Tips: 503235



Dumont #5

Material: Dumostar, non-magnetic,

non-corrosive stainless steel

Length: 11 cm long (4.75 in.) 0.025 x 0.015 mm Tips:

Tip Profile: 500085

1:1

Dumont #5

Stainless steel Material: 11 cm (4.3 in.) Length: Tips: 0.025 x 0.005 mm 501985

available 1:1

Dumont #5

Material: Dumostar, non-magnetic, non-corrosive

Length: 11 cm (4.3 in.) Tips: 0.1 x 0.06 mm Tips

Tip Profile:

Dumont #5

Material: Dumoxel, non-magnetic stainless steel

Length: 11 cm (4.3 in.) 0.1 x 0.06 mm tips Tips:

Tip Profile:



Dumont #5

Material: Stainless steel 11 cm (4.3 in.) Length: 0.10 x 0.06 mm Tips:

500342 Tip Profile:



Dumont #5

Material: Stainless steel, Biology Length: 11 cm (4.3 in.) 0.05 x 0.01 mm tips Tips:

Tip Profile:



Dumont #5

Material: Stainless steel, Medical Biology

11 cm (4.3 in.) Length: 0.05 x 0.01 mm Tips: 14095.....

Tip Profile:



Dumont #5B

Material: Stainless steel, Biology, bent at 45°

11 cm (44.3 in.) Length: Tips: 0.05 x 0.01 mm

Tip Profile:



Dumont #7

Material: Stainless steel 12 cm (4.75 in.) Length: 0.17 x 0.1 mm Tips:

14097.....

Tip Profile: 1:1

Forceps

Iris Forceps, serrated Stainless steel 10 cm (4 in.) long, straight 0.8 mm tips

Standard German

15914-G ...

Iris Forceps, serrated

Stainless steel 10 cm (4 in.) long, curved 0.8 mm tips

Tip Profile: 1:1

1:1

1:1

www.wpiinc.com

Adson Forceps, serrated

Stainless steel 12 cm (4.75 in.) long Tip Profile:

Standard 14226 . . . German 14226-G



German

Tip Profile:




Kelly Hemostatic Forceps

Stainless steel

tandard	German

501241 Straight, 14 cm (5.5 in.) long..... 501241-G **501288** Curved, 14 cm (5.5 in.) long 501288-G **501714** Straight, 15.5 cm (6 in.) long.....

501714-G 501715-G



Hundreds of additional surgical instruments and related items are available. See the complete catalog.

501715 Curved, 15.5 cm (6 in.) long

Tip Profile:

Titanium Forceps

- 100% non-corrosive (great for seawater procedures)
- 100% non-magnetic (MRI compatible)
- 40% lighter than stainless steel (reduces hand fatigue)
- Forcep tips coated with tungsten carbide for increased gripping power
- Anodized, non-glare blue finish

Straight Forceps

Length: 11.5 cm (4.5 in.)

Tips: 1 x 2, 0.12 mm teeth, 5 mm tying platform

Titanium 555047FT.....



Straight Forceps

Length: 11.5 cm (4.5 in.)



Calibri Tying Forceps

Length: 11.5 cm (4.5 in.)

Tips: 3 mm, 4 mm tying platform **Titanium 555008FT**.....



Kelman-McPherson Forceps

Length: 8.5 cm (3.3 in.)

Tips: 45°angle, 7.5 mm smooth jaw **Titanium 555190FT**.....



Kelman-McPherson Forceps

Length: 11.5 cm (4.5 in.)

Tips: 45°angle, 7.5 mm smooth jaw



McPherson Forceps

Length: 11.5 cm (4.5 in.)

Tips: 45°angle, 6 mm tying platform

Titanium 555009FT.....



McPherson Forceps

Length: 8.5 cm (3.3 in.)

Tips: angled, 4.5 mm tying platforms

Titanium 555005FT.....



Forceps

Length: 8.5 cm (3.3 in.)

Tips: 12 mm, straight, 4.5 mm tying platforms



Forceps

Length: 8.5 cm (3.3 in.)

ips: straight, 0.12 mm teeth (1x2), 5 mm tying platform

Titanium 555041FT



Tip Profile:

Vannas Scissors, Super Fine, stainless steel

Length: 8 cm

Blades: straight 3 mm Tips: 0.015 x 0.015 mm

501778

Length: 8 cm Blades: curved 3 mm Tips: 0.015 x 0.015 mm

501839



Vannas Scissors, stainless steel

Length: 8.5 cm (3.3 in.)

Blades: straight 7mm 0.025 x 0.015 mm Tips:

Length: 8.5 cm (3.3 in.) curved 7mm Blades:

0.025 x 0.015 mm Tips:

500086.....

501232



Vannas Scissors

Length: 8 cm (3 in.) straight 5 mm Blades: 0.1 mm tips Tips:

Standard

German

14003 14003-G



Vannas Scissors

Length: 8 cm (3 in.) curved, 5 mm Blades: 0.1 mm tips Tips:

Standard

German

14122..... 14122-G



Vannas Scissors

Length: 8 cm (3 in.)

Blades: 45° angled to side, 5 mm

0.1 mm tips Tips:

German Standard

500260 500260-G.....



McPherson-Vannas Scissors

Length: 8 cm (3 in.) Blades: straight 5 mm 0.1mm Tips:

Length: 8 cm (3 in.) curved 5 mm Blades: Tips:

0.1mm

Standard German

14124 14124-G.....

501234 501234-G

Spring Scissors

12 cm (4.75 in.) Length:

straight 12 mm extra fine and long Blades:

Standard 14125..... German 14125-G

Length: 12 cm (4.75 in.)

Blades: curved 12 mm extra fine and long

14126.....



Tip Profile:

Iris Scissors

Length: 9 cm (3.5 in.)

sharp, 11 mm straight Tips: Titanium 555560S.....



Vannas Scissors

Length: 9.5 cm (3.7 in.) angled sharp, 10.5 mm Tips:

<u>Titanium</u> 555584S



555583S

Vannas Capsulotomy Scissors

Length: 9.5 cm (3.7 in.) Tips: curved, sharp, 12.5 mm

Length: Tips:



Castroviejo Straight Scissors

Length: 10.5 cm (4.1 in.)

Tips: 11 mm



Castroviejo Curved Scissors

10.5 cm (4.1 in.) Length:

small Tips:



Mcpherson-Westcott Stitch Scissors

Length: 11.5 cm (4.5 in.) 10.5 mm Tips:

Titanium 555540S.....



Mcpherson-Westcott Conjunctiva Scissors

10.5 cm (4.1 in.) Length: curved, blunt Tips: Titanium 555500S.....



Iris Scissors

10.5 cm (4.1 in.) Length: Tips: sharp, 11 mm curved


Titanium 555562S



21



12.5 cm (5 in.) curved......



Dissecting Scissors

	8		10
Standard		German	
14393	10 cm (4 in.) straight	14393-G	Tip Profile:
14394	10 cm (4 in.) curved	14394-G	
15922	12.5 cm (5 in.) straight	15922-G	

15923-G

Iris Scissors, stainless steel

Length: 11.5 cm (4.5 in.)


Blades: straight, Tungsten Carbide Standard German

500216 500216-G

Length: 11.5 cm (4.5 in.)

curved, Tungsten Carbide Blades:

500217 500217-G



Iris Scissors, stainless steel

Length: 11 cm (4.3 in.) Blades: straight Standard German


501758 501758-G Length: 11 cm (4.3 in.) Blades: curved

501759 501759-G

Mini-Iris Scissors

8 cm (3.1 in.) straight, sharp tips. . 503670

503671 8 cm (3.1 in.) curved, sharp tips



Iris Scissors, SuperCut, stainless steel

one edge micro serrated, one edge honed to the sharpness of a knife edge 10 cm (4 in.)

Standard German Straight 14218.... 14218-G Curved 14219..... 14219-G

11.5 cm (4.5 in.)

Straight 503259

Curved 503260




504519 WPI Swiss Scissor, Stainless Steel, 9cm, Blade Extra Fine WPI Swiss Scissor, Stainless Steel, 9cm, Blade Fine Sharp...... 504520

504521 Wpi Swiss Scissor, Stainless Steel, 9cm, Heavy Duty, Rounded.....

504522



504519



	orial profial p	301/34	301/34-G		
	Blunt/Blunt	501743	501743-G		
1.5 cm (4.5 in.) Curved	Sharp/Sharp	501755	501755-G	Tip Pro	files:
	Sharp/Blunt	501756	501756-G		
	Blunt/Blunt	501757	501757-G		
)
4 cm (5.5 in.) Curved	Sharp/Sharp	501220	501220-G		
	Sharp/Blunt	501221	501221-G		
	Blunt/Blunt	501222	501222-G		\oslash
14 cm (5.5 in.) Straight	Sharp/Blunt	14192	14192-G		
	Sharp/Sharp	501218	501218-G		Ø
	Blunt/Blunt	501219	501219-G		

Surgical Scissors, stainless steel

One tip sharp, one tip blunt

16 cm (6.3 in.) Straight 501223 501223-G Curved 501226 501226-G	
Curved 501226 501226-G	
	1:1
	-

501223

Needle Holders

Needle Holder 12.5 cm (5 in.) long straight serrated jaws extra delicate 14109 Standard Tip Profile: German 14109-G 1:1 1:1 Mayo-Hegar Needle Holder, Stainless Steel 15 cm long (6 in.) V503382.....









Physiology Kit I

SuperCut Tenotomy Scissors curved (#14396)

Rongeur 3mm jaw (#14091)

Utility Scissors (#501322)

Operating Scissors straight Sharp/Blunt (#14192)

Dumont Tweezer #5 (#14098) Flat Jaw Tweezers (#501303)

Probe 1.0 mm diameter, blunt (#501313)

KIT-PHYSIO-I

Mouse Kit

Dumont #5 (#14098)

Vannas Scissors (#14003)

Iris Forceps, curved, serrated (#15915)

Dissecting Scissors, straight 10cm (#14393)

Wire Retractor (#14130) Needle Holder (#14109)

Blunt Probe, *1.0 mm diameter (#501313)*

MOUSEKIT

Mouse Kit.....

Physiology Kit II

Vannas Scissors (#14003)

SuperCut Iris Scissors *straight* (#14218)

Rongeur 1.3mm jaw (#14292) Utility Scissors (#501322)

Probe 1.0 mm diameter, blunt (#501313)

Operating Scissors straight Sharp/Blunt (#14192)

Dumont Tweezer #5 (#14098)) Stevenson Retractor (#14131)

Iris Forceps curved (#15915)

Adson Forceps 1x2 teeth (#500092)

Olsen-Hegar Needle holder (#500227)

PHYSIO-II

Physiology Kit II.....

Rat Kit

Dumont #5 (#14098)

Vannas Scissors (#14124)

Iris Forceps, serrated (#15915)

SuperCut Iris Scissors, straight 10cm (#14218)

Alm Retractor (#14240)

Needle Holder (#14110)

Blunt Probe, 1.0 mm diameter (#501313)

Student **Dissecting Kit**

Includes:

Dissecting Scissors, 4.5" Dissecting Scissors, 4" Dressing Forceps, 5.5" Micro Dressing Forceps, 4" #4 Knife Handle Straight Teasing Needle Angled Teasing Needle #22 Surgical Blade Canvas Roll-up





501336

501838

Canvas Roll-up only

Clips and Clamps



for 7 mm clips only 500343



Reflex Clip 9 mm

for use with #500345 100/box, Stainless Steel, non-sterile 500346

Reflex Clip 7 mm

for use with #500343 100/box, Stainless Steel, non-sterile

500344

Micro Bulldog Clamps



Curved serrated jaws

14119

32mm long, 12mm jaw

German

38mm long, 9mm jaw

14119-G

Retractors

Agricola Retractor

Length: 4 cm (1.6 in.)
Blades: 3 x 3 sharp prongs

Small, self-retaining retractor perfect for small animal surgery and dissection. Maximum spread, 2.5 cm.

501846.....



Reflex Clip Applier

11 cm (4.3 in.) long

for 9 mm clips only

500345

Scalpels, Knives, Blades and Handles

Scalpel Handle #3

13 cm (5 in.) long, stainless steel

500236.....

German 500236-G.....



For Blades #10 ~ #15

Sterile stainless steel blades All stainless steel scalpel blades are made by Feather®, using a precise beveling technique to create the edge's micron sharpness. They are the finest blades available.

500239	Scalpel Blades #10, stainless steel, sterile (100 / box)
500240	Scalpel Blades #11, stainless steel, sterile (100 / box)

Animal Lancets



504539	Animal Lancet, 3mm, Mice,1-2 drops, 200/box
504540	Animal Lancet, 4mm, Mice under 2 mo, 200/box
504550	Animal Lancet, 5mm, Mice 2-6 mo,200/box
504551	Animal Lancet, 5.5mm, Mice over 6 mo, 200/box
504552	Animal Lancet, 5mm, Rats under 3 mo, 200/box
504553	Animal Lancet, 6mm, Rats under 3mo, 200/box
504554	Animal Lancet, 7mm, Rats 3-4mo, 200/box
504555	Animal Lancet, 3mm, Primates, 200/box
504556	Tatoo Identification Lancet, 3mm, 200/box

Disposable Scalpels





DCAP Guillotine for Rodents and other small animals (opening 1.5 x 1.5 in.) DCAP-M Guillotine for large rodents and other medium animals (opening 2.5 x 2.5 in.) DCAP-L Guillotine for larger animals (opening 4" x 4")

- Large, stable base
- Hardened blades for long service
 - Ambidextrous configuration

The small animal guillotine has been completely redesigned for ease of use and extra added safety features. The blades are drawn together by magnetic force to ensure a clean and precise cut through very strong bones and skin.

There is a large base for stability, long handle for extra leverage, spring action so the blades can not fall down unexpectedly, hardened stainless blades for endurance, simplified construction

for easy maintenance. The fluoropolymer coated surface on the base makes cleaning easy.

The guillotine is considered one of the most humane methods to sacrifice a subject.

OmniDril[|]

This line-powered micro drill will make easy work of grinding, finishing, cutting, and drilling bone, teeth, and other material. The high-torque 35,000 rpm (maximum) motor is quiet and has minimal vibration which reduces wear on the motor and provides greater comfort for the user. It also features a forward and reverse switch, "E Type" handpiece, and handpiece holder. The handpiece has a removable nose cone that can be cleaned and sterilized. It accepts 3/32" and 2.33 mm bur shanks. Unlike battery-powered drills, this unit will maintain consistent power for the duration of use. The wide range of speeds allows the user to control the amount of heat generation.

The following accessories are included with the Micro Drill System:

Qty Description

Abrading Tip, Rubber Abrading Tip, Stone

Accessory Stand

Ball Mill, Carbide, #1, .031" Diameter Ball Mill, Carbide, #2, .039" Diameter

Ball Mill, Carbide, #3, .047" Diameter

Ball Mill, Carbide, #4, .055" Diameter


Ball Mill, Carbide, #5, .063" Diameter

Ball Mill, Carbide, #6, .071" Diameter Ball Mill, Carbide, #7, .083" Diameter

Ball Mill, Carbide, #1/4, .019" Diameter Ball Mill, Carbide, #1/2, .027" Diameter

Cutoff Disk

Mandrel, Screw



OMNIDRILL35 SPECIFICATIONS

110V. 50/60 Hz **INPUT** OUTPUT 0-32 Vdc **FUSE** 1 amp **OPERATING SPEED RANGE** 0-35,000RPM

DIMENSIONS 178 x 114 x 89mm (7 x 4.5 x 3.5in.)

WFIGHT 1.7kg (3.75lbs)

	normas	1, 11
	SE STANSON STA	0
	((6)	
4	math	
4		

503598	OmniDrill35 Micro Drill System, 110 V
503599	OmniDrill35 Micro Drill System, 220 V
REPLACEN	MENT ACCESSORIES
501850	Abrading Tip, Rubber, pk of 20
501851	Abrading Tip, Stone, pk of 5
501852	Accessory Stand
501853	Ball Mill, Carbide, #1, .031" Diameter, pk of 5
501854	Ball Mill, Carbide, #2, .039" Diameter, pk of 5
501855	Ball Mill, Carbide, #3, .047" Diameter, pk of 5
501856	Ball Mill, Carbide, #4, .055" Diameter, pk of 5
501857	Ball Mill, Carbide, #5, .063" Diameter, pk of 5
501858	Ball Mill, Carbide, #6, .071" Diameter, pk of 5
501842	Ball Mill, Carbide, #7, .083" Diameter, pk of 5
501860	Ball Mill, Carbide, #1/4, .019" Diameter, pk of 5
501861	Ball Mill, Carbide, #1/2, .027" Diameter, pk of 5
501862	Cutoff Disk, pk of 20
501863	Mandrel, Screw, pk of 5
501864	Mandrel, Threaded, pk of 5
504459	Footswitch
502237	Stereotaxic Holder for OmniDrill35 Microdrill

BioClave Mini Research Autoclave

- Fully automatic-just press START
- Rapid, complete sterilization
- Extremely compact
- Built in water tank does not require connection to an external water supply
- Benchtop model

The BioClave Mini comes pre-programmed for basic sterilization. Choose either the 120°C or 134°C temperature mode, and the cycle time is fixed. This is the perfect autoclave for small spaces.

The stainless steel sterilization chamber accommodates a variety of liquids, media,



instruments, glassware and plasticware. A mechanical and electrical safety interlock prevents the door from being opened until the pressure is released (0 PSI).

Simply press the Start button, and the entire sequence runs-fill, sterilize, exhaust and dry.

504187	BioClave	Mini	(115)	V)
504188	BioClave	Mini	(230	V)


BIOCLAVE MINI SPECIFICATIONS

CHAMBER VOLUME...... 8L MAX. PRESSURE 29 PSI/2 Bar

CHAMBER DIMENSIONS 6.7 x 12.5"/17 x 31cm TRAY DIMENSIONS. 4.75 x 9.75"/12 x 24.8cm EXTERIOR DIMENSIONS 20 x 13.5 x 12.9"/51 x 24 x 33cm

WEIGHT 74 Lb. (33kg)

Aluminum Sterilizing Trays



These sterilization containers are made of hard anodized aluminum and can withstand autoclaving and ethylene oxide sterilizing. The nonperforated base has two locking clasps that secure the lid in place. The perforated lid with disposable filter allows steam and gas to move freely in and out of the container during sterilization. A silicone mat may be purchased separately and cut to fit the base of the tray.

501913	Small Aluminum Tray, 300 x 140 x 40 mm
501914	Large Aluminum Tray, 300 x 140 x 70 mm
501915	Single-use Filters, 100 per package
500255	Silicone Mat, 38 x 25.5 cm

Dry Sterilizer

Sterilize your microdissecting and tissue culture instruments, thoroughly and conveniently, in seconds. No chemicals. No flames. No risk of burns. No disinfectant fluids. Glass beads heated to 260°C kills all viruses, aerobic and anaerobic bacteria, yeasts and spores. (1.5 mm lead-free glass beads included.)



GERMINATOR SPECIFICATIONS

I.D..... 5.1 x 10.2 cm (2 x 4 in.) WEIGHT 2.3 kg (5 lb.)

Dry Sterilizer (110 V)
Dry Sterilizer (220 V)
Extra Glass Beads, 300 g (11 oz)

SweepZone® Ultrasonic Cleaning System

SweepZone operates with a special cleaning wave of ±2KHz energy that literally sweeps back and forth throughout the ultrasonic tank. Pulse Width Modulation (PWM) enables the device to sense cleaning loads, solution levels and voltage fluctuations.

It automatically adjusts to different conditions in order to create more uniform cleaning



 $C \in$

power by eliminating hot spots and weak areas in the tank. This results in unsurpassed cleaning performance with every application.

User-friendly digital touch pad controls make it easy and convenient to set any of the five preset timer settings from a range of six to sixty minutes. Durable stainless steel housing and tank are corrosion resistant and long lasting. Includes timer, drain and stainless steel cover Dimensions: 32.2 x 1.75 x 22.2 cm (12.68 x 6.88 x 8.75 in.)

Shipping Weight: 5.4 kg (12 lb) Tank Capacity: 4.0 L (1.06 gal.)

Tank Internal Dimensions: 29.8 x 15.2 x 10.2 cm (11.75 x 6 x 4 in.)

504217 SweepZone Ultrasonic Cleaning System

Quantrex® **Ultrasonic Cleaning Systems**

Provides super-strength cleaning every time. Versatile enough for a variety of cleaning applications. When used with L&R's specialty formulated solutions, the self-contained, compact unit offers efficient troublefree cleaning. Each Quantrex machine comes standard with increased power—strength you can see as soon as you turn the unit on.



60 Minute Timer.

Constructed with vinyl-clad steel and stainless steel.

Fourteen quality inspection steps for strength and durability.

Stainless steel drain with multi-positional outlet for easy removal of solution.

Dimensions: 26.0 x 16.5 x 21.0 cm (10.25 x 6.5 x 8.25 in.)

Shipping Weight: 4.5 kg (10 lb)

Tank Capacity: 3.2 L (0.85 gal.)

Internal Dimensions: 23.8 x 13.7 x 10.2 cm (9.38 x 5.38 x 4.0 in.)

Includes Timer, Drain and Cover. Heater Optional.

504216 Quantrex Ultrasonic Cleaning System

Ultrasonic Cleaner

Half liter stainless steel tank. Durable and compact, robust all-metal construction allows for continuous duty. Includes lid.







14342

UBATH SPECIFICATIONS

INPUT POWER PEAK OUTPUT 70W, 55 kHz TANK CAPACITY 0.53L (18 oz.) TANK I.D. 12.1 x 8.6 x 6.6 cm (4³/4 x 3³/8 x 2⁵/8 in.)

13.7 x 10.5 x 12.1 cm TANK O D

(53/8 x 41/8 x 43/4 in.) 1.8 kg (4 lb.)

503737

UBATH-Y Ultrasonic Cleaner (110 V) **UBATH-Z** Ultrasonic Cleaner (220 V, Euro plug) **UBATH-B** Ultrasonic Cleaner (220 V, UK plug) Specify line voltage

	specijy iine voitage				
13740 Ultrasonic Detergent (4 lb)					
14342	Ultrasonic Detergent (1 gal)				
503737	Mesh Draining Basket				
503738	Bur Sterilizing Tray				
503739	Positioning Cover				
504766	Sterilization mesh casette 40x40x20mm				
504767	Sterilization mesh casette.80x80x34mm				
504768	Sterilization mesh casette.105x70x25mm				

SHIPPING WEIGHT



The PI-MC2010 is a micro-cauterizing unit that embraces thermodynamic technology to perform micro-surgical tasks with efficiency and cleanliness. Cells just beyond the cutting line are left healthy and functional.

How It Works

A current passes through the platinum cutting wire and super-heats it at a preprogrammed rate for the span of a few microseconds, then abruptly stops. A thin, gaseous-tube of solution surrounding the wire is created because the rate of heating is much faster than the velocity of the solution molecules due to their average thermal energy at room temperature. This creates an intense but finite shockwave as the superheated tube of vapor expands from the wire inducing turbulence. The quick heating and cooling of the platinum wire causes a flushing effect as cool water rushes in after the brief shockwave and turbulence. Debris from the lysed cells directly adjacent to the wire is flushed away immediately. Cellular debris does not stick to the wire and, even at the highest power settings, only an extremely small quantity of heat is generated and then quickly cooled so as not to increase the temperature of the surrounding environment.

Features

- 13 & 18 micron diameter cautery tips with 1 & 1.5mm loop
- 15 Power settings
- Light weight pencil handle
- Foot switch activated
- Simple operation

PI-MC2010	Micro-Cautery Instrument
PI-IVICZU IU	Micro-Cautery instrume

Comes with foot switch, pencil, tip, holder, power cord, operation manual, and choice of one box of tips.

MICRO-CAUTERY ELECTRODE TIPS

504045	18 micron wire tip cautery electrode, red-black cover,10mm loop, 5/pkg
504046	18 micron wire tip cautery electrode, red cover,1mm loop, 5/pkg
504047	18 micron wire tip cautery electrode, red cover,1.5mm loop, 5/pkg
504048	13 micron wire tip cautery electrode, yellow,1mm loop, 5/pkg
504049	13 micron wire tip cautery electrode, yellow,1.5mm loop, 5/pkg

Economy Electrosurgical Unit

Electrosurgery utilizes alternating current at radio frequencies to cut and coagulate. Using this method, the current enters the patient's body and the patient becomes part of the electrical circuit. This requires the use of a return, or indifference plate.

This economically priced electrosurgical unit has 10 levels of output intensity, three operational modes (cut, coagulate, and cut/coagulate), and various choices of electrodes.

The unit comes complete and ready-to-use with a handpiece, indifference plate, footswitch, and one of each electrodes. All accessories can also be ordered separately.

ELECTROSURGICAL UNIT SPECIFICATIONS

OPERATION FREQUENCY 1.5 MHz STABLE & FINE POWER SETTING 10 Steps

POWER SUPPLY

115V ± 10% - 50/60Hz 1.8A, 210VA

230V ± 10% - 50/60Hz 0.9A, 210VA

OUTPUT POWER 70 WATTS ± 5% WORKING FROUENCY 1.5-1.7 MHZ ± 5%

24cm x 22cm x 8.5cm (lxwxh) DIMENSIONS

SHIPPING WEIGHT 10 lb (4.5kg)



501274	Electrosurgical Unit, 110V					
501285	Electrosurgical Unit, 220V					
501273	Handpiece for electrodes, Ø 1.6mm shaft					
501275	Indifference Plate					
501277	Footswitch					
501278	Diamond Shape Loop Electrode, Ø 1.6mm shaft					
501279	Small Loop electrode, Ø 1.6mm shaft					
501280	Large Loop Electrode, Ø 1.6mm shaft					
501281	Fine Wire electrode, Ø 1.6mm shaft					
501282	Wire Electrode, Ø 1.6mm shaft					
501283	Small Oval Loop Electrode, Ø 1.6mm shaft					
501284	Ball Electrode, Ø 1.6mm shaft					
501286	Set of 7 Electrodes					



Features

- Integrated turnkey solution for small animal anesthesia
- Safe for surgical personnel, 90% below OSHA Isoflurane limit
- Designed by veterinarians
- Compact and portable
- Time efficient and cost effective
- Virtually stress-free for the animals
- Easy to setup and use, simplifying the training of new staff and reducing the threat of human error
- Speedy recovery time

EZ Anesthesia is the system of choice for anesthetizing small animals, and it comes with a variety of choices. Animals to be anesthetized are placed in the acrylic induction chamber and the system delivers a precisely blended mixture of oxygen and isoflurane. An activated charcoal air filter canister at the top of the chamber releases safe, filtered air back into the room. A water-heated cage warmer or warming plate (ATC2000) is used to retain the animal body temperature while in the induction chamber. After the initial anesthetizing, the animal can be moved to the heated surgical water bed and positioned properly in the snugly fitted nosecone. A highly sensitive valve regulated by the animal's breathing works with the nosecone to ensure non-rebreathing efficiency, and allow safe anesthesia for up to several hours.

The breathing device also includes an air filter that releases safe, filtered air back into the room.

EZ Systems include all necessary hardware components and connections. Oxygen and liquid Isoflurane are not supplied. These are required to operate the system.

Each system comes with all necessary components and connections for immediate use, including::

- Oxygen regulator
- Vaporizer Unit
- One Induction Chamber (standard EZ-178)
- One Breathing Device (standard **EZ-103**)
- Case of charcoal filters

Other available options include:

- Heated Beds and Heating Pump
- Additional Breathing Devices

EZ-FF9000 Fixed Flow System

The Fixed Flow System is our most advanced system and provides preset, fixed flow rates with no need for adjustment. The system offers five gas outlets, each with an individual ON/OFF switch. The system ensures consistent, precise gas flow when connected to any pressurized gas source higher than 7 psi.

- Induction Chamber output is fixed at 1 Lpm
- Four Breathing Circuit outputs are fixed at 0.5 Lpm

This unit features an oxygen flush button that purges the induction chamber with pure oxygen, thereby protecting personnel from exposure to anesthesia gas when opening the chamber. A handle on the vaporizer makes the unit easy to safely transport.

EZ-FF9000 Fixed Flow System

EZ-7000 Classic System

The Classic System is the upgraded version of our original EZ-Anesthesia system. The new version has a reconfigured flow meter that is more user-friendly. The system offers five gas outlets that can supply a single induction chamber and four breather circuits simultaneously.

The upgraded system also features an oxygen flush button that purges the induction chamber with pure oxygen, thereby protecting personnel from exposure to anesthesia gas when opening the chamber. A handle on the vaporizer makes the unit easy to safely transport.

EZ-7000 Classic System

EZ-B800 Basic System

The Basic System is designed for use with a single animal. It utilizes one output directly from the vaporizer into a Y-splitter which creates a dual feed that can be directed to the induction chamber or the breather circuit.

This unit incorporates an oxygen flush system that purges the induction chamber with pure oxygen, thereby eliminating personnel exposure to anesthesia gas when opening the chamber.

EZ-B800 Basic System

> For the full line of anesthesia products, see www.wpiinc.com/anesthesia

EZ-1130 Connection Kit

This kit is used for connecting additional components to the EZ-Anesthesia Systems. It includes 6 ft of clear PVC tubing with a quickdisconnect fitting.

EZ-1130 Connection Kit

EZ-830 Ventilator Connection Kit

This kit is used to connect the SAR-830 ventilators with EZ-Anesthesia Systems. It includes all required components, pre-assembled for simple connection between the ventilator and the anesthesia system.

EZ-830 Ventilator Connection Kit



heavy duty stainless steel lids have a stem fitting for connection to the guick-disconnect fitting on the hose from the gas source. A foam lid gasket ensures a good seal on the cage. Multiple lids may

EZ-177 and EZ-1785 Induction Chambers

These chambers incorporate a positive seal O-ring gasket for containment during use.

EZ-177	Sure-Seal Mouse Chamber (5"L × 4.75"W × 4.38"H)
EZ-178	Sure Seal Mouse/Rat Chamber (9.75"L × 4.75"W × 4.38"H)
EZ-1785	Large Mouse/Rat Induction Chamber (7"W \times 11"D \times 6"H inside)
EZ-179	Rabbit/Guinea Pig Induction Chamber (18.75"× 12.75" × 12.75")

EZ-320 and EZ-330 Oxygen Regulators

Regulators are preset to 50 psi.

- EZ-320 utilizes a CGA-540 connection for large "H" tanks.
- EZ-330 utilizes a CGA-870 connection for small "E" tanks.

EZ-320	Oxygen	Regulator	for	large tanks
EZ-330	Oxygen	Regulator	for	small tanks

EZ-340-XXXX Custom Hose Assemblies

Custom built hose assemblies for wall or ceiling outputs specific to facility needs: Chemetron, Ohmeda, Schrader or DISS.

EZ-20027	Small Lid (13" x 9") Fits old-style mouse cages
EZ-20028	Small Lid (16" x 10") Fits new-style mouse cages
EZ-20030	Square Lid (13" x 13") Fits Thorn cages
EZ-20032	Medium Lid (20.5" x 11") Fits rat cages
EZ-20034	Large Lid (23" x 16.5") Fits guinea pig cages
EZ-20029	Lid Storage Bracket (wall-mounted,holds up to four lids)
EZ-20027G	Small Lid Gasket (13"x9")
EZ-20028G	Small Lid Gasket (16"x9")
EZ-20030G	Square Lid Gasket (13"x13")
EZ-20032G	Medium Lid Gasket (20.5"x11)
EZ-20034G	Large Lid Gasket (23"x16.5")

Vaporizer Pole Mount

be used to treat several cages at once.

The **EZ-E28000** is a mounting option for the EZ-Anesthesia system, combining system portability with a small footprint. Constructed of Aluminum and Stainless Steel, it features a five leg base for maximum stability and an oxygen "E" tank mounting bracket. Anesthesia system must be purchased separately

E-28000 Vaporizer Pole Mount



EZ-104A	Versaflex Non-Rebreathing Unit
EZ-103A	Microflex Non-Rebreathing Unit
EZ-107A	Rat Stereotaxic Non-Rebreathing Unit
EZ-109	Multi-Animal Non-Rebreathing Unit
EZ-211	Mouse/Rat Thin-Line Heated Waterbed
EZ-212	Mouse/Rat Standard Heated Waterbed

Rodent Isoflurane Facemask Kit

Mouse and rat facemasks, medical grade PVC tubing (12 feet), 1/4-inch adapters (2), 3/8-inch adapters (2).



OC-3FW-KII	Siliali Rodelit i aceillask Rit (olie iliask, tubilig, adapters)
OC-MFM-KIT	Medium Rodent Facemask Kit (one mask, tubing, adapters)
OC-LFM-KIT	Large Rodent Facemask Kit (one mask, tubing, adapters)
OC-XLFM-KIT	XLarge Rodent Facemask Kit (one mask, tubing, adapters)
OC-MOUSE-KIT	Mouse Facemask Kit (sm & med masks, tubing, adapters)
OC-RAT-KIT	Rat Facemask Kit (lg & xlg masks, tubing, adapters)
OC-ALLFM-KIT	Complete Rodent Facemask Kit
	(sm, med, lg & xlg masks, tubing, adapters)

Mobile Workstations

Two mobile workstations, constructed of heavy-duty stainless steel with locking casters, integrate all your EZ-Anesthesia components into one portable unit. Open side shelves accommodate 20 lb. cylinders, and convenient 2" port holes allow for easy rigging of gas and electrical lines. Below the work surface of each mobile workstation is an open shelf and a locking cabinet. t provides a 42"x24" work surface and holds up to four cylinders. EZ-E27000 has a 22"x21" work surface and holds up to two cylinders. These systems are easy to set up and provide maximum flexibility and mobility.

EZ-E25000	Mobile Workstation, 42" x 24" Top
EZ-E27000	Mobile Workstation, 22" x 21" Top

Analgesia

Electronic von Frey Anesthesiometer

- Plug up to three probes into a single unit
- LCD readout (Floating or last maximum/minimum)
- Rigid tips up to 800 gm
- "Supertips™" 15 up to 65 gm
- 1,000 gm probe available
- Independence from temperature
- Optional analog output cable for chart recorder
- Pipette tips can be customized to any specification
- Microprocessor electronics 0.1 gm plug-in probes

11-2390	Electronic von Frey Anesthesiometer, rigid tips, 90 gram range
II-2391	Electronic von Frey Anesthesiometer, rigid tips, 800 gram range
II-2392	Electronic von Frey Anesthesiometer, rigid & 15 super tips, 90 gram range
II-2393	Electronic von Frey Anesthesiometer, rigid & 15 super tips, 800 gram range
II-23931	Electronic von Frey Anesthesiometer, custom rigid tips, 1000 gram range
II-2394	von Frey Probe, 90 gram range
II-2395	von Frey Probe, 800 gram range
II-2396	von Frey Probe, 1000 gram range
II-2397	MRI Probe Option (add to price of probe above)
II-2400	Analog Output Cable

To assess mechanical allodynia, which is a painful response to a light touch or pressure from a stimulus that is not normally painful, the Electronic von Frey Anesthesiometer was developed. The Electronic von Frey meter uses one of 15 different flexible von Frey hairs called "SuperTips™" (or rigid tips up to 800 grams). Each hair, regardless of model chosen, is exactly 0.8 mm in diameter. This uniformity of design eliminates false readings and allows for comparison of test results. The Electronic von Frey can be used with chart recorders and analog/digital converters, and it never needs calibrated. This system includes either a 90, 800 or 1,000 gram probe. Mesh stands are available in a variety of sizes for large group studies.



Trio

Get three test systems in one package with the Trio, featuring the Electronic von Frey, Plethysmometer and Randall Selitto Meters. Just like the Quattro package, the modular design allows these three test systems to communicate with the same electronic controller. The stand and sling for the Randall Selitto test are sold separately.

II-2888 Trio 3-in-1 System

Quattro

This special package offers four tests, including Electronic von Frey, Plethysmometer, Randall Selitto and the Grip Strength Meter. You get all four test modules and the electronic controller that is interchangeable with all four systems. The electronic controller has up to three inputs, so you can perform up to four unique tests with only one electronic system. If you prefer, you may build your system as you grow. Because of the modular design of these four systems, you need to order only one complete system. Then, the modules for the other three tests, which integrate into the system, can be purchased separately, as needed. The stand and sling for the Randall Selitto test

II-2889 Quattro 4-in-1 System



Plethysmometer (Paw Volume) Meter



- No wetting solution needed
- One calibration/year
- Battery-operated or line powered controller
- One-year warranty

Test the effectiveness of anti-inflammatory agents and agents to reduce edemic conditions with the Plethysmometer. Simply insert the rat or mouse paw into water. A sensor in the water notes a pressure change when the paw is immersed. Pressure is calibrated in 0.1 milliliters and displays on the battery-powered monitor. The acrylic stand offers full visibility of the subject throughout the testing.

II-520MR Paw Volume Meter for Mouse & Rat

Randall Selitto Paw Pressure Meter

- Hands free operation with footswitch
- Visible force limit indicator
- Portable electronic display
- Battery-operated or line powered controller
- 3-oz. handheld probe with accuracy of 0.5%
- No calibration required

The Randall Selitto Paw Pressure Meter for analgesia testing is digitally controlled. Use the handheld instrument to apply force to an animal's extremity and get instantaneous, live readings. You can even view the last maximum force applied during the test. The new limit

indicator lets you select the maximum force limit, and then indicates with a warning light when the system reaches that limit. This unit comes with an acrylic stand to allow for easy viewing of the display. Stand and sling are sold separately.

Randall Selitto Paw Pressure Meter, 800 g pressure applicator II-2500

Grip Strength Meter

- Mouse and rat wire mesh grid plates included
- Maximum force range is 2,000 gm (10% over range allowed). Higher ranges are available by special order.
- 1 gram increment readings
- Suction feet on the heavy, anodized base plate resist even large pulling forces
- Battery-operated or line powered controller

One year warranty
Measure muscle hyperalgesia in rats and mice


with the Grip Strength Meter, which gauges the forelimb grip force using a digital force transducer. Simply hold the animal by the tail and gently dangle it over the wire mesh plate until the animal grasps the plate with its forepaws. The force transducer, connected with the wire mesh plate, measures the strength of the animal at the time of the test. The battery-operated, electronic control device calculates the average of three measurements, and it holds the last maximum force in a "peak and hold" type readout until you reset it.

Grip Strength Meter for Mouse & Rat **II-2200**









Analgesia

Incremental Hot Cold Plate Analgesia Meter

Heat or cool, from 0-70 °C

Ramping temperatures for threshold & latency results

• Rapid increase or decrease in temperature

Precise programmable digital control

Printout of data

Temperature stability is 0.1°C

Includes clear animal enclosure

Plate size 4 x 8"

Two-year warranty

This safe, humane device for rats and mice is used for latency and threshold-based nocicetption, ramping temperatures for 0-70°C. Because this hot cold plate is incremental, it measures latencies of much more than just the strong narcotic agents, broadening dramatically the range of analgesia research with devices of this type.

Microprocessor-controlled, the Incremental Hot Cold Plate can heat or cool in increments of 0.1°C, at a rate of 1-10°C per minute. With uniform heating and cooling and upper/ lower cut-off limits, this device is predictable

TTE-WE

and safe. It can also function as a constant temperature plate with great stability (0.1 °C). As soon as a reaction is observed from the chosen paw, the unit reverses to the standby temperature.

II-PE34

Incremental Hot Cold Plate Analgesia Meter for Mouse & Rat

Incapacitance Meter for Mouse & Rat

Precise programmable digital control

Start, stop and reset test from controller's front panel

180-270 gram holder included (other sizes available)

Reaction detected automatically

Manual override of all timer functions

Alphanumeric readout

5 to 999 seconds test period

Alphanumeric readout

All functions and parameters entered via key-pad

Two-year warranty

II-600MR Incapacitance Meter for Mouse & Rat



Test pain and inflammation in the hind limbs of mice, rats or birds with the Incapacitance Meter. It uses a technique called dual channel weight averaging, which tests both hind limbs. This gives you a clean, stressfree correlation of the paw pressure test. Conduct control and testing of the animal at the same time. Place the animal in the holder with its hind limbs resting on the two weight-averaging platform pads. The controller records the average weight (grams) over the test period as the animal shifts its weight from each pad

Hot Plate Analgesia Meter

- Includes both mouse and rat enclosure
- Temperature indicated in 0.1°C increments
- Holding accuracy is +/- 0.1 °C
- Digitally controlled
- Two-year warranty

Use the Hot Plate Analgesia meters for latency testing in rats and mice. Simply place the animals on a black anodized, aluminum plate (11 \times 10.5 \times 0.75", 275 \times 263 \times 15 mm) and set the plate's surface temperature to the desired setpoint (up to 75°C). The plate maintains a consistent temperature throughout the test.

II-39

Hot Plate Analgesia Meter for Mouse & Rat



Plantar Test Apparatus/Tail Flick Test Analgesia Meter



- Includes three acrylic animal enclosures that each hold two rats or four mice
- Precise programmable digital control
- User-defined humane cutoff feature
- Adjustable beam intensity in 1% increments up to 250°C
- Reaction is detected automatically
 - automatic response in conditioned animals. You can set a humane cutoff timer that

its reaction time may be slower. This unique system offers a heated glass option that prevents the glass enclosure from acting as

with the Tail Flick meter) a heat sink, giving a more accurate reading. An optional tail temperature monitor can also be selected for use with the Tail Flick meter. This option actually preheats the tail

before experimentation. Once the preset tail

available in two sizes for large group studies.

temperature is reached, the test and timer

automatically begin. A glass stand is also

automatically shuts off the heat if no response is observed during the designated time frame. When an animal is placed on cold glass,

testing with a single unit. Either testing system
is also available individually. In plantar mode,
the visible light/heat source is directed at the
paw or other desired body part, and in tail
flick mode it is directed at the subjects' tails.
Test up to 12 mice or 6 rats simultaneously.
If desired, other animals like cats and rabbits
may also be used. Tests are simple to setup.
The focused, radiant heat/light source creates
a 4 x 6 mm intense spot. Because the light
is visible, you know when the test starts and
ends. The equipment is silent (no whining

or clicking sounds) to avoid triggering an

narcotics and strong non-narcotic drugs, offers

both Plantar (Hargreaves Method) and Tail Flick

This unit, which is designed for testing

II-336T	Combination Plantar/Tail Flick Meter, non-heated glass and tail temperature for mouse and rat
II-336TG	Combination Plantar/Tail Flick Meter, tail temperature and heated glass for mouse and rat
II-390	Plantar Test Analgesia Meter, non-heated glass for mouse and rat
II-390G	Plantar Test Analgesia Meter, heated glass for mouse and rat

Alphanumeric readout

Heated glass option

via key-pad

Manual override of all timer functions

All functions and parameters entered

Tail temperature monitor option (for use



Animal Temperature Controller

ATC2000 is a low noise heating system for maintaining animal body temperature during experimental procedures. The DC heater is extremely quiet in terms of electromagnetic radiation. This is essential in electrophysiological recordings which are very sensitive to electromagnetic interference.

The controller uses proportional, integral, and derivative (PID) technology in adjusting the DC voltage output. Compared with switched on/ off type controllers, PID controllers provide a much more precise and stable control of temperature. The PID approach is also more immune to the variation of the experimental conditions such as change in animal size and unexpected disturbances.

The controller has three temperature sensing inputs.

- The resistive temperature device (RTD) Probe input can be used to monitor an RTD rectal probe to control the animal temperature or to monitor ambient temperature, induction chamber temperature or any other temperature.
- When using a thermocouple probe, the thermocouple (TC) probe input can be used just like the RTD input. (A T type thermocouple must be used.)



- Low Noise DC Heater
- Adaptive mode-auto adjust PID regardless of animal size
- PID control for maximum temperature stability
- Three temperature inputs-RTD, TC, plate
- Automatic shutdown if the plate reaches 45°C
- The heater plate also has an internal temperature sensor. The ATC2000 monitors this sensor to prevent the localized hot spots under animal.

The controller has three operational modes:

- Normal mode uses the configured sensor (RTD or TC) or the plate sensor to drive the control loop.
- Adaptive mode uses the temperature of the heated plate and the temperature of the subject to control. This approach is less prone to overshoot, but somewhat slower the normal mode, depending on the sampling rate used.

Shutdown is a fail safe mode used if the plate temperature ever exceeds 45°C.

The ATC2000 is tuned at the factory. However, the PID parameters may also be set manually. The temperature resolution of the controller is 0.1°C. The rectal temperature probe has a 6-ft ultra-flexible, shielded cable and an RTD

The metal heating plates (available separately) have built-in temperature sensors. Compatible with stereotaxic systems, the rigid, flat surface fits under the U-frame. Plates are washable with water and detergent.



ATC2000 SPECIFICATIONS

TEMPERATURE RANGE.....up to 45°C

WEIGHT...... 5 kg (11 lb)

RESOLUTION 0.1°C	-
ACCURACY	°C
RAT SENSORRTD,	OD 2.0mm tube with 3.5mm ball
head	(Optional mouse sensor is
available)	
MAXIMUM DC OUTPUT 10V,	3A
POWER	ersal AC Adapter 90–264V Input,
1200	94.5A maximum output
DIMENSIONS 8.9 x	20.8 x 27.6 cm (3.5 x 8.2 x 10.8 in.)

ATC2000 Animal Temperature Controller **NEEDLE MICROPROBES** (see page 49 for details)

MT-29/1 29 ga 1 cm Needle Microprobe, 5-ft cable MT-23/3 23 ga 3 cm Needle Microprobe, 5-ft cable MT-D Needle Microprobe, 5-ft cable

FLEXIBLE IMPLANTABLE PROBES (see page 49 for details)

IT-18 0.025-inch diam Flexible Implantable Probe, 3-ft cable IT-23 0.009-inch diam Flexible Implantable Probe, 3-ft cable IT-1E 0.025-inch diam Flexible Implantable Probe, 3-ft cable

300443 RTD Extension Cable, 3 m

REQUIRED ACCESSORIES (select one sensor and one plate) 61800 Heating Plate with built-in RTD sensor, 15 x 25 cm

61830 Heating Plate with built-in RTD sensor, 15 x 10 cm 61840 Heating Plate with built-in RTD sensor for stereotaxic frame, 15x4 cm

REPLACEMENTS

61824 RTD Rectal Temp Probe, 1.25 mm shaft diam., 2.5mm ball tip TC Rat Rectal Temp Probe, 1 mm shaft diam., 3.2mm ball tip RET-2 RET-3 TC Mouse Rectal Temp Probe, 1 mm shaft diam., 1.6mm ball tip 503573 Silicone Pad for ATC2000 (10 x 15 cm)

Blood Pressure Measurement for Rats and Mice



- Test up to 200 animals at a time
- Sensor is MRI compatible
- Quick and accurate blood pressure measurement at temperatures as low as 32°C

This revolutionary design brings non-invasive blood pressure testing to a new level — a true turn-key system for accurate, consistent blood pressure measurement on mice, rats or any other laboratory animal

It is a compact, simple yet versatile system that can test from one to 200 animals at a time with independent control of each channel. Simple daisy-chaining allows expansion of up to 200 independently controlled systems.

All components are built into one small unit — controls, inflation of tail cuff, warming environment with whisper-quiet fans — providing an ideal system for teaching facilities and for the pharmaceutical industry when high throughput is a must.

Single animal systems are controlled from the touch screen, which allows keying in all necessary test setups. Touch screen control allows ease of operation, supplying automatic evaluation of test results systolic, diastolic, mean and heart rate.

- Highly sensitive photoelectric sensor for blood pressure detection
- Monitor, record, store or export real time systolic, diastolic, mean and heart rate

Data is collected, stored, displayed and can be transferred to the supplied memory stick. The USB interface allows for software control of multichannel systems. For single animal systems, built-in software lets you view and export data. Reports are in an Excel-style format and may be easily exported.

No computer is required. However, the analog output may be interfaced with your own data acquisition software.

System is easily cleaned. Removable trays are included with each system.

In addition to the standard one-year warranty on the system, tail cuff sensors have a lifetime warranty..

Includes all accessories necessary to run system.

Mouse Blood Pressure System II-MRBP-M II-MRBP-R Rat Blood Pressure System

Call for pricing on multiple channels.

Blood Pressure Monitor and Transducer



- Monitors animal arterial or venous blood pressure
- Displays systolic, diastolic or average blood pressure

BP1 accepts WPI's BLPR blood pressure transducer (below) as well as other blood pressure transducers.

An audio monitor provides a signal with variable pitch and amplitude, allowing you to hear changes in blood pressure. Digital LCD display provides average or peak signal values from 0 to 1999 mV. With an optional pressure gauge (not provided — see PM015D, page 185), the user may calibrate the display to read mm Hg. Recorder output connector allows direct connection to a pen recorder, oscilloscope or computer via a data acquisition system.



BLPR2 can be used for the direct arterial and venous pressure measurement in animal blood vessels. Supplied sterile, BLPR2 is accurate, linear and stable with temperature. May be sterilized cold with Cidex or a similar bactericide.

BLPR2 is equipped with a twelve-foot cable and connector compatible with WPI's four-channel signal conditioning unit, TBM4M Transbridge, and the single-channel BP1 blood pressure monitor. Cable has moisture-resistant locking connector. A continuous, uniform lumen



eliminates places for bubbles to form and lodge. The clear fluid path is easy to inspect. Easy to mount — slotted transducer body accepts Velcro strap.

To facilitate setup and operation, a four-way stopcock that allows easy filling, flushing, and zeroing of the transducer is included. Typically, the stopcock is located between the transducer and the animal catheter where it can be used to quickly zero, flush, or de-bubble the transducer.

SYS-BP1	Pressure Monitor (transducer & cable not included)		
	Specify line voltage		
ACCESSOR	IES		
BLPR2	NEW Blood Pressure Transducer & Cable		
BPCABLE2	NEW Cable (12 ft) with DIN connector for BLPR2		
503067	BLPR2 Transducer without cable		
13024	Single Rack Mount Kit		
13025	Dual Rack Mount Kit		
3491	Extension Cable, 5 ft		
500184	BNC-to-BNC Cable, 10 ft		
14036-15	4-Way Luer Stopcock, Blue Tint (package of 15)		
KZ1101	Micro Cannula, 3 inch		

Note: BLPR2 is intended for animal research only and may not be used for human blood pressure measurement.

BP1 SPECIFICATIONS

AMPLIFICATION x1, x10, x100, variable (x5 to x1000)

OUTPUT VOLTAGE SWING + 5 volts MAXIMUM OUTPUT CURRENT

INPUT IMPEDANCE, EACH INPUT

100 kΩ | | 0.01 μF

TRANSDUCER APPLIED VOLTS

10 V nominal, varies with load. 25 mA,

maximum

95-135 V or 220-240 V, 50/60 Hz

8.5 x 5.12 x 10 in. (21.6 x 13 x 25.44 cm)

SHIPPING WEIGHT 11 lb (5 kg)

BLPR2 SPECIFICATIONS

WORKING PRESSURE -50 to + 300 mm Hg OVERPRESSURE -400 to +4000 mm Hg 1-10 VDC or RMS to 5 kHz **EXCITATION VOLTAGE**

SENSITIVITY 5 μV/V/mm Hg DYNAMIC RESPONSE 100 Hz

EIGHT HOUR DRIFT 1 mm Hg after 5 minute warm-up

MAXIMUM FRROR Total combined effects of Sensitivity, Linearity,

Hysteresis (at 25°C and 5 µV/V/mm Hg) do not exceed ±2% or 1 mm Hg, whichever is greater.


SHIPPING WEIGHT

www.wpiinc.com

POWER

DIMENSIONS





3-dimensional placement of any probe or accessory within the working space of a stereotaxic frame. No computer is required. The MTM-3 supplied with WPI stereotaxic frames is also compatible with standard stereotaxic frames and can be adapted to existing frames of other manufacturers. Single and dual manipulator arm motorized systems are available.



Increased Precision and Repeatability of motion over traditional manual Stereotaxic frames

- Accurate microstepping motor drive for high resolution placement
- Set the "final approach" speed between 2mm/sec and .02mm/sec

Increased Convenience, decreased error measurement

- No more error resulting from reading Vernier scales
- Brain atlas coordinates can be input into the controller, no computer required
- Coordinate distances are automatically calculated
- Touch screen for ease of control
- Graphic controller display for instant operational feedback

MTM-3 Operation

Manual mode

- Move the actuator using the touch panel or the intuitive manual 3 axis wheel controller
- Individual axes may be easily disabled/enabled to ensure motion on only the desired axis
- Controller allows for three different speed sensitivity levels

Coordinates mode

- Specify retracted position
- Store three origin definitions (references)
- Probe position may be displayed with respect to any of the references, as well as absolute coordinates

MTM-3	3-Axis Stereotaxtic Controller (Left or Right)	
MTM-3B	3-Axis Stereotaxtic Controller with Bluetooth (L or R)	Ī
MTM-6	6-Axis Stereotaxtic Controller (Left and Right)	
MTM-6B	6-Axis Stereotaxtic Controller with Bluetooth	
MTM	3-Axis Stereotaxtic Arm and Controller only (L or R)	
MTMB	3-Axis Stereotactic Arm and Controller only, Bluetooth	

- Position plot cursor graphically represents the actuator arm position with respect to any of the stored references
- All functions are accomplished without the use of a computer (computer interface available, if desired)

Optional Computer Control

- Remote computer control through a USB port
- Wireless Bluetooth computer control for the MTM-3 models equipped with a wireless option
- Text based command set provided
- Use a terminal program to create simple scripts for repetitive operations



Precision Stereotaxic Instruments for Small Animals


WPI's Precision Stereotaxic Instrument is built around the time-proven U-frame design concept, providing stability, and adaptability to most species. Precision alignment ensures accurate placement of electrodes, micropipettes, and other devices. It is ideal for researchers in need of a versatile, reliable instrument for stereotaxic procedures with small animals

Versatility of positioning

The manipulator arm controls medio-lateral and vertical positioning via lead screws with 80 mm of travel. This allows the fastest positioning possible, consistent with lining up the scales easily at a given coordinate. The antero-posterior movement is controlled via a dovetail slide movement, with 80 mm of travel possible in each direction. A universal joint allows the investigator to change the angle of the probe up to 90° in either the antero-posterior or medio-lateral planes. The locking mechanism will hold any angle position without drift or creep. It also provides an absolute lock at 90° vertical.



Above: 502600 Precision Stereotaxic Frame. At right: 502603 Dual Manipulator Stereotaxic Frame.






Easily read scales

All scales are oriented to be read easily from the open end of the "U". This is the position from which most scientists prefer to work. The numerals on the scales are clear and easy to read. Precise alignment with facing vernier scales gives accurate resolution to 0.1 mm.

Convenient for electrophysiology

The entire Stereotaxic frame including the dovetails, manipulator arms and base are electrically continuous. Grounding of the entire frame including the base plate can be accomplished by connecting the provided grounding stud to earth. This is ideal for electrophysiological studies where the animal and surrounding structures need to be grounded to reduce electrical noise.

502600	Precision Stereotaxic Frame with 18°Ear Bars
502650	Precision Stereotaxic Frame with 45°Ear Bars
502603	Dual Manipulator Stereotaxic Frame with 18°Ear Bars
502653	Dual Manipulator Stereotaxic Frame with 45°Ear Bars
TAXIC-600	Stereotaxic Frame with 18°Ear Bars plus UMP3-1 System
TAXIC-650	Stereotaxic Frame with 45°Ear Bars plus UMP3-1 System
TAXIC-603	Dual Manipulator Stereotaxic Frame with 18°Ear Bars plus UMP3-1 System
TAXIC-653	Dual Manipulator Stereotaxic Frame with 45°Ear Bars plus UMP3-1 System



Digital Stereotaxic Frame with LCD Display



This new Digital Stereotaxic Frame features sealed electronic sensors and an easy-to-read LCD display with 10-micron resolution in all three axes. A zeroing function aids in targeting specific coordinates. The battery-powered display is electronically quiet, making it useful in electrophysiology experiments as well as keeping the workbench tidy.

Features

- Adaptors available for use with rats, mice, birds, cats, geckos, guinea pigs and other
- 80 mm of vertical, lateral and anteriorposterior travel
- 180° rotation and lock at any vertical angle

- 360° rotation and lock at any level angle
- Syringe pump, microcamera and drill can be attached directly
- Dual-lead screws provide stable and accurate movement
- Maintains accuracy and flexibility under different temperatures
- Zeroing function for targeting specific
- 10µm precision of manipulator in all directions
- Desktop display is separate from sensors, making it easy to read
- Battery-powered sensors, without electronic noise, are suitable for electrophysiology experiment

- Extended base plate 400mmX255mm suitable for various animals
- Vertical lock and fixing knob are separated to ensure accurate position at any angle
- The precisely-designed rotary knob and U frame allow sufficient space for the anterior-posterior operation
- Laser engraved scales and darkened rod make numerals easy to read
- Curved design of nose clamp fixes the head of the animals more securely

502300 Digital Stereotaxic Frame with 18°Ear Bars 502350 Digital Stereotaxic Frame with 45°Ear Bars 502303 Dual Manipulator Digital Stereotaxic Frame with 18°Ear Bars 502353 Dual Manipulator Digital Stereotaxic Frame with 45°Ear Bars **TAXIC-300** Digital Stereotaxic Frame with 18° Ear Bars and UMP3-1 TAXIC-350 Digital Stereotaxic Frame with 45° Ear Bars with UMP3-1

Also available with UMP3-1 UltraMicroPump special pricing!



502227	Stereotaxic Frame System with one manipulator for Cat and Monkey
502228	Stereotaxic Frame System with two manipulators for Cat and Monkey
502229	Stereotaxic Frame System with three manipulators for Cat and Monkey
502230	Stereotaxic Frame System with four manipulators for Cat and Monkey
502231	Stereotaxic Frame System with one manipulator for Dog
502232	Stereotaxic Frame System with two manipulators for Dog
502233	Stereotaxic Frame System with three manipulators for Dog
502234	Stereotaxic Frame System with four manipulators for Dog

For more options and accessories, visit the website — www.wpiinc.com

Stereotaxic Accessories

ADAPTORS

502063	Mouse and Neonatal Rat Adaptor
502213	Platform, Gas Anesthesia, with mouse Mask (use with 502063)
502062	Mouse Adaptor
502204	Rat Adaptor with a pair of ear bars, 18°
502226	Cat/Monkey Adaptor for 502600 series
502238	Spinal Adaptor for Rat
502060	Guinea Pig Adaptor for 502600 series
502241	Dog/Monkey Adaptor for Parallel Rail Stereotaxic instruments, with a pair of ear bars, 18°

The WPI Mouse and Neonatal Rat Adaptor (502063) employs light, Delrin® adjustable ear bars with tapered points on one end and thumbscrew on the other to facilitate surgery on mice and rat pups. Adjustable ear bars may be independently adjusted in height to level the skull. Laser engraved scales show the vertical positions of the ear bars. A tooth bar and nose clamp secures the nose. A well in the thick aluminum body may be filled with dry ice and alcohol for hypothermic anesthesia of neonatal animals. The adaptor clamps securely on the right side of the "U" frame of the stereotaxic instrument. When gas anesthesia is needed, a gas mask assembly (502213) may be mounted on the adaptor.
















502235	502056	502225	502242
Bullium and bulliu	infinite a justinified	2 2 D	a minimum o
		1	

OTHER A	ACCESSORIES
502053	Mask, Gas Anesthesia, Mouse
502054	Mask, Gas Anesthesia, Rat
502201	V-Clamp, 10/32 screw
504608	V-Clamp screw
502213	Platform, Gas Anesthesia, with mouse Mask (use with 502063)
502243	Adjustable Stage Platform for 502600 series, 2cm high
503598	Micro-Drill, 35K rpm, 110/220VAC, w/ a set of bits
503599	Micro-Drill, 35K rpm, 240VAC, w/ a set of bits
61840	Heating Plate for 502063, 4X15cm, 5mm thick (use with ATC2000)

		502244
EAR BARS		
502055	Ear Bars, Rat, 18°, (pair)	
502056	Ear Bars, Rat, 45°, (pair)	
502224	Ear Bars, Cat, 18°, (pair)	
502225	Ear Bars, Cat, 45°, (pair)	
502235	Ear Bars, Mouse, 60°. Non-rupture, (pair)	
502242	Ear Bars, Rat, Hollow. 1.5mm hole for auditory stimulation	



502259—Spinal Adaptor for Rat

Microprobe Thermometers



A Microprobe Thermometer is the instrument of choice for biological and laboratory temperature measurements. These thermometers are very versatile, providing fast response, high accuracy and stability with digital display and analog signal for connection to a computer or recorder. With the wide selection of probes, the instruments can be used in almost any application.

BAT-12 This thermometer has a sealed construction making it water, dust and fume resistant. The BAT-12 has a single microprobe input and a single range with the same high accuracy as the BAT-10. Comes complete with carrying case.

The thermometers can be used with any "Type T" thermocouple. Select a temperature microprobe on the following page for your specific application.

Tripod Stand for BAT-12

BAT-10 This is the most versatile thermometer available. The instrument has a wide temperature range and fast response with most microprobes. The BAT-10

- Super Accuracy Fast Response
- Analog output signal
 Multiple inputs
- Differential Temperature Measurement



accuracy is NIST traceable and in each of the two temperature ranges, the accuracy is the same as the resolution. There are three microprobe inputs, 1 and 2 can be selected as separate inputs while 2 and 3 will read the differential temperature measurement between the two. The instrument has automatic warnings for low battery or faulty probes on the digital display. The linearized analog output (LOP) signal allows ease of connection to a data acquisition system

BAT-10R/LOP	Multiple Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 110 VAC adapter (microprobes ordered separately)	
BAT-10R/LOP220	Multiple Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 220 VAC adapter (microprobes ordered separately)	
BAT-12R	Single Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 110 VAC adapter (microprobes ordered separately)	
BAT-12R-220	Single Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 220 VAC adapter (microprobes ordered separately)	
OPTIONAL ACCESSORIES		
EXT-6	Probe Extension Lead 180 cm long	

	BAT-10	BAT-12
TEMPERATURE RANGE & RESOLUTION	-200°C to +400°C, 1°C resolution -100°C to +199.9°C, 0.1°C resolution	-100°C to +199.9°C, 0.1°C resolution
DIFFERENTIAL TEMP. RANGE	-19.99°C to +19.99°C Linearization centered at 40°C 0.01°C resolution	N/A
ACCURACY 1° Range 0.1° Range Diff. Range	1°C ±1 least significant digit 0.1°C ± 1 least significant digit 0.01°C ± 1 least significant digit	0.1°C ± 1 digit between 0-50°C 0.1% ± 1 digit over full range
REPEATABILITY	± 1 least significant digit	
CALIBRATION CONFORMITY	Conforms to NIST tables	Follows NIST thermocouple tables within 1 digit
DISPLAY	3½ Digit LCD	3½ Digit LCD
INPUT SOCKET	Miniature, quick disconnect, copper-constanta	nMiniature, quick disconnect, copper-constanta
ANALOG OUTPUT	Non-linearized set at 1.6 V, corresponding to temperature of 401°C	≈ 10 mV per degree C
POWER SUPPLY / BATTERIES	BAT-10: 4 alkaline "C" cells (life: 1000 hr) BAT-10R: 4 Ni-Cad "C" cells (rechargeable unit)	BAT-12: 9V cell BAT-12R: 9V Ni-Cad with charger
SENSORS	Three Type T thermocouple inputs	One Type T thermocouple input
AMBIENT OPERATING RANGE	15-45°C	Auto-compensated to 0.1°C from 0°C to 50°C
DIMENSIONS	21.6 x 22.9 x 8.9 cm (8.5 x 9 x 3.5 in.)	12.7 x 6.4 x 15.2 cm (5 x 2.5 x 6 in.)
WEIGHT	1.6 kg (3.5 lb), including carrying case	1 kg (2 lb), including carrying case

Temperature Probes


- Flexible Teflon microprobes are used for implantation in tissue, in spectrophotometer cuvettes, rectally in neonatal mice, in water baths, PCR thermal cyclers, etc.
- Animal rectal temperatures during surgical procedures and pyrogen testing.
- Skin temperature measurement during exercise physiology studies.

When precise temperature measurements are required, WPI can provide you with a very accurate monitor and thermocouple microprobes. WPI monitors have both resolution and accuracy of 0.1°C in the 0-50°C range and are traceable to NIST standards, whereas, other competitive electronic

thermometers have an accuracy that is usually to 0.5°C or worse. Furthermore, all our type T clinical probes are guaranteed accurate to 0.1°C, due to our stringent wire standards. These are five times more accurate than competitive probes made with regular "Special Limits" wire.



Probe Type	Size	Style	Time Constant	Isolated	Max. Temp.	Lead Length	Description
NEEDLE M	ICROPROBI	ES					
ast-response r	needle probes for	instant read	dings in tissue, se	misolids, liquids	, very small s	pecimens, powder	s and materials. Needle tip is sealer
ensure only:	stainless steel con	itacts specin	nen.				
/IT-29/1	29 ga / 1 cm	А	0.125 sec	No	200°c	5 ft	29 gauge approximately 0.013-in
/IT-29/2	29 ga / 2 cm	A	0.125 sec	No	200°c	5 ft	-
/IT-29/3	29 ga / 3 cm	A	0.125 sec	No	200°c	5 ft	-
/IT-29/5	29 ga / 5 cm	A	0.125 sec	No	200°c	5 ft	-
IT-26/2	26 ga / 2 cm	A	0.1 sec	No	200°c	5 ft	26 gauge approximately 0.018-in
MT-26/4	26 ga / 4 cm	A	0.1 sec	No No	200°c	5 ft	-
/T-26/6	26 ga / 6 cm	A	0.1 sec	No	200°c	5 ft	22
/IT-23/3	23 ga / 3 cm	A	0.15 sec	No No	200°c	5 ft	23 gauge approximately 0.125-in
<u>/IT-23/5</u> /IT-23/8	23 ga / 5 cm	A A	0.15 sec 0.15 sec	No No	200°c 200°c	5 ft 5 ft	
/11-23/6 /IT-4	23 ga / 8 cm	A	0.15 sec 0.025 sec	No No	200°c	5 ft	Similar to MT-29/1 but has a blunt tip.
11-4	29 ga / 1 cm	Α	0.023 Sec	INU	200 C		istant skin and surface temperatures, liquids
/IT-D		C	0.025 sec	No	200°c	5 ft	Fast response surface probe (stainless steel for locating
מ-ווו		C	0.023 300	INU	200 (310	inflammation, arteries, etc. Also for dental use.
							illiallillation, arteries, etc. Also for delital use.
	IMPLANTAI						
esigned for hi		tremely sm	all specimens suc	n as insects, see	ds, etc. Maxir	num insertion dep	th 1/8". Totally sheathed in chemical resistant Teflon.
	Sensor						
	Lead Diameter		0.2	V	1500-	2.0	
Γ-14 Γ-10	0.050" dia	D	0.3 sec	Yes	150°c	3 ft	-
Γ-18 Γ-19ΕΥΙΙΟΝΙΟ	0.025" dia	D	0.1 sec	Yes	150°c	3 ft	
I-18EXLUNG	0.025" dia.	D	-	Yes	150°c	5 ft	-
Г-21	0.016" dia	D	0.08 sec	Yes	150°c	1 ft	
Γ-23	0.010 dia	E	0.005 sec	Yes	150°c	3 ft	For ultra fast measurements and for use on micro-size
1-25	0.000 uld	L	0.000 300	103	150 €		Tissue implantable with 239a. Needle
							Rather fragile. Teflon coated.
T-1E	0.025" dia	F	0.005 sec	Yes	150°c	3 ft	As IT-18 sensor except bead exposed. Combines ultra-
1-16	0.025 dia	'	0.005 300	103	150 €		e of IT-23 with sheath strength of IT-18.
ECTAL D	2005					iust repons	e of the 25 which shedder still right of the 10.
ECTAL PI	KORF2	G	0.0	Ma	125°c	г.	Dartel and a far arte torically for fact interestate at
RET-2	-	G	0.8 sec	No	125-0	5 ft	Rectal probe for rats typically for fast intermittent measurements. Smooth ball tip (0.125-in. dia.) with
						111 1 (0.0	
ET-3		G	0.5 sec	No	125°c	5 ft	59-in. dia) stainless steel shaft. Rectal probe for mice similar to RET-2 except tip diam.
E1-3	-	u	0.5 260	INU	123 (nd shaft 3/4-in. long (0.028-in. diam.)
						0.005-III. di	10 Stratt 3/4-111. 1011g (0.026-111. utatit.)
	PURPOSE						
IT-1		Н	0.5 sec	No	400°c	5 ft	"Workhorse" probe for liquids, gases, semi-solids.
							Plastic handle with straight stainless steel shaft. Not
			2.5		1000		rface temperatures.
IT-2		H	0.5 sec	No	400°c	5 ft	Like HT-1 except shaft length is 9-in.
BT-1		I	0.15 sec	No	240°c	5 ft	Plastic handle with welded stainless steel, immersible
							for surface temperatures of solids, liquids,
					0		o is 0.02-in. diam., at right
				angle to	nrohe to fac	ilitate surface mea	curament



See page 36

Neuroscience Cannulas

For in vivo investigation of rodents



Cannula assembly: fixing screw, internal cannula with attached tubing, and guide cannula.



Internal cannula secured with fixing screw.

Exceptional Quality Best Prices Rapid Order Response

This cannula system for neuroscience study and pre-clinical research includes an entire range of cannula options. The three primary components include the Guide Cannula, the Internal Cannula and the Dummy Cannula

- WPI cannulas are beveled inside and out and then polished to remove any burrs and ensure that the inside diameter is perfectly cylindrical. This limits undesired trauma to tissue and ensures smooth
- WPI offers exceptional quality at the best
- Quantities of each item are kept in stock so you can order as needed.

Guide Cannulas	Gauge	OD	ID
	22	0.64	0.46
A STATE OF THE STA	24	0.56	0.40
-	26	0.48	0.32
944444		-0	

The Guide Cannula is a surgical grade, stainless steel tube that is implanted into a rodent's skull and cemented into place using dental cement and screws. It guides the Internal Cannula to the specific injection site.

Internal Cannula	Gauge	OD	ID
1	22	0.41	0.26
• 0.7mm	24	0.36	0.21
	26	0.30	0.16
A			

The Internal Cannula is inserted into the Guide Cannula to sample or inject fluid.

Dummy Cannula



The Dummy Cannula has a stainless steel wire

core, and it is placed in the Guide Cannula when the Internal Cannula is removed. It seals the opening and prevents tissue from entering the bottom of the Guide Cannula. The Dummy Cannula is threaded to securely tighten it so that the animal will not unscrew it while grooming.

Ordering

Order the Guide Cannulas based on the gauge and length from the base and the Internal or Dummy Cannulas based on the length of the Guide Cannula and the projection from the Guide Cannula tip.

Understanding Part Numbers

Cannula Type

GC22-37 Cannula Type Gauge of Guide Cannula Length of Cannula or Projection

GC-Guide Cannula

INC-Internal Cannula **DUMC**-Dummy Cannula

Gauge-Choose the gauge of the Guide Cannula that will be used, even if you are specifying an

Internal or Dummy Cannula. Choices include 22, 24 and 26 gauge.


Length–Guide Cannulas can be ordered in a range of lengths from 1.0-9.9mm. Length can be specified to 0.1mm, with a tolerance of ±0.07mm. For example, a **GC24-60** is a 24ga. cannula that is 6.0mm long.

For other cannulas, the length is determined by the desired projection length and the guide cannula length. The projection can extend beyond the tip of the Guide Cannula up to 1.0mm.

DUMC Length = GC Length + Projection

INC Length = GC Length + Projection

For example, for a Dummy Cannula flush with the end of a 26ga Guide Cannula that is 6.0mm long, order a **DUMC26-60**. For an inernal cannula that projects 0.5mm beyond the 6.0mm Guide Cannula (6.5mm), order an INC26-65.



* DUMC Length=GC Lenth + Projection. If the cap is screwed on too tightly, the projection will be longer than expected.
** INC Length=GC Lenth + Projection. Internal cannula mounts flush and does not screw into place.

Cannulas		Qty 1-9	Qty 10-19	Qty 20-49	Qty 50-99	Qty 100+
GC22-X	Guide Cannula, 22 Gauge, X.0mm					
INC22-X	Internal Cannula, 22 Gauge, X.0mm					
DUMC22-X	Dummy Cap Cannula, 22 Gauge, X.0mm					
GC24-X	Guide Cannula, 24 Gauge, X.0mm					
INC24-X	Internal Cannula, 24 Gauge, X.0mm					
GC26-X	Guide Cannula, 26 Gauge, X.0mm					
INC26-X	Internal Cannula, 26 Gauge, X.0mm					
DUMC26-X	Dummy Cap Cannula, 26 Gauge, X.0mm					

Flexible	PE Tubing	Qty 1-9	Qty 10-19	Qty 20-49	Qty 50-99	Qty 100+
504278	0.25mm ID, 0.5mm OD, 1m long					
504279	0.42mm ID, 0.85mm OD, 1m long					
504280	0.6 mm ID, 1.1mm OD, 1m long					
	Flexible PE tubing are recommended	d to conn	ect to Intern	al Cannula	(INC) only,	

NOT Guide Cannula (GC) or Dummy Cap (DUMC). 504278 matches INC26-XX 504279 matches INC24-XX 504280 matches INC22-XX

Misc. Accessories Qty 1-9 Qty 10-19 Qty 20-49 Qty 50-99 Qty 100+ Fixing Screw for connection of INC & GC 504281 504282 Instant Cyanoacrylate Adhesive, 3g 504283 Instant Cyanoacrylate Adhesive, 20g

Sensors

MACRO SENSORS


SPECIES	Carbon Monoxide	Nitric Oxide	Hydrogen Peroxide	Oxygen	Hydrogen Sulfide
Order Number	ISO-COP-2	ISO-NOP	ISO-HPO-2	ISO-OXY-2	ISO-H2S-2
Price					
Available Diameters	2 mm	2 mm	2 mm	2 mm	2 mm
Response Time	< 10 sec	< 5 sec	< 5 sec	< 10 sec	< 5 sec
Detection Limit/Range	10nM to 10μM	1 NM to 40μM*	< 100nM to 100μM	0.1%-100%	< 5nM-100μM
Sensitivity	~0.5 pA/nM	≤2 pA/nM	8 pA/μM	0.3-0.6nA/%	2 pA/nM
Drift	<1pA/min	<1pA/min	0.1pA/min	< 1%/min	
Temperature Dependent	Yes	Yes	Yes	Yes	Yes
Physiological Interference	nitric oxide	NaNO ₂ (10 ⁻⁶ or better) None	None	None
Replacement Sleeves (pkg of 4)	#95620	#5436	#600012	#5378	#600016
Filling Solution	#95611	#7325	#100042	#7326	#100084
Start-up Kit	#95699	#5435	#600011	#5377	#600015

^{*} Higher detection limit available on request — call for custom pricing.

MICRO SENSORS												
	ISO-NOPF200	ISO-NOPF200-Lxx ³	ISO-NOPF100 ISO-NOPF100-L ³	ISO-NOP70Lxx³	ISO-NOPF500-Cxx	ISO-NOP3005	ISO-NOP3020	ISO-NOP30L³	ISO-NOP007	ISO-NOPNM	ISO-HPO-100	ISO-H2S-100
Species					Nitric	Oxide					H ₂ O ₂	H ₂ S
Price	(pkg of 2)	(pkg of 3)	(pkg of 2)	(pkg of 2)	(pkg of 2)	(pkg of 3)	(pkg of 3)	(pkg of 3)	(pkg of 3)	(pkg of 3)	(pkg of 3)	(pkg of 2)
Fiber Diameter (µm)	200	200	100	70	500	30	30	30	7	7 Conical tip: 100nm	100	100
Tip Length ² (mm)	1-5 ¹	1-10 ¹	1-5 ¹	3	5–10	0.5	2	3	2	2	1-41	2-5 ¹
Response Time (sec.)	< 5	< 5	< 5	< 3	< 10	< 3	< 3	< 3	< 3	< 3	< 5	~5
Lowest Detection Limit/Range (nM)	0.2	0.2	0.2	1	0.2	1	1	1	0.5	0.5	1	<5
Nominal Sensitivity-New Sensor ² (pA/nM)	≥20	≥50	≥10	≥10	≥ 20	≥1	≥1.5	≥1	≥1	≥0.5	≥1	1-4
Baseline Drift (pA/min)	none	none	none	none	none	none	none	none	none	none	<2.0	<2

¹Sensor available in 1mm length increments (for example, 1mm, 2mm, 3mm...).

Some nitric oxide sensors are available in custom lengths. When ordering custom lengths, use the part numbers ISO-NOPF100-Cxx or ISO-NOPF200-Cxx and replace the xx with the desired length. For example, for a 1mm flexible sensor tip, the part number should be ISO-NOPF200-C01. Sensors can be ordered in the following custom lengths: 1mm, 2mm, 3mm, 4mm or 5mm.



²Sensor sensitivity varies with length and diameter.

Any 100µm sensor can be purchased with a hypodermic sheath. Add a -H to the end of the part number (for example, ISO-HPO-100-H).

³L-shaped sensor for use with a tissue bath.

Nitric Oxide Sensor Guide

All WPI NO sensors are 100% compatible with ISO-NO Mark II (NOMK2), APOLLO 4000, APOLLO 1000, and TBR4100 Free Radical Analyzer.

APPLICATIONS	in vivo	Cell Cultures, NO ₂ , NO ₃ ,	Tissue Bath	Tissue Bath	Microvessels	Microvessels	Single Cell
SENSOR	ISO-NOPF	ISO-NOP	ISO-NOP30L	ISO-NOP70L	ISO-NOP30	ISO-NOP007	ISO-NOPNM
SENSOR DIAMETER	100, 200, or 500 μm	2 mm	30 micron	70 micron	30 micron	7 micron	100 nm
RESPONSE TIME (with NOMK2)	< 5 sec	< 5 sec	< 3 sec	< 3 sec	< 3 sec	< 3 sec	< 3 sec
LOWEST DETECTION LIMIT	0.2 nM	1 nM	1 nM	1 nM	1 nM	0.5 nM	0.5 nM
TEMPERATURE SENSITIVITY	some	yes	yes	yes	yes	yes	some
DRIFT	none	none	none	none	none	none	none
SENSITIVITY	10 pA/nM	2 pA/nM	1.4 pA/nM	1.4 pA/nM	1.4 pA/nM	1 pA/nM	0.5 pA/nM
PHYSIOLOGICAL INTERFERENCE	none	none	none	none	none	none	none

Selectivity of WPI's NO sensors

The ideal NO sensor should be insensitive to other reactive species likely to be present within the measurement environment. The conventional Nafion coated carbon fiber NO sensor exhibits a large response to such

species. WPI's unique NO sensor technology utilizes a novel surface membrane which amplifies the response to NO while eliminating responses to a vast range of reactive species, including nitrite, absorbic acid, hydrogen peroxide, catecolamines, and much more.

SO-NOP - The original nitric oxide probe - ideal for cell cultures, cell suspensions and many other applications



The ISO-NOP is a popular, robust and high performance sensor encased within a 2 mm diameter disposable stainless steel protective sleeve. The tip of the sleeve is covered with a NO-selective membrane. Replacement membrane sleeves can be purchased separately (WPI #5436) and require an internal electrolyte (WPI #7325).

NO3 and NO2 Detection


A simple change in experimental protocol will enable the ISO-NOP to be conveniently used for indirect rapid accurate determination of nitrite (NO₂) and nitrate (NO₃) concentration in samples. Using this method a detection limit for NO₂ or NO₃ as low as 1 nM is routinely possible.

> Abdominal X-ray showing the appratus consisting of two customized ISO-NOP nitric oxide probes, 4-channel pH catheter, and Teflon nasogastric tube. (Courtesy Prof. K.E.L. McColl, University Department of Medicine and Therapeutics, Western Infirmary, Glasgow, Scotland.) lijima, K., et al. Gastroenterology 2002: **122**: 1248-1257.



ISO-NOP	Replacement 2 mm shielded sensor and cable
5435	ISONOP Startup Kit (recommended with first purchase)
5436	Replacement Sleeve Kit for 2 mm sensor, pkg of 4
7325	ISO-NO Electrolyte (10 mL)
7521	ISO-NO Electrolyte, CO2-insensitive (10 mL)
5399	T-Adapter Kit (pkg of 3) for ISO-NOP
7357	Nitrite Standard Solution, 1 gram/liter (100mL)
	, ,


ISO-NOPNM - The world's smallest nitric oxide NanoSensor, designed for measurement of NO at the cellular level.



Schematic drawing of the new integrated NO nanosensor. (US Patent Pending)

The ISO-NOPNM NanoSensor has a tip diameter of just 100 nm (0.1 µm) and a detection limit for NO of less than 0.5 nM — making it indisputably the smallest and most sensitive NO sensor in the world!

The ISO-NOPNM is based on a novel design in which an electrochemically "activated" composite graphite nanofibre is used as the NO-sensing element. The surface of the NanoSensor is then modified using a unique multi-layered NO-selective membrane. Figure at right illustrates the response of the ISO-NOPNM following successive additions of nanomolar concentrations of NO. The ultra-low noise of the ISO-NOPNM (0.5 pA) enables a detection limit of just 0.5 nM NO. The response time of ISO-NOPNM is less than 3 seconds.




Amperometric response of the NO nanosensor (ISO-NOPNM) to the successive additions of 2nM, 4 nM, 8 nM NO into 0.1 M PBS (pH=7.4).

ISO-NOPNM	100 nm NanoSensor, pkg of 3 (requires cable #91580)
91580	Microsensor Adapter Cable
SNAP50	SNAP, 50 mg vial

ISO-NOPF - unique flexible NO sensor! Designed for arteries, microvessels, in vivo applications, and similar applications.

ISO-NOPF electrodes are the newest addition to WPI's nitric oxide sensor family and are available in 100 µm and 200 µm diameters. Utilizing the latest advances in nanotechnology and material science, scientists at WPI's Sensor Laboratory have created these completely flexible and virtually unbreakable NO sensors. The new sensors are based on a composite graphite NO-sensing element combined with a reference electrode. The surface of the sensor is then coated with a unique multi-layered NO-selective membrane.





Sensor Body

 $60\text{mm x} < 2\text{mm}\varnothing$

Response of ISO-NOPF to NO.

Schematic		

Ag/AgCl

flexible

100, 200


or 500µm

ISO-NOPF100	100 μm Flexible NO Sensor, pkg of 2
ISO-NOPF200	200 μm Flexible NO Sensor, pkg of 2
ISO-NOPF500	500 μm Flexible NO Sensor, pkg of 2
91580	Microsensor Adapter Cable
SNAP50	SNAP, 50 mg vial


ISO-NOP007 ISO-NOP30 - 7 and 30 micron sensors with exceptional performance — ideal for tissues and microvessels

The ISO-NOP007 and ISO-NOP30 have recently been improved in design and performance. The ISO-NOP007 has a tip diameter of just 7 microns and a length of 2 mm. The ISO-NOP30 has a tip diameter

of 30 microns and is available in two different tip lengths (i.e., ISO-NOP3020 has tip length of 2 mm, ISO-NO3005 has tip length of 0.5 mm). The response of the ISO-NOP007 and ISO-NOP30 is linear over a wide dynamic concentration range of NO. The design of both electrodes is based on a single carbon fiber coated with WPI's NOselective membrane. A detection limit of approximately 1 nM NO makes these electrodes ideal for use in tissues and microvessels.




Schematic drawing of ISO-NOP007 and ISO-NOP30.



The response of a 7 µm NO sensor (ISO-NOP007) to successive additions of NO (100 nM). Inset shows the linearity of the resulting calibration plot.

ISO-NOP007	7 μm Nitric Oxide Sensor (pkg of 3)	
ISO-NOP3020	30 µm Sensor Tips (2 mm length), pkg of 3 (requires #91580)	
ISO-NOP3005	30 µm Sensor Tips (0.5 mm fiber), pkg of 3 (requires #91580)	
91580	Microsensor Adapter Cable	
SNAP50	SNAP, 50 mg vial	



Simultaneous measurement of force (top trace) and changes of NO concentration (lower trace) in (A) the rat superior mesenteric artery relaxed with ACh and (B) a small human artery relaxed with ACh and SNAP. In this artery oxyhaemoglobin(oxyHb) partly reversed the increase in NO concentration, with only a small change in force. [U. Simonsen, et al., J. Physiol., 1999, **516**: 271-282.1


ISO-NOP30L ISO-NOP70L ISO-HPO-100-L ISO-NOPF200-L10

L-shaped sensors for tissue bath & cell culture studies

The ISO-NOP30-L is a unique L-shaped nitric oxide sensor designed specifically for use in tissue bath studies and similar applications (e.g., see WPI's MYOBATH). The shape of the sensor has been engineered to facilitate

placement of the electrode within the lumen of the tissue vessel under study. The ISO-NOP70-L is similar in construction to the ISO-NOP30 but with the advantage of having a flexible tip (70 µm diameter). The ISO-NOPF200-L10 is designed specifically for cell culture studies.

ISO-NOP30-L	NO Sensor, L-Shaped 30-micron (pkg of 3)		
ISO-NOP-70-L	NO Sensor, L-Shaped 70-micron (pkg of 2)		
ISO-NOPF200-L10	NO Sensor, 200 μm Flexible L-shaped (pkg of 2)		
ISO-HPO-100-L	HPO Sensor, L-Shaped 100-micron (pkg of 2)		
91580	Microsensor Adapter Cable		



HYDROGEN PEROXIDE & OXYGEN SENSOR GUIDE

	ISO-HPO-2	ISO-HPO-100	ISO-HPO-100H	ISO-HPO-100-L	ISO-OXY-2
APPLICATION	Cell Cultures, etc.	Tissue/ Microvessels	Hypodermic Sheath	Tissue Bath	Cell Cultures, etc.
SENSOR DIAMETER	2.0 mm	100 micron	100 micron	100 micron	2.0 mm
RESPONSE TIME	<5 SEC (90%)	<5 SEC (90%)	<5 SEC (90%)	<5 SEC (90%)	<10 SEC (90%)
DETECTION LIMIT	<100 nM to 100 mM	1 nM to 1 mM	<10 nM to 1 mM	1 nM to 1 mM	0.1 % to 100%
DRIFT	<0.1 pA/min	<1.0 pA/min	1.0 pA/min	<1.0 pA/min	<1%/min
SENSITIVITY	8 pA/μM	1 pA/nM	1 pA/nM	1 pA/nM	N/A
PHYSIOLOGICAL INTERFERENCE	none	Contact WPI	Contact WPI	Contact WPI	none

Hydrogen Peroxide Sensors

Hydrogen Peroxide is produced in biological systems by controlled pathways at low concentrations that impact on cell signaling. At higher concentrations inflammatory cells produce local intense amounts of this oxidant to kill pathogens. In the progress of human disease, uncontrolled formation of hydrogen peroxide from the mitochondrial respiratory chain and enzymes. such as xanthine oxidase, can occur (Prof. Victor Darley-Usmar, Univ. of Alabama, personal communication). Despite the recognized importance of this oxidant in biology real-time measurements at low concentration have been difficult. The hydrogen peroxide sensors developed by WPI are designed to compliment existing high sensitivity fluorescent approaches with direct quantitative measurement in biological samples in the low nM range.

Four hydrogen peroxide sensors are currently available. The ISO-HPO-2 is a 2.0 mm stainless steel sensor, with replaceable membrane sleeves (#600012) and an internal refillable electrolyte (#100042). The sensor is designed for use in cell cultures and similar applications.

The ISO-HPO-100 is a 100 micron tip diameter hydrogen peroxide micro sensor designed for use in tissues and similar applications. The sensor design is based on a flexible "activated" carbon fiber sensing electrode coated with a proprietary membrane that enhances hydrogen peroxide detection.

ISO-HPO-2

Both sensors incorporate WPI's proprietary combination electrode technology whereby the hydrogen peroxide sensing element and separate



reference electrode are encased within a single Faraday-shielded probe design. This design has been shown to enhance performance during measurements and minimizes overall sensor size.




600011	ISO-HPO Startup Kit (recommended with first purchase)	
ISO-HPO-2	2mm Shielded HPO Sensor & Cable	
ISO-HPO-100	100 μm HPO Sensor*, pkg of 3	
ISO-HPO-100-L	100 μm HPO Sensor, L-shaped*, pkg of 3	
ISO-HPO-100H	100 μm HPO Sensor in hypodermic sheath*, pkg of 3	
600012	Replacement Sleeve Kit for ISO-HPO-2, pkg of 4	
100042	ISO-HPO-2 Electrolyte (10 mL)	
91580	Microsensor Adapter Cable	
	* Requires 91580 Microsensor Adapter Cable	

Oxygen Sensors

This sensor incorporates WPI's proprietary combination electrode technology whereby the oxygen-sensing element and separate reference electrode are encased within a single shielded sensor design. A gas-permeable polymer membrane is fitted over the end of the sleeve, which allows oxygen to pass while blocking liquids, ions and particulate matter.

Oxygen diffuses through the membrane. The voltage applied to the sensor is held at -0.7V when the monitoring device is on and



the sensor is properly connected. The magnitude of the generated electrical current is determined by the rate of diffusion through the membrane. The rate is proportional to the partial pressure of oxygen outside the membrane. The current serves as a measure of the partial pressure of O2.

The ISO-OXY-2 is a 2.0 mm stainless steel sensor, with replaceable membrane sleeves (#5378) and an internal refillable electrolyte (#7326). The sensor is similar in design to WPI's popular OXELP oxygen sensor (see page 68).

ISO-OXY-2	2 mm Shielded Oxygen Sensor & Cable	
5377	ISO-OXY Startup Kit (recommended with first purchase)	
5378	Replacement Electrode Sleeve Kit, pkg of 4	
7326	ISO2 Filling Solution (electrolyte)	

Temperature Sensor



The temperature sensor (#ISO-TEMP-2) is based on a 2.0 mm tip diameter high quality miniature platinum RTD (Resistance Temperature Detector) electrode. This design has been shown to provide greater accuracy, stability and interchangeability during temperature measurements than traditional thermistor and thermocouple sensors. The ISO-TEMP-2 is included with the purchase of a system.


ISO-TEMP-2	2 mm Platinum RTD Temperature Sensor (requires #91580)	
91580	Microsensor Adapter Cable	

ISO-H2S-2 ISO-H2S-100-Cxx


Although hydrogen sulfide (H₂S) is generally thought of as a poisonous gas, it is endogenously produced in many mammalian tissues. It has been detected in micromolar amounts in blood and brain tissue. Hydrogen sulfide is reported as having a broad range of biological functions and although its potential to participate in cell signaling is clear, this biological role is not well understood. H2S is strongly anagolous to nitric oxide (NO) because they share several physical and metabolic properties.

Like NO, H₂S is a potent vascular signal that can mediate vasoconstriction or vasorelaxation depending on the O2 level and tissue. In the rat aorta, H₂S concentrations that mediate rapid constriction at one O₂ level will cause rapid relaxation at lower O2 levels.

The ISO-H2S sensor is a low detection limit sensor to record H2S in vitro. It is the only sensor available that measures H₂S.



The stepped response to increasing concentrations of H2S are linear (see inset, R=.9963).



The sensor is insensitive to competing species such as ascorbic acid, dopamine, nitrate, nitrite, and glutathials.

ISO-H2S-2	2 mm Shielded Hydrogen Sulfide Macro Sensor	
ISO-H2S-100-Cxx	Hydrogen Sulfide Micro Sensor (pkg of 2)	
600016	Replacement Sleeves for ISO-H2S-2 (pkg of 4)	
100084	Filling Solution for ISO-H2S-2	
600015	Start-up Kit for ISO-H2S-2	



ISO-COP-2 SPECIFICATIONS

Sensor diameter	2mm
Sensor sensitivity	~0.5 pA/nM
Detection limit	~10nM
Linear range	10nM-10µM
Response time	<10 seconds

ISO-COP-2

Electrochemical CO Sensor for In Vivo Measurements for CO dissolved in solution

Carbon monoxide (CO) is a versatile mediator of physiological processes. Carbon monoxide (CO) formed by internal mechanisms (endogenous) is measured in a variety of ways, but standard measurement methods are of limited utility in most biological systems. WPI's ingenious ISO-COP-2 CO sensor measures CO in vivo or in vitro in real time!

This CO sensor is comprised of a 2.0mm stainless steel body with a replaceable membrane-covered sleeve. The sleeve is filled with electrolyte. It is an amperometric sensor designed for use in cell culture and similar applications.

In principle, CO diffuses through the gas-permeable membrane and

is then oxidized to CO₂ on the working electrode of the sensor. This oxidation creates a current with a magnitude directly related to the concentration of CO in solution.

It is designed for use with WPI's **TRB4100** (4-Channel) or **TBR1025** (1-Channel) Free Radical Analyzers.

References

Motterlini, M., Sawle, P., Bains, S., Hammad, J., Alberto, R., Foresti, R. and Green, C, "CORM-A1: A new pharmacologically active carbon monoxide-releasing molecule," FASEB Journal, November 19, 2004, express article 10.1096/fj.04-2169fje.

ISO-COP-2	2 mm Shielded Carbon Monoxide Sensor	
95620	Replacement Sleeves (pkg of 4)	
95611	Filling Solution	
95699	Start-up Kit	
	•	



IGS100 – Implantable glucose sensor

Measuring glucose in vivo over the long term is challenging and difficult. Previous measurement systems were limited to acute studies or a few days at best. WPI introduces a new kind of implantable glucose sensor based on a patented technology. This sensor provides a tool for researchers to directly detect glucose in chronic studies in vitro or in vivo. The sensor is fully compatible with WPI's Apollo system.

IGS100	Implantable Glucose Sensor (pkg of 2)
91580	Microsensor Adapter Cable

GLUCOSE SENSOR SPECIFICATIONS

IN VITRO PRECISION Coefficient of Variation (CV) ≤5% **GLUCOSE RANGE** 36 - 450 mg/dl (or 2-25 mM/L) RESPONSE TIME (sec) 100 - 300s

IN VIVO CALIBRATION In vivo calibration

Acetaminophen, ascorbic acid, uric acid INTERFERENCE SPECIES

LENGTH 5 cm SENSOR SIZE $0.6 \times 0.7 \text{ mm}$ REFERENCE ELECTRODE Ag/AgCl

POLARIZATION VOLTAGE (V) 0.65 -0.7V vs. Ag/AgCl

SENSOR LIFE 3-4 months in solutions at room temperature under continuous polarization; 15-30 days in vivo

SHELF LIFE 6 months

OPERATION CONDITIONS 20° to 40° C (68° to 104° F) STORAGE CONDITIONS 10° to 25° C (50° to 77° F)



Four-Channel Free Radical Analyzer

- Real-time detection using electrochemical microsensors.
- Integrated system includes one temperature sensor, your choice of two additional sensors, and a start-up kit.
- Measure up to four different species in the same preparation or simultaneous measurement in four different preparations.
- Current measurement range from 300 fA to 10 μA (four ranges) permits wide dynamic range for detection.
- Wide bandwidth allows recording of fast events.

- Measure nitric oxide from < 0.3 nM to 100 μM.
- Measure hydrogen peroxide < 10 nM to</p> 100 mM.
- Measure hydrogen sulfide
- Measure glucose
- Measure oxygen from 0.1% to 100%.
- Isolated architecture allows Lab-Trax interface to simultaneously measure free radical and independent analog data (i.e., ECG, BP, etc.) data on any channel.

Real-time detection and measurement of a variety of redox-reactive species is fast and easy using the electrochemical (amperometric) detection principle employed in the new TBR4100. This optically isolated four-channel free radical analyzer has ultra low noise and independently operated channels.


For use with WPI's wide range of nitric oxide, hydrogen peroxide, hydrogen sulfide and oxygen sensors, the TBR4100 can measure four different species simultaneously in the same preparation. Simply plug a sensor into any one of the input channels on the front panel and select the current range. Poise voltage can be selected from a range of values tuned for optimal response from WPI sensors. An independent output for real-time monitoring of temperature is also included.

The TBR4100 analyzer utilizes PC-based data acquisition via our Lab-Trax interface; data traces are displayed and recorded in real-time. The Data-Trax software comes preconfigured for single or multiple electrode recording; filters, gains, and smoothing are all set for optimal results. Data can be viewed making adjustments to smoothing and filter settings without affecting the original stored raw data. Electrode calibration from multiple concentration readings can be input into the software's Multipoint Calibration utility quickly provides a plot and slope calculation for electrode sensitivity determination. Alternately, the **Lab-Trax** data interface can be used for providing simultaneous acquisition of Free Radical data along with other physiological data (ECG, HR, BP, etc.) as each of the four input channels has its own independent input, filters, and 24-bit

See www.wpiinc.com/TBR4100 for more information on Lab-Trax data acquisition.



Don't need four channels? The single-channel TBR1025 packs the power of its big brother in a small, economical package.




Multipoint electrode calibration and slope determination can be quickly derived from recorded calibration data.

TBR4100	SPECIFICATIONS
Power	Analog Output
Operating Temperature (ambient)0 - 50°C (32 - 122°F)	Amperometric Input
Operating Humidity (ambient)15 – 70% RH non-condensing	Number of Amperometric Channels4
Warm up Time<5 minutes	Signal Bandwidth0-3 Hz
Dimensions	Polarization Voltage (selectable via rotary switch) Nitric Oxide
Analog Output Range	Current measurement Performance Range Analog Output Noise @ 3Hz * Noise @ 0.3 Hz * +/- 10 nA 1 mV / 1 pA < 1 pA < 0.3 pA +/- 100 nA 1 mV / 10pA < 7 pA < 3 pA +/- 1 μA 1 mV / 100pA < 70 pA < 30 pA +/- 10 μA 1 mV / 1μA < 700 pA < 300 pA *Instrument performance is measured as the (max-min) over 20 seconds period with open input. Typical values are given at 3 Hz and 0.3 Hz bandwidth. Typical sensor performance with TBR4100 ISO-NOPF100 noise

TBR4100-416	Four-Channel Free Radical Analyzer with Lab-Trax 4/16 Data Acquisition System
	Includes TBR4100 analyzer & power cord, Lab-Trax-4/16 data acquisition system & USB cable, 4 BNC cables,
	1 electrode adapter cable, 1 temperature probe, 2 sensors of your choice, and sensor start-up kit(s), if applicable.
TBR4100-424T	Four-Channel Free Radical Analyzer with Lab-Trax 4/24T Data Acquisition System
	Includes TBR4100 analyzer & power cord, Lab-Trax-4/24T data acquisition system & USB cable, 4 BNC cables,
	1 electrode adapter cable, 1 temperature probe, 2 sensors of your choice, and sensor start-up kit(s), if applicable
TBR1025	Single-Channel Free Radical Analyzer — <i>Includes 1 sensor of your choice & 1 start-up kit</i>

RECOMMENDED ACCESSORIES

SNAP50 SNAP S-Nitroso-N-acetyl-D-penicillamine, 50 mg vial



Micro C Advanced

Neurotransmitter **Detection**

Measures oxidizable compounds such as catecholamines (epinephrine, norepinephrine, dopamine), indolamines (serotonin, melatonin), ascorbic acid and Fe (II) with exquisite sensitivity, low noise and site specificity

1-picoampere to 2 microamperes. The built-in carbon electrode activation feature allows the easy renewal of electrode sensitivity. In addition, MicroC features a low-pass filter and the option of applying DC potential externally. A wide range of compounds can be detected: dopamine, epinephrine, norepinephrine, serotonin, etc. Other compounds, such as glutamate, glucose, acetylcholine and alcohol, can also be detected with MicroC using enzyme-modified biosensors.

See Application Notes, "Carbon Fiber Microelectrodes", available as a PDF file which may be downloaded from WPI"s web site.

The MicroC Potentiostat is supplied with a carbon electrode probe, with 5 feet triax cable, which accepts 0.79 mm connector pin, and a reference electrode with a 4 mm Ag/AgCl half cell (see page 121). For applications where smaller half cells are needed, please call WPI for more information.

MICROC SPECIFICATIONS

MicroC, WPI's low cost and elegant instrument for electrochemical

detection using carbon microelectrodes, will record the presence and

concentration of oxidizable biological compounds in vivo and in vitro.

It also features inherently low noise and a sensitivity of 1 millivolt per

cholamine release is less than 1 millisecond. When used with carbon

fiber microelectrodes, redox current can be recorded over a range of

picoampere of oxidation current. Response time to quantal cate-

METHOD 2 electrode, DC potentiostat APPLIED POTENTIAL 0.65 V, variable ± 2.5 V **CURRENT RANGES** 2000 pA, 20 nA, 200 nA, 2 μA BANDWIDTH 1.67 Hz, 167-1000 Hz

NOISE

DISPLAY 3½-digit LCD display, ±2 V

RECORDER OUTPUT + 4.5 volts RISE TIME < 1 millisecond ELECTRODE PROBE/CABLE LENGTH Triax shielded, 5 feet

POWER Six 1.5 volt alkaline batteries (included)

BATTERY LIFE > 1000 hours, est. SHIPPING WEIGHT 4 lb (1.8 kg)

References

G. A. Gerhardt, "Nafion-coated electrodes with high selectivity for CNS electrochemistry" Brain Research, 290: 390-395 (1984).

R. M. Wightman, et al., "Temporally resolved catecholamine spikes correspond to single vesicle release from individual chromaffin cells." Pro. Nat'l Acad. of Sci. 88: 10754-58, (1991).

Z. Zhou and S. Misler, "Action Potential-induced Quantal Secretion of Catecholamines from Rat Adrenal Chromaffin Cells", J. Biol. Chem. 270; 3498-3505, (1995).

SYS-MICROC Potentiostat Replacement Probe for MicroC **MICROCP OPTIONAL ACCESSORIES** 300305 ATP Adapter (0.031" pin to 1 mm socket)

(nack of 5)

Carbon Fiber Microelectrodes

Sensitive, renewable/durable and economical carbon fiber electrodes for electrochemical detection of oxidizable compounds

Carbon fiber microelectrodes have been used in both the detection of oxidizable compounds (Gonon, et al., 1978; Cahill and Wightman, 1995) and extracellular single-unit recording (Armstrong-James and Millar, 1979). WPI's ultra-sensitive and low-noise carbon fiber (CF) electrodes can be applied, with our Micro-C Potentiostat or similar instruments, in the electrochemical detection of catecholamines (epinephrine, norepinephrine and dopamine), indolamines (serotonin, 5-HT or melatonin), ascorbic acid, Fe (II), and other oxidizable compounds.

CF electrodes (diameter of 10 or 30 µm) respond with an excellent linearity to the oxidizable compounds (Fig. 1) and can detect the compounds as low as 0.2 nM. While the shorter (25-100 μm) CF electrodes are suitable for in vivo amperometric and voltammetric measurements, the longer CF electrodes provide higher sensitivity and are especially useful for the *in vitro* studies (amperometric or differential pulse voltammetry). When used with the Micro-C Potentiostat, these CF electrodes can be activated and renewed in sensitivity for multiple use. The selective detection of catecholamines can be achieved with our Nafion-coated CF electrodes. For selective detection of 5-HT and ascorbic acid, please contact WPI for more information




Fig. 1 — Excellent linearity in the response of carbon fiber electrode (CF30-500) to dopamine recorded on Micro-C. Courtesy: Drs. D. Yeomans and X.-T. Wang, University of Illinois at Chicago.




Fig. 2 — Extracellular recording using a carbon electrode in CA1 region of the hippocampus in an anesthetized rat shows ultra**low noise (-<5 μV-).** Courtesy: Dr. Carolyn Harley of Memorial University, Newfoundland, Canada.

References

P. S. Cahill, R. M. Whightman, Anal. Chem., 67, 2599-2605 (1995). F. Gonon, et al., Hebd Seances Acad. Sci. Ser. 286, 1203 (1978). M. Armstrong-James, J. Millar, J. Neurosci. Methods, 1, 279 (1979).


CARBON FIBER MICROELECTRODES, UNCOATED

	Diameter	Length	(pack of 5)
CF10-100	10 μm	100 μm	
CF10-250	10 µm	250 μm	
CF10-500 *	10 µm	500 µm	
CF30-50 *	30 µm	50 μm	
CF30-100	30 µm	100 μm	
CF30-500 *	30 µm	500 μm	
CF30-1000 *	30 µm	1000 μm	

CARBON FIBER MICROELECTRODES, NAFION-COATED Diameter Langth

	Diameter	Length	(pack of 3)
CFN10-50 *	10 μm	50 µm	
CFN10-100 *	10 µm	100 µm	
CFN10-250 *	10 µm	250 µm	
CFN30-50 *	30 µm	50 µm	
CFN30-100 *	30 µm	100 µm	
CFN30-250 *	30 µm	250 µm	
CFN30-500 *	30 µm	500 μm	
CFN30-1000 *	30 µm	1000 µm	

* Built to order — allow up to 4 weeks manufacturing time.



FLOX is a device for measuring fluorescence lifetime, phase and intensity. It uses LED excitation and photodiode detection with filterbased wavelength selection for easy experimental set-up and control. Because the unit is self-contained, it is invariant to fiber bending and stray light, and has a wide dynamic range of optical intensity as well as low optical and electronic crosstalk, and low drift and phase noise. **FLOX** is especially useful for oxygen sensing applications where stability and sensitivity to drift is important and where sample set-ups must be left undisturbed for long periods of time.

- Fluorescence-sensing detector for optical sensors, a viable alternative to traditional chemical sensing devices
- Self-contained, benchtop system-invariant to fiber bending or
- Affordable-Half the cost of comparable phase measurement
- Excellent stability, extremely low drift and phase noise
- Simple calibration, setup and control

The new oxygen sensing system measures fluorescence lifetime, phase and intensity, using LED excitation and photodiode detection with filter-based wavelength selection. The system is simple to set-up and control. The compact, self-contained unit makes it invariant to fiber bending and stray light. It also has a wide dynamic range of optical intensity, as well as low optical and electronic crosstalk, and low drift and phase noise.

When stability and sensitivity to drift are important in your oxygen sensing experiment, this unit is ideal. It is perfect for applications where sample set-ups must be left undisturbed for long periods of time.

Three coatings available for probes and patches

The optical sensors consist of transducer materials, applied to the tips of optical fibers or to substrates such as patches or cuvettes, which change optical properties in response to specific analytes in their immediate environment.

OXY-The standard oxygen sensor designed for monitoring oxygen partial pressure in gas and aqueous solutions is a fiber optic fluorescence probe with a proprietary oxygen sensing coated tip.

HIOXY–Designed for monitoring oxygen partial pressure in non-aqueous vapors and solutions. The sensor coating chemistry is compatible with oils, alcohols, and hydrocarbon-based vapors and liquids.

FOSPOR-A new generation of highly sensitive sensor coating can be used for monitoring traces of oxygen in gas and liquids.

Fluorescence-

based optical

sensor system

The oxygen sensor probes are low-power and offer high sensitivity, reversibility and stability, ideal for remote monitoring. The thin coating on the probe tips consumes no oxygen, allowing for continuous contact with the sample. They are ideal for viscous samples and are immune to interference caused by pH, ionic strength or salinity fluctuations or biofouling.

AL300

A 500 µm OD (300 µm core diameter) aluminum-jacketed optical fiber probe for applications that require fine spatial resolution. 300 µm alµminµm-jacketed fiber assembly; 500 µm OD, 1 m length

The OR125 is a 1/8" OD optical fiber probe

100 µm optical fiber, stainless steel ferrule;

OR125

used as a direct replacement for 1/8" OD O2 electrodes.

3.175 mm OD, 63.5 mm length

INTRODUCTORY PRICE

FLOX-PATCH Non-Invasive Oxygen Monitoring Kit, including phase measurement system, temperature probe; select sensor patches when ordering

R Sensor Probe

OR125G and OR125GT

3.175 mm OD, 63.5 mm length

1.587 mm OD, 152.4 mm length

electrodes.

These probes are 1/8" OD optical fiber probes

used as direct replacement for 1/8" OD O2

1000 µm optical fiber, stainless steel ferrule;

1000 µm fiber in a stainless steel 1/16" ferrule;

FLOX-PROBE In Situ Oxygen Monitoring Kit, including phase measurement system, temperature probe; select sensor when ordering



The PI600 is a silicone-jacketed, polyimidecoated optical fiber probe used in environments where a non-metallic probe is

200 or 600 µm optical fiber with silicone jacketing; 710 µm OD, 2 m length

Dissolved oxygen meter & oxygen electrode

Electrically isolated, the electrode can function in recording situations where ground-referred electrodes would fail



- Measures gas and dissolved oxygen
- Low O₂ consumption
- 2 mm tip

Isolation prevents adverse interaction with other electrodes and instruments. Highly accurate and stable ISO2 measures oxygen concentrations in aqueous solutions and in gas mixtures. Measurement modes are percent oxygen, parts per million, and oxygen reduction current in nanoamperes.

The lower detection limit (DL) for ISO2 (in gaseous or in liquid phase) is 0.1 ppm or 0.1%. Oxygen concentration around 0.5 ppm or 0.5% and up can be routinely measured in the gaseous or in the liquid phases.

DL for a particular sensor tends to be different in gaseous and liquid phases when the baseline noise level is also different. In the case of ISO2, however, the baseline noise level is very close whether the electrode is immersed in solution or used in the gaseous phase, so DL for the ISO2 remains relatively close in either phase.

The small tip size (2 mm diameter) and low oxygen consumption of the **OXELP** electrode make it ideal for measurements in vivo or in vitro. With a T-Adapter (#5399), the sensor probe can be used also for continuous-flow monitoring of oxygen in small fluid volumes. OXELP also features a fast response time, typically 10 seconds. Optional BNC-to-double banana adaptor (#13347) and BNC cable (#500184) allow ISO₂ to be connected directly to your chart recorder.

SYS-ISO2 Dissolved Oxygen Meter & Electrode Replacement Oxygen Electrode for ISO2 5378 Replacement Electrode Sleeve Kit (pkg of 4) Four sleeves with membranes, plus 10 mL refill solution. 7326 ISO2 Filling Solution (10 mL) 5377 Replacement ISO2 Start-up Kit Includes Calibration Bottle, 10mL Refill Solution, 1 cc Syringe, 2 Replacement Membranes Sleeves, MicroFil (28 ga.) 5399 T-Adapter Flow-Through Kit Includes 3 female luer T's, 3 luer lock fittings, 3 2mm gaskets, 6

male luer to 1/8-in. tubing, 3 luer lock fittings

13347 Chart Recorder Adapter (requires BNC cable) 500184 **BNC Cable**



Also see WPI's temperature stabilzed chamber, page 65.

ISO₂ SPECIFICATIONS

MEASUREMENT MODES

% O2: 0-100% 0-20 ppm: 0-200 nA Current: **RESOLUTION** 0.1 ppm

ACCURACY ± 1.5%

RECORDER OUTPUT 1000 Ω resistance for chart recorders

DISPLAY

POWFR 2 nine-volt alkaline batteries, supplied

BATTERY LIFE 1000 hours, estimated DIMENSIONS

8 x 4 x 2 inches (20 x 10 x 5 cm)

5 lb (2.3 kg) SHIPPING WEIGHT

OXELP OXYGEN ELECTRODE

2 mm TIP DIAMETER STAINLESS STEEL SLEEVE 66 mm OVERALL LENGTH

CABLE LENGTH

122 cm (4 feet), including BNC connec-

RESPONSE TIME 10 seconds, 90% response typical

in well-stirred solution

DRIFT

OxyMini & OxyMicro

A new generation of fiber optic oxygen sensors based on luminescence lifetime

The OxyMini system is optimized for process control and biotechnology applications. The OxyMicro is designed for biological research applications including implanting into tissues, cell cultures, profiling of biofilms and sediment related bioassays. The measurement principle of the sensor system is based on the detection of oxygen concentration as a function of luminescence lifetime either in dissolved or gaseous phase environments.

OxyMini and OxyMicro Benefits

- Oxygen is not consumed during the experiment
- Immune to electrical and magnetic interference
- Excellent long-term stability
- No lengthy polarization necessary (e.g., as Clark-type oxygen electrodes require)
- Fast response time < 0.5 s for MicroTip sensors
- Probe size of MicroTip sensors as small as 50 μm
- Measurement is feasible in dry gas
- Optical isolation of sensor tip available for fluorescent or photosynthetically active samples

Measurement Principle

Conventional fiber-optic oxygen sensor systems based on intensity measurements are limited in their accuracy by light source stability and ambient light fluctuations. Using a luminescence lifetime detection, measurements are not affected by light source stability, intensity fluctuations caused by fiber bending or changes of the optical properties of the sample (turbidity, refractive index, coloration, etc.). These advantages make WPI's OxyMini and OxyMicro the most advanced and reliable fiber-optic oxygen system available.

Calibration: The sensors can be calibrated by a simple two point calibration, 100% air-saturation and 0% air saturation.

OxyMini and OxyMicro oxygen meters: The OxiMini and OxyMicro fiber optic oxygen meters are compact, easy to transport. The instruments are designed for in/outdoor use and can be connected to a PC via a RS232 interface. Data can be visualized, analyzed and stored with the supplied software. A full range of sensors covering most biomedical applications are available.

OXY MICRO

OxyMicro systems

The OxyMicro is a single channel fiber optic oxygen meter for WPI's fiber optic oxygen microsensors. Applications include:

- Oxygen profiles of marine sediment, soils, or tissue
- Implantation into living tissue (e.g., heart or muscle tissue)
- Control of cell culture media in Biotechnology.

MicroTip

The MicroTip (WPI **#501656**) is a needle-type (27 ga.) oxygen micro sensor designed for applications where a small tip size of 50 μ m and a fast response time (t₉₀) of 1 s are necessary. The oxygen sensitive sensor tip consists of 140 μ m fiber

tapered to a 50 μ m tip. The sensor is housed inside a stainless steel needle of 22 mm length and 0.4 mm diameter. This allows penetration through a septum rubber or similar material. These sensors are ideal for oxygen profiling in sediments and biofilms.



MicroFlow

The MicroFlow fiber optic oxygen sensor (WPI #501657) is a miniaturized fiber optic



chemical sensor optimized for fast response time (t₉₀ < 1 sec in gases, < 5 sec in liquids). The tiny probe has a tip size of 50 µm and is integrated in a T-shape flow cell for easy connection via Luer-Lock adapters to external tubings. Liquids (like water, blood, etc.) can be pumped through the cell.

MicroImplant

The MicroImplant fiber optic oxygen sensor (WPI **#501658**) is an implantable probe (IMP) with a tiny probe tip size $50 \mu m$, an exposed fiber length of 5-mm and a jacket diameter of $900 \mu m$. The IMP sensor was successfully implanted in crabs, fishes and soil.

OxyMini systems

The OxyMini is a single-channel fiber optic oxygen meter for WPI's fiber optic oxygen minisensors. These sensors are based on 2 mm polymer optical fibers and have a length of 2.5 m. A wide range of applications is possible with these sensors.

- Process control: bottling plant in breweries and quality control of packages
- Biotechnology: Control of cell culture media and non-invasive control of bioreactors
- Implantation of oxygen sensors into soil and trees.

MiniTip

This oxygen dipping probe (WPI **#501641**) has a tip diameter of 4 mm and consists of a polymer optical fiber, with an oxygen sensitive coating. The MiniTip's range is 0 to 100%. This robust sensor has a response time (t₉₀) of approximately 40 s.

MiniFlow

The MiniFlow oxygen probe (WPI **#501642**) is a miniaturized fiber optic chemical sensor integrated in a T-shape flow through cell. The standard T-shape



flow cell can be easily connected via Luer-Lock adapters to external tubings. Liquids (e.g., water, blood, etc.) can be pumped through the cell. The sensor has a response time (t₉₀) of approximately 40-s and an excellent long-term stability.

MiniFoil

WPI offers the sensor material on a 1 cm² support disk made of polyester. This material can be glued, for example, inside glass vials and the oxygen concentration can be measured noninvasively and non-destructively from outside through the wall of the flask. A plastic fiber optic cable (WPI #501644, WPI #501645)



is used to illuminate the sensor foil. The wall of the flask must be transparent/ non-fluorescent. Response time (t_{90}) of approximately 50 s. The material can be implanted into animal tissues or custom-made housings.

MINISENSOR SYSTEM

AOT Fiber-optic Oxygen Meter for Minisensors *
ORS (not interchangeable with Microsensors)
MiniTip, fiber-optic oxygen sensor
MiniFlow, flow-through cell with integrated planar oxygen sensor
MiniSpot, planar oxygen-sensitive spot, 5 mm diam. (includes 10)
Requires #501644
Polymer optical fiber with 1 SMA connector

MICROSENSOR SYSTEM

OXY-MICRO	-AOT Fiber-optic Oxygen Meter for Microsensors *
MICROSENS	ORS (not interchangeable with Microsensors)
501656	MicroTip, needle-type housing fiber-optic oxygen sensor, 50 μm tip
501656-C	MicroTip, needle-type housing, 50 µm tip, optical isolation
501656-F	MicroTip, needle-type housing, 140 μm flat tip
501657	MicroFlow, flow-through housed oxygen microsensor
501658	MicroImplant, implantable oxygen microsensor, 50 µm tip
501658-F	MicroImplant, 140 μm flat tip
*	Meter contains two analog outputs and one trigger input

	MiniTip	MiniFlow	MiniSpot	MicroTip	MicroFlow	MicroImplant
Measurement Range dissolved/gaseous	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg
Response Time [t ₉₀] dissolved/gaseous	40 s 10 s	40 s 10 s	40 s 10 s	< 2 s < 0.5 s	< 2 s < 0.5 s	< 2 s < 0.5 s
Sterilization (EtOH, H ₂ O ₂) autoclavable (130°C, 1.5 atm)	Y N	Y Y	Y Y	Y N	Y Y	Y Y
Drift (100,000 datapoints, 20°C)	< 0.1%	< 0.1%	< 0.1%	< 0.3%	< 0.3%	< 0.3%
Accuracy (20°C)	0.2%					
Resolution (20°C)	2.75 ± 0.01 ppm, 9.00 ± 0.05 ppm, 220 ± 0.15 ppm, 45.0 ± 0.25 mmHg, 150 ± 0.75 mmHg, 375 ± 2.6 mmHg					
Temperature Range	-10°C to 50°C					
Probe Assembly Length	World Precision Instruments www.wpiinc.com					

Optical sensor

ment or bag

with transparent window

Optical fiber

pH*Optica*™

A novel fiber optic pH system

pH*Optica*[™] is a pH measuring system which uses fiber optic sensors and patented DLR technology. This method allows referenced measurements with single excitation to be implemented.

Features of pHOptica Meter

- Single-channel, compact, easy to transport fiber-optic meter for pH measurements with miniature sensors.
- Two 12-bit, programmable analog outputs, with electrical isolation.
- One external trigger input, with electrical isolation.
- Computer with RS232 interface required for operation.
- User-friendly software saves and visualizes measured values.
- Several **pH***Optica* meters can be connected to one computer.
- Temperature variation is recorded using a temperature sensor.
- pH optical sensors No reference electrode is needed
- Immune to electrical interferences and magnetic fields
- Low drift

Features of

PHOptica

- High spatial resolution due to small tip size
- Measurement in very small sample volumes
- Additional optical isolation of the sensor tip is available for measurements in colored or photosynthetically active samples.

pHOptica Micro system

The pH Optica micro system is a single channel pH system for use with fiber optic micro sensors. The applications include:

- Penetration or implantation into living tissue (heart, muscle or animal blood vessels).
- Soil implantation for pH measurement.
- Implantation into customer-made housing.




Needle-Type Housing Sensor the glass-fiber with its pH-sensitive tip is protected inside a stainless steel needle (18 ga.); fiber has to be extended during measurement; penetration through septum.

To protect the small glass fiber tip against breaking, suitable housings and tubings around it, depending on the respective application, were designed.



pHOptica Mini system

The pH Optica mini system is a single channel pH system for use with fiber optic mini sensors, foil and spot surface sensors. The applications include:





Dipping probe pH measurements

pH mini sensors

● OD of the dipping sensor is 4 mm. ● Sterilization of the pH sensor spots via gamma radiation. • The pH mini sensor meter is based on 2 mm PMMA waveguides. • Drift of 0.1 pH units for 10,000 measurements (4 days measurement in the 30 sec data update mode).

Two different housings and sensor spots (sensorfoils) are offered:



POF Coated with a pH-Sensitive Foil—Small and robust pH dipping sensor: no reference electrode needed.



Flow-Through Cell with Integrated pH Sensor—On-line monitoring; can be easily connected via Luer-Lock adapters.



Planar pH Sensitive Foils and spots—non-invasive and non-destructive measurement from outside through the wall of the flask; online monitoring.

pH micro sensors

- Tip size 140 micrometer.
- Drift of 0.1 pH units for 2000 measurements (16 hours measurement in the 30 sec data update mode).

PHOPTICA SPECIFICATIONS

DATA INTERFACE SAMPLE RATE MEASURING pH RANGE RESOLUTION (at 20 °C)

RESPONSE TIME

RS232 1 sample per second

< 1 min

± 0.03 (microsensors); ± 0.01 (minisensors)

DIMENSIONS WEIGHT POWER SUPPLY

185 x 110 x 45 mm 630 g 100 - 220 V AC

MINISENSOR SYSTEM (cannot be used with microsensors)

PH-OPTICA-MINI Fiber Optic pH Meter for minisensors, foils and spots

503538* pH MiniTip, fiber optic pH sensor dipping probe, disposable (4 mm OD), pkg of 3

502120* pH MiniFlow, fiber optic pH flow sensor, pkg of 3

*Requires #503110

502122** pH MiniSpot, fiber optic pH spot sensors, pkg of 10, OD 5 mm **Requires #501644

501644 Polymer Optical Fiber with 1 SMA connector **503110** Fiber Optic Cable with 1 SMA connector

MICROSENSOR SYSTEM (cannot be used with minisensors)

PH-OPTICA-MICRO Fiber optic pH meter for microsensors

502123 pH MicroTip (needle-type), fiber optic pH sensors (140 μm OD), pkg of 3 **502124** pH MicroImplant, fiber optic pH implantable sensor (140 µm OD), pkg of 3

Beetrode® Micro pH Electrodes

with 100-micron sensing tips!

- NEW—IMPROVED—super-miniature, superfast coated wire pH electrodes
- 100-micron diameter ideal for monitoring fast pH changes in very small places



Beetrode® is a solid state pH sensor with ideal characteristics over a wide pH range. Exhibits a larger E₀ than conventional glass electrodes. Requires a separate reference electrode, such as WPI's Dri-Ref Series. Beetrodes (except the NMPH2) connect to your pH meter via a BNC-terminated cable. Beetrodes generate mV readings on standard pH meters. To obtain pH-scale readings on standard pH meters, use **BEE-CAL™**, a small, battery-operated compensator (AA battery included) that adjusts the electrode offset potential so that Beetrodes will produce standard pH-scale readings. (NMPH2B requires ZBEECAL.)

BEETRODE SPECIFICATIONS

TIP DIAMETER TIP I FNGTH

100 μ (0.1 mm) ~2 mm, except

~2.5 mm on NMPH3L

~5 mm on NMPH5 0.187 in. (5 mm)

BODY DIAMETER (NMPH 1, 3, 3L, 5) BODY LENGTH (NMPH 1, 3, 5) BODY LENGTH (NMPH 3L)

1.875 in. (48 mm) 0.75 in. (19 mm) 1 s (90%) typical

RESPONSE TIME pH RANGE SI OPF

0 - 14Nernstian 100 k Ω (max)

RESISTANCE **SELECTIVITY**

No significant interference by K⁺, Na+, Ca++ in 0.1 to 1 M solutions

DRIFT

< 2.5 mV / 5 min.

BEE-CAL BEE CAL	ZBEECAL
	SEE-CAL-
1358	3508
	Q.

NMPH1*	Beetrode—2 mm receptacle
NMPH2	Beetrode with 1 m cable (unterminated)
NMPH2B	Beetrode w/ BNC cable, 1 m cable
NMPH3*	Dental Beetrode, 45° Bend, 2 mm receptacle
NMPH3L*	Dental Beetrode w/ 2 mm loop, 2 mm receptacle
NMPH5*	Beetrode, 2 mm receptacle
SYS-BEECAL	Beetrode Offset & Cable (BNC, 2 mm pin)
ZBEECAL	Bee-Cal Level Shifting Device
3508	BNC-to-US Standard Adapter
1358	BNC-to-2 mm Pin Adapter

*Requires BEE-CAL or 1358 BNC-to-2 mm Pin Cable (4-ft) for connection to a pH meter.

Dri-Ref™ Reference Electrodes



Dri-Ref™ reference electrodes were developed by WPI to have extremely low electrolyte leakage properties, hence the name "Dri-Ref". In addition to this key feature, these electrodes exhibit stable and reproducible potential and low resistance. Stored in KCl when not in use, they have a long life expectancy.

Although the internal filling solution contains KCI, the low fluid leakage means Dri-Ref may be used in combination with ion selective electrodes, including those for K⁺ and Cl⁻, without significant contamination from the reference electrode.

The Dri-Ref electrodes are chemically resistant to strong acids and alkalines. Dri-Ref electrodes are not suitable for use in organic solvents. In addition, the long, thin FLEXREF may be easily manipulated to accommodate a difficult experimental setup.

SUPER-Dri-Ref

With a diameter of 2-mm, SUPER-Dri-Ref does not leak electrolyte at all. Exhibiting the electrical stability of a classic flowing junction reference cell, this electrode exhibits low resistance and a stable half-cell potential essentially independent of sample electrolyte concentration. SUPER-Dri-Ref is ideal for small volume and low salt concentration measurements.

Micro-Reference Electrode

Only 450 µm in diameter and 1 inch long, WPI's new **DRIREF-450** reference electrode can be used along with other sensors in space-restricted areas and very small sample volumes.

Luer-Tip Reference

The male luer fitting at the front of the DRIREF-L allows it to be easily connected to a female luer port (see WPI's luer fittings kit, page 138) to form a tight seal — a very convenient installation for a flow-through system.

FLEXREF	Flexible Dri-Ref, 1.5 mm diam.
DRIREF-2	Dri-Ref, 2 mm diam.
DRIREF-2SH	Dri-Ref, 2 mm diam. (Short)
DRIREF-5	Dri-Ref, 4.7 mm diam.
DRIREF-5SH	Dri-Ref, 4.7 mm diam. (Short)
SDR2	SUPER-Dri-Ref, 2 mm diam.
DRIREF-450	Micro-Dri-Ref, 450 μm diam.
DRIREF-L	Reference Electrode with Luer Tip

			DRI-REF SI	PECIFICATION	ONS			
	DRIREF-450	DRIREF-5	DRIREF-2	FLEXREF	SDR2	DRIREF-L	DRIREF-5SH	DRIREF-2SH
LENGTH	2.54 cm	9 cm	13 cm	13 cm	9 cm	7.5 cm	3.5 cm	2 cm
DIAMETER	450 µm	4.7 mm	2 mm	1.5 mm	2 mm	Standard Luer	4.7 mm	2 mm
CONSTRUCTION	Coated Glass	Epoxy	Isoplast™	Teflon™	PVC	Polypropylene	Epoxy	Isoplast™
LEAD LENGTH	30 in. (76 cm)	30 in. (76 cm)	30 in. (76 cm)	30 in. (76 cm)	30 in. (76 cm)	30 in. (76 cm)	30 in. (76 cm)	30 in. (76 cm)
CONNECTOR	2 mm pin	2 mm pin	2 mm pin	2 mm pin	2 mm pin	2 mm pin	2 mm pin	2 mm pin
RESISTANCE (typical)	< 5 K Ω	~500 Ω	~2.7 KΩ	~2.7 KΩ	< 5 KΩ	~500 Ω	~500 Ω	~2.7 KΩ
FILLING SOLUTION	KCI	KCI	KCI	KCI	KCl	KCl	KCl	KCI
ELECTROLYTE LEAKAGE (mL/hr)	_	~7.4×10 ⁻⁷	~5.7×10 ⁻⁸	~5.7×10 ⁻⁸	_	~7.4×10 ⁻⁷	~7.4×10 ⁻⁷	~5.7×10 ⁻⁸
Isoplast is a trade mark of Dow Ch	hemical Teflon is a	trade mark of DuPor	nt					

Calcium Calibration Solutions

CALBUF-1

For use with calcium electrodes

A set of eight calcium buffers covering the range of concentration from 10⁻¹ to 10⁻⁸ M Ca⁺⁺. Each buffer contains 20 mL of solution and enough potassium chloride to set the ionic strength

to 0.1 M. Limited shelf life; use within 30 days. Concentration: 1x10⁻¹, 1x10⁻², 1x10⁻³, 1x10⁻⁴, 1x10⁻⁵, 1x10⁻⁶, 1x10⁻⁷, 1x10-8 M at 20°C. Limited shelf life; use within 30 days.

CALBUF-1 Kit of 8 Calcium Buffer Solutions



CALBUF-2

For use with calcium fluorescent indicators

CALBUF-2 is especially suitable for calibrating fluorescent Ca⁺⁺ indicators. It provides eleven buffer standards in the 10-4 to 10-8 M Ca++ range, whereas other commonly used fluorescent Ca⁺⁺ indicators have the apparent K_d in the range of 100 to 300 nM. As with any ionic sensitive indicator, the sensitivity range of these indicators is about 1.0 log unit above and below the K_d. CALBUF-2 provides seven calibration points in this sensitivity range. It has an osmolarity of 0.305, which is isotonic with most mammalian cells.

Concentration: 1x10⁻⁸, 4x10⁻⁸, 1x10⁻⁷, 2.5x10⁻⁷, 5x10⁻⁷, 7.5x10⁻⁷, 1x10⁻⁶, 4x10⁻⁶, 1x10⁻⁵, 4x10⁻⁵, and 1x10⁻⁴ M at 20°C. Ionic strength: 0.150 M. 11 bottles, 20 mL each. Limited shelf life; use within 30 days.

CALBUF-2 Kit of 11 Calcium Buffer Solutions

Kwik-Tip"

Ion selective electrodes

- Superior, stable PVC membrane
- Fast response
- 2mm diameter tips
- Interchangeable tip holder

These highly stable electrodes accurately measure calcium, potassium, hydrogen and TPP ion activity. Tips consist of 2 mm diameter plastic tubes sealed at one end with an ion-sensitive

membrane. After filling with electrolyte solution, the user inserts the tube into the holder and connects it to a pH meter. Tips and holders are interchangeable, so one tip may be replaced with another sensitive to a different ion. Replacing a tip takes less than a minute. Electrode tips normally last several months, when stored properly in saline solution. When replacement is necessary, only the tip need be replaced.

Kwik-Tip electrodes are available separately and as kits. Each "KWIK" Electrode Holder kit includes a reusable holder and three removable tips. In addition to a 4-foot BNC cable and an electrolyte filling syringe; "TIP" Electrode Kits contain three electrode tips for a specific ion. A

separate reference electrode, such as WPI's Dri-Ref™, is also required.

If your pH meter requires a US Standard connector, also order Part #3508 (BNC-to-US pH Standard adapter).



COMPLETE	KITS

KWIKCAL-2	Holder & 3 Calcium Electrodes
KWIKH-2	Holder & 3 Hydrogen Electrodes
KWIKPOT-2	Holder & 3 Potassium Electrodes
KWIKTPP-2	Holder & 3 TPP (Tetraphenylphosphonium) Electrodes

HOLDERS AND REPLACEMENT TIPS

KWIK-2	Electrode Holder with BNC cable
TIPCA	Calcium Electrode Tips (3)
TIPH	Hydrogen Electrode Tips (3)
TIPK	Potassium Electrode Tips (3)
TIPTPP	TPP+ (Tetraphenylphosphonium) Electrode Tips (3)

ACCESSORIES

BNC-to-US pH Adapter 3508

Also see Dri-Ref Reference Electrodes

TIP ELECTRODE SPECIFICATIONS

Part No.	Electrode	Color Code	Recommended Filling Solution	Min. Slope / Decade	Concentration Range	Selectivity Coefficients (-log)
TIPCA	Calcium	Green	0.1 M CaCl ₂	28 mV	0.1 M - 10 ^{-6.75} M	Na+ 5.5, K+ 5.4, Mg++ 4.9
TIPH	Hydrogen	Orange	1 M Citric Acid, 0.01 M NaCl, pH 5.6	54 mV	pH 5.0 - 12	Na ⁺ 10.4, K ⁺ 9.8, Ca ⁺⁺ 11.1
TIPK	Potassium	Yellow	0.1 M KCl	54 mV	0.1 M - 10 ^{-4.5} M	Na ⁺ 4.0, Ca ⁺⁺ 3.9, Mg ⁺⁺ 3.0
TIPTPP	TPP+	Purple	10 mM TPP+	54 mV	0.001 M - 10 ⁻⁴ M	K+ 6.0

Liquid Ion Exchangers

FFICIENTS 12.7	1.97	Ca ⁺⁺ IE 200 5.5 4.9
E 010 FFICIENTS 12.7	IE 190 S* 1.97	IE 200 5.5
FFICIENTS 12.7	5* 1.97	5.5
12.7	1.97	
_	2.95	4.9
_	_	5.4
_	2.7	_
2-10	4-10	4-10
56 mV	58 mV	28 mV
oH 4-12	pK 0-3	pCa 1-7
_	Corning	ETH1001
	477317	
	56 mV bH 4-12	2-10 4-10 56 mV 58 mV DH 4-12 pK 0-3 — Corning



When used in micropipettes to record cellular ion concentrations. consider using WPI's Duo 773 electrometer (channel A).

IE010	Hydrogen Ion Exchanger (0.1 mL)
IE190	Potassium Ion Exchanger (1.0 mL)
IE200	Calcium Neutral Ion Exchanger (0.1 mL)

Multi-Port Measurement Chamber



Temperature stabilized four-port closed chamber for measurements of nitric oxide, oxygen hydrogen peroxide and other species in cell culture

- Four port (WPI #NOCHM-4) chamber accommodates WPI's 2-mm sensors for nitric oxide (ISO-NOP), oxygen (OXELP), hydrogen peroxide, and WPI's KWIK-TIP ion selective electrodes in combination with WPI's 2 mm **Dri-Ref™** reference electrodes.
- Two additional top ports for injection of reagents using WPI's MicroFil™ syringe
- Closed chamber design greatly reduces the surface area of the solution exposed to air
- One top port and up to three side ports configuration provides adequate space for convenient sample and electrode manipulation
- Temperature control through an external circulating bath
- The chamber can be used for nitric oxide and other species calibration at temperatures from 4-40 °C

NOCHM-4 Four-Port Closed Chamber, for use with WPI's 2.0 mm electrodes (e.g., ISO-NOP and OXELP, etc.)


Spare Plug-adapter for ISO-NOP nitric oxide electrode NOCHM-P

800100-5 Spare Center Chamber Gasket (package of 5)

 $C \in$

Pre-polarizer

Keep extra NO sensors ready to use



Achieve a stable background current quickly. This small battery-powered device applies a potential to the NO electrode equivalent to the potential applied by the ISO-NO meter. Consequently, a sensor which has been connected to the activator may be transferred to the meter for immediate use. For use with all WPI NO electrodes.

ISO-NO Activator NSA-3

ISO-NOP Rejuvenator

After an ISO-NOP 2-mm sensor is used for long periods, sensitivity may become reduced and response time may increase. This little device can restore ISO-NOP performance to original levels by applying an electric waveform for a few seconds. 9v alkaline battery included.

JUV ISO-NOP Rejuvenator

GSNO

S-nitrosoglutathione

GSNO has been identified in vivo as a potential storage and transport vehicle for NO in the body. GSNO has been used in clinical trials to treat a form of preeclampsia and to prevent platelet aggregation. It also has considerable potential as NO donor in medicine.

M.W.336.3 • $C_{10}H_{16}N_4O_7S$ • Purity > 98% • Soluble in water or DMSO · Storage: -20°C

GSNO-50 GSNO 50 mg vial GSNO-100 GSNO 100 mg vial

S-Nitroso-N-acetyl-Dpenicillamine

SNAP is a stable green crystalline S-nitrosothiol compound that mimics the



action of nitric oxide in vivo. It has vasodilatory properties and has been shown to relax isolated bovine

coronary artery rings by activating soluble granulate cyclase. This reagent also actuates apoptosis in mouse thymocytes and has been accounted for reversible inactivation of protein Kinase C. SNAP can be used for calibration of all WPI NO sensors.

M.W. 220.2 • Purity > 98% by NMR or TLC

SNAP25	SNAP, 25 mg vial
SNAP50	SNAP, 50 mg vial
SNAP100	SNAP, 100 mg vial

Dual Channel Differential Electrometer



• High input impedance (10 15 Ω) • Differential (A-B) output • Low noise and wide bandwidth Electrode resistance test circuitry
 Probe test circuitry
 Driven guard shield

The **FD223a** is a dual channel differential, high impedance amplifier/ electrometer designed specifically for electrochemical measurements using ion specific (K⁺, Na⁺, C1⁻, etc.) or pH glass microelectrodes.

The instrument is very stable, drift free, and features a built in provision for measuring and adjusting input leakage current. DC levels may be independently adjusted for each probe channel.

The ability to locate the sensing probes directly at the measurement site overcomes the noise introduced by the long cables usually needed to bring the measured potential to the instrument. Signal-driven guards at the probe input maintains the specified high resistance and reduces the stray capacitance of the probes.

Careful design, coupled with quality component selection, particularly in the headstage, results in an excellent amplifier with low noise and wide bandwidth. The FD223a will faithfully reproduce the measured signal.

To reduce the noise and stray capacity even farther the probe housing includes a signal driven guard. A portion of this inner driven shell is exposed at the probe tip allowing a spring shield to be extended over the electrode holder and microelectrode.

The amplifier features a probe test port that permits testing of the electrode test feature and setting of the probe leakage current, (IG). A standby mode is included and should be used when attaching glass microelectrodes or electrode holders to the probe input. While in the standby mode the voltage at the probe input is clamped near zero volts thus protecting the input.

#2547 Driven Guard Shield

FD223A SPECIFICATIONS

INPUT IMPEDANCE $> 10^{15} \Omega$, shunted by 0.5 pF

INPUT CAPACITANCE 1 pF, nominal LEAKAGE CURRENT 75 fA max **GAIN** $1.000 \pm 0.1\%$

 50Ω **OUTPUT RESISTANCE** INPUT SWING VOLTAGE ±10 V

RISE TIME (10 TO 90%) 5 µs, small signal

NOISE (0.1 HZ TO 10 KHZ) <100 μV p-p, input shorted

BASELINE STABILITY ±0.1 mV/day POSITION CONTROLS RANGE ±600 mV

PHYSICAL DIMENSIONS Case: 8.8 x 21.0 x 17.5 cm (H x W x D)

Probe: 12.7 x 65 mm (D x L), 1.8 m cable

POWFR 90-265 VAC, 50/60 Hz, 10 VA

PROBE HANDLE 6.5 x 65 mm (D x L)

SHIPPING WEIGHT 2.5 kg

OPERATING CONDITIONS Equipment is intended to be operated in a controlled laboratory environment. Temperature: 0-40 °C; altitude: sea level to 2000 m; relative humidity: 0-95%.

FD223a Dual Channel Differential Electrometer FD223A

2 probes, driven guard shields and micropipette holder MEH1SF included for all glass microelectrodes O.D. 1.0 mm, 1.2 mm, 1.5 mm, or 2.0 mm.

OPTIONAL ACCESSORIES

M3301L	Micromanipulator (specify left- or right-handed)
M-3	80° Tilting base
RC1T	Reference cell (Ag/AgCl)
2547	Driven guard shield for FD223AP Probe
MEH1SF	Microelectrode holder
FD223AP	Replacement probe (includes calibration)

See cables and connectors, page 98 See microelectrode holders, page 106 See capillary glass, page 110

WPI's Quietest Intracellular Amplifier!



amplifier well-suited for the student lab Photo shows two units arranged for differential recording. Manipulators not included.

- Remote Headstage Easily mounted in any manipulator, this small probe, containing the first stage of amplification, includes a microelectrode holder, which plugs directly into the probe input.
- Battery Power Four 9V alkaline batteries (included) power the Electro 705 for approximately 500 hours giving a super clean low noise source of power making the Electro 705 the quietest amplifier available. Batteries can be easily tested by the press of a button.
- Capacitance Compensation Corrects for loss of rise time caused by the presence of electrode capacity. Up to 50 pF of electrode shunt capacity may be neutralized.
- Driven Guard Shield Stray capacitance can be further reduced by placing the driven guard shield (included) over the microelectrode holder at the input end of the probe.

- Tickler Circuit A momentary oscillation that helps achieve cell penetration
- Electrode Resistance Test The 705 provides a 1 nA electrode test current. Electrode resistance is monitored at the 1X output as a voltage (1 mV/M).
- Probe Test Port Allows the convenience of testing the amplifier's intrinsic noise and gain without cumbersome external test hookups. Gate leakage current can also be adjusted with minimum effort.
- Baseline Position Control Adds or subtracts up to 300 mV to the headstage output, allowing artifact voltages such as liquid junction potentials to be nulled prior to recording.
- **Differential Output** Two Electro 705s can be connected in tandem to create an optional differential amplifier probe system.

ELECTRO 705 SPECIFICATIONS

INPUT IMPEDANCE 1012 Ohms, shunted by 1 pF

100 Ohms, both outputs **OUTPUT IMPEDANCE**

GAIN X1: ±0.1% INPUT VOLTAGE RANGE RISETIME 15 µs, 10-90%

NOISE LEVEL 500 μV peak-to-peak*

INPUT CAPACITANCE COMPENSATION 0-50 pF

GATE LEAKAGE CURRENT ±10 pA, adjustable to zero

ELECTRODE RESISTANCE TEST 1 mV/ M Ohms DC POSITIONING + 300 mV

COMMON MODE REJECTION >104 (in differential mode) **POWER** Four 9V alkaline batteries, supplied

DIMENSIONS 8.5 x 3.5 x 2.2 in. (22 x 9 x 6 cm)

SHIPPING WEIGHT 5 lb (2.3 kg) * Full band width, with 20 M Ohms source

Electro 705 Electrometer SYS-705

Probe, driven guard shield and micropipette holder MEH1SF included for glass microelectrodes O.D. 1.0 mm, 1.2 mm, 1.5 mm, or 2.0 mm.

OPTIONAL ACCESSORIES

·	7.0000000		
M3301L	Micromanipulator (specify left- or right-handed)		
M-3	80° Tilting base		
RC1T	Reference cell (Ag/AgCl)		
2541	Driven guard shield for 705PF Probe		
MEH1SF	Microelectrode holder		
705PF	Replacement probe (includes calibration)*		

*Instrument must be returned to WPI for free calibration with new probe. See cables and connectors, page 98

See microelectrode holders, page 106 See capillary glass, page 110

Reference

Koch, U. (2000) "Interdependence of spatial and temporal coding in the auditory midbrain." Journal of Neurophysiology 83, 4, 2300-2314

Duo 773 Dual Microprobe System

2-channel intracellular amplifier



For intracellular dual or differential studies, the Duo773 has separate negative capacity controls and built-in active filtering that allows the precise balancing of time constants for artifact-free differential measurement. Comes complete with two probe headstages, 1015 Ohms & 1011 Ohms probes to monitor signals from ion-specific micro-electrodes as well as KCI-filled electrodes.

HE/

ACT

GAI

OU.

OU MAX PRC DC FIF

INP NO

* Although injected currents are "constant," the maximum current in a given situation will always be limited by the system compliance of 10 V.

**The 712P headstage may be used on either A or B channels, however Current Injection specifications do not apply when used on channel A. The 715P headstage may not be used on the R channel

References

L. Pluja (2000) "Electrical and mechanical effects of vasoactive intestinal peptide and pituitary adenylate cyclase-activating peptide in the rat colon involve different mechanisms." European Journal of Pharmacology 389, 217-

G. X. Wang, X. B. Zhou, et al. (2000) "Effects of mitoxantrone on excitation-contraction coupling in guinea pig ventricular myocytes." Journal of Pharmacology and Experimental Therapeutics 293, 2, 501-508.

S. Tsuruoka (2000) "Acute effect of cadmiummetallothionein on glucose and amino acid transport across the apical membrane of the rabbit proximal tubule perfused in vitro." Journal of Pharmacology and Experimental Therapeutics 292, 2, 769-777.

DUO 773 SPECIFICATIONS

ADSTAGE (PROBE)	712P (red, port "B")	715P (blue, port "A")
TIVE PROBE INPUT IMPEDANCE	>10 ¹¹ Ω	$10^{15}\Omega$
AIN .	x1, x10	x1
JTPUT RESISTANCE	100 Ω	100 Ω
JTPUT VOLTAGE RANGE	±10 V	± 10V
AXIMUM INPUT VOLTAGE	±15 V	±15 V
OBE LEAKAGE CURRENT	5 X 10 ⁻¹² A	10 ⁻¹⁴ A
POSITION ADJUST RANGE	± 300 mV	± 300 mV
ECTRODE RESISTANCE TEST CURRENT	1 nA	1 pA, 1 nA selectable
PUT CAPACITY COMPENSATION	+10 to -50 pF	0 to -10 pF
DISE		
Input shorted 20 M Ω carbon resistor	<50 μV p-p 10kHz bandwidth <200 μV p-p 10kHz bandwidth	<50μV p-p 10kHz bandwidth <200μV p-p 10kHz bandwidth

1 μs, typical 25 μs, typical (712P only)** ± 50 nA low range, ± 500 nA high range ± 500nA low range, ± 5 μA high range 20 mV/nA low range, 2 mV/nA high range 100 mV/nA low range, 10 mV/nA high range 3V low range, 10V high range

0-100 M Ω , 0-1000 M Ω

40 dB/decade, continuously variable 1-30 kHz

3.5-digit LED 200 mV, 2000 mV, 20 V, 200 nA, 2000 nA

17 x 5.25 x 10 in. (43 x 13 x 25 cm) Diameter: 12 mm Length: 34 mm 95-135 V or 220-240 V, 50/60 Hz 15 lb (7 kg) CF. CSA

10-90% direct input small signal

10-90% through 20 MΩ (-C "on")

Externally commanded Current

External current command factor

CURRENT INIECTION

Internal DC Current

Current monitor Compliance

Bridge balance

Accuracy and resolution **DIMENSIONS**

I OW PASS FILTER

METER SECTION

Display

Ranges

Probe

POWER

Instrument

SHIPPING WEIGHT

CERTIFICATION

Bridge amplifier gain

Headstage—Two gold-plated, epoxy sealed miniature active probes can be positioned directly to the measurement site. Microelectrode holders containing an Ag/AgCl electrochemical half-cells plug directly into the probes. Stray capacitance can be reduced by placing the included driven guard shield over the microelectrode holder at the end of the probe.

Capacity Compensation — Channel A can compensate up to 10 pF of electrode shunt capacity and Channel B can compensate up to 50 pF.

Tickler Circuit — Assists in cell penetration. The frequency and amplitude of the oscillations may be varied for differences in membrane thickness or cell size. The duration of tickle can be controlled either by using the momentary switch, a foot switch, or by applying a signal to the remote tickler input.

Active Filters — Low pass settings on a -40 dB/decade active filter vary the cutoff from 1 to 30 kHz. Either probe or bridge outputs may be selected for filtering.

Current Injection — Channel B can eject current through the microelectrode by applying a command signal to the stimulus input connector; the resulting output from the probe will then be a constant current replica of the input signal. Two ranges of current delivery are provided: 50 nA and 500 nA or by an external source. This source can be useful for delivering hyperpolarizing currents to stabilize the cell membrane potential and as a holding current for microiontophoresis.

Compliance Alarm — When the electrode voltage exceeds the probe input maximum allowed voltage, an audible over-compliance alarm will

Bridge Balance—Subtracts the excess electrode voltage associated with delivering current through the recording micropipette. Electrode resistances up to 1000 M Ω can be balanced in two ranges. The balanced signal is available from x10 or x50 front panel output connectors.

Independent Outputs — The Duo773 has an output for each probe independent of gain filtering or balancing. In addition the Duo773 has a 10x and a 50x output for easy integration to most data aguisition

Digital Meter — The Duo773 comes complete with a 3½-digit display for monitoring injection current or the voltages for either probe (single ended or differential).

SYS-773 Duo 773 Electrometer

Specify line voltage

Includes two probes (712P and 715P or two 712P) with driven guard shields and eight MEH1SF microelectrode holders for 1.0 mm, 1.2 mm, 1.5 mm, or 2.0 mm glass electrodes.

OPTIONAL ACCESSORIES


712P	Replacement probe (includes calibration)*
715P	Replacement probe (includes calibration)*
*Instrui	ment should be returned to WPI for free calibration with new probe.
2022	Pack Mount Kit 51/ in high

2933	Rack Mount Kit, 5¼-in. high
2547	Driven Guard Shield for 712P & 715P Probes
15790	Replacement Probe Handle
TW100F-4	Glass capillary with filament
TW150F-4	Glass capillary with filament


See Dri-Ref, page 63. See cables and connectors, page 98.



2547 Driven Guard Shield



Optional Holders for Intracellular **Amplifiers**



See Microelectrode Holders, page 106

WPI's Low-Noise Amplifiers Outperform Cheap Imitations



n amplifier, in simplest terms, is an electronic device that magnifies an input signal. However, the way an amplifier is designed to handle noise and bandwidth limitations greatly affects the quality and sustainability of the final output signal.

Defining terms

To knowledgeably discuss amplifiers, let's define a few terms.

- Gain The gain is the multiplier defining how much the amplitude of an input signal is increased. A signal with an ×1 gain is not amplified. An ×10 gain produces an output signal ten times greater than the input signal.
- Noise Any unwanted signal fluctuations are called noise. While noise can also result from external sources, for the purpose of this discussion, we are primarily concerned with the noise resulting from the inner workings of the electronic device, our amplifier. This intrinsic noise is called shot (or schott) noise.
- Signal to Noise Ratio (SNR) The ratio of the output signal to the noise of the amplifier is called the signal to noise ratio. The smaller the shot noise signal in an amplifier in comparison with the output signal, the easier the desired signal is to discriminate. When engineering an amplifier, the SNR may be improved by boosting the first stage gain to yield a larger output signal or by using quality components to minimize the shot noise level of the amplifier.
- Output Range The output range determines the maximum output signal that can be generated with the amplifier. It is

by the maximum voltage of the power supply. If

the amplitude of the output signal is too large for the output range, part of the signal is cut off (clipped).

- Rail The upper or lower limit of the amplifier range is called a rail. Signals that exceed the rail cannot be faithfully
- **DC Offset** DC offsets can appear in biological preparations. This offset is the amount the output signal is displaced away from a zero reference point, and it is usually a result of the potential difference at the electrode's tip.

How does an amplifier work?

Power Supply Rails Limits the Range

In a perfect world an input signal can be infinitely multiplied by the gain factor to determine the output signal. For example:

Input Signal	Gain	Output Signa
2mV	×1	2mV
2mV	×2	4mV
2mV	×10	20mV
2mV	×100	200mV
2m\/	×10 000	201/

In the real world, however, the power supply rails limit the possible output range of the amplifier. For example, a bio-amplifier could have a range of ±5.0V. In order for the output signal to be faithfully reproduced, the input signal times the gain factor must fall within the voltage window set by the power rails. Otherwise, the output signal will go off scale, and the input signal will not be faithfully

reproduced. This is called "hitting the rail."

In our example, a 1.0µV input signal at an ×10⁶ gain would generate a 1.0V output signal. Since the power supply is rated up to +5.0V, this output signal is clearly visible. If the input signal in this example is greater than $5.0\mu V$, the output signal would be greater than +5.0V. Since 5.0V is the top of the range that the power supply is capable of producing, the output signal hits the upper rail and gets cut off. This amplifier will give a +5.0V DC output signal for all input signals greater than or equal to 5.0µV. In this instance, a smaller gain factor should be used to bring the output signal back into the dynamic output range of the amplifier.

Noise Limits Amplifier Useability

All electronic devices produce their own internal electronic noise, an unavoidable signal that can mask the output signal. For example, if the input signal is 2mV and the noise is 1mV, the signal to noise ratio is two to one (2:1), and the output signal would be undetectable. In this case, it is nearly impossible to discern which part of the output is generated by noise and which part is the desired signal. (Fig. 1)




Fig. 1-The higher the signal to noise ratio, the more discernable the desired signal.

Ideally, the signal to noise ratio should be at least 50 to 1 to produce a quality output signal. A good signal to noise ratio can be achieved in one of two ways:

- · Boost the output signal by increasing the
- Reduce the noise.

While increasing the gain is the simplest solution, too much gain can impose a limitation on the dynamic range of the amplifier. Reducing noise is a more complicated solution, but it offers a greater range and more stability in the end.

Two-Stage Amplifiers

FD223A


DC

Bio-amplifiers usually involve multiple stages of amplification.

Stage One – The unadulterated signal coming into the amplifier is unaffected by the intrinsic noise of the amplifier. Then, it runs through the critical first stage of amplification where the signal is boosted by the primary gain factor to produce an output signal with the desired signal to noise ratio. The intrinsic noise is not amplified in the first stage. Higher gain factors used in the first stage of amplification can seriously limit the dynamic range available at output stage. Large stage one gains also limit the gain factor available in the second stage of amplification.

Stage 2 – The stage one output signal enters the second stage of amplification where both the signal and the noise from the first stage are amplified together by the second stage gain factor so that the signal is large enough to be seen on a chart recorder or data acquisition system. The second stage amplification is the gain the user controls. It does not change the signal to noise ratio.

Result of Amplifier Using Gain to Control Signal to Noise Ratio






Fig. 2-As the offset naturally increases over time, a poorly constructed amplifier will not be able to faithfully reproduce the signal. This offset can also be a result of gain drift which can occur as the temperature rises.

2 mm nin

Instead of using high gains in the first stage of amplification, a well constructed bio-amplifier that uses high quality components, like WPI's DAM series amplifiers, minimizes the noise in the first stage of amplification so that the dynamic range is retained throughout the amplification process. A poorly designed amplifier will simply increase the gain of the first stage amplification until the desired signal to noise ratio is reached.

Why not boost the power rails?

Theoretically, increasing the voltage rails powering the amplifier will increase the available dynamic output range. It would seem natural to increase the power supply rails coming into the amplifier in order to provide the capability for greater first stage gains. However, most data acquisition systems are limited to a maximum input signal ranging between ±10.0V. Therefore, it is not practical to increase the power rails of bio-amplifier

2

beyond ±10.0V. Since the industry standard limits us to ± 10.0 V power supply rails, the only way to improve the signal to noise ratio is to minimize the shot noise in the first stage of amplification. This is why high quality amplifier components are imperative.

Why does my signal flatline?

Regardless of the amplifier used, biological potentials are often accompanied by a DC offset, because the electrodes polarize over time. The DC offset naturally increases over time. Since the poorly constructed amplifier that utilizes greater first stage gain has restricted its dynamic range, it has limited ability to handle this offset. As the offset continues to increase, the output signal may eventually be forced by the offset into the rail causing the flat line (clipping the signal). (See Fig. 2.)

The amplifier that minimizes the noise in the

first stage amplification offers a larger dynamic output range and handles a much greater offset

Amplifier	AC/DC	Differential	Head- stage	EMG EKG	Stimu- lation	Isolated	Multi- channel	Battery Powered	Connectors
Intracellular Bioamplifiers									

IDZZJA	DC	•	•				2		Ζ ΙΙΙΙΙΙ ΡΙΙΙ
Electro 705	DC		*					*	2 mm pin
Duo773	DC	•	*		•		2		2 mm pin
Extracellula	ır Bioampl	ifiers							
ISODAM8A	DC	•	opt	*		•	4 - 8		Mini Banana or 8-pin DIN
ISO80	AC	•	•	•	•	•		•	Mini Banana
DAM50	AC/DC	*		*				*	RJ-11
DAM80	AC	*	*	•	•			*	Mini Banana
Transducer	Amplifiers	s							
BRIDGE8	DC	•					4 - 8		8-pin DIN WPI transducers
TBM4M	DC	•					4		8-pin DIN WPI transducers
Epithelial V	oltage/Cu	rrent Clan	p Bio Amp	lifier					
EVC4000	DC				•		1 - 4		Ussing 2 mm

WPI's amplifiers

The purchase of a low-noise amplifier pays dividends in the end. WPI's amplifiers were engineered for the bio-medical researcher. While 20-30µV of noise is common in bio-amplifiers, WPI's DAM series amplifiers generate 0.4µV RMS (root mean squared) at 0.1-100Hz. (That's equal to 2µV peak-to-peak.) The chart at left compares WPI's bio-amplifiers.

See www.wpiinc.com/amplifiers

Ш

DAM Series Bioamplifiers

A family of very low noise battery-operated amplifiers

· Gated or manual current generation for histological marking, iontophoresis, or cell stimulation.

> A very low noise remote active headstage (DAM80 only) is useful for very high impedance amplification utilizing glass or metal electrodes.

> > DAM series amplifiers can be used as standalone units on any tabletop, or use optional clampmounting hardware to locate them conveniently within the work area. Alternatively, a pair of amplifiers can be mounted into a standard equipment rack with a rack mount kit (#3484). A variety of hook up accessories are available to configure your application.

> > > Optional probe #5489 for use with DAM50 also includes microelectrode

adapter #5469.



WPI's DAM series amplifier's are well known as a standard of the industry for extracellular potential amplification. These battery powered bio-amplifiers are designed with a compact chassis profile that enables the user to locate the unit closer to the preparation and thereby minimize long lead lengths which contribute to noise. Each amplifier is equipped with selectable high and low filters, and a position control to offset galvanic potentials which may develop during recording. A choice of models offer additional features that are useful for certain applications:

#3294 Grounding Clip

#5469 Metal

Microelectrode

Extracellular recording

#5489 Probe

#2033 Adapter

Metal Microelectrode

(i.e., TM33B01)

Application Examples

Central Ground

DAM50

DAM50 #5447 Adapter #5371 Cable #5371 Cable Light Stimulation #300102 electrode extension, held by micromanipulators #3578 Cable Metal Microelectrodes Central Ground Immobilized specimen - EP2 Ag/AgCl Half Cell Ag/AgCl pellet ERG recording of fly eyes using metal microelectrode

Micromanipulator



DAM80, an AC amplifier only, features a very low noise headstage probe which can be mounted in micromanipulators for up-close cortical recording, for extracellular recording from high impedance glass or metal microelectrodes. Also provides a gated current for tissue marking. Microelectrode holder MEH3SB is recommended.

Included with the DAM-80 is a Startup Kit containing the following accessories needed for basic metal electrode electrophysiology research:

CBL102 Cable, BNC-to-3.5mm plug, 6 ft (2m) (two)

5469 Adapter, mini-banana to 0.031 skt. (two)

13388 Adapter, mini-banana to 2mm skt. (two)

3294 Cable, ground clip to wire, 3 ft

2033 Mini-banana plug, black

2034 Mini-banana plug, red

2035 Mini-banana plug solderable turrent (two) EP1 Ag/AgCl pellet (70 mm wire) 1mm diam x 2.5 mm long M3301EH Electrode Holder, 14cm (two) **5470** 0.031-inch jack on 12-inch wire (package of 4)

DAM SERIES SPECIFICATIONS

INPUT IMPEDANCE INPUT LEAKAGE CURRENT MAX. DC DIFFERENTIAL SIGNAL

COMMON MODE REJECTION RATIO

INPUT CAPACITANCE AC MODE NOISE AC MODE NOISE

DC MODE NOISE (DAM50) BANDWIDTH FILTER SETTINGS

AC Mode

AC Mode (DAM80) DC Mode (DAM50)

OUTPUT CONNECTORS

OUTPUT VOLTAGE SWING OUTPUT IMPEDANCE BATTERY TEST CALIBRATOR SIGNAL POSITION

CURRENT SOURCE

DAM80: DC Generator EXTERNAL COMMAND

AC OR DC CURRENT WAVEFORM

BATTERIES **DIMENSIONS** DAM50 DAM80

SHIPPING WEIGHT

 $10^{12} \Omega$, common mode and differential

50 pA (typical) ± 2.5 V (DAM 50)

AC: 100x, 1000x, 10000x DC: 10x, 100x, 1000x (DAM50)

100 dB @ 50/60 Hz

20 pF

 $0.4~\mu V$ RMS (2 μV p-p) 0.1-100 Hz 2.6 μV RMS (10 μV p-p) 1 Hz-10 kH 7.5 μV RMS (30 μV p-p) 3-10 kHz

Low frequency, 0.1, 1, 10, 300 Hz High frequency, 0.1, 1, 3, 10 kHz High frequency, 0.1, 1, 3, 10 kHz BNC on DAM50; 3.5 mm MiniPhone

connector on DAM80

±8 V 470Ω Audible tone 10 Hz square wave Approximately 250 mV

0 to ±50 μA, variable Input Voltage ±10 V commands ±50 μA max. amplitude @ 200 KΩ 2 x 9V alkaline (included)

8 x 4 x 1.75 in. (20.3 x 10.2 x 4.4 cm) $7 \times 4 \times 1.75$ in. $(17.8 \times 10.2 \times 4.4 \text{ cm})$

3.5 lb (1.6 kg)

DAM80 AC

Yes Yes

differential 100-10K (AC) yes yes

3.5 mm mini phone mini banana

(2) nine volt alkaline batteries

317 0 0.031 men jack on 12 men vine (package or 1)		
#3294 Grounding Clip (included) DAM80 or ISO-80 DAM80 or ISO-80	FEATURE	DAM50
#13388 Adapter #5371 Cable #300102 electrode extension, held by micromanipulators	Input Mode Input configuration Gain Range High / Low Filters Offset position control Current Generator Remote Active headstage Output connection	AC/DC differential/single ended 100-10K (AC), 10-1K (DC) yes yes No No BNC
FMG recording of a rat Metal Microelectrodes Ag/AgGl skin electrode pad (#EL203)	Standard input connection* Power supply *see optional accessories for	unterminated wire (2) nine volt alkaline batteries additional alternatives

SYS-DAM50	Bio-amplifier
SYS-DAM80	Bio-amplifier with active probe (DAM80P)
OPTIONAL AC	CESSORIES
DAM80P	Replacement Probe
3072	6 Replacement Modular Cables (DAM50)
3517	2 Optional Shielded Modular Cables (DAM50)
CBL102	3.5 mm Phone plug-to-BNC Cable
2851	BNC-to-BNC Cable
2033	Black Insulated Mini-Banana Plug
2034	Red Insulated Mini-Banana Plug
2035	Uninsulated Mini-Banana Plug
2101	9V Alkaline Battery, each (2 required)

Rack Mount Kit (for 1 or 2 DAM preamps)
Ringstand Mounting Kit
Electrode Adapter (DAM50)
Metal Microelectrode Adapter for DAM80
(mini-banana plug to 0.031 in. (0.79 mm) socket)
Adapter for Metal Microelectrode (DAM50)
Adapter, mini-banana plug to 2mm socket
Cable, Low Noise (2 mm pin to 2 mm pin)
Adapter Cable for Ag/AgCl pellets (2 mm pin)
Electrode Extension, 4-inch
9V NiMH Battery

Also see cables and connectors, metal microelectrodes, carbon-filled micropipettes.

Isolated Differential Amplifier

The improved ISO-80 provides low noise AC coupled amplification and offers excellent recording performance for monitoring extracellular nerve action potentials in vitro and in living animals. The ISO-80 is provided with a remote headstage (1 m cable) which incorporates an electrode impedance test function and a constant current stimulator.

The constant current stimulator can be used for cell marking, stimulation or electrode cleaning. Typical applications include measuring EMG, EEG, extracellular and action potentials in vitro or in vivo. The ISO-80 system is DC isolated from the subject ground and employs state of the art electro-magnetic shielding for improved noise rejection. The amplifier employs both high pass and low pass filtering with gain from 100 to 10,000. The lowest low-pass setting is 5Hz and

the upper passband is 10 kHz.



INPUT RESISTANCE >1011 Ohms, Common Mode and

differential

INPUT LEAKAGE CURRENT 50 picoamperes, max. AMPLIFICATION $\times 10^{2}$, $\times 10^{3}$, $\times 10^{4}$ COMMON MODE REJECTION RATIO 100 dB typ. @ 50/60 Hz EOUIVALENT NOISE SIGNAL INPUT 0.4 microvolts rms (0.1-100 Hz) 2.0 microvolts rms (1 Hz - 10 kHz)

FILTER SETTINGS

5, 10, 100, 300 Hz Low frequency High frequency 100 Hz, 1, 3, 10 kHz

MAX. OUTPUT VOLTAGE SWING

ELECTRODE IMPEDANCE RANGE 100 kOhm - 10 MOhm @ 300 Hz

0 to ±20 micro amperes (constant cur-STIMULATION CURRENT

MAXIMUM STIMULATION VOLTAGE ±15 volts MAXIMUM ELECTRODE VOLTAGE ±40 volts DISPLAY 3½-digit LCD **BATTERY TEST** Low battery display

POWER Two 9-volt NiCad batteries & charger,

supplied

SHIPPING WEIGHT 4 lb (1.8 kg) Included with the ISO-80 is a Startup Kit containing the following accessories needed for basic metal electrode electrophysiology research:

 ϵ

CBL102 Cable, BNC-to-3.5mm plug, 6 ft (2m) (two)

5469 Adapter, mini-banana to 0.031 skt. (two)

150.80

ISOLATED BIO-AMP

13388 Adapter, mini-banana to 2mm skt. (two)

3294 Cable, ground clip to wire, 3 ft

2033 Mini-banana plug, black

2034 Mini-banana plug, red


2035 Mini-banana plug solderable turrent (two)

EP1 Ag/AgCl pellet (70 mm wire) 1mm diam x 2.5 mm long

M3301EH Electrode Holder, 14cm (two)

5470 0.031-inch jack on 12-inch wire (package of 4)

ISO-80	Isolated Bioamplifier w/ active probe (ISO80P)
	Specify line voltage
OPTIONA	L ACCESSORIES
ISO-80P	Replacement ISO-80 Probe
CBL102	3.5 mm phone plug-to-BNC cable



4-Channel Transducer Amplifier

Transbridge (TBM4M) is a four-channel analog transducer manifold, specifically designed to amplify output voltage signals from pressure, force, displacement, and temperature transducers as well as a wide variety of other signal sources. Analog output signals are available from



Transducers available separately

each channel for input to a data acquisition system for digital signal processing in a computer.

Each channel contains a regulated 10-volt power supply (+5 and -5 volts with respect to signal ground) to provide DC power to transducers, and a precision differential amplifier with selectable voltage amplification and variable position adjustment control.

Transducers can be connected to Transbridge via any of the 8-pin connectors on the front panel. Four spare 8-pin DIN plugs are provided with each instrument to allow you to rewire cables of other manufacturers' transducers and connect them to Transbridge. Each Transbridge channel may be used in either Full Bridge or the Half Bridge mode independently. For transducer types other than resistive bridges, such as active transistor circuits, magnetic, photocell or piezoelectric devices, the instrument's differential amplifiers may still be used effectively for signal amplification in differential (full bridge) and single-ended (half bridge) modes.

SYS-TBM4M	Transbridge Transducer Amplifier	\$2,085
	Specify line voltage	
OPTIONAL AC	CCESSORIES	
13024	Single Rack Mount Kit	\$11 4
13025	Dual Rack Mount Kit	\$ 114
500184	BNC-to-BNC cable, 10 ft	§31
3161	8-pin DIN plug	§28
3718	Package of 4, 8-pin DIN (startup kit)	\$98

Dual Microiontophoresis Current Generator



The Dual Microiontophoresis Current Generator (Model 260) is an electrically isolated, battery-operated instrument designed for the electro-iontophoresis of dyes, drugs and charged substances from micropipettes. Two identical battery operated current generators are available. In ordinary use, the two current generators are operated in parallel providing two distinct currents; one for preventing substances in the micropipette from outward diffusion (the retain or hold current) and the second for actively ejecting charged material. For pipettes with submicron tips, a hold current may not be necessary if there is little outward diffusion of pipette material. Model 260 is powered by two 9-volt alkaline batteries per side (four, in total); unique circuitry converts the ±9 V to ±100 V without a transformer, yielding an exceedingly quiet output.

SYS-260	Dual Microiontophoresis Current Generator					
OPTIONAL	ACCESSORIES					
2933	Rack Mount Kit, 5¼-in. high					

Ωmega-Tip-Z™

Millivolt and megohm meter measures impedance of metal or glass capillary microelectrodes

Omega-Tip-Z was created especially for measuring impedance in etched tungsten, platinum-iridium* and steel microelectrodes, as well as electrolyte-filled micropipettes. The meter's AC impedance-

measuring circuit is unaffected by electrode offset or tip junction potentials. The gold-plated miniature probe lets you conveniently monitor micro-electrode impedance in electrolytes, and an electrode tip cleaning feature lets you remove buildup quickly. Omega-Tip-Z can also measure DC electrode tip potentials up to 2000 millivolts. The instrument operates for hundreds of hours without battery failure.

NOTE: Metal microelectrodes which have been precalibrated at 1 kHz should be baselined for use with Omega-Tip-Z.



*See Metal Microelectrodes, page 102.

SYS-OMEGAZ	Omega-Tip-Z with Probe & Holder			
711P	Replacement Probe			
5468	Adapter to connect metal microelectrodes to probe, 2 mm socket to .031 in. receptacle			
OPTIONAL AC	CESSORIES			
Z-LITE	Fiber Optic Illuminator (115v, 60Hz, beige case)			
Z-LITE-Z	Fiber Optic Illuminator (230v, 80Hz, black case)			
500186	Bifurcated Light Guide with lenses			
Z-LITE-186	Z-Lite Illuminator and bifurcated light guide			



FORCE TRANSDUCERS

These rigid-lever force transducers transform applied force into proportional voltage. Using balanced strain gauges, FORT transducers produce linear output voltage vs. applied force input with very little

To use, clamp the handle of the FORT transducer in a horizontal position and apply the forces to be measured to a rivet or hook mounted in the hole at the end of the flat sensing leaf.

FORT100	Force Transducer (100 g)
FORT250	Force Transducer (250 g)
FORT1000	Force Transducer (1000 g)
FORT5000	Force Transducer (5000 g)

FORT SPECIFICATIONS

FORCE RANGES, FULL SCALE OUTPUT SENSITIVITY (± 10%) INPUT & OUTPUT RESISTANCE **RESOLUTION** RESONANT FREQUENCY LINEARITY ERROR MAX. OPERATING VOLTAGE MAXIMUM APPLIED FORCE DRIFT

DIMENSIONS

WEIGHT (excluding cable)

FORT100 100 grams 7 μV/V/g 350 Q 300 Hz 10 V AC or DC

0.01% of full scale force Less than 0.1% of full scale 3× rated full scale force thermally compensated 0.3 inch diam $\times 4$ in. $(7.6 \text{ mm diam} \times 10.2 \text{ mm})$ 0.3 oz (8 g)

FORT250 250 grams 3 μV/V/g 350 Q 0.01% of full scale force 300 Hz Less than 0.1% of full scale 10 V AC or DC

3× rated full scale force thermally compensated 0.3 inch diam $\times 4$ in. $(7.6 \text{ mm diam} \times 10.2 \text{ mm})$ 0.3 oz (8 g)

FORT1000 1000 grams

0.84 µV/V/g 350 Q 0.01% of full scale force

300 Hz

Less than 0.1% of full scale 10 V AC or DC 3× rated full scale force thermally compensated

0.3 inch diam $\times 4$ in. $(7.6 \text{ mm diam} \times 10.2 \text{ mm})$ 0.3 oz (8 g)

FORT5000

5000 grams 0.38 µV/V/g 350 Q

0.1% of full scale force

60 Hz

Less than 0.1% of full scale

10 V AC or DC

3× rated full scale force thermally compensated

0.3 inch diam $\times 4$ in. $(7.6 \text{ mm diam} \times 10.2 \text{ mm})$

0.3 oz (8 g)



10g & 25 Force **Transducers**

FORCE RANGE, FULL SCALE **OUTPUT SENSITIVITY** INPUT & OUTPUT RESISTANCE **RESOLUTION** RESONANT FREQUENCY LINEARITY ERROR MAXIMUM APPLIED FORCE

MAXIMUM OPERATING VOLTAGE DRIFT **DIMENSIONS**

FORT10g 0-10 grams

10 mV/gm, nominal

 1500Ω < 1mg

450 Hz

<0.2% of full scale

10 V DC (-5V ~ +5V or 0 ~ 10V)

2× rated full scale force

<30 mg/hr

40 x 22 x 19 mm Handle 88 mm

100 gram

FORT25

0-25 grams

3 mV/gm, nominal

1500 Ω

< 2mg

450 Hz

<0.2% of full scale

10 V DC (-5V ~ +5V or 0 ~ 10V)

3× rated full scale force

<50 mg/hr


40 x 22 x 19 mm Handle 109 mm

100 gram

These 10-gram and 25-gram force transducers are reliable tools for high precision force measurement. Using balanced semiconductor strain gauges, both produce linear output voltage vs. applied force input with very little deflection. The rigid lever force transducer transforms the applied force into a proportional voltage. Featuring a temperaturecompensated, full-bridge configuration with four high sensitivity semiconductor strain gauges, these transducers have broad dynamic measuring range and very high sensitivity.

To use, clamp the handle of the FORT10 or FORT25 transducer in a horizontal position and apply the forces to be measured to a rivet or hook mounted in the hole at the end of the flat sensing leaf.

WEIGHT



- Simultaneously measures electric potential and pressure
 - Preset internal microelectrode pressure
 - Air-filled system no debubbling

Model 900A is designed to measure hydrostatic pressure in small vessels and cells. Pressure ranges of -200 to +400 mm Hg can be measured with stability and accuracy. The system's sensing element is an electrolyte-filled glass microelectrode with a tip diameter range of 2 to 5 microns.

Pressures of electrolyte solutions are measured by maintaining a salt concentration gradient at the tip of the sensing electrode in dynamic equilibrium by applying an equal air pressure inside the microelectrode. The pressure reading appears on the front panel display and via the BNC recorder

Because the piezoelectric pressure controller uses external pressure and vacuum sources, pressures lower than -200 to greater than +400 mm Hg can be quickly and accurately measured at the microelectrode tip.

The open pressure chamber is almost immune to vibrations and movements and, unless they are extremely large, the open system is unaffected by leaks. The pressure controller is contained in a small, lightweight enclosure that can easily be mounted near the micropipette to help reduce dead space. It includes an amplifier, a piezoelectric valve and a pressure transducer. The user supplies fluid-filled microelectrodes, +500 mm Hg pressure source and a -300 mm Hg vacuum source.

Measuring electric potential and pressure simultaneously lets you use potential recording as an additional cue for locating the electrode where visibility is limited, or correlate pressure and potential when this is meaningful.

The unique "Set Pressure" mode lets you preset

the internal pressure of the microelectrode select a positive pressure for flushing the tip, or a negative pressure for pulling solution into the tip. By disconnecting the microelectrode holder and attaching the tubing to a manometer, you can check the calibration against a standard.

A built-in alarm sounds to indicate maximum pressure. The alarm also sounds when the tip is blocked or electrical continuity is broken (e.g., the microelectrode comes out of the solution, too little filling solution to cover the Ag/AgCl pellet, disconnected ground reference, etc.).

The piezoelectric pressure controller regulates internal pipette pressure by controlling air flow into and out of a small pressure chamber. A vacuum source is connected on the outlet side of the chamber, and a piezoelectric valve meters air entering the pressurized chamber. The residual volume of the pressure chamber includes the micropipette, the

connecting tubing and the pressure transducer on the outlet side of the piezoelectric valve. The 900A accurately controls and adjusts the pressure in the chamber to match pressures applied externally to the microelectrode tip.

The response time of the piezoelectric valve is 0.5 ms from fully closed to fully open. Overall system response time depends largely on the amount of residual volume in the tubing. When this volume is small, the system responds very rapidly (typically 10 milliseconds).

The lightweight pressure controller pod may be mounted close to the microelectrode using smallbore tubing, to minimize system dead space.

Microelectrode holders MEH6RF and MEH6SF for 1.0 mm O.D. capillary glass included. (1.2, 1.5 and 2.0 mm also available — please specify O.D. requirement when ordering.)

SYS-900A Micropressure System

System price includes a one-day technical training session at WPI in Sarasota, Florida.

Specify line voltage

OPTIONAL ACCESSORIES 900AP Replacement Probe Pressure Calibration Chamber CAL900A Probe Extension Cable 3491 Rack Mount Kit 2933 5332 Replacement Liquid Trap **MEH6RF** Micropipette Holder (1.0, 1.2, 1.5 or 2.0 mm — Specify O.D.) **MEH6SF** Micropipette Holder (1.0, 1.2, 1.5 or 2.0 mm — Specify O.D.) TIPTW900A Prepulled Micropipette for 900A (1 mm thin-wall, 2 μTip) (pkg of 10) Replacement Pressure Pod 900APP SYS-PM015D Pressure Manometer (15 psi)

Automated TEER Measuring **System**



The **REMS AutoSampler** automates measurements of electrical resistance of transepithelial, transendothelial or Caco-2 cell membranes being grown to confluence on microporous filters of high throughput screening (HTS) 24- and 96-well microplates. It is a PC-controlled tissue resistance measurement system that offers reproducibility, accuracy, flexibility, and ease-of-operation for this kind of measurement. Automated measurement of tissue resistance in cell culture microplates provides the important advantages of speed, precision, decreased opportunity for contamination and the instant availability of measured resistance data on a computer. These measurements are useful in applications such as drug bioavailability studies and studies on the mechanisms of drug transport.

The main components of the REMS AutoSampler include: the robotic sampler that moves the electrode over each well of the microplate, the electrode which is located on the robotic arm, a base plate for the 24- and 96-well tray, a Windows-based data acquisition card, the REMS interface unit and the REMS software to operate the system on a Windows-based

The REMS AutoSampler automates TEER measurements previously made with WPI's **EVOM** Epithelial Voltohmmeter. Automated tissue resistance measurements up to 20 k Ω can be performed on 24- or 96-well HTS microplates. Microplates presently supported include the Corning Costar HTS Transwell-24, Falcon HTS Multiwell insert systems, and Millipore Multiscreen™ CaCo 96-well plate.

The REMS AutoSampler is designed to facilitate integration with other robotic systems. Special locating bars are installed on the REMS base platform that allow other system robots to place an HTS tray into a precise location on the REMS base.

The REMS AutoSampler will automatically measure and record tissue resistance from a user-specified matrix of culture wells on the microplate. According to the specified sequence, the robotic arm moves over the identified wells taking TEER measurements. By means of a x-y-z locating system, the electrode-containing arm is positioned precisely

SYS-REMS Automated Tissue Resistance Measuring System

Includes robot sampler, base plate, data acquisition board; computer, display, keyboard, mouse; software for Windows XP or Vista; and electrode for either 24-well plate (Corning Costar HTS Transwell-24 or Falcon HTS Multiwell) or 96-well plate (Millipore Multiscreen CaCo) — SPECIFY WHEN ORDERING.

ACCESSORIES

ACCESSORIE	
REMS-24	Replacement REMS STX Electrode for 24-well HTS Plate
REMS-96	Replacement REMS STX Electrode for MilliporeTM 96-well Plate
REMS-96C	Replacement STX Electrode for Corning 96-well plate

Contact WPI for detailed information.

REMS AUTOSAMPLER SPECIFICATIONS

MEMBRANE RESISTANCE RANGE 0 to 2000 Ω and 0 to 20 k Ω AC SQUARE WAVE CURRENT +/- 20 µA @ 12.5 Hz **ELECTRODE POSITIONING** Resolution in X, Y and Z: +/- 1 mm ELECTRODE PERFORMANCE Repeatability in X, Y and Z: +/- 0.25 mm ELECTRODE ARM SPEED X- and Y-axis: 250 mm/sec Z-axis: 247.3 mm/sec

TYPICAL MEASUREMENT TIME

SCAN PATTERN LINE VOLTAGE

DIMENSIONS

WEIGHT

1 min, 10 sec

Choice of any well pattern sampling User specified: 100/120 V or 220/240 V

 $53.5 \times 43.7 \times 37.1$ cm $(21\%2 \times 17\%6 \times 14\% \text{ in.})$

24 kg (52 lb)

and reproducibly over each well. The ability of the REMS AutoSampler to reproducibly and precisely locate the electrode results in highly reproducible TEER measurements. TEER measurements are stored in the computer as the electrode moves from one well to the next. The Windowsbased software provides user-friendly features to acquire, display and store the tissue resistance measurements.

The REMS electrode is very compact and robust in design. Each of two rod-shaped probes, 1.5 mm in diameter, consists of a pair of electrodes: one electrode for injecting current and the other for measuring the voltage. The use of two pairs of electrodes eliminates the error caused by the electrode-liquid interface. To take a measurement, the robot inserts one probe into the center of the filter well and the other into the opening slot of the 24- or 96-well plate. The use of AC current to measure resistance provides several advantages over DC current, including:

- Absence of offset voltages on measurements;
 - There is a zero net current being passed through the membrane and, therefore, it is not adversely affected by a current charge;
 - No electrochemical deposition of electrode metal.

The REMS AutoSampler also features a rinse and calibration check station. If occasional rinsing of the REMS electrode is required it may be sent to a rinse station by pressing the rinse station button on the menu bar.

Trans Epithelial Electric Resistance (TEER) Measurements

During the last two decades TEER measurements have become universally established as the most convenient. reliable and non-destructive method to evaluate and monitor the growth of epithelial tissue cultures in vitro. The confluence of the cellular monolayer is quickly

determined by a sharp increase in TEER.

TEER measurement technology, which was first introduced by WPI in the mid-1980's, has since been perfected and expanded to include a range of TEER related manual and automatic instrumentation.

EVOM

Epithelial Voltohmmeter

Manual TEER measurements of epithelial cells in 6-, 12-, and 24-well

Electrically isolated meter that plugs into a standard outlet for continual readout without push buttons

Compatible with Endohm chambers

STX2 manual electrodes and test electrode included with every meter

Free standing with tilt bail, making viewing results easy

The EVOM was the first instrument designed specifically to perform routine Trans Epithelial Electrical Resistance (TEER) measurement in tissue culture research. EVOM2 is the next generation, redesigned for ease of use. The EVOM2 not only qualitatively measures cell monolayer health, but also quantitatively measures cellular confluence. The unique electronic circuit of the EVOM2 and the included STX2 electrode detect the

confluence of the cellular monolayer. When combined with WPI's Endohm chamber, the EVOM2 can also be used to perform more accurate quantitative measurements or lower resistance measurements like transendothelial electrical resistance measurements.

The isolated power source of the EVOM2 was specifically designed to avoid adverse effects on tissue and the formation of electrode metal deposits, even when it is plugged into a standard wall outlet. Now, the EVOM2 is always on when you need it. In addition, its rechargeable battery allows up to 10 hours of mobile use. The four and a half digit readout provides a range of 1-9,999 Ω . The included



EVOM2	Epithelial Tissu	ie Voltohmmeter	(includes STX2	electrode set)
-------	------------------	-----------------	----------------	----------------

REPLACEMENTS AND ACCESSORIES

STX2 Replacement "Chopstick" Electrode Set

STX3 Adjustable Tip Spread "Chopstick" Electrode Set 3993 Electrode Adapter (for electrodes with 2 mm pins)

Replacement Battery, Rechargeable NiMH 91736

91750 EVOM2 Test Resistor test electrode lets you calibrate the resistance measurements for an accurate reading every time, and the voltage meter never needs calibration. An analog BNC output is standard with the EVOM2, providing an output port for recording data or remote display of the EVOM2 output.

EVOM2 comes complete with the popular STX2 "chopstick" electrodes, 4 mm wide and 1mm thick. Each stick of the electrode pair contains a silver/silver-chloride pellet for measuring voltage and a silver electrode for passing current. The small size of each electrode is designed to facilitate placement of the electrodes into a variety of standard cell culture wells.

EVOM2 SPECIFICATIONS

MEMBRANE VOLTAGE RANGE RESOLUTION RESISTANCE RANGE

RESISTANCE RESOLUTION 1 Ω AC SQUARE WAVE CURRENT

POWER

12 VDC

NOMINAL BATTERY RUN TIME **BNC OUTPUT**

DIMENSIONS WEIGHT

ELECTRODE CONNECTION

TEST RESISTOR **ENVIRONMENTAL RANGE** ±200 mV 0.1 mV 0 to 9999 Ω

±10 µA nominal at 12.5 Hz Internal rechargeable 6V NiMH 2700 mAH battery with external

OUTPUT

supply for recharging

1-10 V (1 mV/ohm)

19 x 11 x 6 cm (7.25" x 4.25" x 2.30")

1.4 kg (3 lb)

RJ-11connector (telephone style)

External, 1000 Ω 10-38°C (50-100°F)

0-90% non-condensing relative humidity

TEER measurements in High Throughput

STX100

Series Electrodes

- Designed for 24-well HTS plate (Corning Costar and BD Falcon) and with 96-well plates (Millipore and BD Falcon)
- Improved accuracy down to 5 Ohm
- Sterilized with EtO, alcohol or bactericide

With the development of a High Throughput Screening (HTS) protocol for faster drug discovery, a new line of cell culture filter plates have been introduced by several major cell culture insert manufacturers. These HTS plates normally have either 24 or 96 individual cell culture inserts "bonded" together as one plate so that it can be handled by a robot apparatus. In response to these developments, WPI has developed an automatic REMS system and a manual electrode, STX100, for TEER measurements using HTS plates.

STX100's design is based on the same reliable design principle as the universally used STX2 electrode, with several important modifications. The size of the electrode tip has been reduced to 1.5 mm to facilitate positioning through the narrower slit of the HTS plate. The STX100 electrode itself is constructed using a stronger material for higher durability and maximum usage applications. The bottom section of the electrode is shaped to fit neatly into the "keyhole" shaped filter well. This enables the STX100 electrode to produce increased accuracy and reproducibility of TEER readings ($\pm 5\Omega$) compared to the standard STX2. Several versions of STX100 are available, designed to fit the Corning Costar 24-well HTS plate, the Falcon 24 well HTS plate,

and the Millipore Multiscreen CaCo 96-well plate. Measurement can be directly performed when the HTS plate is in either a common or divided tray, reducing the possibility of contamination as well as mechanical damage to the cultured cells.

STX100C	STX100 for Corning Costar HTS Transwell-24
STX100F	STX100 for BD Falcon HTS Multiwell Insert System
STX100M	STX100 for Millipore Multiscreen™ HTS 96-Well Plate
STX100C96	STX100 for Corning HTS 96-Well Plate

OPTIONAL ACCESSORIES

OFIIOIVAL	ACCESSORIES
13685	Modular Cable, 7 ft
13347	Chart Recorder Adapter
2851	Standard BNC Cable, 5'2"
500184	Standard BNC Cable, 10 ft (3m)

CaliCell[™]

Cell culture cups with synthetic membrane for testing STX electrodes, Endohm and Ussing chambers

It takes a long time and a lot of work to grow a batch of cells, so you will want to make certain that your test apparatus is functioning properly. The CaliCell™ provides a quick and positive way to test STX electrodes, EVOMs, Endohm, and Ussing chamber.

The CaliCell™ is a major improvement in TEER electrode calibration. Its membrane makes use of our unique electric current constriction technology to produce resistance readings comparable to those obtained with real cell cultures. The CaliCell™ does not have to be refrigerated, and can be cleaned and sterilized with alcohol. Readings will not drift over time as long as the unit is kept in good physical condition.

CALICELL-12

12 mm Calibration Cell for Endohm-6/Endohm-12



Screening (HTS) cell culture filter plates



Endohm™

For TEER measurement of endothelial cell cultures in individual cups

- Compatible with EVOM2
- Improved accuracy of 1-2 Ohm
- Accommodates 6mm, 12mm, 24mm cups and Costar Snapwell cup
- Sterilized with EtO, alcohol or a bactericide

Using WPI's EVOM2 resistance meter, Endohm chambers provide reproducible resistance measurements of endothelial tissue in culture cups. Transfer cups from their culture wells to the Endohm chamber for measurement rather than using hand-held electrodes. The chamber and the cap each contain a pair of concentric electrodes: a voltage-sensing silver/silver chloride pellet in the center plus an annular current electrode. The height of the top electrode can be adjusted to fit cell culture cups of different manufacture. Endohm's symmetrically apposing circular disc electrodes, situated above and beneath the membrane, allow a more uniform

current density to flow across the membrane than with STX2 electrodes. The background resistance of a blank insert is re-

duced from 150 Ω (when using WPI's hand-held STX2 electrodes) to less than 5 Ω . With Endohm's fixed electrode geometry, variation of readings on a given sample is reduced from 10-30 Ω with STX2 electrodes (depending on the experience of the user) to 1-2. Compared with other resistance measurement methods, Endohm with EVOM2 offers a much more convenient and economic solution to "leaky tissue" measurement. Because of the uniform density of the AC square wave current from EVOM2, errors owing to electrode polarization or membrane capacitance are largely eliminated. Endohm together with EVOM2 offers the most accurate and economical endothelial ohmmeter now avail-able. To date, cups from Costar, Millipore, ICN Biomedicals, and Falcon have been tested. Endohm chambers may be sterilized with EtO, alcohol or a bactericide (also see: Cidex, Microsurgery section); not autoclavable.

ENDOHM-6	Endohm for 6 mm culture cup (24 wells per plate)						
ENDOHM-12	Endohm for 12 mm culture cup (12 wells per plate)						
ENDOHM-24SNAP	Endohm for 24 mm & Costar Snapwell™ cup (6 wells per plate)						
Requires EVOM2, EVOM, EVOMX or Millicell ERS-2							
53330-01 Replacement Endohm Cable							

For electrophysiological investigation of epithelial transport

- Direct connect low-resistance electrodes
- Simplified operation, easy to control temperature and clean after use
- Luer type leak-free attachment of tubing and electrodes
- Recessed electrode ports to avoid air bubble formation
- Secure membrane holding by sharp stainless steel pins or O-ring
- Specialized chamber adapts cell culture insert (Costar Snapwell) for monolayer cell culture
- Chambers with rectangular openings for tubular tissues from small animals

WPI's Ussing System offers researchers a quick, effective means of making low-resistance electrical connections to the Ussing chamber without need of long agar bridges or Calomel half-cells. Ag/AgCl half-cells screw into short tubes which plug firmly into place in the chamber's luer ports. These direct-connect electrodes eliminate the inconvenience and expense of Calomel half-cells in open liquids. The system includes one Ussing Chamber (eight sizes available), support stand, electrode kit, glass circulation reservoir (two sizes available), and a tubing start-up kit (25 feet of 0.375-in. tubing, 10 feet of 0.156-in. tubing, plus four male luer fittings, two compressor clamps, one Y-connector, and one clip). Sixteen possible system configurations are listed at right. Components are also available separately. (Preamplifier in photo not included.)

Ussing Chambers

WPI's classical Ussing Chambers are well established perfusion chambers that are easy to operate, easy to control temperature, and easy to clean after use. Hundreds of them are used daily by scientists in the field. Ussing Chambers are machined from solid acrylic with eight entry ports for fluid lines, electrodes, or agar bridges. For easy, leak-free attachment of tubing and electrodes, all eight ports are luer type. The four ports for voltage and current electrodes are recessed to prevent formation of air bubbles in the chamber. The fluid compartments in each side of the chamber are separated by the epithelial membrane being studied. Sharp stainless steel pins on one side of the chamber hold the membrane in position and mate with holes in the opposite chamber interface. (In the CHM4, tissue is held by an O-ring instead of pins.)

The CHM5 chamber adapts the Costar Snapwell, a cell culture insert for monolayer cell culture, into WPI's "classical" epithelial voltage clamp system. Until now, classical Ussing Chambers have not been widely used for monolayer cell culture inserts because most inserts have a very deep profile, limiting good fluid perfusion at the surface of the membrane and limiting voltage electrodes from measuring the potential close to the surface of the membrane. CHM5 solves these problems: Perfusion fluid is introduced into the chamber at an angle so that it will flow directly to the surface of the membrane. The voltage electrode is also inserted into the chamber at an angle so as to reduce the distance between the surface of the membrane and the electrode.



Complete Ussing System includes stand, glass reservoir, electrodes, Ussing chamber and tubing (EVC3 preamp and post mounting kit not included—see page 18).


Two small chambers with rectangular openings are designed for tubular tissue from small animals such as the mouse intestinal tract membrane (CHM6) and rat intestinal tract membrane (CHM7). The rectangular opening more closely matches the shape of the tissue than would a circular opening, significantly increasing the membrane area available for testing. The larger membrane area increases the transport rate of low permeability chemicals; it also reduces the electrical resistance of the system for easier current clamping.

Optional Drains

Drains may be added to Ussing chambers to allow quick and complete evacuation of radioactive or toxic substances. To have drains added at the time of order, add a "D" to the part number (such as "USS1LD"); cost of the drain (\$105 to \$140) will be added to the cost of the chamber or system ordered.

Cartridge Electrodes

The Electrode Kit contains four voltage/current electrodes, plus four luer-tipped cartridges. Electrodes are threaded and screw securely into the end of each cartridge. The luer tip then plugs securely into the luer openings of the chamber. The cable from each electrode terminates with a 2 mm pin which may be plugged into voltage/current clamps such as WPI's DVC1000 or EVC-4000.



Cartridge Electrodes

The miniature electrode-gel cartridge is a small plastic tube with a male luer tip identical to those at the tip of hypodermic syringes. The tube may be filled with different gel materials; agar is commonly used but other gel materials may also be satisfactory.


CHM2

СНМ5

СНМ4

СНМ7

- 30 mm



Assembled chambers are 101.6 mm (4 in.) long.

	CHM1 (Medium)	CHM2 (Small)	CHM3 (Large)	CHM4 (Extra Small)	CHM5 (Snap)	CHM6 (Rect., Small)	CHM7 (Rect., Large)	CHM8 (Extra Small)
Reservoir Opening	12 mm	9 mm	13.5 mm	4 mm	12 mm	5 x 14.5 mm	7 x 30 mm	4 mm
Half-Chamber Volume	1.0 mL	0.75 mL	1.2 mL	0.5 mL	1.7 mL	0.8 mL	5.5 mL	0.5 mL
Pin Circle Diameter	17 mm	12 mm	18.5 mm	6 mm*	6 mm* N/A 7		9 x 32 mm	5.5 mm
				*O-ring diam.				



Circulation Reservoirs

Hand-blown borosilicate glass, with jacketed chambers for temperature control. Available in two sizes — #5210 holds 20 mL per side, and #5362 (at left) holds 10 mL per side (useful when expensive chemicals are involved). Reservoir condenser caps prevent air bubbles and turbulence in fluid reservoirs.

JSSING SYSTEMS, LARGE RESERVOIR

EVC4000-4

USS1L	Medium Chamber, Stand, Reservoir, Electrodes, Tubing						
USS2L	Small Chamber, Stand, Reservoir, Electrodes, Tubing						
USS3L	Large Chamber, Stand, Reservoir, Electrodes, Tubing						
USS4L	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing						
USS5L	Snap Chamber, Stand, Reservoir, Electrodes, Tubing						
USS6L	Small Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing						
USS7L	Large Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing						
USS8L	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing						
USSING SYSTEMS, SMALL RESERVOIR							
USS1S	Medium Chamber, Stand, Reservoir, Electrodes, Tubing						
USS2S	Small Chamber, Stand, Reservoir, Electrodes, Tubing						
USS3S	Large Chamber, Stand, Reservoir, Electrodes, Tubing						
USS4S	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing						
USS5S	Snap Chamber, Stand, Reservoir, Electrodes, Tubing						
USS6S	Small Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing						
USS7S	Large Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing						
USS8S	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing						
*	* Add EVC4000 at reduced price when buying Ussing System with equivalent number of channels						
EVC4000-1	1-Channel Voltage Clamp & Preamps						
EVC4000-2	2-Channel Voltage Clamp & Preamps						
EVC4000-3	3-Channel Voltage Clamp & Preamps						

4-Channel Voltage Clamp & Preamps

Water Bath

The Julabo circulating bath (see page 85) is ideal for controlling temperatures of external systems. With a powerful 15L/min flow rate, the pump provides optimum heat exchange. The tap water cooling feature is standard with a range of 20-100°C. The bath opening is 15cm x 15cm x 15cm and can hold 3-4.5L of liquid.

System components also available separately:

System com	ponents also avallable separately:
xxxxD	Drain option (add "D" to part number of chamber or system)
CHM1	Medium Chamber
CHM2	Small Chamber
СНМ3	Large Chamber
CHM4	Extra Small Chamber with O-Ring Seal
CHM5	Snap Chamber (fits Costar Snapwell cups)
СНМ6	Small Rectangular Chamber
СНМ7	Large Rectangular Chamber
CHM8	Extra Small Chamber with Mounting Pins
EK1	Ussing Electrode Kit (2 voltage, 2 current)
EKC	Extra Ussing Current Electrode (red) (each)
EKV	Extra Ussing Voltage Electrode (blue) (each)
5210	Large Glass Circulation Reservoir, (20 mL per side)
5233	Replacement Condenser for 5210
5362	Small Glass Circulation Reservoir, (10 mL per side)
5361	Replacement Condenser for 5362
3955	EKV Cartridges, 35 mm (pkg of 12)
3960	EKC Cartridges, 58 mm (pkg of 12)
3669	Tubing Kit (flexible hose and luer fittings)
3579-20	Replacement luer fittings for tubing connections (pkg of 20)
5153	Support Stand
3485	Post Mounting Kit for Preamp
J. F	

Julobo

^{*} Drain option, \$105 to \$140, varies according to chamber.

Multi-Channel Voltage / Current Clamp



More channels and a wider range of voltage clamp commands than WPI's classic DVC-1000. The superior design of the cartridge electrodes makes 100-volt current excursion unnecessary, so this safe, low-voltage system is easier to adjust and use.

EVC4000 employs the voltage clamp technique to monitor membrane permeability as a function of membrane voltage or applied chemicals. When combined with WPI's patented EKC and EKV cartridge electrodes, EVC4000 can efficiently voltage or current clamp up to four sample membranes simultaneously using safe moderate voltages on the current wire leads. The superior design of the cartridge electrodes makes 100-volt current excursion unnecessary, so this safe, lowvoltage system is easier to adjust and use. Extremely stable and accurate, each module, with its companion preamplifier, can operate independently in one of three different modes: Voltage Clamp (VC), Current Clamp (CC), or Open Circuit Potential (PD) measurement. EVC4000 can be controlled from the front panel of the instrument or from computer generated commands applied at the rear panel of the instrument. A feature unique to EVC4000 is an electronic potentiostat in the preamplifier box that maintains the serosal electrode invariant potential at zero relative to system ground. The preamplifier apparatus actively maintains one surface of the test membrane close to ground potential under all operating conditions.

References

W. K. MacNaughton (2000) "Role of constitutive cyclooxygenase-2 in prostaglandin-dependent secretion in mouse colon in vitro." Journal of Pharmacology and experimental Therapeutics 293, 2, 539-544

EVC4000-4	4-Channel Voltage Clamp & preamps (shown above)
EVC4000-3	3-Channel Voltage Clamp & preamps
EVC4000-2	2-Channel Voltage Clamp & preamps
EVC4000-1	1-Channel Voltage Clamp & preamp

Specify line voltage

EVC4000 SPECIFICATIONS

가	<	L	A	N	1	۲	L	I۲	1	E	h	ć		
								_						

Input Resistance Input Leakage Current Maximum Input Voltage

VOI TAGE CLAMP

Panel Display Clamp Voltage / External Input

Range of Voltage Electrodes Max. Clamp Voltage

Fluid Resistance Compensation

CURRENT CLAMP

Panel Display

Maximum Clamp Current Current Clamp Output

DISPLAY RESOLUTION

Voltage Current

DIMENSIONS

SHIPPING WEIGHT (EVC4000-4)

1012 Ohms 100 pA, max.

±15 volts

±200 mV ±0.1 mV 100 mV per Volt

±32 Volts ±100 mV 0 to 1000 Ohms

±999 μA ±1 μA ±1 milliampere 1 uA / mV

0.1 mV 1 uA

> 18.25 x 7.2 x 9.6 in. (46 x 18 x 24 cm)

26 lb (12 kg)

OPTIONAL ACCESSORIES

SYS-EVC4000	Replacement Voltage Clamp & EVC3 Preamplifier
EVC3	Replacement Preamplifier Module
EK1	Ussing Electrode Kit (2 voltage, 2 current)
EKV	Extra Ussing Voltage Electrode (each)
EKC	Extra Ussing Current Electrode (each)
2851	BNC Cable
3485	Post Mounting Kit for Preamp

Circulating Bath

The Julabo circulating bath is ideal for controlling temperatures of external systems. With a powerful 15L/min flow rate, the pump provides optimum heat exchange. The tap water cooling feature is standard with a range of 20-100°C. The bath opening is 15cm x 15cm x 15cm and can hold 3-4.5L of liquid.

- LED temperature display (0.1°C resolution)
- Stainless steel bath tank
- Adjustable high temperature cut out and dry running protection.
- PID temperature control
- Large capacity for temperature applications with larger external systems and open systems
- Internal bath for simultaneous applications with smaller objects
- Built in cooling coil for tap water connection when you require a temperature less than the ambient temperature




JULABO BATH SPECIFICATIONS
Temperature Selection digital
Temperature Indication LED
Resolution
Temperature Control PID1
Heater Wattage 230V–2.0kW
115V–1.0kW
Circulating Pump
Discharge, max a 0 bar15L/min
Pressure, max. at 0 I 0.35 bar
Ambient Temperature5-40°C
Mains Power Connection (230V/50Hz) 190-253V/50Hz
Current Input (230V)
Mains Power Connection (115V/60Hz) 103-127V/60Hz
Current Input (115V)
Bath Tank Stainless steel
Working Temperature Range 20-100°C
Temperature Stability ±0.03
Bath Opening (WxL)
Bath Depth15cm
Filling Volume
Dimensions
Weight
Recommended Bath Fluid soft/decalcified water

IIII ADO DATH CDECIEICATIONS

503843	Julabo Circulating Bath, 115V
503844	Julabo Circulating Bath, 230V
504122	Julabo Circulating Bath, 100V (Japan)
504142	Julabo Circulating Bath, 13L volume, 110V
504141	Julabo Circulating Bath, 13L volume, 220V

Rodent Brain Matrices





WPI offers one of the largest selections of brain matrices available. Made of acrylic, or stainless steel, these matrices are sturdy and can be heated, chilled, autoclaved (stainless steel only), scrubbed — and stand up to rigorous daily use. Coronal matrices have the additional feature of a mid-line sagittal cut to facilitate splitting of the left and right hemispheres. Sections can be as fine as 1-mm. The olfactory/ spinal/notch is cut into each matrix.

Order #	Subject	Material	Section	A	В	С	D	E	Cavity Depth	Weight
RBMA-200C	Adult Mouse, 40-75g	Acrylic	Coronal	3.18	11.1	8.73	19.1	12.2	7.4	0.5 lb
RBMA-200S	Adult Mouse, 40-75g	Acrylic	Sagittal	3.18	11.1	8.73	19.1	12.2	7.4	0.5 lb
RBMA-300C	Rat, 175-300g	Acrylic	Coronal	4.7	15.9	12.7	36.6	23.8	7.61	0.5 lb
RBMA-300S	Rat, 175-300g	Acrylic	Sagittal	4.76	15.9	12.7	36.6	23.8	10.91	0.5 lb
RBMA-600C	Rat, 300g-600g	Acrylic	Coronal	4.76	19.8	14.7	36.6	24.7	10.91	0.5 lb
RBMA-600S	Rat, 300g-600g	Acrylic	Sagittal	4.76	19.8	14.7	36.6	24.7	10.91	0.5 lb
RBMS-200C	Adult Mouse	Stainless Steel	Coronal	3.18	11.1	8.73	19.1	12.2	7.4	1.0 lb
RBMS-200S	Adult Mouse	Stainless Steel	Sagittal	3.18	11.1	8.73	19.1	12.2	7.4	1.0 lb
RBMS-300C	Rat, 175-300g	Stainless Steel	Coronal	4.76	15.9	12.7	36.6	23.8	7.61	1.0 lb.
RBMS-300S	Rat, 175-300g	Stainless Steel	Sagittal	4.76	15.9	12.7	36.6	23.8	7.61	1.0 lb
RBMS-600C	Rat, 300g-600g	Stainless Steel	Coronal	4.76	19.8	14.7	36.6	24.7	10.91	1.0 lb
RBMS-600S	Rat, 300g-600g	Stainless Steel	Sagittal	4.76	19.8	14.7	36.6	24.7	10.91	1.0 lb



Combining the accuracy of digital electronics with the convenience of analog controls

A pulse generator/stimulator combining the reproducibility and accuracy of digital electronics with the fine resolution and continuous adjustment possible with analog circuitry. All timing parameters are entered with tenturn readable potentiometers and six-position range switches. Outputs are accurate to within 1% of the set value.

Pulses can be created in continuous run, single-shot, or train/burst modes. Duration of the train/burst is easily controlled using the onboard envelope generator or by using either of two external gating inputs. Used in conjunction with the A360, A365, A385, or A395, bipolar pulses or trains may be easily produced. Output stimulus can be fed through the Duo 773 for iontophoresis. Footswitch allows hands-free operation.

Three separate outputs are available on the front panel. A Monitor output provides 10-15 V signals (up to 50 mA) for viewing the output on an oscilloscope or for controlling other devices. The stimulator's signal, simultaneously available at the Isolator output, is sufficient to drive any WPI A300 Series stimulus isolator (A360, A365, or A385) and is also TTL and CMOS compatible. The Variable output can provide signals varying between ±10 V with a resolution of 1 mV. Separate variable outputs are provided for positive and negative signals.

Accupulser™ Signal Generator SYS-A310

Specify line voltage

OPTIONAL ACCESSORIES

3259 Footswitch for A310 2933 Rack Mount Kit, 51/4 in. high

A310 ACCUPULSER™ SPECIFICATIONS

TIMING PARAMETERS

EVENT INTERVAL 100 μs to1000 s* EVENT DELAY 10 µs to 100 s * **PULSE WIDTH** 10 μs to 100 s * TRAIN DURATION (ENVELOPE) 100 µs to 1000 s* PULSE INTERVAL 20 µs to 100 s*

OUTPUTS

SYNC 5 μs, TTL, and 5 V CMOS compatible,

20 mA max.

MONITOR 10-15 V, 50 mA max.

TTL & 5 V CMOS compatible, 20mA max. ISOI ATOR

VARIABLE (Pos or Neg)

PUI SED/DC LOW RANGE HIGH RANGE 0 to ± 1 V 0 to ±10 V Range Resolution 1 mV 10 mV

NOISE

Pulsed at 100 kHz bandwidth <500 µV DC Wide Band <500 uV OUTPUT IMPEDANCE <1.O

INPUTS

EXTERNAL SYNC Accepts 1-us minimum pulses

TTL, CMOS compatible

EXTERNAL GATE Accepts 1-µs pulse to continuous

TTL, CMOS compatible

95-130 V or 190-260 V, switch selectable **POWER**

single phase, 50/60 Hz

DIMENSIONS 17 x 5.25 x 10 in. (43 x 13 x 25 cm)

SHIPPING WEIGHT 14 lb (6.4 kg)

*Continuously variable in six ranges. All accuracies better than 1% of set value. 50kHz maximum pulse frequency.

A300 Pulsemaster™ **Multi-Channel Stimulator**



The Pulsemaster™ (Model A300) is WPI's third generation, multichannel, pulse/train generator/stimulator that combines the superb accuracy of digital electronics with the "you-seewhat-you-get" displays only available on single-channel products. In one compact rack mountable enclosure, the Pulsemaster contains an event interval generator, five pulse train channels, two mixing channels and a very quiet variable voltage output channel. System timing is accurate to 100 ppm; output timing is continuously variable in 0.1% of full scale increments over a range of eight orders of magnitude. Bright, threedigit LED displays continuously and simultaneously show all the variable timing parameters.

An integrated five-channel pulse generator/stimulator including one interval generator, five pulse or train channels, two mixer channels and one very quiet variable voltage output stage

The Pulsemaster is designed for ease of use and flexibility. Each channel can be operated synchronized with the onboard event interval generator, triggered manually from any other channel or external source, and as an independent asynchronous pulse generator. Except for the external source, all channel interconnections are accomplished on the panel, without the use of cables. The output from each channel is compatible with standard digital circuitry and is also designed to drive WPI's A300 series stimulus isolators. If desired, any channel's output may be internally connected to the variable channel, whose amplitude can be continuously adjusted from millivolts to ten volts.

Pulsemaster™ Multi-Channel Stimulator **SYS-A300**

Specify line voltage

A300 PULSEMASTER SPECIFICATIONS

EVENT INTERVAL CHANNEL

Operating Modes EXTernal SYNC, SINGLE EVENT, CONTINUOUS ON

EXT SYNC accepts ≥ 1µs pulses; TTL, CMOS, Input

RS232C compatible

Timing EVENT INTERVAL 10 µs to 999 s (100 kHz - 0.001 Hz),

±0.1% of full scale, continuously variable in 0.1% of full scale increments, through three orders of magnitude,

in six ranges

SYNC OUT pulse of ≈6 µs, TTL, 5 V CMOS compatible Output

PULSE TRAIN CHANNEL (5 provided)

Operating Modes EXTernal SYNC, SELF SYNC, manual SINGLE event,

sync from Event Interval, sync from any of other four

Pulse Trains, sync from one of the MIXers, off,

TRAIN/PULSE

EXT SYNC accepts ≥ 1µs pulses; TTL, CMOS Input

Timing DELAY and WIDTH 10 µs to 999 s, ±0.1% of full scale

Output OUTPUT PULSE/TRAIN of preset timing, TTL, 5 V

CMOS compatible, 4 mA sink and source

MIXER CHANNEL (2 provided)

Any combination of an EXTernal pulse, the outputs of Inputs

the five Pulse Train channels, and DC continuous ON/ DC MOMentary EXT INPUT accepts ≥ 1µs pulses

OUTPUT, TTL, 5V CMOS compatible, 4 mA sink and Output

source

VARIABLE CHANNEL

Output from any one PULSE TRAIN channel Inputs

or one of the two MIXER channels or DC

Output 0 to +1 V low range, 1 mV resolution

0 to +10 V high range, 10 mV resolution

5 mA max sink and source

Output Impedance <1 ohm

<500 μV peak @ 100 kHz bandwidth, PULSED mode Noise

<500 μV, wide band, DC mode

Signal Ground Floating, i.e., not connected to chassis

POWER 95-135 V or 220-240 V, 50/60 Hz

BATTERIES Three 1.2 V DC, size AA, NiMH batteries

DIMENSIONS 8.5 x 19 x 8.75 in. (22 x 45 x 22 cm)

SHIPPING WEIGHT 21 lb (9.5 kg)

89

Isostim™ Stimulator/Isolator

Combining the ease of use and accuracy of WPI's 300 Series stimulators with the power output of a stimulus isolator



Timing

Pulse interval and width are set with single-turn continuously variable

controls from 5 ms to 5.5 s in three ranges.

Pulse width is continuously variable from 50 µs to 550 ms in four ranges.

Modes of operation

In FREE RUN, Isostim™ generates continuous square waves. In EXT GATE or EXT SYNC modes, externally applied pulses can generate trains or single events. Single pulses of finite duration can be produced using a pushbutton on the instrument's front panel. EXT/DC mode converts Isostim to a passive stimulus isolator.

Dual tone audible alarm

A tone sounds when an open circuit is detected or when system compliance is reached. A second tone, which sounds when a signal is applied to the input, can only be heard if the batteries have sufficient charge to operate the isolator. A violation light advises when pulse width exceeds the interval.

ISOSTIM™ SPECIFICATIONS

TIMING PARAMETERS

Interval

Pulse width

INPUT

External sync External gate External command voltage

Trigger threshold OUTPUT

Current ranges

Waveform

Load voltage excursion (compliance) 100 V nom., 150 V max. Output polarity Current rise time and delay Current fall time and delay

Optocoupler

POWER

Dry Cell (Version D) Rechargeable (Version R) 16 alkaline 9V batteries included 16 rechargeable NiMH 9V batteries

(included)

DIMENSIONS 8.5 x 3.5 x 4.9 in (22 x 9 x 12 cm) SHIPPING WEIGHT 4 lb (1.8 kg)

5 ms to 5.5 s continuously variable in three ranges (0.18 to 200 Hz)

50 µs to 550 ms continuously variable in four ranges

Accepts 1 µs minimum pulses

Accepts 1 µs pulse to continuous 5.0 V at 3.0 mA (TTL level), 10 V max. 2.0 V at 0.5 mA

DC, pulse from internal timing or externally generated pulse

0-1 mA, 0-10 mA Reversible, manual switch

8 μ s, typical (1 K Ω load) 10 μ s, typical (1 K Ω load) 10¹² Ohms

Leakage resistance, output to ground 2500 V rated min. breakdown voltage

Current delivery

Stimulus currents up to 10 mA can be set on the front panel with a control knob and a two-position range switch. Output current is load-independent.

Isostim model A320D is powered by readily obtainable 9-volt alkaline batteries (included). Under average use these will last several months before replacement is required. The rechargeable A320R is supplied with a nickel metal hydride battery stack which provides 10-12 hours of operation before recharge is required. The A362 Battery Charger must be used with the A320R.



A362 Battery Charger

Required for A320R, A365R and A395R

Recharges the high-voltage nickel-cadmium or NiMH battery stack in the A320R, A365R or A395R. LED lamp indicates charging status. Full charge overnight. Dimensions: 2.8 x 4.1 x 5 in. (7 x 10 x 13 cm). Shipping weight: 4 lb (1.8 kg).

SYS-A362	Battery Charger for A320R, A365R, A395R
A320RC	A320R with Charger (A362)
SYS-A320D	Isostim™ Stimulator/Isolator
SYS-A320R	lsostim™ Stimulator/Isolator (rechargeable)

Specify line voltage

OPTIONAL ACCESSORIES

• · · · · · · · · · · · · · · · · · · ·	O. 1101012 ACC2220 ACC22				
DRL	Dummy Load Resistor Kit (set of 3)				
13347	BNC-to-Double Banana Adapter				

Stimulus Isolator / Precision Current Source



Activated by conventional logic-level commands, Model A365 can be gated by any pulse generator, stimulator, or computer output; automated bipolar pulsing for zero net charge on biological preparations.

Dual tone audible alarms — A tone sounds when an open electrode circuit is detected or when system compliance is reached. A second optional tone sounds when a signal is applied to the input. A test switch is also provided to check battery charge.

Current delivery — Stimulus currents are set using a three-digit control

knob and a three-position range switch. Output current tracks control settings to better than 1%. Output current is load independent; voltage sufficient to push the desired current through the load is automatically developed, subject only to compliance limits. Model A360LA produces up to 10 milliampere current, in three ranges, at more than 100 volts compliance.

Polarity — Output polarity is determined by a push switch on the front panel. Bipolar current is toggled by the command waveform, setting alternating pulses as positive or negative.

SYS-A365D	High Voltage Isolator, Bipolar, alkaline batteries
A365RC	A365R with charger (A362)
SYS-A365R	High Voltage Isolator, Bipolar, rechargeable
SYS-A362	Battery Charger for A320R, A365R, A395R
	Specify line voltage
OPTIONAL A	CCESSORIES
DRL	Dummy Load Resistor Kit (set of 3)

BNC-to-Double Banana Adapter



DRL — **Dummy Load Resistor Kit** Converts current output to precise voltages.

A365 SPECIFICATIONS

OUTPUT WAVEFORM DC or current pulse **OUTPUT CURRENT RANGES** 0.1, 1.0, and 10 mA **CURRENT AMPLITUDE ERROR** 0.5% of full scale, max. 0.1% of full scale, typical CURRENT RESOLUTION **OUTPUT LOAD VOLTAGE EXCURSION (COMPLIANCE)** 100 V EXTERNAL COMMAND VOLTAGE 5.0 V at 3.0 mA (TTL level), 10 V max. TRIGGER THRESHOLD 2.0 V at 0.5 mA Reversible, manual switch or automatic OUTPUT POLARITY **CURRENT RISE TIME & DELAY** 6 μ s, typical (1 K Ω load) **CURRENT FALL TIME & DELAY** 10 μ s, typical (1 K Ω load) **OUTPUT TO GROUND RESISTANCE OPTOCOUPLER** 2500 V, rated min. breakdown voltage

Model A365D (dry cell)

Model A365R (rechargeable) incl.

DIMENSIONS SHIPPING WEIGHT

POWFR

16 alkaline 9 V batteries, included 16 rechargeable NiMH 9 V batteries

 $8.5 \times 3.5 \times 5$ in $(22 \times 9 \times 12$ cm)

4 lb (1.8 kg)

91

13347



A385 High **Current Stimulus** Isolator

Delivers positive, negative, or bipolar currents. For bipolar delivery, polarity of the output is toggled to the opposite state with each pulse presented to the input. Pulse duration is controlled by an externally applied voltage. Input connector is a standard BNC, allowing signals from any source — such as computer D/A or I/O lines — to be used.

Output amplitude is set on a 3-digit, ten-turn dial as a percentage of the range selected: for example, a setting of 45.6 in the 0-10 mA range translates to 4.56 mA at the output. Accuracy and repeatability are excellent. Designed for subcutaneous stimulation, maximum output voltage at the stimulating electrodes is 36 volts, reducing the possibility of serious accidental transcutaneous shocks. A compliance/output alarm sounds

when the 36-volt limit is reached. Internal circuitry maintains electrodes short-circuited during inactive periods ("electrode exhauster" feature). A385 is not appropriate for transcutaneous stimulation.

The 1.2 amp-hour rating of the six heavy-duty lead-acid rechargeable batteries ensures that experiments will not be interrupted by dead batteries — even at peak currents. Indicator lights and audible alarms keep the user constantly apprised of battery charge status. These batteries must be recharged by the A382 System Charger designed especially for the A385.

A385RC	A385R with A382 Charger
SYS-A385R	High Current Isolator, rechargeable
SYS-A382	Battery Charger for A385 (see below)

Specify line voltage

A385 SPECIFICATIONS

OUTPUT WAVEFORM **OUTPUT CURRENT RANGES CURRENT AMPLITUDE ERROR** CURRENT RESOLUTION REPEATABILITY

OUTPUT LOAD VOLTAGE EXCURSION (COMPLIANCE)

EXTERNAL COMMAND VOLTAGE TRIGGER THRESHOLD

OUTPUT POLARITY

CURRENT RISE TIME AND DELAY CURRENT FALL TIME AND DELAY **OUTPUT TO GROUND RESISTANCE** OPTOCOUPLER

POWER

DIMENSIONS SHIPPING WEIGHT DC or current pulse 1, 10, and 100 mA 0.5% of full scale, max

0.1% of full scale, typical

5.0 V at 3.0 mA (TTL level), 10 V max.

2.0 V at 0.5 mA

Reversible, manual switch, or electronically switched bipolar delivery

6 μ s, typical (1 $K\Omega$ load)

10 μ s, typical (1 $K\Omega$ load)

10¹² O

2500 V, rated minimum breakdown

Six rechargeable lead-acid batteries (Requires companion charger A382)


8.5 x 3.5 x 5 in. (22 x 9 x 12 cm)

5 lb (2.3 kg)

A382 Battery Charger

An innovative threestep charger, A382

employs fast, medium, and



A382 SPECIFICATIONS

POWER 95-135 V or 220-240 V, 50/60 Hz **DIMENSIONS** $8.5 \times 3.5 \times 5$ in $(22 \times 9 \times 12$ cm) SHIPPING WEIGHT 5 lb (2.3 kg)

A395 Linear Stimulus Isolator


Replicates a programmed waveform of any shape or polarity



All WPI stimulus isolators are designed to supply constant current because current threshold (not voltage) is the most quantitatively reproducible parameter for stimulation of nerve and muscle. Model A395 dispenses current reproducibly from its Output terminals; the amplitude being determined by the selected current RANGE and the input voltage. Current amplitude is "constant", that is, load resistance independent, provided that the I x R (load) product does not exceed the available battery supply voltage. A visual indicator (the compliance LEDs) displays if I x R reaches this limit. When the unit is out of compliance, one of the two LEDs (labeled - and +) illuminate, depending in which direction the current is flowing. Model A395 D can generate a voltage of 70 volts or more across its OUTPUT terminals. Thus, the user can be sure that the amplitude of the current will be as dialed as long as the voltage drop across the load (stimulus electrode path) does not reach the magnitude of the supply voltage. The compliance LEDs will then be visible. The user would then know that (a) too much current was dialed for a given load or (b) inter-electrode resistance was too high or the electrode circuit path was open.

Model A395 generates an output current of arbitrary (user-defined) wave shape; DC, AC, pulse, and combinations thereof. Battery operated, and photoelectrically-isolated from the input voltage drive, the instrument regenerates output currents which are linearly proportional to the analog voltage waveforms provided by your D/A converter or signal generator (see diagram below).

The A395 is ideally suited for data acquisition and stimulator generators. It can be easily daisy-chained for mutiple channel requirements.



Current Delivery —A 10 V input produces the maximum output current for the current range selected, i.e., 100 µA, 1 mA, or 10 mA. Front panel controls allow DC current to be generated. Externally applied signals can be superimposed simultaneously (DC offset). Warning lamps indicate open circuit or excessive current conditions.

Digital Meter — Measures DC or average output current.

Overload Lamps — Indicate when output voltage has reached positive or negative compliance voltage limit.

A395RC	A395R with Charger (A362)
SYS-A395D	Linear Stimulus Isolator
SYS-A395R	Linear Stimulus Isolator, Rechargeable
SYS-A362	Battery Charger

Specify line voltage

A395 SPECIFICATIONS

OUTPUT CURRENT, Imax 3 ranges: 100µA, 1 mA, and 10 mA

OUTPUT VOLTAGE RANGE ± 70 V

OUTPUT BANDWIDTH 10 kHz (measured across 1K Ω load

INPUT RESISTANCE $> 20 \text{ M}\Omega$ INPUT VOLTAGE @ Imax ± 10 volts INPUT/OUTPUT LINEARITY ERROR < 0.5% RISE, FALL TIME 26 μs @ 10 ΚΩ

POWER

17 alkaline 9 V batteries Model A395D Model A395R 17 rechargeable NiMH 9 V batteries 6.5 x 4 x 3.5 in. (16 x 10 x 9 cm) DIMENSIONS

SHIPPING WEIGHT 4 lb (1.8 kg)

A362 Battery Charger

Required for A320R, A365R and A395R

Recharges the high-voltage nickel-cadmium or NiMH

indicates charging status. Full

charge overnight.

Dimensions: 2.8 x 4.1 x 5 in. (7

x 10.5 x 12.7 cm).

Shipping weight: 4 lb (1.8 kg).



Adhesives Application Guide

Part No.	Description	Curing time	Useful Applications and Characteristics
Epoxies	Form strong bonding. Used in wire bonding	ng applications.	
4898	Silver filled conductive Epoxy	12 hrs@50C; 5 min @150C	Connecting conductors that can't be soldered. Constructing or connecting Ag/AgCl pellets.
7335	Carbon filled conductive Epoxy	48 hrs@25C; 5 min @150C	Constructing carbon electrode.
4886	High performance Structural Epoxy	12 hrs@25C.	Forms a strong and slightly flexible bond on plastic, metal, & glass. Bonds some low surface.
Hot melt (EVA)	Easy to use for bonding, needs large gap f	illing	
13316	Mini Glue Gun with glue sticks	As soon as it cools down	Bonds wood, glass, metals, and many plastics.
Silicone Adhesiv	es/Sealants/Primers	Good moisture resistant and ela	astic. Low toxic.
1571	Room temperature vulcanizing (RTV) adhesive. Acyloxy/moisture cure system. Acetic acid is cure by-product.	24hrs@25C	Has the best adhesion property in this silicone family. Will bond to many materials.
7128	RTV sealant. Alkoxy/Moisture cure system. Methanol as cure by-product.	72hrs@25C	Good for bonding or sealing electronics circuits (metal).
SYLG184	Sylgard, Two parts, vinyl/platinum cure sealant. Hydrogen as cure by-products. Very low toxic	24hrs@25C, 15 min.@150C	Coating Patch Clamp electrodes, Cell culture dish, making dissection pads.
Kwik-Sil	Two part, adhesive. Vinyl/platinum system, Hydrogen as cure by-products. Very low toxic.	< 5 min@25C	Live tissue and nerve studies. Medium strength adhesion.
Kwik-Cast	Two part sealant. Vinyl/platinum cure system. Hydrogen as cure by-products. Very low toxic.	< 5 min@25C	Sealant for live tissues. Embedding peripheral nerves with electrodes.
6820	Primer for silicone	N/A	Enhances adhesion of silicone adhesives for difficult to bond plastic surfaces
Cyanoacrylate	Forms an instantaneous bonding.		
7341	Ethyl Cyanoacrylate, low viscosity 90-120 cps	<10 seconds	Mounting rat/mouse brain slices. Ideal for relatively small gaps and smooth surfaces.Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.
7342	Ethyl Cyanoacrylate, high viscosity 1100-1600 cps	<30 seconds	Use on brain slice exp. Ideal for larger gaps, allows slightly longer bonding time. Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.
Vetbond	Butyl Cyanoacrylate, Low toxic	<10 seconds	Bonds tissues, alternative to suture, helps small wound healing. Antimicrobial effect. Used in forensic science.
503763	Octyl Cyanoacrylate, Low toxic	<15 seconds	Suitable for surface wound bonding, protection, holding a sensor or other device on the tissue.

S Kwik-Gard™

Kwik-Gard is specially packaged Sylgard 184 silicone for quicker and easier application, eliminating the messy procedure of preparing the mixture before application. Its special cartridge controls the precise mixing ratio to ensure proper curing. The disposable tip

mixes resin and hardener as they are dispensed. Since no air is introduced during mixing, the resin does not need degassing for most applications. The mixed silicone is applied directly to the site, reducing preparation time and material waste.

Each Kwik-Gard cartridge contains 37 mL of resin and hardener. The dispensing tip has a dead volume of 0.75 mL.



KWIKGARD	Kwik-Gard Start-up Kit
(incl.	dispenser, 1 cartridge, 5 tips)
KWIKGLUE	Kwik-Gard Refill (2 cartridges,
10 d	ispensing tips)
KWIKMIX	Dispensing Tips (pkg of 10)
KWIKGUN	Kwik-Gard Dispenser

Scotch-Weld 2216 **Structural Epoxy**



Probably still the best epoxy for bonding plastic, often used as the benchmark for testing the binding strength of other adhesives. The slightly rubbery texture also makes it less easy to break off. It is the only epoxy known that can bond PEEK. Color: gray. Cures at room temperature.

Shipping weight: 1 lb. (0.5 kg)

Scotch-Weld 2216 (2 oz.)

Low toxicity 5-minute adhesives for live tissues!

- Specially formulated for WPI
- Extremely low toxicity
- Excellent moisture resistance
- Cures at room temperature

Kwik-Sil and **Kwik-Cast** silicones are based on technology with vinyl terminated siloxane and platinum complex catalysts. In order to gain enough speed to cure at room temperature, special cross linkers and high catalyst concentration is used. Although the high concentration of catalysts makes these products more costly than traditional RTV silicones and less attractive for general usage, they provide an excellent value in applications for the biological research field.

Both **Kwik-Sil** and **Kwik-Cast** have very low toxicity before, during and after curing. In traditional RTV silicone systems, a by-product of the condensation (curing) is either acetic acid or alcohol, which are toxic to living cells. In vinyl systems, the by-product of condensation is a small amount of hydrogen gas, which is much less toxic to the cell.


Kwik-Sil and **Kwik-Cast** curing speed is hundreds of times faster than traditional RTV silicones. A curing time of a few minutes at room temperature is especially useful for encapsulation of live tissue or implanting into a live animal.

Unlike many vinyl-based silicones in which the platinum complex catalysts are easily poisoned by contamination from amines and animal tissue, **Kwik-Sil** and **Kwik-Cast** are not sensitive to contamination from animal tissue.



Kwik-Sil™ is a translucent, medium-viscosity silicone adhesive, developed for chronic peripheral nerve studies such as anterograde tracing with fluorescent indicators or electrode recording. Good adhesion and mechanical properties (tear strength and elongation) allow days of study without breaking of the bonding. Curing speed is very reproducible.

Kwik-Cast[™] is a very low viscosity silicone sealant developed to embed peripheral nerves with electrodes for acute multifiber recordings. It flows easily, filling the small spaces around the nerve and leaving no channels through which peritoneal



fluid can travel and thus short the nerve/electrode contact. Equally important is the ability of the material to flow into itself and create one continuous mass from underneath the nerve all the way to the top of the nerve/electrode contact to ensure long-term recording stability. **Kwik-Cast** is color-coded to make the mixing foolproof. The catalyst is yellow and the base is blue. When uniformly mixed, it is green. **Kwik-Cast** can be applied and cured underneath mineral oil. After recording, electrodes are easily recovered due to the low tear strength.

KWIK-CAST & KWIK-SIL SPECIFICATIONS

		Kwik-Sil	Kwik-Cast
	Mix Ratio	1 to 1	1 to 1
	Working time	< 5 minutes*	4 minutes
	Setting time (room temp., 1:1 ratio)	5-10 minutes**	<10 minutes
	Cure time	~15minutes	
	Viscosity, cps	15,000	10,000
	Shelf life at 23 °C	1 year	1 year
	Volume	5 mL	5 mL
	Number of mixing tip	10	10
	Dead volume of the mixing tip	<0.12 mL	<0.12 mL
A	AFTER CURING 24 HOURS:		
	Tear Strength, ppi	90	44
	Elongation %	650	60
	Durometer (shore A-2)	30	36
	Color	translucent	green
	Volume Resistivity, W/cm	1x10 ¹⁵	1x10 ¹⁵

- * 3 minutes average with about 90 seconds of liquidity
- ** no longer mixable at this point

KWIK-SIL	Silicone Adhesive Compound (two 5-mL syringes)
KWIK-CAST	Silicone Casting Compound (two 5-mL syringes)
600009	Replacement KWIK-CAST Mixing Tips (pkg of 10)
600022	Replacement KWIK-SIL Mixing Tips (pkg of 10)

PRICE BREAK

PRICE BREAK								
	1-5 pkg	6-9 pkg	10 or More					
KWIK-CAST								
KWIK-SIL								

"Super" Adhesives for Life Science Research

Cyanoacrylate adhesives have been on the market since 1958. Most industrial or household grade cyanoacrylate is made of shorter alkyl chain derivatives such as methyl or ethyl cyanoacrylate (WPI's #7341 and **#7342**). They are very useful for temporarily holding tissues such as mounting specimens for microtome sectioning. However, they are not suitable for bonding wounds on live animals. The difficulties of using cyanoacrylate for bonding live animals are: (1) a strong, irritating odor;



(2) quick loss of bonding strength due to breakdown of the bonding by hydration; (3) the breakdown products, cyanoacetate and formaldehyde, are toxic and can cause inflammatory reactions; and (4) they have low flexibility and tend to be brittle.

To overcome these problems, several longer alkyl chain cyanoacrylates have been developed especially for veterinary and human use. The first longer alkyl chain product is butyl cyanoacrylate. This product has been used for animal and human applications outside the USA since 1970. It is much less toxic and has a lower odor than the methyl or ethyl cyanoacrylate. The butyl cyanoacrylate offered by WPI is **Vetbond**™.

A family of adhesives containing

octyl cyanoacrylate, a plasticizer and stabilizer, was developed In the 1990's (one of them approved by FDA for human use). When bonding to tissue, these new adhesives are four times stronger and less toxic than butyl cyanoacrylate. Compared with the

traditional suture, the new super adhesive has several advantages. On average, it takes only one-tenth of the time to close an incision. The **bonding strength** is equal to 5-0 monofilament suture. It also has a mysterious antimicrobial effect that can decrease infection rates in contaminated wounds. Bonding will slough off naturally in 5 to 7 days. Cosmetic appearance of the healed incision is also better.



Gluture Topical Tissue

Adhesive #503763 forms a strong and flexible film and is thus more suitable for surface wound bonding, protection, and holding a sensor or other device on the tissue. Setting time is about 10 seconds, which gives ample time for application. It can also be used for temporarily holding a live tissue. For example, there is a report of using it to hold nematodes on a glass slide for patch-clamp neurons recording.

All of the products offered by WPI are veterinary grade (not suitable for human application). Though very similar to the grade for human use, they are not sterile and do not have FDA approval.

503763	Gluture Topical Tissue Adhesive (10 tips), 1.5 mL
7341	Cyanoacrylate Adhesive, low viscosity—90-120 cps (package of 10 vials, each approx 1.5mL)
7342	Cyanoacrylate Adhesive, high viscosity—1100-1600 cps (pkg of 10 vials, each approx 1.5mL)
VETBOND	3M Vetbond™ Adhesive (3 mL)

Sylgard



A two-part silicone elastomer, ideal for potting and encapsulating applications. Very low dielectric constant sealing compound used in patch clamping and many other lab applications. After cure, will withstand -55° to 200 °C. Shipping weight: 2 lb. (1 kg)

Sylgard (1.1 lb) **SYLG184**

Silicone Dissecting Pad Kit



Make your own silicone dissecting pads easily and quickly. Mix the 2-part silicone right in the plastic petri dishes and allow to cure 24 hours at room temperature. Kit includes enough silicone to prepare 20 dishes.

Kit Includes:

2-Part Sylgard silicone elastomer 20 plastic petri dishes with lids, 65mm Pins

501986 **Silicone Dissecting Pad Kit**



501601	Digital Caliper
502157	Replacement Battery (package of 10)



501321 Roger Wirecutting Scissors

Screwdriver Set



This production grade precision 5-piece screwdriver set is the highest quality tool you can find on the market. Made by German craftsmen, the chrome-vanadium tips will fit any screw securely without leaving marks. The set contains an ESD safe handle with 8 interchangeable blades. Phillips Sizes: 000, 00, 0, 1. Slotted Sizes: 1.5, 2.0, 3.0, 3.5 mm. The 000 size Philips blade is the smallest you can find anywhere — it can fit smallest screw on a 35 mm camera. Weight: 0.24 lb

501635 Professional Quality Screwdriver Set

Electrically Conductive Silver Epoxy

Two-component silver-filled epoxy for electrical connections which cannot be soldered, such as Ag/AgCl pellets. This widely used silver-filled epoxy features low viscosity and smooth flowing character. Pure silver is dispersed in both resin and hardener. Cures in 15 minutes at 120 °C. Mix ratio 1:1. May be premixed and frozen for later use.

Shipping weight: 1 lb. (0.5 kg)

4898 Silver Epoxy (1 oz.)

Electrically Conductive Carbon Epoxy

Two-component carbon-epoxy, curable at room and elevated temperatures. Ideal for electrostatic discharge protection and EMI/RFI shielding. 1:1 mix ratio. May be premixed and frozen for later use. Shipping weight: 1 lb. (0.5 kg)

Carbon Epoxy (2 oz.)



Silicone RTV adhesive (non-acidic)

Because it is non-corrosive, this material is ideal for use on metal, for encapsulating small circuits on connectors. After cure, will withstand -55° to 200 °C. No mixing required.

Shipping weight: 1 lb. (0.5 kg)

RTV Coating (3 oz.)

RTV Prime Coat

Enhances adhesion of silicone adhesives to many difficult-to-bond plastic surfaces.

Shipping weight: 1 lb. (0.5 kg)

RTV Prime Coat, 400 ml (13.5 oz.)



Silicone RTV adhesive

Clear silicone sealant provides good bonding to plastic. After cure, will withstand -55 to 200 °C. No mixing required. A handy, general purpose laboratory sealant. (Releases acetic acid during curing.) Shipping weight: 1 lb. (0.5 kg)

RTV Sealant (4.7 oz.)

Mini Glue Gun

Comes with three sticks of special formula hot melt glue. UL approved. 110V 60Hz only. Shipping weight: 1 lb. (0.5 kg)

Mini Glue Gun

Luer Valve Assortment Kit



A useful kit (above) for building your own liquid flow experiment. It provides the means to start, stop, add, divide and control a flow of liquid or gas. Included in the kit are over 200 assorted parts such as one-way and three-way stopcocks, manifolds, Y-connectors, injection sites, male and female luer caps, check valves, syringe-activated check valves, slide clamps, roller clamps, and pinch clamps. All (except clamps) have a luer fitting for quick and easy connecting and disconnecting. Includes assorted luer fittings for use with flexible tubing.

Luer Valve Assortment Kit



Barb-to-Tubing Coupler Assortment Kit

Barb-to-Tubing Assortment Kit (at left) includes three different sizes of tubing and two boxes with different fittings, T-connectors, elbow connectors, check valves and plugs.

Barb-to-Tubing Assortment Kit (polypropylene)

Includes 25ft. each of three tubing sizes: 1/16" ID, 1/8" ID, 1/4" ID

Luer-to-Tubing Coupler Assortment Kit

Assemble guick-disconnect luer fittings for use with flexible tubing with internal diameters of \%", \%2", \%" and \%2". A variety of quick-disconnect connectors can be quickly made for connecting small diameter flexible tubing; 3-way connections can be made with the use of the 3-way luer tee; luer plugs, tees, connectors, bulk-head mounts, color coding rings, locking nuts, male and female luers—are all included to enhance the versatility of this kit. The kit has 253 assorted parts and is offered in two different types of materials. Polypropylene fittings are chemically inert and resistant to most organic and inorganic solvents. Nylon fittings are strong and can be bonded with adhesive.

14012 Luer-to-Tubing Coupler Assortment Kit (Polypropylene) 500895 Luer-to-Tubing Coupler Assortment Kit (Nylon)





Parts now sold individually

Buy parts from Luer Valve Assortment Kit 14011 individually



Parts in kit may differ slightly in appearance from those pictured.

Cables & Connectors











#3508





















































PART #	APPLICATION/DESCRIPTION	CONNECTOR A	CONNECTOR B	CABLE LENGTH	
1358	Beetrodes	BNC (male)	2 mm pin	3 ft (0.9 m)	
2026-10	2 mm socket, unwired (pkg of 10) (Not Shown)	2 mm socket	unwired	none	
2851	Standard BNC cable	BNC (male)	BNC (male)	6 ft (1.8 m)	
3142	Mini-Banana Adapter	Screw Terminals	Dual Mini-Banana	none	
3161	Connector for input to TBM4M and BP-1	DIN (male)	unwired	none	
3294	Ground wire for DAM80 probe	Clip	none	3 ft (0.9 m)	
3417-10	2 mm plug, unwired (pkg of 10)	2 mm pin	unwired	none	
3491	Extension for any 8-pin DIN	DIN (male)	DIN (female)	5 ft (1.5 m)	
3492	Connector, adapts WPI transducers to non-WPI equipment	DIN (female)	unwired	none	
3508	Adapts BNC pH electrode to pH meter with	(/			
	"U.S. Standard" input	BNC (male)	US Standard	none	
3517	DAM50, DAM60, DAM70, shielded (two cables/pkg)	Modular phone plug, 4 wire	none	3 ft (0.9 m)	
3578	Adapter cable for Ag/AgCl pellets	2 mm pin	none	5 ft (1.5 m)	
3670	Double banana plug with solder turret terminals	Dual Banana (male)	Dual Banana (female)	none	
5371	Low-noise cable for microelectrode holders	2 mm gold pin	2 mm gold pin	2 ft (0.6 m)	
5372	Low-noise cable for microelectrode holders	2 mm gold jack	2 mm gold jack	2 ft (0.6 m)	
5373	Low-noise cable for microelectrode holders	2 mm gold pin/jack	2 mm gold pin/jack	2 ft (0.6 m)	
5374	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold pin	4 ft (1.2 m)	
5375	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold jack	4 ft (1.2 m)	
5385	Cable, shielded transducer stock	none	none	25 ft (7.6 m)	
13324	Adapter	Double-banana (female)	BNC (male)	none	
13347	ISO2 (chart recorder adapter)	Double-banana (male)	BNC (female)	none	
13388	Electrode adapter for DAM probes	Miniature banana (male)	2 mm jack	none	
13451	Adapter: Iso-DAM, Iso-DAM8	BNC (female)	two 2 mm pins	6 in. (15 cm)	
13555	Serial Cable (not shown)	DB9 (male)	DB9 (female)	6 ft (1.8 m)	
13620	Low-noise cable for microelectrode holders	2 mm gold pin	2 mm gold jack	2 ft (0.6 m)	
13685	SP Series pump-to-pump linking cable	Modular phone plug	Modular phone plug	7 ft (2.1 m)	
13776	Adapts reference electrode to VF4 ground jack	Banana (male)	2 mm jack	none	
13854	BNC T-connector, male to:	BNC (female)	BNC (female)	none	
14254	BNC Straight Adapter	BNC (female)	BNC (female)	none	
15623	Serial cable and adapter, SP Series pump	SP Pump	IBM 9-pin "D" connector	5 ft (1.5 m)	
15975	Adapter	2 mm socket	1 mm pin	none	
15976	Adapter	1 mm socket	2 mm pin	none	
300040	Adapter Extension	2 mm socket	2 mm socket	4 in. (10 cm)	
500040	Standard BNC Cable	BNC (male)	BNC (male)	10 ft (3 m)	
500184	BNC Right Angle Adapter	BNC (male)	BNC (female)	none	
500257	Standard BNC Cable	BNC (male)	BNC (male)	6 in. (15 cm)	
500257	Standard BNC Cable	BNC (male)	BNC (male)	12 in. (30 cm)	
500258	Standard BNC Cable Standard BNC Cable	BNC (male)	BNC (male)	18 in. (46 cm)	
				10 ft (3 m)	
503301	Cable, Extension	8-pin miniDIN (male)	8-pin miniDIN (female)		
503536 504713	Cable, USB	USB (male) Banana (male)	USB (female)	10 ft (3 m)	
	Cable (red and black pair)		Banana (male)	36 in. (91 cm)	
504714	Cable (red and black pair)	Banana (male)	Mini-Gator	36 in. (91 cm)	
504715	Cable (red and black pair)	Banana (male)	Mini-Clip	36 in. (91 cm)	
504716	Cable (red and black pair)	Banana (male)	Micro-Clip	36 in. (91 cm)	
CBL100	MiniPhone Patch Cable	3.5 mm MiniPhone plug	3.5 mm MiniPhone plug	6 ft (1.8 m)	
CBL102	DAM Series, PM Series	3.5 mmMiniPhone plug	BNC (male)	6 ft (1.8 m)	
















FrameWorks

Non-Magnetic Bases, Stainless Steel Rods, & Clamps

These high quality components are made of stainless steel and polymer that resist organic solvents and corrosion. They can be easily assembled to make a stand-alone setup for student labs or to make a complicated frame for research labs.



503041

Large Clamp with Rod (157 mm), opens up to 85 mm



Medium Clamp with rod (157 mm), opens up to 45 mm



503086

Small Clamp with rod (157 mm), opens up to 16 mm



14073-4

Open-sided Frame Clamp (pkg of 4)



503082-4

Board Frame Clamp, opens to 8.5mm (pkg of 4)



503080-4

Frame Clamp with Parallel Surface Mount (includes mounting screws) (pkg of 4)



502193-4

Parallel Frame Clamp (pkg of 4)




503078-4

T-joint Frame Clamp (pkg of 4)



503079-4

In-line Frame Clamp (pkg of 4)





502190

Heavy Rectangular Base (with M8 thread mount and thumbscrew mount), 23×15.6 cm, 4 lb

Light Rectangular Base (with M8 **503083** thread mount and thumbscrew mount), 23×15.6 cm, 0.5 lb



503085

Large 10-in. V-base with M8 Thread Mount



503084

Small 6-in. V-base with M8 Thread Mount



503081-4

Vertical Surface Mount, M8 Threaded

g-SPIN™ Microcentrifuge

• Small size with rubber feet keeps the centrifuge stable

 Supplied with interchangeable rotors and adapters for 0.5mL-2.0mL microtubes and PCR strips

On/off switch lets you start and stop in seconds

Safety switch stops rotor without cover in place

6000 rpm, fixed

6 x 0.5/2.0 mL tubes 2 x 0.2 mL strip tubes

> 6-place 0.5—2.0-mL tubes 2 strips PCR 0.2 mL tubes

2000× g

Fixed angle

6"W x 5"H x 7"D

110V 60 Hz or 220V 50 Hz

g-SPIN SPECIFICATIONS

SPEED RANGE MAX. RCF

TUBE CAPACITY

DIMENSIONS POWER

ROTOR

Also available:
10,000 RPM
units — call for details.
ecuis.





G-SPIN6	Microcentrifuge, 6000 RPM, 110V 60 Hz
G-SPIN6-220	Microcentrifuge, 6000 RPM, 220V 50 Hz, CE
503529	Microcentrifuge tube, 0.5 mL, natural, bag/1,000
503530	Microcentrifuge tube, 1.5 mL, natural, bag/1,000
503531	Microcentrifuge tube, 2.0 mL, natural, bag/1,000
503532	Microcentrifuge tube strips, 0.2 mL, and domed cap strips, bag of 250
503537	Microcentrifuge tube strips, 0.2 mL, and flat cap strips, bag of 250

Rackmounting Hardware

Many instruments may be mounted in standard 19-inch instrument racks with the appropriate rackmount kit, as noted on the page featuring the instrument.

Dual rackmount kits allow many smaller instruments to be joined by bolting the chassis together and mounting the pair into a standard rack.

Brackets for instruments which are less than 17.5 inches wide have "wings" which extend to the standard rack width.



3484 Rackmount Kit for DAM Series amplifiers

Rackmount Kit, 3½-in. high (121)





DAM Series amplifiers and Iso-DAM amplifiers may be mounted to the rack panels above by fastening bolts (included) through holes in the panel.

Rackmount Kit, 8 ¾-in. high




Single Rackmount Kit (TBM4M)


Metal Microelectrodes

Superior microelectrodes for outstanding extracellular recording — tungsten, iridium, platinum-iridium, and Elgiloy®

EXPOSED TIP DIMENSIONS (nominal)					
Nominal Impedance	Tungsten	Elgiloy	Platinum Iridium	Pure Iridium	
10 k Ω 50 k Ω	250 μ 200 μ	_	_	_	
0.1 MΩ	100 μ	1 20 μ	6 0 μ	45 μ	
$0.5~\mathrm{M}\Omega$	55 µ	66 µ	18 μ	14 μ	
1.0 MΩ	30 μ	36 µ	10 μ	10 μ	
2.0 MΩ	12 μ	15 μ	6 μ	5 μ	
5.0 Meg Ω	5 μ	6 μ	3 μ	2.5 μ	




Kapton* tubing, indicated by "KT" in the part number, extends from the connector to within 5 mm of the tip, providing stiffness and additional insulation to the electrode shaft. Kapton-clad electrodes are recommended when the electrode is to be inserted through a cannula for extra deep penetration.



Note: Electrode diagrams not shown to scale.




Type B



Concentric Bipolar Electrodes

Excellent for shielded macro recording as well as evoked potentials — especially well suited for bipolar stimulation

The tungsten electrode is sharpened to a point and is 75 microns in diameter. The outer stainless steel conductor is insulated with Polyimide tubing to within 0.2 mm of the end of the stainless steel tube. Also available without the outer Polyimide insulation.



Heat Treated Tip (above) is ideal for penetrating tough membranes (not

recommended for chronic

implantation). This process is performed using a microforge in which the heating element is positioned in close proximity to the tip in order to melt the Parylene-C distal to the exposed metal. It provides a smooth transition and produces better adherence of the Parylene-C to the metal.


To have your electrodes heat treated, just add the suffix "H" to any of the "KT" numbers on the facing page. Cost of the treatment is \$10 per package of 10 electrodes (\$20 per package for Type B).

* Parylene is a trade mark of Union Carbide. Kapton is a trade mark of DuPont. Elgiloy is a trade


ACCESSORIES

300102	Micromanipulator holder, 4 in., 2mm to 0.031 socket			
5468	2 mm receptacle to 0.031-inch jack (for Omega-TipZ)			
5469	Adapts mini banana plug (DAM80) to 0.031-inch			
	receptacle (metal microelectrode)			
5470	0.031-inch jack, 28 ga. wire, 12 inch (pkg of 4)			
5482*	Pins, 0.031-inch, gold-plated (pkg of 50)			
5483*	Sockets, 0.031-inch gold-plated (pkg of 50)			

*Gold-plated pins (#5482) and sockets (#5483) may be attached to 24-, 26-, or 28-gauge wire.



Gold-plated pins (#5482) and sockets (#5483) may be attached to 24-, 26-, or 28-gauge wire.



Introductory Assortments

Each of these assortment kits includes electrodes with different impedance within each style. Use an assortment kit to determine which electrode you need for your experiment. Ten electrodes per box, no mixing.

Item	Contains the following electrode impedances by quantity (pkg of 10)	Price
TM31/33Axx	TM33A05 (2), TM33A10 (3), TM33A20 (3), TM31A50 (2)	
TM31/33AxxKT	TM33A05KT (2), TM33A10KT (3), TM33A20KT (3), TM31A50KT (2)	
TM33BxxKT	TM33B01KT (3), TM33B05KT (2), TM33B10KT (3), TM33B20KT (2)	
TST33AxKT	TST33A05KT (3), TST33A10KT (4), TST33A20KT (3)	

Concentric Electrodes*

Item	Metal Core	Length	Imp	Probe Outer Diameter (total)	Tip Diam.	Core diam.	Y dim.	X dim. w/ polyimide	Price (pkg of 5)
TM33CCNON	Tungsten	3" (76)	10-15K	0.013" uninsulated (325 µm)	3-4 μ	.003" (76 μm)	0.4 mm	.005" (127 μm)	
TM33CCINS	Tungsten	3" (76)	10-15K	0.016" insulated (400 μm)	3-4 µ	003" (76 μm)	0.4 mm	.005" (127 μm)	
TM53CCINS	Tungsten	5" (127)	10-15K	0.018" insulated (450 μm)	3-4 µ	005" (127 μm)	0.4 mm	.008" (203 µm)	
PTM23CC001NON	Pt/Ir	2" (51)	10K	0.020" uninsulated (525 µm)	3-4 µ	0.01" (254 µm)	0.4 mm	.014" (356 µm)	
PTM3CC02INS	Pt/Ir NS fine	3" (76)	200K	0.013" insulated (325 μm)	2-4 µ	0.002" (50.8 μm)	.25 mm	.004" (114 µm)	

^{*}All have a stainless steel outer shaft

Selection Guide for Metal Electrodes


Item	Length	Insul. Thick.	Shaft Diam.	Nominal Impedance (± 20%)		Typical Use
Tungsten — P	rofile A					Package of 10
TM31A10	76 mm	1 μ	0.127 mm	1.0 MΩ	1 μ	Multi unit and single unit recording and microstimulation
TM31A20	76 mm	1 μ	0.127 mm	$2.0~{ m M}\Omega$	1 μ	Multi unit and single unit recording and microstimulation
TM31C05	76 mm	1 µ	0.085 mm	$0.5~{ m M}\Omega$	1 µ	Recording from small tightly packed cells
TM33A05	76 mm	3 µ	0.127 mm	$0.5~{ m M}\Omega$	1 µ	Multi unit and single unit recording and microstimulation
TM33A10	76 mm	3 μ	0.127 mm	$1.0~{ m M}\Omega$	1 μ	Multi unit and single unit recording and microstimulation
TM33A20	76 mm	3 µ	0.127 mm	$2.0~{ m M}\Omega$	1 µ	Multi unit and single unit recording and microstimulation
TM33B01	76 mm	3 µ	0.254 mm	$0.1~{ m M}\Omega$	1-2 µ	Single and multi unit recording and microstimulation
TM33B05	76 mm	3 µ	0.254 mm	$0.5~{ m M}\Omega$	1-2 µ	Single and multi unit recording and microstimulation
TM33B10	76 mm	3 μ	0.254 mm	$1.0~{ m M}\Omega$	1-2 µ	Single and multi unit recording and microstimulation
TM33B20	76 mm	3 μ	0.254 mm	$2.0~{ m M}\Omega$	1-2 µ	Single and multi unit recording and microstimulation
TM33C05	76 mm	1 μ	0.085 mm	$0.5~{ m M}\Omega$	1 μ	Single unit and stim / chronic use
TM33C10	76 mm	1 µ	0.085 mm	$1.0~{ m M}\Omega$	1 μ	Single unit and stim / chronic use
Tungsten — P	rofile C	•				Package of 10
TM31A10KT	76 mm	1 μ	0.216 mm	1.0 MΩ	1 μ	Multi unit and single unit recording and microstimulation
TM33A10KT	76 mm	3 µ	0.216 mm	1.0 MΩ	1 µ	Multi unit and single unit recording and microstimulation
TM33B01KT	76 mm	3 µ	0.356 mm	0.1 MΩ	1-2 µ	Single and multi unit recording and microstimulation
TM33B05KT	76 mm	3 μ	0.356 mm	$0.5~{ m M}\Omega$	1-2 µ	Single and multi unit recording and microstimulation
TM33B10KT	76 mm	3 µ	0.356 mm	1.0 MΩ	1-2 µ	Single and multi unit recording and microstimulation
Elgiloy®/Stainl	less — Pro	file A				Package of 10
SSM33A70	76 mm	3 u	0.229 mm	7.0 MΩ	1-2 µ	Recording and Stimulating (Prussian blue staining)
SSM33A120	76 mm	3 µ	0.229 mm	12.0 MΩ	1-2 µ	Recording and Stimulating (Prussian blue staining)
Elgiloy®/Stainl	less — Pro	file C				Package of 10
SSM33A20KT	76 mm	3 μ	0.356 mm	2.0 MΩ	1-2 µ	Recording and Stimulating (Prussian blue staining)
Tungsten — P	rofile B					Package of 10
TST33A001KT	76 mm	3 u	0.356 mm	10 kΩ	1 μ	Tissue slice stimulation
TST33A05KT	76 mm	3 µ	0.356 mm	0.5 MΩ	1 µ	Stereotrode / Bipolar, differential measurements
TST33A10KT	76 mm	3 µ	0.356 mm	1.0 MΩ	1 µ	Stereotrode / Bipolar, differential measurements
TST33A20KT	76 mm	3 µ	0.356 mm	2.0 MΩ	1 µ	Stereotrode / Bipolar, differential measurements
TST33C05KT	76 mm	3 µ	0.216 mm	0.5 M Ω	1 u	Stereotrode / Bipolar, differential meas. — extra fine (75 µm separation)
TST53A10KT	127 mm	3 μ	0.356 mm	1.0 MΩ	1-2 µ	Stereotrode / Bipolar, differential measurements
Pure Iridium -	– Profile A					Package of 10
IRM23E10	50 mm	3 μ	0.106 mm	1.0 MΩ	1-2 µ	Single and multiunit recording and stimulation
IRM23E15	50 mm	3 µ	0.106 mm	1.5 M Ω	1-2 µ	Single and multiunit recording and stimulation
IRM23E25	50 mm	3 µ	0.106 mm	$2.5~\mathrm{M}\Omega$	1-2 µ	Greater selectivity - small cells
IRM23E30	50 mm	3 µ	0.106 mm	3.0 M Ω	1-2 µ	Greater selectivity - small cells
Pure Iridium -	– Profile C					Package of 10
IRM23E01KT	50 mm	3 μ	0.180 mm	0.1 M Ω	2-3 μ	Multiunit & ERP recording & stimulation
IRM23E20KT	50 mm	3 µ	0.180 mm	$2.0~{ m M}\Omega$	1-2 µ	Greater selectivity & microstimulation
IRM23E25KT	50 mm	3 µ	0.180 mm	$2.5~{ m M}\Omega$	1-2 µ	Greater selectivity - small cells
IRM23E30KT	50 mm	3 μ	0.180 mm	$3.0~{ m M}\Omega$	1-2 µ	Greater selectivity - small cells

Eligoy Steel *Cobalt/chromium/nickel alloy. The KT suffix refers to Kapton™ cladding.

All Metal Microelectrodes are available in custom lengths, blunt or heat treaded (extra charge).

Add the $\vec{\bf B}$ suffix where blunt electrodes are desired. (For example, an IRM123A10KT ordered as a blunt will be IRM123A10KTB.) Add the **H** suffix where heat treated electrodes are desired. (For example, an IRM123A10KT ordered as with heat treatment will be IRM123A10KTH.)


Additional metal microelectrodes available on website www.wpiinc.com





Ag/AgCl Half-Cells


New, improved sintered pellets with lower resistance and high strength. Stable and well balanced in the presence of current, these small and inexpensive half-cells are easy to work with as bath electrodes.


RC1	Reference Cell with 1.5 m lead
RC1T	Reference Cell, 1.5 m lead, 2 mm pin
RC2	Reference Cell with 2.0 mm pin
RC2F	Reference Cell with female connector
RC3	Reference Cell with epoxy body, 4.5 mm diam x 50 mm
RC6	Reference Cell with glass body, 1.5 mm diam x 50 mm
EP05	Ag/AgCl Electrode 0.5 mm diam x 20 mm
EP08	Ag/AgCl Electrode 0.8 mm diam x 20 mm
EP1	Ag/AgCl Electrode 1.0 mm diam x 3 mm
EP2	Ag/AgCl Electrode 2.0 mm diam x 4 mm
EP4	Ag/AgCl Electrode 4.0 mm diam x 1 mm
EP8	Ag/AgCl Electrode 8.0 mm diam x 1 mm
EP12	Ag/AgCl Electrode 12.0 mm diam x 1 mm
3578	Adapter Cable for Ag/AgCl Pellets











Precious Metal and Specialty Wire



New! Micro coaxial cables (**MAXxxxx**) are ideal for microelectrode fabrication and construction of similar research tools. The dual shielding eliminates electrical interference caused by radio frequencies (RF), electrostatic and microphonics (e.g., bending and vibration. Available with single or dual (twin) conductors.

Teflon-coated stainless steel (type 304) wire (SSTxxxx) is available in 25-ft and 50-ft lengths. The Teflon coating is 150 micro-in.

Carbon wire (C3005) is a single 30-micron fiber of electrochemically activated carbon. This fiber is especially useful in micro-electrochemical experiments.

Platinum/iridium wire — uncoated (**PTxxxx**) and Teflon-coated (PTTxxxx) — is an alloy of 90% platinum and 10% iridium, giving excellent tensile strength and corrosion resistance. Uncoated pure platinum wire (PTPxxx) is 99.95% pure.

Indium wire (IN1003) is 99.99% pure, with a melting point of 156.4°C.

Annealed silver wire (AGWxxxx), 99.99% pure, is available in five diameters; three of those sizes are also available with a Teflon coating (AGTxxxx).

Tungsten wire (TGWxxxx), available in three diameters, is 99.95% pure.

Gold wire (AUWxxxx) is 99.99% pure. Titanium wire (Tlxxxx) is 98.9% pure, annealed, in two diameters.

Stainless steel wire (SSxxxxx) is type 316.

Catalog No.	Metal	Coating	AWG*	Diameter	Precut Length
AGT0510	Silver	Teflon	36	0.005 in. (0.125 mm) ¹	10 ft (3 m)
AGT0525	Silver	Teflon	36	0.005 in. (0.125 mm) ¹	25 ft (7.6 m)
AGT05100	Silver	Teflon	36	0.005 in. (0.125 mm) ¹	100 ft (30 m)
AGT1010	Silver	Teflon	30	0.010 in. (0.25 mm) ¹	10 ft (3 m)
AGT1025	Silver	Teflon	30	0.010 in. (0.25 mm) ¹	25 ft (7.6 m)
AGT10100	Silver	Teflon	30	0.010 in. (0.25 mm) ¹	100 ft (30 m)
AGT1510	Silver	Teflon	26-27	0.015 in. (0.38 mm) ¹	10 ft (3 m)
AGT1530	Silver	Teflon	26-27	0.015 in. (0.38 mm) ¹	30 ft (9.1 m)
AGW0510	Silver	_	36	0.005 in. (0.125 mm)	10 ft (3 m)
AGW0530	Silver	_	36	0.005 in. (0.125 mm)	30 ft (9.1 m)
AGW1010	Silver	_	30	0.010 in. (0.25 mm)	10 ft (3 m)
AGW1030	Silver	_	30	0.010 in. (0.25 mm)	30 ft (9.1 m)
AGW1510	Silver	_	26-27	0.015 in. (0.38 mm)	10 ft (3 m)
AGW1530	Silver	_	26-27	0.015 in. (0.38 mm)	30 ft (9.1 m)
AGW2010	Silver	_	24	0.020 in. (0.5 mm)	10 ft (3 m)
AGW2030	Silver	_	24	0.020 in. (0.5 mm)	30 ft (9.1 m)
AGW4010	Silver	_	18	0.040 in. (1.0 mm)	10 ft (3 m)
AUW0170	Gold	_	50	0.001 in. (0.025 mm)	70 ft (21 m)
AUW201	Gold	_	24	0.020 in. (0.5 mm)	1 ft (30 cm)
C3005	Carbon	_	49	0.0012 in. (30 μm)	5 ft (1.5 m)
PT1002	Platinum / Iridium	_	30	0.010 in. (0.25 mm)	2 ft (61 cm)
PT0402	Platinum / Iridium	_	38	0.004 in. (0.102 mm)	2 ft (61 cm)
PT0203	Platinum / Iridium	_	44	0.002 in. (0.051 mm)	3 ft (91 cm)
PT0110	Platinum / Iridium	_	50	0.001 in. (0.025 mm)	10 ft (3 m)
PTP101	Platinum	_	30	0.010 in. (0.25 mm)	1 ft (30 cm)
PTP201	Platinum	_	24	0.020 in. (0.5 mm)	1 ft (30 cm)
PTP401	Platinum	_	18	0.039 in. (1.0 mm)	1 ft (30 cm)
PTP406	Platinum	_	18	0.039 in. (1.0 mm)	0.5 ft (15.2 cm)
PTT0502	Platinum / Iridium	Teflon	36	0.005 in. (0.125 mm) ¹	2 ft (61 cm)
PTT0203	Platinum / Iridium	Teflon	44	0.002 in. (0.051 mm) ¹	3 ft (91 cm)
PTT0110	Platinum / Iridium	Teflon	50	0.001 in. (0.025 mm) ¹	10 ft (3 m)
SS31605	Stainless Steel	_	36	0.005 in. (0.125 mm)	50 ft (15.2 m)
SS31614	Stainless Steel	_	27	0.014 in. (0.36 mm)	30 ft (9.1 m)
SST30407-25	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) ³	25 ft (7.6 m)
SST30407-50	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) ³	50 ft (15.2 m)
TGW0325	Tungsten	_	40	0.003 in. (0.075 mm)	25 ft (7.6 m)
TGW0515	Tungsten	_	36	0.005 in. (0.125 mm)	15 ft (4.6 m)
TGW1510	Tungsten	_	26-27	0.015 in. (0.38 mm)	10 ft (3 m)
Microcoaxial	Cables				
MAX3820	Tinned Cu Alloy	Coaxial		0.0173 in. (0.44 mm)	20 ft (6 m) ⁴
MAX4020	Tinned Cu Alloy	Twin Coaxial	0.015	8x0.024 in. (0.4x0.61mm)	20 ft (6 m) ⁵

^{*}Brown & Sharpe

¹ Plus 0.002 in. for Teflon coating

³ Teflon adds 0.00015 in. (4 µm) to diameter

⁴ Impedance: 50 ohm; capacitance: 95 pF/m; resistance: 5 ohm/m

⁵ Impedance: 100 ohm; capacitance: 54 pF/m; resistance: 1.9 ohm/m

Micropipette Holders & Half-Cells

We offer a large variety of micropipette holders. Our popular ones are stock items. Custom holders (designated by ♦) can be manufactured on demand but require an additional setup fee. Call for a quote. See all the options at

www.wpiinc.com/MEH


WPI's microelectrode holder-half-cells couple fluid-filled glass micropipettes to high input impedance amplifiers. A Ag/AgCl pellet (or a silver wire) molded into the holder body provides stable potential. Electrical connection is made via male 2 mm pins or female 2 mm sockets. The pipette may be mounted axially or at right angles to the holder. Pipettes are held with screw-caps or rubber gaskets (without caps). Filling WPI microelectrode holders with electrolytes containing chloride results in stable electrode potential. Suitable electrolytes include KCl, NaCl and CaCl2. Holders are supplied for standard WPI single capillary tubing of 1.0, 1.2, 1.5 and 2.0 mm outside diameters. (Call WPI regarding custom designs for other glass diameters.) The holder style you select will depend on your experimental application, space, and


Hints for selecting and ordering micropipette holders


- 1. Determine the required electrical connection on the holder: for example, if you wish to connect the holder to a 2 mm pin you should select a holder equipped with a 2 mm jack. Most WPI probes require a holder equipped with a 2 mm jack.
- 2. Decide on the required alignment of the electrical connection: either in-line with the glass pipette, or at a right angle to it. Space considerations in your experimental setup and requirements imposed by other pieces of equipment typically determine which alignment is appropriate.
- 3. Determine if you want to hold the glass pipette by a rubber gasket (e.g., MEH1S) or a screw-cap (e.g., MEH3S). Rubber gaskets offer easier insertion and removal of glass pipettes whereas screw-caps provide more secure mounts for micropipetters.
- 4. Choose a holder with either a silver wire or a silver/silver chloride pellet for the metal/liquid coupling. Silver/silver chloride pellets provide a more stable low-noise baseline which is important for low-noise DC recording. Pellets require the glass pipette and holder to be free of air bubbles to achieve a good connection. Silver wire holders are durable and are easier to use when the holder is equipped with a pressure port because the fluid in the pipette does not have to be filled to the top of the pipette to achieve a good electrical connection
- 5. Choose a holder equipped with a pressure port only when you want to pressure inject liquid from the pipette. Two types of ports are available: 2.0 mm O.D. and standard "syringe-style" luer. The luer port is often recommended because it makes assembly and disassembly much easier. Quick-connect luer fittings for four common sizes of tubing (1/16", ³/₃₂", ¹/₈", ⁵/₃₂" I.D.) are included with each luer-equipped holder.
- **6.** Some non-WPI preamplifiers or headstages cannot be mounted on micromanipulators. In such cases, a holder equipped with a rod (e.g., MEH8) permits the holder to be conveniently mounted on a micromanipulator.
 - 7. Finally, remember to specify the O.D. of the glass you will be using when you place your order.


MEH6RF/SF is designed primarily for use with the Model 900A Micropressure System; EHB1 for use in electrode beveling; and **MEH3SW** for microtitration of chloride with a silver wire as the electrode and a solution of silver nitrate filling the holder. MPH models do not contain Ag/AgCl half-cells and are used for pressure injection of substances through microelectrodes. PicoNozzle, used for pressure injection with PV800 Series PicoPumps, includes an MPH6S holder which may also be used to couple a micropipette to a syringe. APEH models are also designed for use in pressurized injection procedures.


 denotes holder sizes manufactured for you as custom orders. Call for price.


























Additional holders shown on pages 108-109.

MICROELECTRODE HOLDERS

= custom order (call for price)

Order Number Replace XX with		ass C amet			Electric Connection	Connector	Half-Cell	Pressure	Screw	Designed for
glass diameter 🛨		1.2	•	•	Angle	,		Port	Сар	WPI Products
APEH1xx	•	•	•	•	Right	Male	Pellet	No Port	2 Caps	
APEH2xx					Right	Male	Wire	No Port	2 Caps	
EHB1					Straight	Male	Wire	No Port	N/A	MBS, 48000
EHBF					Straight	Female	Wire	No Port	N/A	MBS, 48000
FOIMPH					Straight	Fiber Optic	None	No Port	w/Cap	MBS, 48000
FOIMPH-LF					Straight	Fiber Optic	None	Male Luer	w/Cap	MBS, 48000
MBMPH5-7						None	None	Female Luer	w/Cap	For P-5 or P-7 glass only
MEH145xx	•	•	•	•	45°	Male	Pellet	No Port	No Cap	
MEH1F45xx	•	•	•	•	45°	Female	Pellet	No Port	No Cap	705, 773, 767, 721, FD223
MEH1Rxx	•	•	•		Right	Male	Pellet	No Port	No Cap	
MEH1RFxx	•	•	•	•	Right	Female	Pellet	No Port	No Cap	705, 773, 767, 721, FD223
MEH1Sxx	•	•	•	•	Straight	Male	Pellet	No Port	No Cap	, , , ,
MEH1SFxx					Straight	Female	Pellet	No Port	No Cap	705, 773, 767, 721, FD223
MEH2Rxx	•	•	•	•	Right	Male	Pellet	Male Luer	w/Cap	
MEH2RFxx	•	•	•	•	Right	Female	Pellet	Male Luer	w/Cap	705, 773, 767, 721, FD223
MEH2RFWxx	•	•	•	•	Right	Female	Wire	Male Luer	w/Cap	705, 773, 767, 721, FD223
MEH2RWxx	•	•	•	•	Right	Male	Wire	Male Luer	w/Cap w/Cap	. 55, 7, 5, 757, 721, 15225
MEH2Sxx		*	*	•	Straight	Male	Pellet	Male Luer	w/Cap w/Cap	
MEH2SFxx					Straight	Female	Pellet	Male Luer	w/Cap w/Cap	705, 773, 767, 721, FD223
MEH2SFWxx		•		•	Straight	Female	Wire	Male Luer	w/Cap w/Cap	705, 773, 767, 721, FD223
MEH2SWxx		*		•	Straight	Male	Wire	Male Luer	w/Cap w/Cap	. 33, 7, 3, 7, 7, 7, 7, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1,
WEH345xx					45°	Male	Pellet	No Port	w/Cap w/Cap	
MEH3F45xx	•	•	•	•	45°	Female	Pellet	No Port	w/Cap	705, 773, 767, 721, FD223
MEH3FW45xx	•	Ť	<u> </u>	<u> </u>	45°	Female	Wire	Port	w/Cap	703,773,707,721,10223
MEH3Rxx	•			•	Right	Male	Pellet	No Port	w/Cap	
MEH3RFxx					Right	Female	Pellet	No Port	w/Cap	705, 773, 767, 721, FD223
MEH3RFWxx	•	•	•	•	Right	Female	Wire	No Port	w/Cap	705, 773, 767, 721, FD223
MEH3RWxx	Ť	Ť	<u> </u>		Right	Male	Wire	No Port	w/Cap	703,773,707,721,10223
MEH3Sxx	•	•	•	•	Straight	Male	Pellet	No Port	w/Cap	
MEH3SBWxx	•	•		•	Straight	Banana	Wire	No Port	w/Cap w/Cap	ISO-80, ISO-DAM8A
MEH3SFxx	•	*		•	Straight	Female	Pellet	No Port	w/Cap w/Cap	705, 773, 767, 721, FD223
MEH3SFWxx				•	Straight	Female	Wire	No Port	w/Cap w/Cap	705, 773, 767, 721, FD223
MEH3SWxx				•	Straight	Male	Wire	No Port	w/Cap w/Cap	100, 110, 101, 121, 1022
MEH3W45xx	•	•	•	•	45°	Male	Wire	No Port	w/Cap w/Cap	705, 773, 767, 721, FD223
MEH6RFxx	•	▼	•	•	Right	Female	Pellet	2.0-mm Port	w/Cap w/Cap	705, 773, 767, 721, FD223
MEH6RFWxx		•	•	•	Right	Female	Wire	2.0-mm Port	w/Cap w/Cap	705, 773, 767, 721, FD223
MEH6SFxx		•	•	•	Straight	Female	Pellet	2.0-mm Port	w/Cap w/Cap	705, 773, 767, 721, FD223
MEH6SFWxx		•		•	Straight	Female	Wire	2.0-mm Port	w/Cap w/Cap	705, 773, 767, 721, FD223
MEH7xx		▼		•	Right	Male	Pellet	2.0-mm Port	w/Cap w/Cap	103, 113, 101, 121, FD223
MEH7Wxx		•			Right	Male	Wire	2.0-mm Port	w/Cap w/Cap	
VIEH8XX		▼			Right	Male	Pellet	No Port	w/Cap w/Cap	
MEH900Rxx					Right	Male	Pellet	2.0-mm Port	w/Cap w/Cap	900A
MEH900SXX					Straight	Male	Pellet	2.0-mm Port	w/Cap w/Cap	900A 900A
								Female Luer		SUUA
MPH1xx						None	None		w/Cap w/Cap	
MPH3xx						None	None	Male Luer		
MPH4xx	_	•	_	_		None	None Pellet	2.0-mm Port	w/Cap	Dicopozzlo Vi+ /E 420 VV
MPH6Pxx	•	•	•	•	Right	Male		Female Luer	w/Cap	Piconozzle Kit (5430-XX)
MPH6Rxx	•	•	•	•	Right	Male	Wire	Female Luer	w/Cap	Piconozzle Kit (5430-XX)
MPH6Sxx					_	None	None	Female Luer	w/Cap	Piconozzle Kit (5430-XX)


[★] Specify O.D. of glass (1.0, 1.2, 1.5 or 2.0 mm) by replacing XX in the Order Number with 10, 12, 15 or 20.


Handles and Accessories (not included)

Handle #2505 is for use with WPI manipulators. The smaller diameter handle #5444 is required for use with Narishige and Zeiss manipulators.


2505	1/4-in (6.3 mm) diameter handle
5444	³/16-in (4.8 mm) diameter handle
GO1-100	Replacement gasket 1.0 mm, Package of 100
GO2-100	Replacement gasket 1.2 mm, Package of 100
GO3-100	Replacement gasket 1.5 mm, Package of 100
GO4-100	Replacement gasket 2.0 mm, Package of 100
1571	Clear Silicone Rubber Sealant (-4.7 oz-)

MEH1RF 20.7 mm Gasket 2 mm Jack to Ag/AgCl Pellet √ 6.3 mm DIA. 1.0, 1.2, 1.5, 2.0


























































































П

M











111

Glass Capillaries

Clean, high quality glass for making micropipette electrodes and other research implements




Filament



Filament

WPI offers a wide spectrum of high-quality glass capillaries. We take pride in our ability to ship your glass order within 48 hours. If you need a special glass that does not appear in our catalog, please call us. We will make every effort to provide it for



Fire-Polished glass capillaries are easier to insert into microelectrode holders without damaging the gasket. More importantly, fire-polished glass won't scratch the chloridized wire used in a recording electrode. Fire-polishing does not affect the glass's mechanical or electrical properties.

Borosilicate glass capillaries: Close dimensional tolerances assure microelectrode uniformity and reproducibility. Available in one-, two-, three-, five- and seven-barrel configurations; a complete range of single barrel thin-wall sizes; and a variety of special configurations. Capillaries with filaments contain a solid filament fused to the inner wall, which speeds filling of electrodes. Capillaries with or without inner filaments are available for making microelectrodes in a wide range of diameters.

Single Barrel standard wall thickness capillaries are offered either with or without inner filaments for quick filling in a variety of lengths and diameters. Two usable electrodes can be made from one 6-inch length. Borosilicate glass is Corning N51A.

Thin Wall single barrel capillaries are offered both with or without inner filaments. The concentricity of this material provides excellent strength. Micropipettes made from thin wall capillaries have fine tips with a short taper.

> Note: Because electrode tips erode when left filled with saline solutions for long periods, electrodes should be made and filled immediately prior to use.

Single-	Barre	el Sta	ndard E	Borosili	icate Gl	ass Tubing
Length	OD (mm)	ID (mm)	Filament	Fire- Polished	Quantity	Item
3 in. (76 mm)	1.0	0.58	V		500	1B100F-3
3 in. (76 mm)	1.0	0.58			500	1B100-3
3 in. (76 mm)	1.2	0.68	V		350	1B120F-3
3 in. (76 mm)	1.2	0.68			350	1B120-3
3 in. (76 mm)	1.5	0.84	V		225	1B150F-3
3 in. (76 mm)	1.5	0.84		V	300	1B150-3
4 in. (100 mm)	1.0	0.58	V	V	500	1B100F-4
4 in. (100 mm)	1.0	0.58		V	500	1B100-4
4 in. (100 mm)	1.2	0.68	V	V	400	1B120F-4
4 in. (100 mm)	1.2	0.68			350	1B120-4
4 in. (100 mm)	1.5	0.84	V	V	300	1B150F-4
4 in. (100 mm)	1.5	0.84		V	300	1B150-4
4 in. (100 mm)	2.0	1.12	V		125	1B200F-4
4 in. (100 mm)	2.0	1.12		V	200	1B200-4
6 in. (152 mm)	1.0	0.58	V		500	1B100F-6
6 in. (152 mm)	1.0	0.58			500	1B100-6
6 in. (152 mm)	1.2	0.68	V		350	1B120F-6
6 in. (152 mm)	1.2	0.68			350	1B120-6
6 in. (152 mm)	1.5	0.84	V		225	1B150F-6
6 in. (152 mm)	1.5	0.84			225	1B150-6
6 in. (152 mm)	2.0	1.12	V		125	1B200F-6
6 in. (152 mm)	2.0	1.12			125	1B200-6

(Schott Duran) Glass Tubing								
OD (mm)	ID (mm)	FIL	Fire- Polished	Length	Quantity	Item		
1.0	0.75	~		3 in. (76 mm)	500	TW100F-3		
1.0	0.75			3 in. (76 mm)	500	TW100-3		
1.2	0.90	~	✓	3 in. (76 mm)	400	TW120F-3		
1.2	0.90			3 in. (76 mm)	350	TW120-3		
1.5	1.12	~		3 in. (76 mm)	225	TW150F-3		
1.5	1.12		V	3 in. (76 mm)	300	TW150-3		
1.0	0.75	~		4 in. (100 mm)	500	TW100F-4		
1.0	0.75		V	4 in. (100 mm)	500	TW100-4		
1.2	0.90	~		4 in. (100 mm)	350	TW120F-4		
1.2	0.90			4 in. (100 mm)	350	TW120-4		
1.5	1.12	~		4 in. (100 mm)	225	TW150F-4		
1.5	1.12		V	4 in. (100 mm)	300	TW150-4		
1.0	0.75	~		6 in. (152 mm)	500	TW100F-6		
1.0	0.75		V	6 in. (152 mm)	500	TW100-6		
1.2	0.90	~	V	6 in. (152 mm)	400	TW120F-6		
1.2	0.90			6 in. (152 mm)	350	TW120-6		
1.5	1.12	V		6 in. (152 mm)	225	TW150F-6		
1.5	1.12		V	6 in. (152 mm)	300	TW150-6		

Thin-Wall Single-Barrel Standard Borosilicate

Patch Clamp Capillary Glass

To select the best patch clamp glass for your experiments, we have evaluated available glass types in terms of the four properties most crucial to successful patch clamp studies:

The softening temperature determines how easily each glass type can be pulled to the desired shape and the extent to which it can be heat polished. Glass with a high softening temperature is difficult to pull and causes unnecessary wear on the heating element of the puller. This makes it very hard to make electrodes that are reproducible and of consistent quality. Patch clamp glass with a low softening temperature is preferred; however, higher softening temperature glass is stronger.

Electrical properties determine how much noise the glass is likely to produce in recording situations. The lower the product of dielectric constant times the loss factor, the smaller the equivalent noise current the glass will produce (Rae and Levis, Methods in Enzymology, 207, p67, 1992). Patch clamp glass with good electrical properties is critical especially in single-channel recording.

Sealability: It is not clear what factors determine the sealing ability of the patch to the glass. Almost any glass can form a gigohm seal under the right conditions. Different glass types vary, however, in how easily they form a seal. It is important to select a patch clamp glass that seals easily. Good fire polish is critical for seal (see DMF1000).

Leachable components: Substances leached from glass can alter channel behavior. Since different channels are sensitive to different glass components, it is best to record one type of channel with several different kinds of pipette glass to eliminate any artifact due to the glass.

WPI offers capillary tubing made from two glass

types widely used in constructing patch clamp electrodes. The significant characteristics of each are as follows:

PG52151-4, PG52152-4 and PG52165-4 are prepared from Schott #8250 glass (equivalent to Corning #7052), one of the most widely used patch clamping glasses. This is a specially formulated borosilicate glass with a softening temperature that is 110°C lower than regular borosilicate glass (Corning 7740, or Pyrex). It has excellent sealing properties for most cells. Electrical properties are also very good.

PG10150-4 and **PG10165-4** are composed of Corning #0010 glass, a high lead content (22% PbO) glass. Its thermal and electrical performance is between the Schott #8250 and Corning #8161 glasses described above. It is much more economical than Corning #8161 glass. It has been found that this glass causes much less alteration in channel behavior than Corning #8161 and Schott #8250 glass (Furman and Tanaka, Biophys. J. 53, p287,

Patch clamp capillaries do not have microfilaments.



PATCH CLAMP CAPILLARY GLASS								
Catalog#	Glass Type	OD/ID (mm)	Dielectric Constant	Softening Point °C	Quantity			
PG52151-4	#8250	1.5/1.0	4.9	720°	100			
PG52165-4	#8250	1.65/1.1	4.9	720°	100			
PG10150-4	#0010	1.5/0.75	6.7	625°	100			
PG10165-4	#0010	1.65/1.1	6.7	625°	100			


Glass Handling Forceps

Ever had difficulty picking up a glass capillary? Special tips on these forceps solve the problem, holding glass firmly without risk of breakage. They also keep the glass clean and avoid contamination from skin oils.







77020 Glass Handling Forceps



Three-Barrel



Five-Barrel



Seven-Barrel

Multi-barrel configurations are designed especially for microiontophoresis. Because the capillaries are fused together during manufacture, you will not need to twist them while pulling to seal the tips together. An inner filament in each barrel makes filling easy and fast.

Also see PolyFil for a novel way to connect multi-barrel pipettes

Multi-Barrel Borosilicate Glass Tubing with Filaments									
Length	Description	OD/ID (mm)	Filament	Quantity	Item				
4 in. (102 mm)	Two-Barrel	1.5/0.84	V	100	2B150F-4				
4 in. (102 mm)	Three-Barrel	1.2/0.68	V	100	3B120F-4				
4 in. (102 mm)	Five-Barrel	1.2/0.68	V	65	5B120F-4				
4 in. (102 mm)	Seven-Barrel	1.2/0.58	V	60	7B100F-4				
4 in. (102 mm)	Seven-Barrel	1.2/0.68	V	75	7B120F-4				
6 in. (152 mm)	Two-Barrel	1.5/0.84	V	100	2B150F-6				
6 in. (152 mm)	Three-Barrel	1.2/0.68	V	100	3B120F-6				
6 in. (152 mm)	Five-Barrel	1.2/0.68	V	65	5B120F-6				
6 in. (152 mm)	Seven-Barrel	1.0/0.58	V	60	7B100F-6				

Special Configuration Borosilicate Capillaries


Septum Theta offers superior cell impalement. The natural bevel resulting from the prominent spear-like projection of the septum gives microelectrodes a sharp, spear-point tip. This style has low resistance for use as a single microelectrode, and it can be used to make superior double-tipped microelectrodes with low trans-tip coupling. The natural bevel of Septum Theta also significantly increases the effective tip

cross-section. As supplied, the width of the septum is approximately 0.2 mm; wall thickness is approximately 0.2 mm.

Piggyback glass consists of a pair of borosilicate capillaries fused together during manufacture. One barrel is larger than the other, and both have inner filaments for quick filling. Piggyback glass makes it simple to fabricate two-barrel electrodes with a significant tip diameter differential.



Septum Theta



Piggyback

Special Configuration Borosilicate Glass Tubing

					_
Description	OD/ID (mm)	Length	Quantity	Item	
Septum Theta	1.5/1.02	6 in. (152 mm)	100	TST150-6	
Piggyback	1.51/0.84 0.75/0.35	4 in. (102 mm)	50	PB150F-4	
Piggyback	1.51/0.84 0.75/0.35	6 in. (152 mm)	50	PB150F-6	

Borosilicate glass rod

1.0 mm diameter — for making tools, probes, tips

Borosilicate Glass Rod						
Description	OD (mm)	Length	Quantity	Item		
Glass Rod	1.0	4 in. (102 mm)	500	GR100-4		
Glass Rod	1.0	6 in. (152 mm)	500	GR100-6		

Micropipette Storage Jar

Stores up to 30 micropipettes, filled or unfilled, up to three inches in length. A gentle sliding action inserts or removes pipettes without damage to the delicate tips.



E210	Storage Jar for 1.0 mm OD Micropipettes
E212	Storage Jar for 1.2 mm OD Micropipettes

Storage Jar for 1.5 mm OD Micropipettes **E220** Storage Jar for 2.0 mm OD Micropipettes

DEDLACEMENT DADTO

KLFLF	ACLIVILITY FARTS
1965	Foam Ring for 0.75 - 1.0 mm glass
1966	Foam Ring for 1.2 - 1.5 mm glass
1967	Foam Ring for 2.0 mm glass



WPI's MicroFil™ fills micropipettes easily and reliably. Its long and fine tip allows you to start the filling very close to the pipette tip, eliminating both air bubble formation and clogging due to the washing down of dust particles. The transparent amber MicroFil needle is constructed from a combination of plastic and fused silica — no metal components are used. The MicroFil needle can be stored for days with the filling solution inside without clogging.


The MicroFil's tip elasticity is sturdy and very flexible though not unbreakable. Since it is more flexible than stainless steel needles, moderate bending will not block or damage the MicroFil needle. The combination of plastic and fused silica in the MicroFil tip is sturdier than plastic tips, allowing easy and repeated insertions into micropipettes. MicroFil's luer fitting allows easy coupling to syringes and syringe filters.

		1-5 pkgs	6-10 pkgs
MF34G-5	MicroFil, 34 ga., 67 mm long (pkg of 5)		
MF28G-5	MicroFil, 28 ga., 97 mm long (pkg of 5)		
MF28G67-5	MicroFil, 28 ga., 67 mm long (pkg of 5)		

CUSTOM MICROFIL

All MicroFil products, including custom orders, can be shipped immediately. Custom orders for special needs can be made using nine sizes of MicroFil tubing in lengths up to 50° cm — except for CMF90UxxL which has a maximum length of 10 cm because of its high resistance to flow. Quantity discounts available. Specify length when ordering by inserting the length (in centimeter increments) into the catalog number in place of the XX's.

CMF20GxxL	MicroFil, 20 Gauge, 700 μm ID, 850 μm OD (pkg of 4)
CMF22GxxL	MicroFil, 22 Gauge, 530 μm ID, 700 μm OD (pkg of 4)
CMF23GxxL	MicroFil, 23 Gauge, 530 μm ID, 665 μm OD (pkg of 4)
CMF26GxxL	MicroFil, 26 Gauge, 320 μm ID, 430 μm OD (pkg of 4)
CMF28GxxL	MicroFil, 28 Gauge, 250 μm ID, 350 μm OD (pkg of 4)
CMF31GxxL	MicroFil, 31 Gauge, 100 μm ID, 238 μm OD (pkg of 4)
CMF34GxxL	MicroFil, 34 Gauge, 100 μm ID, 164 μm OD (pkg of 4)
CMF35GxxL	MicroFil, 35 Gauge, 75 µm ID, 144 µm OD (pkg of 4)
CMF90UxxL	MicroFil, approx. 36 Gauge, 20 µm ID, 90 µm OD (pkg of 4)





PolyFil allows easy and secure coupling of a multi-barrel micropipette to a pressure source. Coupling is achieved by bonding temperature-resistant and flexible MicroFil to the capillary tube with hot melt adhesive. The luer end of each MicroFil is connected to PVC tubing (200 PSI rated). Kits also include a five-port manifold that allows use of a single PV800 Series PicoPump to drive up to six micropipette barrels independently by switching on only the barrels to be injected. All connections are locking luers — pressure safe and convenient.

Kit includes: 1 pipette holder/handle, plastic; 7 pieces MF28G MicroFil; 7-pieces tubing with male luer lock fittings; 1 flow-thru manifold with five luer lock ports; 1 hot melt glue gun(110V only); 3 glue sticks.

5440	PolyFil Multi-Barrel Micropipette Coupling Kit
13316	Mini Glue Gun and (3) glue sticks



U ip Borosilicate glass micropipettes



Eliminate the cost and trouble of making your own micropipettes — WPI can quickly supply your need for consistently sized pre-pulled glass micropipettes for injection of dyes or proteins into cells, oocytes and for many other biomedical laboratory applications. Tip diameters (ID) range from 0.1 to 10 micrometers.

- Schott Duran borosilicate glass
- 0.5 micrometer and smaller ID micropipettes include an internal glass fiber for easy filling
- Tip inner diameter tolerance ±20%
- Short taper yields high strength
- Nominal length ≈ 50 mm
- OD:ID = 1.33:1
- Standard capillary outer diameters are 1.0 mm (thin-wall) or 1.14 mm
- Every pipette individually tested and inspected
- Vacuum packed

Silanized Tips (Luer Shank)

Silanization waterproofs the glass to retard water when inserting into cell. This will not let the outside fluid run down the pipette and get inside so easily.

Shank	Tip I.D.	Shank Length	Glass O.D.	Filament	Fire Polished	Catalog #
PLAIN	0.1 µm	_	1.0 mm Thin-Wall	Yes	No	TIP01TW1F
	0.2 µm	_	1.0 mm Thin-Wall	Yes	No	TIP02TW1F
	0.3 µm	_	1.0 mm Thin-Wall	Yes	No	TIP03TW1F
	0.4 µm	_	1.0 mm Thin-Wall	Yes	No	TIP04TW1F
	0.5 µm	_	1.0 mm Thin-Wall	Yes	No	TIP05TW1F
	1 µm	_	1.0 mm Thin-Wall	No	Yes	TIP1TW1
	2 µm	_	1.0 mm Thin-Wall	No	Yes	TIP2TW1
	5 µm	_	1.0 mm Thin-Wall	No	Yes	TIP5TW1
	10 µm	_	1.0 mm Thin-Wall	No	Yes	TIP10TW1
	10 µm	_	1.14 mm A203XV glass *	No	Yes	TIP10XV119
	30 µm	_	1.0 mm Thin-Wall	No	Yes	TIP30TW1
.UER	0.1 µm	_	1.0 mm Thin-Wall	Yes	_	TIP01TW1F-L
	0.2 µm	_	1.0 mm Thin-Wall	Yes	_	TIP02TW1F-L
	0.3 µm	_	1.0 mm Thin-Wall	Yes	_	TIP03TW1F-L
	0.5 µm	_	1.0 mm Thin-Wall	Yes	_	TIP05TW1F-L
	1 µm	_	1.0 mm Thin-Wall	No	_	TIP1TW1-L
	2 µm	_	1.0 mm Thin-Wall	No	_	TIP2TW1-L
	5 μm	_	1.0 mm Thin-Wall	No	_	TIP5TW1-L
	10 µm	_	1.0 mm Thin-Wall	No	_	TIP10TW1-L
	30 µm	_	1.0 mm Thin-Wall	No	_	TIP30TW1-L
.UER/SILANIZED	5 μm	1 inch	1.0 mm Thin-Wall	No	_	TIP5TW1LS01
	5 µm	2 inch	1.0 mm Thin-Wall	No	_	TIP5TW1LS02
	10 µm	1 inch	1.0 mm Thin-Wall	No		TIP10TW1LS01
	10 µm	2 inch	1.0 mm Thin-Wall	No	_	TIP10TW1LS02
	30 µm	1 inch	1.0 mm Thin-Wall	No	_	TIP30TW1LS01
	30 µm	2 inch	1.0 mm Thin-Wall	No	_	TIP30TW1LS02

^{* 10} μ (ID), 1.14 mm capillary pipettes are for use in WPI's **Nanoliter 2000**.

μTIP SAMPLER	μTIP SAMPLER ASSORTMENTS		
TIPMIX01-05	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 µm ID, plain shank		
TIPMIX05-10	Two each, 0.5, 1, 2, 5, 10 µm ID, plain shank		
TIPMIX01-05-L	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 μm ID, Luer		
TIPMIX05-10-L	Two each, 0.5, 1, 2, 5, 10 µm ID, Luer		

Micro Cannula

- 0.4mm O.D., 0.2mm I.D. tubing
- Autoclavable
- Biocompatible Perfluorocarbon tubing material

KZ1101

Micro Cannula, 3-inch

This micro cannula is ideal for placement in the carotid or femoral artery of mice, rats, and other small animal blood vessels. It can be used with a pressure transducer (WPI'S BLPR2) for blood pressure measurement, or in conjunction with a micro-syringe injection system (like WPI's UMPIII or MMP pumps). The incorporated standard female luer fitting makes connecting to existing experimental plumbing quick and easy. The cannula is provided with a contoured-tip stainless steel stylet (trocar) to facilitate placement using established techniques. A movable "shoulder" ring provides a tie-in point to prevent accidental removal. The cannula may be left in place for 2 hours or more, and with proper care and cleaning, may be re-used multiple times. Instructions for use included.

Sometimes the simplest designs work best.

The MF200 Microforge is a versatile instrument designed specifically for the fabrication of glass micropipettes and other related tools. The system was developed in collaboration with Dr. Ming Li of the Department of Pharmacology, University of South Alabama. It is perfect for patch pipette tip polishing, tip size reduction, contact stretching, in vitro fertilization pipette production and a variety of other pipette configurations. The MF200 simple, reliable and is priced economically.

Features of the MF200

The MF200 system includes: An easy to use analog temperature controller, a specially configured WPI model W30S-LED research grade compound microscope, 40x longworking distance objective and 10x eyepiece. 40x magnification is essential when polishing pipettes as small as half a micron (0.5 μm) in diameter. Compared to a conventional 40x objective, the long working distance objective reduces the danger of damage to the pipette and/or objective lens during the polishing process. It is also the only commercial microforge using the Kohler illuminator and Abbe condenser for illumination. This provides less glare and sharper image of the pipette than frosted glass illuminator, which was used on all of the other commercial Microforge.

MF200 Microforge







	S
MF200-1	Complete Microforge System, incl. W30S-LED Microscope (110 v)
MF200-2	Complete Microforge System, incl. W30S-LED Microscope (220 v)
MF200-M1	MF 200 without microscope (110v)
MF200-M2	MF 200 without microscope (220v)
*Above MF20	0 microforges include 40X long working distance objective
OPTIONAL	ACCESSORIES
500292	Optional 15× Eyepieces (pair)
	Note: No reticle available for 15x eyepieces
500329	25× Long-Working Distance Objective
	(fits most microscopes with a 160 mm Focal Length)
13142	Optional foot switch
REPLACEME	NT ACCESSORIES
MF200-H2	Replacement heating filament (large gauge)
MF200-H3	Replacement heating filament (medium gauge)
MF200-H4	Replacement heating filament (small gauge)
75090	Filament Adjustment Assembly for 22mm OD Objectives
75050	Replacement Micropipette Slide
75040	Replacement Filament Cable

MF200 SPECIFICATIONS

100-240 VAC 50/60 Hz AC POWER MODULE FILAMENTS (3) H2. H3. H4

FILAMENT ON Pushbutton Controlled or Optional Foot

Switch Controlled FILAMENT ADJUSTMENT ASSEMBLY For 40x and 25x Long-Working

Distance Objectives: mounts on objective **OBJECTIVE** 40× Long-Working Distance (3 mm)

OPTIONAL 25× Long-Working Distance (5 mm) **EYEPIECE**

RETICLE (10x eyepiece only)

1.25 µm/division (at 40x) 0-90° Angle at 5°/division

OPTIONAL EYEPIECE 15x (pair) GLASS HOLDER

Mounts on Microscope Stage DIMENSIONS: Control Unit $4 \times 7 \times 1\%$ in. $(10.2 \times 17.8 \times 4.8 \text{ cm})$ SHIPPING WEIGHT 3 lb. (1.4 kg)

MICROSCOPE See W30S 16 lb. (7.3 kg) SHIPPING WEIGHT

MICROFORGE FEATURE COMPARISON

	MF200	DMF1000
W30S-LED Microscope	· ·	V
40x Long Working Distance Objective	V	V
Analog Controller	Marian P	Pour Advisoration Continue
Digital Controller	- // =	
Pressurized Air Control	-	
Microinjection Capability	-	
Optional Foot Switch	V	
Memory	-	
Auto-sense of Filament Type	_	V
Digital Temperature Control	_	V

117

US: Tel: 941-371-1003 • sales@wpiinc.com

Digital Signal Processor (DSP) Technology

The DMF1000 is powered by the latest digital signal processor (DSP) technology. A digital timer is used to precisely control the polish heating time. Ten memories can be used to store settings of the heating power and heating duration. All of the settings are controlled and displayed digitally for better accuracy and reproducibility. Two different operating modes are provided: Manual and Auto. In the Manual mode. the DSP will memorize the duration of the time that is used to achieve a desired polishing. In Auto mode, the heat will be applied for the duration of the timer




The DMF1000 system includes a specially configured WPI model W30S-LED research grade compound microscope equipped with a high quality metallurgic 40x long-working distance objective and a pair of 10x evepieces. It is the most powerful long-working distance objective currently available on any commercial microforge. The long working distance objective reduces the danger of damage to the objective lens during the heating process.

Other benefits of the DMF1000 design include the use of a Kohler illuminator and Abbe condenser, which provide the reduced glare and sharper image contrast necessary when polishing pipettes as small as half a micron (0.5 µm) in diameter.


Pressure Polishing


The DMF1000 incorporates a unique digital pneumatic pressure feature that enables pressurized air to be delivered through the pipette during fire polishing. In the fabrication of patch pipettes, the pressurized air can be used to blunt the taper at the pipette tip without changing the size of the tip opening. This reduces electrical resistance of the tip, leading to lower noise during patch-clamp recordings (Goodman & Lockery,




stretching)


Fire Large Tip **Polishing Sharpening** (contact







Tip Bending (holding pipettes)



Sealing in **Plastic Sensor**



Filament Holder mounts directly to objective to provide precise control of heating element position.

Ease of use

The Heating Filament

With a conventional microforge often the most difficult and time-consuming part of using a high magnification objective is being able to move both the heating filament and the pipette into the same viewing area. Finding and moving both the heating filament and the pipette without collision can be a challenge. However, this difficulty is eliminated with the DMF1000 because the heating filament is directly attached to the microscope's objective. Hence it can be easily adjusted to any position within the viewing

The low heat capacity and low thermal coefficient of linear expansion of the filaments are key design features of the DMF1000. The low

heat capacity of the filament allows it to reach fire-polishing temperatures without excessive heat. This permits the user to bring the pipette tip close to the filament during polishing without fear of collapsing the pipette tip. Low heat capacity eliminates the need for an auxiliary aircooling system. The low coefficient of expansion characteristic of the filament ensures minimal displacement of the filament during heating. This feature eliminates much of the guesswork out of tip placement in relation to the filament.

Two different heating filaments are provided with the DMF1000 to accommodate various applications. The **H5** filament is large gauge and can be reformed into a "U" for fabrication of pipettes, air forming of patch pipettes and other applications. The H4 is a smaller gauge filament and is ideal for polishing patch clamp pipettes.

The Pipette and Microscope Stage

The pipette rests on a specially designed holder that sits on top of the microscope stage. The position of the pipette, relative to the heating filament, is controlled by the (X, Y, Z) adjustment of the stage. This unique design makes locating and polishing the pipette extremely easy. The stage of the microscope has a high quality rail that gives precise, smooth and stable control of the pipettes movement. This configuration also eliminates the need and expense of an additional micromanipulator to control pipette

Typical applications of the DMF1000

Polishing the Patch Pipettes

It is well known that the proper fire polishing of patch pipettes is the single most important factor for forming a stable giga-seal in patch clamp recording. This is even more important than the type of glass capillary used. Difficulties often arise in forming giga-seals because the polishing of patch pipettes using a conventional low magnification microforge is inadequate. However, since the DMF1000 uses a 40X long-working distance objective, pipette polishing is much more accurately controlled. Pipettes polished using the DMF1000 achieve excellent stable giga-seals with a wide variety of cells. Both whole cell patch pipettes and single channel patch pipettes can be conveniently polished with the DMF1000 to the highest quality and reproducibility achievable with any microforge.

For the single-channel patch clamp pipettes the pipette needs to be pre-coated with Sylgard 184 before polishing. For this procedure the user can follow a simple and effective coating method described previously (Li, 1999)

Microforging Holding Pipettes

A holding pipette with a large blunt tip and a small opening is used to hold a floating cell in place prior to microinjection by applying suction to the rear of the pipette. The procedure for making holding pipettes involves three steps: squaring off, large bore flame polishing, and tip reducing. These steps are accomplished with a larger heating filament.

Microforging Beveled Injection Pipettes

Occasionally, a beveled large bore pipette is not sharp enough to penetrate a cell without damaging the area around the pipette. With the DMF1000 and the large heating filament, a sharp point can be formed on the beveled tip to assist the penetration of the cell. This process is referred to as contact stretching

Pipette Tip Calibration & Microinjection

The integrated digital pneumatic pressure system can be used to calibrate the precise diameter (I.D.) of a micropipette tip, based on a technique described previously (Hagag & Randolph 1990, Bowman & Ruknudin 1999). The pressure system can also be used separately as a simple but highly accurate controller for microinjection applications.



DMF1000-1	Complete Microforge, incl W30S-LED Microscope (110 v)
DMF1000-2	Complete Microforge, incl W30S-LED Microscope (220 v)
DMF1000-M1	Microforge without microscope (110v)
DMF1000-M2	Microforge without microscope (220v)
*Above DMF100	0 microforges include 40X long working distance objective

75040	Replacement Filament Cable
75050	Replacement Micropipette Slide
MF200-H4	Replacement heating filament (small gauge)
DMF1000-H5	Replacement heating filament (large gauge)
503513	21 mm 10X Eyepiece with 100/10 reticle
800292	40x Long Working Distance Objective, 3 mm 0.25NA
REPLACEMEN [®]	T ACCESSORIES
13142	Optional foot switch
500292	Optional 15x Eyepiece (pair)
500329	25x Long Working Distance Objective, 5 mm 0.50NA
OPTIONAL AC	CESSORIES
^ADOVE DIVIFIUU	0 microforges include 40X long working distance objective

Professional-Grade Microscope

The **W30** professional-grade microscope is a best-seller in universities, medical schools, and reseach laboratories. Equipped for performance, its features include titaniumfinished DIN or Semi-Plan optics and a 30-year anti-fungal coating. The W30 is the choice for superior performance at a great price.



W30S-LED	Binocular Microscope
W30ST-LED	Trinocular Microscope
503513	21 mm 10X Evepiece with 100/10 reticle

DMF1000 SPECIFICATIONS

AC POWER MODULE 100-240 VAC 50/60 Hz TIMER RANGE (for heater & timer) 0.01 to 360 sec NUMBER OF MEMORYS

PRESSURE ADJUSTING RANGE 0.5 - 60 PSI (3.5 - 414 kPa

PRESSURE RESOLUTION 0.1 PSI (0.7 kPa)

FILAMENTS H4 — Small filament for working with 40×

long working distance objective.

H5 — Large filament for working with 10× objective. Filament adjustment assembly

provided for both objectives. Auto or Manual via Pushbutton, TTL,

HEATER AND TIMER CONTROL

or Optional Foot switch.

DIMENSIONS: Control Unit $4 \times 7 \times 178$ in. $(10.2 \times 17.8 \times 4.8 \text{ cm})$

SHIPPING WEIGHT 4 lb. (1.8 kg) W30S-LED (see below) MICROSCOPE SHIPPING WEIGHT 16 lb. (7.3 kg)

W30S-LED SPECIFICATIONS

HEAD Binocular (Seidentopf) Inclined 30°, rotates 360°

Dual diopter adjustment, Interpupillary distance

range 55-75mm

10X/18 wide field eyepieces

NOSEPIECE Quadruple forward-facing nosepiece

OBJECTIVES DIN Plan, anti-fungal 4X, 10X, 40X, 100XR (oil)

Parfocal, parcentric, color-coded

Mechanical stage (140mm x 140mm) **STAGE**

Coaxial drive controls

XY Movement: 73mm x 43mm

FOCUS Coarse adjustment: range of 30mm

Fine adjustment: graduation of 2µm

Tension control knob

ILLUMINATION Moveable Abbe condenser, NA 1.25,

Iris diaphragm

Variable LED light source (3W bulb) 110V/220V switchable electronics

DIMENSIONS AND WEIGHT 15" (38cm) x 9" (23cm) x 7" (17.8cm) (h x l x w)

14 lbs. (6.4kg)

World Precision Instruments

www.wpiinc.com



An optically-flat mirrored glass disk, wetted with an abrasive slurry, spins at 60 rpm (120 V), producing sharply beveled tips on fluid-filled glass microelectrodes of one micron or smaller. This eases cell impalement and improves the electrode's linearity. The microelectrode's resistance can be monitored during beveling with WPI's Ωmega-Tip-Z™ megohm meter (see page 81). The beveler is permanently mounted on a precision magnetic plate that gives stable support for the optional 1350M Micropositioner shown. Start-up kit includes 0.05 µm alumina abrasive powder #3531, wick electrode, and wick support.

SYS-1300M Microelectrode Beveler & Start-Up Kit

OPTIONAL ACCESSORIES

(micropositioner not included)

Specify line voltage.

2478	Replacement Mirrored Disk Alumina Abrasive, 0.05 μm (5 g) fine		
3531			
3551	Alumina Abrasive, 0.30 μm (5 g)		
2479	Replacement "O" Ring		
SYS-OMEGAZ	Ωmega-Tip-Z with Probe & Holder		
1350M	Micropositioner (M3301R) and M10 Magnetic Stand		
711P	Replacement Probe		
5468	Adapter to connect metal microelectrodes to probe,		
	2 mm socket to .031 in. receptacle		
Z-LITE	Z-Lite Fiber Optic Illuminator (115 V, 60 Hz, beige case)		
Z-LITE-Z	Z-Lite Fiber Optic Illuminator (230 V, 50 Hz, black case)		
500186	Bifurcated Light Guide with lenses		
Z-LITE-186	Z-Lite Illuminator & Bifurcated Light Guide		
Z-LITE-Z186	Z-Lite Illuminator & Bifurcated Light Guide		
MES	Microelectrode Beveler System		

Model 1350M Micropositioner — This package (shown with beveler above) includes WPI's M3301R Manipulator and an M10 magnetic stand. The stand-manipulator assembly mounts directly onto the beveler baseplate, allowing convenient positioning of electrodes onto the beveling surface. Three axes of adjustment, including coarse and fine control in the axis of the electrode.

1300M SPECIFICATIONS

BEVELING SURFACE 7.8 cm diameter, optically flat reflective glass

MAXIMUM BEVEL

ALUMINA ABRASIVE POWDER $0.05~\mu$ size supplied (0.3 μ also available) **RPM** 60 rpm at 120 V, 60 Hz; 50 rpm at 240 V, 50 Hz

MOTOR AC synchronous

POWER REQUIREMENTS 95-135 V or 220-240 V, 50/60 Hz

DIMENSIONS

Steel base plate 8.5 x 11 x 0.375 in. (22 x 28 x 1 cm)

Overall height 4 in. (10 cm)

Height of abrasive surface 2.75 in. (7 cm) above base plate

SHIPPING WEIGHT 20 lb (9.1 kg)

MES includes: 1300M Microelectrode Beveler; 1350M Micropositioner & Magnetic Stand; OmegaZ; 5052 Steel Base Plate; 5468 Adapter; 3485 Ringstand Mount. Shipping Weight: 59 lb. (27 kg)



Bevel micropipette tips larger than 1 micron, for applications such as microinjection

- Tool holder on microscope keeps pipette in focus during beveling
- Steel base provides solid support for beveler and other magnetic stands
- Includes stereo zoom microscope **PZMIII** with up to 90x magnification
- Variable speed, reversible
- Abrasive surface can be easily replaced by several types of Diamond and Alumina Lapping Film
- Pipette tip illuminated internally via fiber optic illuminators

The advantage of WPI's MicroBeveler over other types of solidsurface bevelers is that the abrasive surface can be easily refreshed. Instead of using a conventional solid abrasive disk, the MicroBeveler abrasive surface is made of a high quality lapping film, widely used in the fiber optics industry. When the surface is damaged or loaded up with glass particles, the abrasive film can be easily replaced.

The solid polishing surface of WPI's new MicroBeveler, turning at 5,400 rpm, provides sufficient cutting force for a very sharp tip in a very short time. The cutting surface is very flat and turns very smoothly, ensuring an undamaged tip.

48000 SPECIFICATIONS

BEVELING SURFACE 3.5 inch diameter disk ABRASIVE MATERIAL alumina, diamond SPEED OF ROTATION 100 to 5,400 rpm **MOTOR** Reversible Direction POWER REQUIREMENTS 120 volts, 60 Hz or 240 volts, 50 Hz, 20 VA to supplied power supply

DIMENSIONS

 $8.7 \times 11 \times 0.4$ in. (22 × 28 × 1 cm) Base Plate Overall Height 3 in. (8 cm) SHIPPING WEIGHT (48000) 16 lbs. (7.6 kg.) SHIPPING WEIGHT (MBS) 35 lbs. (16 Kg.)

MBS MicroBeveler System

> Includes 48000 MicroBeveler, Z-LITE illuminator, fiber optic cable, PZMIII Stereo Zoom Microscope with tilting base especially adapted for use with MicroBeveler, two clear 20x eyepieces, one 20x eyepiece with reticle, tool holder, and pipette holder FOIMPH.

SYS-48000 MicroBeveler

Specify line voltage

OPTIONAL	ACCESSORIES
48015-03	Lapping Film, Alumina, 0.3 micron (50-pack)
48015-10	Lapping Film, Alumina, 1 micron (50-pack)
48015-30	Lapping Film, Alumina, 3 microns (50-pack)
48014-01	Lapping Film, Diamond, 0.1 micron (3-pack)
48014-05	Lapping Film, Diamond, 0.5 micron (3-pack)
48014-10	Lapping Film, Diamond, 1 micron (3-pack)
48014-30	Lapping Film, Diamond, 3 microns (3-pack)
48025	Fiber Optic Cable for Pipette Illumination
15934	Replacement Beveler Disk Plate
48300	Tilt Base Assembly for PZMIII binocular head
48200	PZM Tool Holder



Store up to 30 programs in memory

Two factory programs installed

Tempered glass cover to reduce the effects of humidity on puller reproducibility

Switchable power supply ensures that line voltage fluctuations don't affect reproducibility

PUL-1000 SPECIFICATIONS

HEATER ELEMENT...... Platinum/Iridium PULLING FORCE Solenoid, adjustable CAPILLARY OD RANGE 1.0-1.5mm MAXIMUM CAPILLARY LENGTH 170mm MINIMUM CAPILLARY LENGTH..... 55mm POWER 90-240VAC, 50/60 Hz

REPLACEMENT FILAMENTS......... 2.5 mm Square Box Filament,

2.5 mm wide (# 13834)

PUL-1000 is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, patch clamp studies, microperfusion or microinjection. The puller was designed with tighter mechanical specifications and precision electronics for complete control of the pulling process and accurate reproducibility. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. This allows graduated cycles for applications like patch clamp recording.

This puller is a reasonably priced, compact, versatile and reliable workhorse. The microprocessor, combined with the LCD display, makes the PUL-1000 easy to use.

Construction

The cover of pulling chamber is made with tempered glass to minimize the humidity effect on the reproducibility of pulled pipettes.

Power Supply

PUL-1000 has a high quality switching power supply for use anywhere in the world without worry about the line voltage differences. Pulling reproducibility is unaffected by line voltage fluctuation. Heating voltage can be controlled to within 0.1% accuracy even when line voltage fluctuates from 90 to 240VAC.

PUL-1000	Micopipette Puller	§3595
13834	Replacement Filaments	[§] 35



Programmable Multipipette Puller

- Produces two identical pipettes every time
- 25 saved programmable sequences
- Optical-digital measurement
- 22 manufacturer preset sequences
- Pull patch clamping pipettes, intracellular electrodes, injection micropipettes and micro-needles
- Pneumatic pressure (instead of gravity or magnetic fields) yields consistent pulling force

The PMP-102 is a microprocessor controlled pipette puller. The PMP-102 is designed to pull a pipette horizontally to produce two identical micropipettes. Different kinds of pipettes can be pulled repeatedly using the preset program sequences. You can create your own patch clamp electrodes, intracellular electrodes, injection micropipette and micro-needles using preset sequences.

Exclusive Optical-Digital Taper Measurement

Instead of the mechanical tip length setting like other pipette pullers, there is an exclusive optical-digital ruler in the PMP-102 to apply precise taper length settings, and for real-time measurement and control. With

this feature, you can handle taper pulling precisely and easily. Equipped with a powerful computerize tip sensing function, the puller can automatically finish the tip pulling. You can pull the ideal tip every time.

Computerize Real-Time Feedback Heater Control

The PMP-102 includes an advanced microcontroller to perform real-time heater monitoring and control. When you select a heating level, the microcontroller measures the actual heating power applied during a pull. The real-time measurement is displayed, and the feedback to the control unit dynamically adjusts the power to match your setpoint. As a result, the puller always provides precise heating power, despite variables like thermal/ electrical changes. When the heating level is set to AUTO, the heater automatically determines the melting point for different glass pipettes. The microcomputer control ensures smart, reliable heating.

Programmable, Savable Sequences

There are 25 manufacturer/user programmable pulling sequences with 18 steps in each sequence. You may easily program different pipette tip sizes, tip lengths and tip shapes in different sequences for a variety of applications. Time and length count, heat level, heat control and action parameters can be individually set in each step.

Pneumatic Pulling Force and Compact Size

Other pipette pullers use gravity or magnetic fields as the pulling force. The PMP-102 applies precisely controlled pneumatic pressure as the pulling force, which gives more control, plus even and consistent dragging characteristics. With double horizontal pulling, the PMP-102 can pull two identical injection tips or microelectrode tips at the same time. The PMP-102 is compact, requiring little bench space. And, it precisely and automatically performs multi-step pulling without manual interruption. A precision micro-linear ball bearing rail and advanced pneumatic components are used to provide no fault pulling movement. A simple 4x4 keypad and a full information display LCD let you easily control and monitor the pulling parameters directly. Parameters include sequences, steps, time, timing, heater level, heater control, tip length and actions. The heater power control and action time count up or down in real-time on the display.

With the versatility of the intelligent PMP-102, pulling an ideal micropipette is no longer an uncertainty of hand skill, but a reproducible, automatic process.

Preset Programs

The PMP-102 pipette puller is well tested and comes with pre-installed programs. The User Manual includes sample pictures along with the characteristics of the pipettes you can expect from each program. Select your sequence and press START. Save time and money by using presets. If you can't find exactly what you need in the presets, you can use the presets as templates for creating new programs.

Micromanipulators

Micromanipulators are used when precision work is conducted under a microscope. A micropipette, electrode or probe can be mounted on a micromanipulator and move as little as a micron at a time. This tool can be used for *in vitro* fertilization, patch clamp experimentation, extracellular recording, microinjection and any application requiring fine mechanical placement (resolution). In addition to micromanipulators, WPI offers tilt bases, piezo translators and a variety of

Micromanipulators can be broken out into three broad categories: Manual, Manual/ Motorized and Motorized.

Micromanipulator	Manual or Motorized	Resolution	Travel	Stands	Tilt Base	Piezo Translator	Notes
SN-PZ-50	Motorized	0.5 µm resolution 30nm min. step					
SM325	Motorized	25nm/step 40,000 steps/rev	25mm (3 axes)	M9, M10, M10L, 501622, 501623	TBS, M-3		Use MCL3 Controller
DC3001*	Manual/ Motorized	Motor–0.5µm Manual–0.1mm	Motor–10mm (3 axes) Manual– X: 37mm Y, Z: 20mm	M9, M10, M10L, 501622, 501623	TBS	MPM10 MPM20 (with STM3 joystick)	Controller is required. Options: • MS314 (DC3314) • MPM10 • MPM20 with STM3 joystick
HS6-3	Manual/ Motorized	10nm/step	25mm (3 axes)				Use MCL3 Controller
M3301*	Manual	0.01mm (X fine) 0.1mm (X,Y, Z)	X(fine): 10mm X: 37mm Y: 20mm Z: 25mm	M9, M10, M10L, 501622, 501623	TBS, M-3	MPM20	
KITE	Manual	0.1mm	X(fine): 10mm X: 35mm Y,Z: 20mm	M9, M10, M10L, 501622, 501623	TBS, M-3		
M325	Manual	10µm	X: 25mm Y,Z: 10mm	M9, M10, M10L, 501622, 501623	TBS, M-3		
ммј	Manual	0.1mm	X: 37mm Y: 20mm Z: 25mm	M9, M10, M10L, 501622, 501623	TBS, M-3		Joystick control
MD4	Manual	10μm (X fine) 100μm (X,Y, Z)	X(fine): 10mm X: 37mm Y: 20mm Z: 25mm	M9, M10, M10L, 501622, 501623	TBS, M-3		Holds two electrodes
HS6	Manual	5µm	25mm (3 axes)				
MM3-3	Manual (mini)	1.5µm	13mm				340g load
MM1-3	Manual (mini)	1.0µm	3mm				225g load

Smallest micromanipulator features full 20 mm of piezo-movement in all axes!

- Ultra-stable and precise electrode positioning
- Smooth movement with minimal vibration
- Penetration mode with high acceleration single steps
- Features true electrode angle adjustment — 4th axis option available soon.
- Back-flip mechanism for easy electrode exchange
- Powered off in standby for zero noise interference
- Ergonomic controller can operate up to 14 manipulators



SN-PZ-50 SPECIFICATIONS

SN-PZ-50 MICROMANIPULATOR

POSITIONING RANGE 20x20x20 mm3 (x-y-z) MINIMUM STEP SIZE 30 nm

CLOSED-LOOP CONTROL 0.5 µm resolution and

3 µm repeatibility for absolute position

~3 mm/s

MAXIMUM LOAD 200 g* TABLE MOUNTING magnet or bolt ELECTRODE DRIVE ANGLE 0 - 35 deg

BACK-FLIP MECHANISM

MAXIMUM SPEED

DIMENSIONS** 32x40x80 mm WEIGHT** 260 g

SN-PCZ CONTROL UNIT

Backlit display

Three rotary knobs with optical encoders Can operate up to 14 micromanipulators*

Two memory positions for each micromanipulator **DIMENSIONS** 190 x 210 x 40 mm

WEIGHT 510 g

110-240 V; 50-60 Hz *1-to-8 connector hubs available for connecting more than

two manipulators

*Load balancing options available for loads exceeding 50 g

**Dimensions and weight without head-stage or electrode holder adapter; axis in initial positions

SN-PCZ-50L	Miniature Piezo Micromanipulator (left-hand version) and Controller
SN-PCZ-50R	Miniature Piezo Micromanipulator (right-hand version) and Controller
SN-PZ-50L	Miniature Piezo Micromanipulator (left-hand version)
SN-PZ-50R	Miniature Piezo Micromanipulator (right-hand version)

Specify line voltage

Battery-operated control unit with rotating knobs and information display. Push buttons allow easy speed adjustment and operation of memory positions. Double-curved profile provides good ergonomics and operation of rotating knobs without moving hand.



Joystick controller also available.

does not rotate as it advances, instruments

Non-rotating Spindle Digital Micrometer Head

Build your own precision micro-positioning device



502102 SPECIFICATIONS

Total Travel Distance 25 mm Resolution 0.001 mm Accuracy \pm 0.003 mm Spindle Diameter Ø8 mm Mounting

Ø 12 mm x 10 mm Total Length 166 mm

Measurement Mode Absolute and incremental

Digital Readout mm or inch Analog Readout mm Data Output RS232 Environmental Protection IP54

Shipping Weight 0.51 kg (1.12 lb)

502102

Non-Rotating Spindle Micrometer Head

Universal Manipulator Stand



Universal Micromanipulator Stands enable scientists to mount their manual and motorized micromanipulators at variable angles and heights. A solid aluminum platform with a grooved tower allows the user to attach any micromanipulator of any size or shape to the post for infinite flexibility. Once mounted, the micromanipulator can be set at any height along the entire length (30 cm or 45 cm). The platform base comes configured with industry standard pre-bored holes (1/4-20 x 1" or M6 x 25 mm), allowing direct mounting to any type vibration-free table (for patch clamp recordings) or optical bench (for laser and optical measurements).

Using additional Rotation Clamps (one included with stand), two or more micromanipulators can be mounted on the stand simultaneously in a space-saving convenient manner. WPI's Universal Stand not only allows the user 360° flexibility in manipulations but also promotes independent angular transitions using a single feather-light tensioner/adjustment screw. This affordable stand is currently the *preferred* choice in micromanipulator stands head and shoulders above set-ups using multiple magnetic-based stands costing hundreds more when combined with antiquated heavy steel base plates.

UMS SPECIFICATIONS

DIMENSIONS

10.0 x 12.5 x 1.5 cm (LxWxH)

Base plate Stand

4.0 x 4.0 x 30 cm (LxWxH) (501622)

4.0 x 4.0 x 45 cm (LxWxH) (501623)

Mounting holes

English 1/4 20 x 1" (2 bolts supplied)

Matrix M6 x 25mm grid (2 bolts supplied)

SHIPPING WEIGHT

5016229 lbs (4 kg)

50162311 lbs (5 kg)

501622	Universal Micromanipulator Stand 30 cm (includes one clamp)
501623	Universal Micromanipular Stand 45 cm (includes one clamp)
501624	Additional Rotation Clamp
VFP	Vibration-Free Platform (24"x30")

Vibration-Free Workstation — see page 137

The M3301. Kite and SM325 (but not M325) can be used with these mounts.

(VFP) not included

High speed penetration and precise control

The piezo-ceramic element in WPI's MPM10 provides high penetration speed over an extremely short distance (0.5 to 10.0 µm). Because the range of travel of a dedicated piezomanipulator is much too limited for it to be useful independently, it must be mounted on a manipulator. The MPM-10 piezo translator combined with a DC3001 motorized micromanipulator (available separately) provides a single electronically controlled system.

When the piezo element is activated, the MPM-10 axis carrying the micropipette shoots forward at a rate which is set on the control panel, then immediately returns (at a slower speed) to its starting position. As soon as the piezo element begins its reverse travel, the motorized manipulator starts to travel forward. The complimentary opposition of these two travel sequences results in the micropipette tip remaining in its advanced absolute position.

The three axes of the DC3001 are controlled by six buttons. Pressing a button for less than 0.3 seconds activates one step, the size of which (0.5 to 10.0 micron) is set on the Step Size control. Pressing longer activates the continuous mode, at the rate set on the Motor Speed control. Pressing the button for the X axis forward direction activates the piezo mode. Advancement speed of the piezo element can be separately adjusted from 1 to 100 mm/sec.

Precise construction and special vibration stabilizers ensure the MPM-10's excellent



puncture characteristics. Lateral deviation from the ideal axis of puncture (measured at the tip of the electrode holder) is \pm 5% of the step size.

Includes controller, piezo translator, electrode holder, cables and mounting bracket. Shipping weight: 6 lb (2.7 kg).

SYS-MPM10	Piezo Controller for DC3001 Motorized Micromanipulator
PM5	Remote Controller for MPM-10
PM6	Replacement Electrode Holder for MPM10
14104	Record/Inject Electrode Holder for MPM10, MPM20
MPH8	Electrode Holder Adapter for MPM10 & MPM20
M-3	80° Tilting Base 6mm x 1mm screw (Shipping Weight: 2 lb)

Specify line voltage



Micromanipulator not included.

Piezo Translator

For use with M3301 and DC3001 micromanipulators

Especially recommended for use with the M3301 micromanipulator, the MPM20 is a very efficient tool for intracellular injection. High penetration speed and precise axial advance allows injection pipette to be brought to its target position with tremendous accuracy. Lateral escape of the cell is almost eliminated, and even tough membranes can be penetrated. Independently selectable reverse speed setting can be used for fast withdrawal, preventing adhesion of the injected cell to pipette tip. Mounts directly onto DC3001 and M3301 micromanipulators. Use with DC3001 requires MS314 controller (for the micromanipulator).

The combination of the MPM20, a micromanipulator, and the PV820 PicoPump (see page 180) constitutes an extremely efficient system for intracellular injection; cell penetration, injection and withdrawal are executed automatically with the press of a button.

Shipping weight: 10 lb (4.5 kg).

SYS-MPM20	Piezo Translator		
PM7	Replacement Electrode Holder for MPM20		
14106	Footswitch for MPM20		
Specify line voltage			

Programmable High Precision

Motorized Micromanipulator

suitable for patch clamp and IVF

WPI introduces a compact high precision motorized micromanipulator (SM325). It features low noise, high stability, a user-friendly software interface and economy that are major concerns in IVF and patch clamp research.

The SM325 is driven in all three axes through high resolution stepping motors, which can achieve 40,000 steps per revolution (25 nm/step) with completely vibration-free motion. In a normal lab environment, it can stay localized overnight without drifting. The 25mm long range of travel makes it unnecessary to have an additional manual coarse adjust-

Its compact construction makes mounting onto the stage plate of a microscope practical. The x-axis can be tilted by 90° that allows for a better positioning of the injection tool. An additional tilting fixture makes it possible to tilt the tool holder for fast and easy cleaning and exchange of the injection tool.

The MCL3 controller features a dynamic micro-step function that makes very quick positioning possible with maximum accuracy. Motor control is achieved with a linear output amplifier, which also drastically reduces electronic noise. Users can control the micromanipulators by joystick, keyboard, mouse or computer. The user-friendly software program can be enabled to remember up to 999 position coordinates from previous procedures and can robotically repeat this same positioning sequence.



SM325 SPECIFICATIONS

CONTROL METHOD Joystick, software, or both TRAVEL DISTANCE 25 mm each axis RESOLUTION 25 nm/step or 40,000 steps/rev MAXIMUM SPEED 4 mm/second POWER SUPPLY 120/240V, 50/60Hz **DIMENSIONS** SM325-M 5x7x5.5 in. (13x18x14 cm) (WxLxH) MCL3 9.8x9x3.7 in. (25x23x9.5 cm) (WxLxH) SHIPPING WEIGHT SM325-M 6 lb. (2.7 kg)

11 lb. (5 kg)

SM325	High Resolution 3-D Motorized Micromanipulator (SM325-M) & Controller (MCL3)		
SM325-M	High Resolution 3-D Motorized Micromanipulator		
MCL3	Controller with Joystick and software for SM325-M		
OPTIONAL ACCESSORIES			
M3301EH	Replacement Electrode Holder, straight, 14cm		
15873	Angled Electrode Holder, 13 cm long		
501622	Universal Micromanipulator Stand, 30 cm high		
501623	Universal Micromanipulator Stand, 45 cm high		
VFP	Vibration-Free Platform		

Vibration-Free Workstation — see page 143

MCL3

Programmable Ultra High Precision Motorized Micromanipulators

M3301EH Electrode Holder included

The HS6-3 is supplied with manual controls and stepper motor drives in all 3 axes. The extremely solid construction eliminates the vibrations and drifts. With the utmost precision and long travel distance in all three

directions, HS6-3 is the ideal tool for patch-clamp or electrophysiological applications. The tilting device is mounted on the base plate serves as coarse height adjustment as well and the tool holder can be swiveled in all directions.

The MCL3 controller features a dynamic micro-step function that makes very quick positioning possible with maximum accuracy and free of vibration. Motor control is achieved with a linear output amplifier, which also drastically reduces electronic noise. Users can control the micromanipulators by joystick,

HS6-3 SPECIFICATIONS

CONTROL METHOD: Joystick, software, or both

TRAVEL (X-Y-Z): 25 mm **RESOLUTION:** 10 nm/step MAXIMUM SPEED: 4.5 mm/sec. STABILITY: 1 nm/hour at 24°C 120/240 V, 50/60 Hz POWER SUPPLY:

DIMENSIONS:

HS6-3: 6.1x9.7x9.9 in (15.5x24.6x25 cm) (WxLxH)

MCI 3 9.8x9x3.7 in (25x23x9.5 cm) (WxLxH)

WFIGHT:

HS6-3: 13.2 lb. (6kg) MCL3: 7.7 lb. (3.5kg) keyboard, mouse or computer. The userfriendly software program can be enabled to remember up to 999 position coordinates from previous procedure and can robotically repeat this same positioning sequence.

HS6-3	High Resolution Motorized HS6 Micromanipulator and
Controller	
	includes HS6-3M and MCL3
HS6-3M	High Resolution Motorized HS6 Micromanipulator
MCL3	Controller with Joystick and software for HS6-3M
OPTIONAL	ACCESSORIES
M3301EH	Replacement Electrode Holder, straight, 14cm
15873	Angled Electrode Holder, 13 cm long
501622	Universal Micromanipulator Stand, 30 cm high
501623	Universal Micromanipulator Stand, 45 cm high
VFP	Vibration-Free Platform

Vibration-Free Workstation — see page 143

The DC3001 is the popular, high-precision manual/ motorized micromanipulator. With 0.5µm resolution, it can be used with an optional joystick or an

MS314 controller. (With the controller, it is sold as SYS-DC3314.) It can be mounted on the TBS tilt base but is too heavy for the M-3 tilt base. When performing intracellular microinjection, the DC3001 can be used with the MPM10 or MPM20 piezo translators.

Manual coarse controls use cross roller bearing slides. Vernier scales allow readings to 0.1 mm. All controls are closely grouped so adjustments can be made in any plane with minimum effort. The DC3001 features DC motor drives and fine control micrometers in all three axes. Left- or right-handed versions of the DC3001 are supplied with a standard 12 mm clamp. Standard accessories provided include one microelectrode holder and a securing bolt and wrench.

The sophisticated MS314 Controller allows control of all three axes. Movements may be continuous through the use of cross switches, or the controller can cause the DC3001 to step in defined increments. Steps as small as 0.5 µm are possible. A popular joystick controller, STM3, allows control of the X, Y and Z axes.

If the MPM10 piezo system is used, it replaces the MS314 controller. This configuration moves the motor forward as the piezo retracts, keeping the micropipette in the penetrated cell. When using the MPM20, the MS314 is necessary as the DC3001 motor controller because the MPM20 controls only the piezo element. The MS314 and MPM20 do not interact with each other.


Manipulator (right-handed) & MS314 Controller SYS-DC3314R SYS-DC3314L Manipulator (left-handed) & MS314 Controller Specify line voltage. System components also available separately: Motorized Manipulator, right-handed DC3001R DC3001L Motorized Manipulator, left-handed SYS-MS314 Controller for DC3001 Joystick Controller for DC3001 STM3 **OPTIONAL ACCESSORIES** Tilt Base with Screw Adjustment TBS Remote controller for MS314 and MPM-10 PM₅ 5-lb Weight for Tilting Base (shipping weight: 7 lb [3 kg]) 5464 M4C Microscope Stage Adapter M3301EH Replacement Electrode Holder (14 cm long) Angled Electrode Holder (13 cm long) 15873 501607 Cable for MS314 and DC3001




MCL3-HS6

	DC3001 SPE	CIFICATION	NS
	TRAVEL RANGE	RESOLUTION	
MANUAL:	X-axis 37 mm	0.1 mm	
	Y-axis 20 mm	0.1 mm	
	Z-axis 20 mm	0.1 mm	
	TRAVEL RANGE	RESOLUTION	MAXIMUM SPEED
MOTORIZED:	X-axis 10 mm	0.5 µm	0.2 mm/sec
	Y-axis 10 mm	0.5 µm	0.2 mm/sec
	Z-axis 10 mm	0.5 μm	0.2 mm/sec
SHIPPING WEI	GHT:		
	DC3001:	3 lbs (1.4 kg)	
	MS314:	1.8 lbs (0.9 kg)	
	STM3:	2.8 lbs (1.3 kg)	

Miniature Micropositioners






Mini Micropositioner

(MM1 and MM1-3)

Single stage measures only $5 \times 11 \times 26$ mm with 3 mm travel. Provides precise and smooth motion with no backlash, positive spring loaded carriage, straight within 1 micron and less than 1 micron maximum wobble. Features fine 80 TPI screw adjustment. 10 mm square mounting surface has a 3.9 mm tapped center hole for transmission and/or mounting. Available in single X (MM1), X-Y, and X-Y-Z (MM1-3) axis configurations.



Micropositioner

(MM3 and MM3-3)

Single stage measures only $7 \times 17 \times 44$ mm with 13 mm travel. Offers precise and smooth motion with no backlash, positive spring-loaded carriage, straight within 1.5 microns, and less than 1.5 microns maximum wobble. Features fine 80 TPI screw adjustment. 13 mm square mounting surface has a 7 mm tapped center hole for transmission and/or mounting. Available in single X (MM3), X-Y, X-Y-Z (MM3-3) axis configurations.



MINI-MICROPOSITIONER SPECIFICATIONS				
	MM1	MM1-3	ММЗ	MM3-3
AXIS	χ	X-Y-Z	Χ	X-Y-Z
STRAIGHT LINE ACCURACY	Within 1 micron over 3 mm travel	Within 1 micron over 3 mm travel	Within 1.5 micron over 13 mm travel	Within 1.5 micron over 13 mm travel
CLEAR APERTURE	3.9 mm tapped hole, 8-32 thread	3.9 mm tapped hole, 8-32 thread	7 mm tapped hole, 5/16-16 thread	7 mm tapped hole, 5/16-16 thread
LOAD CAPACITY	255 g Normal	255 g Normal	340 g Normal	340 g Normal
FINISH	Black Anodized	Black Anodized	Black Anodized	Black Anodized
WEIGHT	3 grams/axis	12 grams/axis	14 grams/axis	48 grams/axis
TYPE	Fine Screw	Fine Screw	Fine Screw	Fine Screw
TRAVEL	3 mm	3 mm	13 mm	13 mm

MM1	Mini Micropositioner, one axis, 3 mm travel
MM1-3	Mini Micropositioner, three axes, 3 mm travel
MM1-A	Mounting Adapter for MM1 and MM1-3
ММ1-С	Clamp for MM1 and MM1-3
ммз	Micropositioner, one axis, 13 mm travel
MM3-3	Micropositioner, three axes, 13 mm travel
ммз-а	Mounting Adapter for MM3 and MM3-3
ММ3-С	Clamp for MM3 and MM3-3
MM3-ALL	Complete 3-Axis Micropositioner & Magnetic Stand
MM1-ALL	Complete 3-Axis Mini Micropositioner & Magnetic Stand



Vernier scales allow readings to 0.1 mm. X-axis fine control allows readings to 10 µm.

Left- or right-handed versions of the KITE micromanipulator are supplied with a standard 12 mm clamp and electrode holder

M3391EH.	Kite Manual Manipulator (right-handed)
KITE-L	Kite Manual Manipulator (left-handed)
KITE-M3-R	Kite (right-handed) + Tilting Base Combo
KITE-M3-L	Kite (left-handed) + Tilting Base Combo

KITE SPECIFICATIONS

	TRAVEL RANGE	RESOLUTION	
X-axis Fine	10 mm	0.01 mm	
X-axis	35 mm	0.1 mm	
Y-axis	20 mm	0.1 mm	
Z-axis	20 mm	0.1 mm	
SHIPPING WEIGHT	3 lbs (1.4 kg)		

OPTIONAL ACCESSORIES

	Also see magnetic stands.			
M4C	Microscope Stage Adapter			
500476	Ball Joint, 4 cm long, for Ø 4 mm Holder			
500475	Ball Joint, 7 cm long, for Ø 8 mm Holder			
	Shipping weight: 7 lb (3 kg)			
5464	5-lb Weight for Tilting Base			
M-3	80° Tilting Base M6 x 1mm screw			
15873	Optional Angled Electrode Holder (13 cm long)			
M3301EH	Replacement Electrode Holder (14 cm long)			
0	7100100011110			

M3301 Manual Micromanipulator

The world's most widely used micromanipulator

Weighing just 550 grams and employing a slim space-saving design, this well-built micromanipulator outsells all others worldwide for high precision experiments where magnification is in the range of up to 250x. Its design allows units to stand tightly grouped — since all control knobs project to the rear. And because control knobs are clustered within an 8 cm area in a single vertical plane, resolution is quick — the hand works blindly while the eye monitors the microscopic image. Vernier scales allow readings to 0.1 mm; x-axis fine control allows readings to 10 microns.

The instrument employs rack-and-pinion drive, V-shaped guideways, and cross roller bearings, so all movement is sure and repeatable, without drift, sideplay, backlash, or sticking. Contact parts are milled of hardened steel for high performance and long life.



M3301R	Manual Manipulator, right-handed
M3301L	Manual Manipulator, left-handed
M3301-M3-R	Manual Manipulator (right handed) & Tilting Base
M3301-M3-L	Manual Manipulator (left handed) & Tilting Base
502105	Axis Adjustment Tool


	M3301 SPECIFICATION	ONS
	TRAVEL RANGE	RESOLUTION
Fine	10 mm	0.01 mm
	37 mm	0.1 mm
	20 mm	0.1 mm
	25 mm	0.1 mm

SHIPPING WEIGHT 3 lbs (1.4 kg)

X-axis I X-axis


Y-axis Z-axis

d new angle. Ball-joint holder attachment 7.9 mm diam. 500475 Ball Joint, 7 cm long, for O.D. 5-9 mm Electrode Holder (shown) Ball Joint, 4 cm long, for O.D. 2.8-4.5 mm Electrode Holder 500476





At last — a tilt base you can operate with one hand!




High Resolution Manual Micromanipulators

HS6 Micromanipulator

Engineered for stability, and built on a twelve-pound steel plate, this instrument is chosen worldwide for high resolution micro-recording. such as patch clamping, and other research requiring solid, drift-free performance. A superior tool in its own right, HS6 serves equally well as a base for other precision microdrives. HS6 can be bolted directly to a lab fixture or vibration-free platform. Resolution is extremely high — each graduation on its large micrometer barrels indicates just 5 micron movements. Rack and pinion drive, V-shaped guideways, and cross roller bearings give sure, repeatable movements without sideplay, slipping, or stick-

ing. All contact parts are milled of hardened steel. A flexible ball-joint assembly allows the electrode to be positioned at any angle relative to the x, y, or z axis. The entire manipulator tilts forward to 25 degrees allowing rapid coarse adjustment of the electrode, and allowing cell penetration along the axis of any of the micrometers. Simple, precise, and durable, the HS6 will provide years of dependable performance.



HS6 SPECIFICATIONS

	TRAVEL RANGE	RESOLUTION
X-axis	25 mm	5 μm
Y-axis	25 mm	5 μm
Z-axis	25 mm	5 µm

SHIPPING WEIGHT 25 lbs (11 kg)

DIMENSIONS 9.9 x 6.6 x 9.9 in. (H x W x L) SYS-HS6 Micromanipulator

M3301EH Replacement Electrode Holder (14 cm x Ø 7.2 mm) 15873 Optional Angled Electrode Holder (13 cm long)

VFP **Vibration-Free Platform**

Vibration-Free Workstation — see page 143

Joystick-Controlled Micromanipulator


Specially adapted for use with the Nanoliter Injector (page 186) for oocyte injection and similar applications, this joystickcontrolled micromanipulator allows an easy "steering" motion that translates normal hand movement into smooth submillimeter shifts. Viewed microscopically, movement of the tooltip corresponds naturally to hand movement, so accurate resolution is intuitive and quick. All fine adjustment can be controlled by the joystick. Pivoting forward, backward, or laterally gives precise x-y adjustment. For added convenience, a separate coarse control lever

is also provided for quick raising and lowering. A stop screw—which is set once resolution is achieved—eliminates refocusing and streamlines repetitive work by guiding the tip to its previous focussing plane. The stop screw also prevents the tool-tip from being broken during sudden lowering and eliminates downward drift—placement is stable enough for even extended use. Because the probe holder tilts a full 90°, the tool-tip pivots easily for precise positioning. Rack-and-pinion drive, V-shaped guideways, and cross roller bearings eliminate backlash, slipping, and sticking. All contact parts milled from hardened steel for precise performance and long life.

Joystick travel: 0.35 mm to 3.5-mm, depending on reduction gear ratio setting (adjustable between 1:15 and 1:150).

MMJR	Joystick Micromanipulator (Right-Handed)
MMJL	Joystick Micromanipulator (Left-Handed)
OPTIONAL	ACCESSORIES
M3301EH	Replacement Electrode Holder (14 cm 🗙 Ø 7.2 mm)
15873	Angled Electrode Holder (13 cm long)
M4C	Microscope Stage Adapter
500475	Ball Joint, 7 cm long, for Ø 8 mm Holder
500476	Ball Joint, 4 cm long, for Ø 4 mm Holder

Also see magnetic stands.



TRAVEL RANGE RESOLUTION 37 mm 0.1 mm 20 mm 0.1 mm 25 mm 0.1mm 0.35~3.5 mm

MMI SPECIFICATIONS

JOYSTICK (X,Y axis) SHIPPING WEIGHT 4 lbs (1.8 kg)

X-axis

Y-axis

Z-axis

Micrometer Slide Micromanipulator



M325 SPECIFICATIONS

TRAVEL RANGE RESOLUTION X-axis 25 mm 10 µm Y-axis 10 mm 10 µm 10 mm 10 μm Z-axis SHIPPING WEIGHT 4 lbs (1.8 kg)

The M325 three-axis manual micromanipulator is built of precision micrometer-actuated linear slides. Each slide is comprised of a large micrometer head and a spring-return linear slide. The micromanipulator has been carefully designed to minimize wear in the moving components to achieve a long operational life without the necessity for frequent maintenance or adjustment. The micrometer head is graduated in 10 micron steps which enable repeatable positioning to an accuracy of ± 2 microns.

A unique spring return mechanism is used to transmit movement of the micrometer spindle to the slide carriage — eliminating backlash, lost motion and reducing thread wear. Each linear slide utilizes ball bearings which enable the M325 to carry loads of up to 1 kg.

The toolholder can clamp onto tools with shaft diameters of 3.0 mm to 12.7 mm and allows rotation around two axes. This provides a wide range of options for incorporating the manipulator into your workstations. The M325 can also be configured very easily in left- or right-handed versions to allow several units to be positioned in close

A quick-release clamp allows easy mounting onto any rod from 10-mm to 12.7 mm diameter.

M325	3-Axis Fine Controlled Manual Micromanipulator		
OPTIONS	AND ACCESSORIES		
M3301EH	Replacement Electrode Holder (14 cm long)		
15873	Optional Angled Electrode Holder (13 cm long)		
500475	Ball Joint, 7 cm long, for Ø 8 mm Holder		
500476	Ball Joint, 4 cm long, for Ø 4 mm Holder		

Also see magnetic stands.

Dual Tool-Holder Micromanipulator

A small and compact micromanipulator for manual manipulation in all three axes (x, y and z), the MD4 is equipped with a mounting bracket for a second tool or electrode holder which can be positioned in the x and y axes independent of the manipulator and may also be tilted and swiveled by two fine-adjust screws. Scales allow readings of coarse adjustment with an accuracy of 100 µm. Additional x-axis fine control is achieved with a micrometer screw with a resolution of 10 µm. Supplied with one M3301EH electrode holder and a 12 mm clamp for mounting on M10 Stand or other 12 mm supports. May also be mounted on optional M-3 Tilting Base.

Travel, standard electrode: x-axis, 37 mm (fine, 10 mm); y-axis, 20 mm; z-axis, 25 mm. Additional electrode: x-axis, 7 mm; y-axis, 10 mm.

MD4R	Double-Holder Micromanipulator (right)
MD4L	Double-Holder Micromanipulator (left)
MD4-M3-R	Double-Holder Micromanipulator (right) + Tilting Base
MD4-M3-L	Double-Holder Micromanipulator (left) + Tilting Base
OPTIONAL A	ACCESSORIES
M3301EH	Replacement Electrode Holder (14 cm long)
15873	Optional Angled Electrode Holder (13 cm long)
M2	Additional Ø 12 mm Clamp
M-3	80° Tilting Base 6mm x 1mm screw
M4C	Microscope Stage Adapter
M5	Additional Ø 10 mm Clamp
M6	Additional Ø ¹ /2-in. Clamp
5464	5-lb Weight for Tilting Base*
500475	Ball Joint, 7 cm long, for Ø 8 mm Holder
500476	Ball Joint, 4 cm long, for Ø 4 mm Holder
	161.

*Shipping weight: 8 lb (3.6 kg)

Also see magnetic stands. Also see Universal Manipulator Stands.



MD4	SPECIFICATION	ONS
	TRAVEL RANGE	RESOLUTION
X-axis Fine	10 mm	10 μm
X-axis	37 mm	100 μm
Y-axis	20 mm	100 μm
Z-axis	25 mm	100 µm
SHIPPING WEIGHT	3 lbs (1.4 kg)	

Magnetic **Holding Devices**



Mechanical clamp tightens three rotatable joints simultaneously with one locking knob. Arm adjusts without distortion. Base exerts a magnetic force of 100 kilos for greatest stability. Equipped with fine adjustment for precise operations.

Magnetic Base:

50 (w) x 60 (l) x 55 (h) mm (2.2 x 2.4 x 2.2 in.)

Vertical Holding Power:

100 kgf (220 lb force)

Arms:

11. 119 mm (4.7 in.) L2: 106 mm (4.2 in.) 25 mm (0.98 in.) L3: Ø 12 mm (0.472 in.)

Clamp Hole:

none

Weight:

1.8 kg (4 lb)

М9

Magnetic Stand

Similar to M1 but with a 12 mm diameter sub pole (fits 12 mm clamp supplied with M3301, DC3001, MD4 and MMJ manipulators).

Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb force)

Main Pole:

diameter: 14 mm (0.55 in.) length: 178 mm (7 in.)

Sub Pole:

diameter: 12 mm (0.47 in.) 165 mm (6.5 in.)

Clamp Hole:

Adjustable from 4.5 mm to 6.5 mm

Weight:

1.8 kg (4 lb)

M10 Magnetic Stand

M10L

Same as M10 but equipped with a taller (14inch) vertical main pole.

Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb force)

Main Pole:

diameter: 14 mm (0.55 in.) 356 mm (14 in.)

Sub Pole:

diameter: 12 mm (0.47 in.) length: 165 mm (6.5 in.)

Clamp Hole:

Adjustable from 4.5 mm to 6.5 mm

Weight:

1.8 kg (4 lb)

M₁₀L Magnetic Stand

Bends freely for maximum flexibility. The connecting arm twists and bends like a snake. Lock the arm in position with a flick of the controlling lever.

Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb force)

Main Pole:

diameter: 16 mm (0.63 in.) 315 mm (12.4 in.) length:

Sub Pole:

none

Clamp Hole:

Adjustable from 6 mm to 8 mm

Weight:

1.4 kg (3 lb)

Magnetic Stand M11

Powerful Ball Joint Rare Earth Magnet

Use in constructing your own holding device for small parts/equipment

- Small but very powerful: holds 2 kilograms (~5 pounds)!
- Steel ball rotates freely 360° on a 180° axis
- M3 mounting screw on ball for attachment to equipment
- Magnet base threaded (M3) for mounting onto a base or equipment




This novel magnetic ball joint has phenomenal holding power for up to 2 kg of attached weight while permitting the ball a full 360° rotation on a 180° axis. You can freely orient your equipment to an infinite number of positions within this rotation. This is made possible by the combination of a steel ball (10 mm diameter) and a powerful rare earth magnet contained in the magnet cylinder (\$\phi\$ 10 x 20mm).

Convenient M3 attachment sites are provided on


both the ball (male) and the magnet base (female).

For use with micromanipulators for the positioning and holding of optical instruments 28 mm including various lighting sources and lasers, pipettes and any small parts that would benefit from the flexibility offered by this new magnetic ball joint.



500871 Magnetic Ball Joint

The base of each stand exerts a powerful magnetic force that holds it solidly on ferrous metal surfaces — even vertically or upside-down



MB2

Mechanical clamping type tightens three joints simultaneously just by on-tough operation. Arm is freely adjustable without distortion. Equipped with fine adjuster and medium size magnet for stabilizing the base. Suitable for performing precision operation.

Magnetic Base:

30 (w) x 30 (l) x 30 (h) mm (1.2 x 1.2 x 1.2 in.)

Vertical Holding Power:

17 kgf (37 lb force)

Arm:

L1: 46 mm (1.8 in.) L2: 46 mm (1.8 in.)

L3: 39 mm (1.5 in.)

Clamp Hole:

Adjustable from 5 to 6.5 mm

Weight:

0.38 kg (0.83 lb)

MB2 Compact Magnetic Stand

A ball joint at the base of the main post allows 360° rotation, offering considerable versatility. The second arm adopts angles up to 75°.

Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb force)

Main Pole:

diameter: 12 mm (0.47 in.) length: 194 mm (7.6 in.)

Sub Pole:

diameter: 10 mm (0.39 in.) 165 mm (6.5 in.) length:

Clamp Hole:

Adjustable from 4.5 mm to 6.5 mm

Weight:

1.8 kg (4 lb)

M8 Magnetic Stand

A precision base providing stable support for such devices as electrodes and manipulators. Adjustable second arm adopts a variety of angles.

Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb force)

Main Pole:

diameter: 12 mm (0.47 in.) 176 mm (6.9 in.) length:

Sub Pole:

diameter: 10 mm (0.39 in.) length: 165 mm (6.5 in.)

Clamp Hole:

diameter: 4.5 mm and 6.5 mm

Weight:

1.8 kg (4 lb)

М1 Magnetic Stand

Same base and support arm as M1, but equipped with a longer (14-inch) vertical post.

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb force)

Main Pole:

diameter: 12 mm (0.47 in.) 356 mm (14 in.) length:

Sub Pole:

diameter: 10 mm (0.39 in.) length: 165 mm (6.5 in.)

Clamp Hole:

diameter: 4.5 mm and 6.5 mm

Weight:

1.8 kg (4 lb)

M₁L Magnetic Stand

Mounting your Micromanipulator

Three of the stands above — M1, M1L, and M8 — have 10 mm diameter mounting rods. The standard mount on several WPI manipulators (DC3001, KITE, M3301, MMJ,

and MD4) accommodates a 12 mm rod. In order to use one of these three stands, you will need to replace the manipulator's standard 12 mm mounting clamp with the




optional M5 clamp. **M5** Ø 10 mm Clamp

World Precision Instruments

www.wpiinc.com

Magnetic **Holding Devices**





Round Base

An ideal accessory for optical tables and vibration-free platform. Reduces experimental set-up time by allowing free positioning and instant clamp down of optical components. Switchable ON/OFF magnetic circuit permits fine adjustment and precise positioning.

- Easy ON/OFF operation using lever
- Thin and powerful magnetic force
- Generous array of tap holes

Holding Power:

20 kgf (44 lb force)

Dimension:

75 (OD) x 20 (h) mm 2.9 (OD) x 0.8 (h) in.

Mounting Hole:

4-M4 x 0.7, depth 6mm * M8 x 1, depth 6mm Span 35mm

Weight:

0.7 kg (1.5 lb)

501651 Magnetic Base, 75mm diameter 503568 Magnetic Base, 50mm diameter

* Posts with M4-threads not available from WPI. See posts with M8 threads on page 142.

Square Base

An ideal accessory for optical tables and vibration-free platform. Reduces experimental set-up time by allowing free positioning and instant clamp down of optical components. Switchable ON/OFF magnetic circuit permits fine adjustment and precise positioning.

- Easy ON/OFF operation using lever
- Thin and powerful magnetic force
- Generous array of tap holes

Holding Power:

20 kgf (44 lb force)

Dimension:

65 (w) x 65 (l) x 20 (h) mm 2.6 (w) x 2.6 (l) x 0.8 (h) in.

Mounting Hole:

8-M4 x 0.7, depth 6mm * M8 x 1, depth 6mm Span 25mm

Weight:

0.6 kg (1.3 lb)

501653	Magnetic Base, 65x65mm
503569	Magnetic Base, 45x45mm
503570	Magnetic Base, 90x90mm
503571	Magnetic Base, 120x120mm

* Posts with M4-threads not available from WPI. See posts with M8 threads on page 142.

MOBITY

MOBITY™ is a new magnetic clamping system. With its ease of use, only one hand is needed to operate the attractive power. The MOBITY™ has a strong 88lbf pull, yet weighs only 1.5 lbs. MOBITY™ meets various applications with 4 tapped holes on the top surface. Requires (1) 9V alkaline battery (included).

Holding Power:

40 kgf (88 lb force)

Dimension:

55 (w) x 73 (l) x 50 (h) mm 2.2 (w) x 2.9 (l) x 2.0 (h) in.

Mounting Hole:

3-M4, depth 20mm * M8, depth 15mm

Weight:

0.7 kg (1.5 lb)

501652 MOBITY Magnetic Clamping System

* Posts with M4-threads not available from WPI. See posts with M8 threads on page 142.

A small holder ideal for use where space is limited. Main post unscrews from base which may then be used alone as a switchable magnetic holder.

Magnetic Base:

30 (w) x 35 (l) x 35 (h) mm 1.2 (w) x 1.4 (l) x 1.4 (h) in.

Vertical Holding Power:

20 kgf (44 lb force)

Main Pole:

Diameter: 7mm (0.28 in.) Length: 52mm (2 in.)

Clamp Hole:

Diameter: 6mm

Weight:

0.36 kg (0.8 lb)

M7 Compact Magnetic Stand

8½"x12" Steel Base Plate #5052

12"x24" Steel Base Plate #5479

ACCESSORIES

5052 Steel base plate, 8½ x 12 in. (10 lb) Steel base plate, 12 x 24 x %-in. (32 lb)

give stability to your experimental setup. Beveled edges make them easy to handle; rubber feet hold them off the benchtop, making them easier to grasp when moving; and the special black coating provides a durable protective finish.

BASE PLATES: A magnetic stand requires a steel mounting surface. WPI's

steel base plates have plenty of mass to

Bench Top Vibration Isolation Platforms

Simple Set-Up and Adjustment



Performance

- · Horizontal frequencies are weight dependent.
- · Horizontal frequency of 1.5 Hz is achieved at or near the upper limits of the payload range.
- At the lower limits of the payload range the horizontal frequency is approximately 2.5 Hz.
- · Vertical frequency is tunable to 0.5 Hz throughout the payload range.

Pictured: MK-BM-8100, 50-105 lb. payload weight range (23 - 48 kg)

These bench top platforms offers 10-100 times better performance than a full size air table in a package only 4.6 inches tall, and without air or electricity! These vibration isolation platforms are extremely easy to use and offer extreme performance -1.5 Hz horizontal natural frequency and 0.5 Hz vertical natural frequency. There are only two adjustments.

This is the thinnest, most portable, and most user-friendly isolator ever offered that is capable of delivering this level of performance.

Weight: Approximately 40 lb (16 kg)

Dimensions: 18" W x 20" D x 4.6" H (457 x 508 x 117 mm)

Model	Payload Range	
MK-BM-825	Bench top Vibration Platform, 10 - 30 lb (4.5 - 14 kg)	Call for Price
MK-BM-850	Bench top Vibration Platform, 25 - 55 lb (11 - 25 kg)	Call for Price
MK-BM-8100	Bench top Vibration Platform, 50 - 105 lb (23 - 48 kg)	Call for Price
MK-BM-8125	Bench top Vibration Platform, 90 - 130 lb (40 - 59 kg)	Call for Price
MK-BM-8150	Bench top Vibration Platform, 125 - 155 lb (57 - 70 kg)	Call for Price
MK-BM-8175*	Bench top Vibration Platform, 150 - 180 lb (68 - 81.5 kg)	Call for Price
MK-BM-8200*	Bench top Vibration Platform, 175 - 205 lb (79.5 - 93 kg)	Call for Price
MK-BM-8225*	Bench top Vibration Platform, 200 - 230 lb (90.5 - 104 kg)	Call for Price
MK-BM-8250*	Bench top Vibration Platform, 225 - 255 lb (102 - 115.5 kg)	Call for Price

* Weight: Approximately 47 lb. (21 kg) / (same dimensions)

Vibration-Free Tables

All buildings vibrate — activities of people, machinery, heating and ventilation systems, and nearby truck or rail traffic cause all types of vibrations. These vibrations, though acceptable to occupants, cannot be tolerated by equipment used in patch clamping, cell injection, analytical balances, and optical microscopes. The short-term effects of such vibrations include inconsistent and unreliable performance. The long-term effects are excessive wear, maintenance, and fatigue failures. In order to protect sensitive instruments and equipment from faulty operation or failure, vibration must be significantly reduced. This can be efficiently accomplished by using Vibration-Free Platforms and Vibration-Free Workstations.



Vibration-Free Workstation

- Vertical and horizontal vibration isolation
- High performance Active-Air Suspension
- Automatic leveling
- VibraDamped Steel
- Class 100 Cleanroom compatible
- Leveling feet

Additional tabletop sizes and finishes are available, as well as optional accessories such as side rails and casters.

Call for more information and prices for the configuration you require.

> For more information, see www.wpiinc.com/vfw

Precision Stereo ZOOM Microscope

- Modular, two parallel beam design
- Planachromatic objectives, no optical distortion
- High-contrast imaging, ideal for observing transparent, low-contrast obiects
- Large zoom ratio: 8:1
- Large zoom range: 0.62x-5.0x
- Long working distance
- Step and continuous zoom
- 5-year warranty

The fourth generation of WPI's precision stereo zoom microscopes uses modular, two parallel beam path design and high quality optical system. The advanced optical design with planachromatic objectives provides sharp and distortion-free contrast image throughout the entire zoom range and comes with an impressive 5-year warranty.

The **PZMIV** is available in a binocular or trinocular version. In addition, an extensive list of optional accessories is available that makes the PZMIV suitable for integrated optical and biological research.

The microscope comes with a track stand, standard 10x eyepieces (wide-field, distortionfree and high eye point), 1x planachromatic distortion-free objective. See the Table on next page for all optical options.





PZMIV	ecision Stereo Zoom Binocular Microscope (Model IV), on Track Stand		
PZMIV-BS	PZMIV Microscope on Boom Stand (tubular)		
PZMTIV	Precision Stereo Zoom Trinocular Microscope (Model IV), on Track Stand		
PZMTIV-BS	PZMTIV Microscope on Boom Stand (tubular)		
PZMTIV-DIG50	PZMTIV Microscope System, including PZMTIV, USBCAM50 USB Computer		
	Camera, 0.5× CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated		
	Optical Fiber Light Guide		
PZMTIV-BS-DIG50	PZMTIV Microscope System, including PZMTIV, USBCAM50 USB Computer		
	Camera, 0.5× CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated		
	Optical Fiber Light Guide, Boom Stand		

PZMTIV-BS-LWD-DIG50 PZMTIV Microscope System, including PZMTIV, USBCAM50 USB Computer Camera, 0.5× CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical Fiber Light Guide, Boom Stand, 0.5× Objective, 20× Eyepieces for Long Working Distance Viewing

See web site for complete configurations.

PZMIV 40 cm





PZMɪv & PZMTɪv Eyepieces and Objectives									
	10x Eyepiece		16x Eyepiece		20x Eyepiece 25x I		yepiece		
Objective	Mag	Field (mm) (Video Field)	Mag	Field (mm) (Video Field)	Mag	Field (mm) (Video Field)	Mag	Field (mm) (Video Field)	Working Distance
0.32x	1.9x - 16x	106 -13.1 <i>(49.8 - 6.1)</i>	3.2x 25.6x	70.6 - 8.8 (49.8 - 6.1)	3.9x - 32x	55.4 - 6.9 (49.8 - 6.1)	5x - 40x	45.4 - 5.6 (49.8 - 6.1)	296 mm
0.5x	3.1x - 25x	67.7 - 8.4 (31.8 - 3.95)	5x - 40x	45.2 - 5.6 <i>(31.8 - 3.95)</i>	6.2x - 50x	35.5 - 4.4 (31.8 - 3.95)	7.8x - 62.5x	29 - 3.6 <i>(31.8 - 3.95)</i>	189 mm
0.63x	3.9x - 31.5x	53.8 - 6.7 (25.3 - 3.15)	6.2x - 50.4x	35.8 - 4.4 (25.3 - 3.15)	7.8x -63x	28.2 - 3.5 (25.3 - 3.15)	9.8x - 78.8x	23 -2.9 <i>(25.3 - 3.15)</i>	149 mm
1.0x (inc)	6.2x - 50x	33.9 - 4.2 <i>(15.9 - 1.97)</i>	9.9x - 90x	22.6 - 2.8 <i>(15.9 - 1.97)</i>	12.4x - 100x	17.7 - 2.2 (15.9 - 1.97)	15.5x - 125x	14.5 - 1.8 <i>(15.9 - 1.97)</i>	80 mm

The Video Field is based on a 1/2-inch CCD (8 mm diagonal) and a 0.5x camera adapter.

502000	PZMIV Binocular Body With 10X Eyepieces, 1× Objective, Eye guards
502001	PZMTIV Trinocular Body With 10X Eyepieces, 1× Objective, Eye guards
502004	Boom Stand (Heavy) W/O Focus Mount (requires 502009 Focus Mount for PZMIV)
502005	Ball Bearing Boom Stand (Heavy) W/O Focus Mount (requires 502009 Focus Mount for PZMIV)
502006	Boom Clamp Stand (Heavy) (requires 502009 Focus Mount for PZMIV)
504123	Extension for Heavy Clamp Stand
502007	Articulated Arm and Table Clamp w/o Focus Mount (requires 502009 Focus Mount for PZMIV)
502009	Universal Focus Mount for 76 mm PZMIV (Required for BS, AAC, BBS, and BCS) (5/8" pin)
504596	76mm Halogen-Halogen Dual Illuminated Track Stand
504597	Replacement Lamp for 504596
502010	10× Wide Field Eyepiece for PZMIV (pair)
502011	16× Wide Field Eyepiece for PZMIV (pair)
502012	20× Wide Field Eyepiece for PZMIV (pair)
502013	25× Wide Field Eyepiece for PZMIV (pair)
500264	10× Eyepiece with Reticle (matches 10× eyepiece #502010)
500266	20× Eyepiece with Reticle (matches 20× eyepiece #502012)
502015	Ring Light Adapter for PZMIV (For R-8-8-WPI01 Ring Light Guide)
502016	0.32×, Planachromatic Objective (Distortion-free) (278 mm WD)
502017	0.50×, Planachromatic Objective (Distortion-free) (174 mm WD)
502018	0.63×, Planachromatic Objective (Distortion-free) (138 mm WD)
502019	1.0×, Planachromatic Objective (Distortion-free) (73mm WD)
500261	0.35× CCD Camera Coupler, C-Mount (Use with USBCAM33)
500262	0.5× CCD Camera Coupler, C-Mount (Use with COLCAM, USBCAM50)
500028	1× CCD Camera Coupler, C-Mount (Use with COLCAM)
502163	Wall Mount Plate for Articulated Arm System
Z-LITE	Z-Lite Fiber Optic Illuminator
500186	Bifurcated Light Guide with Lenses
	Ring Light Guide
Z-LITE-186	Z-Lite Fiber Optic Illuminator with Bifurcated Light Guide and Lenses

LED Ringlight



PZMIV SPECIFIC	ATIONS
EYEPIECES	WFH 10×
AUXILIARY LENSES	1×
ZOOM RANGE	0.62× - 5×
TOTAL MAGNIFICATION	6.2× - 50×
ZOOM RATIO	8:1
FIELD OF VIEW	Ø33.9- Ø4.2 mm
WORKING DISTANCE	80 mm
BINOCULAR TUBE	Inclined 45°
INTERPUPILARY DISTANCE	50 – 75 mm
DIOPTER ADJUSTMENT	± 5 Diopter
MICROSCOPE BODY	Rotatable 360°
OPTIONAL ACCESSORIES	
Eyepieces	16×, 20×, 25×
Auxiliary lenses	0.32×, 0.5×, 0.63×
Total Magnification Field of view	1.9x - 125x Ø106 - Ø1.8 mm
rieiu oi view	ااااا ۱.۵ الا - ۱۱۱۵

Working Distance

SHIPPING WEIGHT

World Precision Instruments

80 – 296 mm


23 lb.

Precision Stereo ZOOM Microscope

Quality and **Precision** to Improve Your Vision

WPI's third-generation stereomicroscope, PZMIII, is an ideal tool for tissue dissection, cell injection, specimen manipulation, electrode inspection, and many other applications that require a magnified, stereo viewing and ample working distance. It offers the leading brand's quality and performance at an affordable price. Advanced optics provide the sharpest image that can only be found among the best of this class. It is superior to many stereomicroscopes costing almost twice as much. Zooming is achieved by a spring-loaded





knob that is smooth and effortless. The compact size and light weight make it more stable and easily manipulated on the boom stand. A specially designed photo/video module is used in the trinocular version of the microscope (PZMTIII) for photo, video,

or digital imaging. In addition, an extensive list of optional accessories is available that can make the PZMIII suitable for almost any bio-research applications requiring a stereomicroscope. See next page for options.

PZMIII-BS	PZMIII Microscope on Boom Stand		
PZMIII-AAC	PZMIII Microscope on Articulated Arm with Table Clamp		
PZMTIII	Precision Stereo Zoom Trinocular Microscope (Model III)		
PZMTIII-DIG50	PZMTIII Microscope System, including PZMTIII, USBCAM50 Computer		
	Camera, 0.5X CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated		
	Optical Fiber Light Guide		
PZMTIII-BS-DIG50 PZMTIII Microscope System, including PZMTIII, USBCAM50 Computer			
	Camera, 0.5X CCD Camera Coupler, Z-LITE Optical Illuminator, Bifurcated		
	Optical Fiber Light Guide, Boom Stand		
PZMTIII-BS-LWD-DIG50 PZMTIII Microscope System, including PZMTIII, USBCAM50			
	Computer Camera, 0.5X CCD Camera Coupler, Z-LITE Optical Illuminator,		
	Bifurcated Optical Fiber Light Guide, 0.5× Objective, 20× Eyepieces for Long		
	Working Distance Viewing		

Precision Stereo Zoom Microscope (Model III), on Post Stand

PZMTIII Microscope on Articulated Arm with Table Clamp All PZMIII and PZMTIII microscopes come with 10x eyepieces and built-in 1x auxiliary lens.

PZMTIII-BS

PZMTIII-AAC

PZMIII

PZMTIII Microscope on Boom Stand



i zimii zyepieces ana objectives									
10× Eyepiece		15× Eyepiece		20× Eyepiece		25× Eyepiece			
Objective	Mag.	Field (mm) (Video Field)*	Mag.	Field (mm) (Video Field)*	Mag.	Field (mm) (Video Field)*	Mag.	Field (mm) (Video Field)*	Working Distance (mm)
0.3×	2x - 13.5x	114 - 17 (53.6 - 8)	3× - 20.3×	84 -13 (53.6 - 8)	4× -27×	69 - 10.3 (53.6 - 8)	5× - 33.8×	44.8 - 6.7 (53.6 - 8)	287 mm
0.5×	3.4× - 22.5×	69 - 10 (32.4 - 4.7)	5× - 33.8×	51 - 7 (32.4 - 4.7)	6.7× - 45×	42 - 6.2 (32.4 - 4.7)	8.4× - 56.3×	26.9 - 4.0 (32.4 - 4.7)	177 mm
0.75×	5x - 33.8x	45 - 7 (21.1 - 3.3)	7.5× - 50.6×	34 - 5 (21.1 - 3.3)	10× - 67.5×	28 - 4.2 (21.1 - 3.3)	12.6× - 84.4×	17.9 - 2.7 (21.1 - 3.3)	117 mm
1.0×	6.7× - 45×	34 - 5 (16 - 4.7)	10× - 67.5×	25 - 3.7 (16 - 4.7)	13.4× - 90×	21 - 3.1 (16 - 4.7)	16.8×-112.5×	13.4 - 2.0 <i>(16 - 4.7)</i>	100 mm
1.5×	10x - 67.5x	23 - 3.4 (10.8 - 1.6)	15× - 101.3×	17 - 2.5 (10.8 - 1.6)	20.1×-135×	14 - 2.1 (10.8 - 1.6)	25.1×-168.8×	9.0 - 1.3 (10.8 - 1.6)	47 mm
2.0×	13.4× - 90×	12 - 2.5 <i>(5.6 - 1.17)</i>	20.1× - 135×	13 - 1.8 <i>(5.6 - 1.17)</i>	26.8×-180×	10 - 1.5 (5.6 - 1.17)	33.5× - 225×	6.7 - 1.0 (5.6 - 1.17)	26 mm


^{*} The video field of view is based on a 1/2-inch (8 mm diagonal) CCD camera and a 0.5× camera adapter.

OPTI	ONAL	ACCESSORIES
VI II	CITAL	ACCESSOIVIES

503102

OPTIONA	AL ACCESSORIES
501352	PZMIII Binocular Body, pair of 10x eyepieces and eyeguards
	13338 Ring Light Adapter NOT included
501353	Fan Post Stand with 76mm Focus Mount
502009	Universal Focus Mount, 76 mm ID for PZMIII Body
502004	Boom Stand (Heavy) without Focus Mount
502005	Ball Bearing Boom Stand (Heavy) without Focus Mount
504123	Post Extension for Heavy Boom Stand
502006	Boom Clamp Stand (Heavy) (requires 502009 Focus Mount for PZMIV)
502007	Articulated Arm with Table Clamp, without Focus Mount (50 cm clamp)
502163	Wall-Mount Plate, 6" × 6" (or 15.24 cm × 15.24 cm)
501369	Wide Field 10× Eyepieces (pair)
501370	Wide Field 15× Eyepieces (pair)
501371	Wide Field 20× Eyepieces (pair)
501372	Wide Field 25× Eyepieces (pair)
504128	10× Eyepiece with Reticle (matches 10× eyepiece on PZMIII)
504129	20× Eyepiece with Reticle (matches 20× eyepiece 501371)
501373	0.3× Long Working Distance Objective Lens
501375	0.5× Long Working Distance Objective Lens
501376	0.75× Long Working Distance Objective Lens
501377	1.5× Long Working Distance Objective Lens
501378	2.0× Long Working Distance Objective Lens
501379	PZMTIII Trinocular Body, pair of 10× eyepieces and eye guards
	13338 Ring Light Adapter NOT included
501381	0.5× C-Mount CCD Camera Coupler
13338	Ring Light Adapter for PZMIII Series
	(included with all microscope configurations on previous page)
503051	Manual Stage for PZMIII

76mm Rectangular Base Post Stand for PZMIII



#503051 Manual Stage — Mounts in the circular opening in the PZMIII base. XY travel distance: 75 x 56 mm. Glass size: 116 x 96 mm. Active diameter: 37.6 mm. Dimensions: 180 x 155 x 27 mm. Fits 503102 base only.

PZMIII SPE	CIFICATIONS
EYEPIECES	WFH 10×
ZOOM RANGE	0.67× - 4.5×
TOTAL MAGNIFICATION	6.77× - 45×
FIELD OF VIEW	Ø 34 MM - Ø 5 MM
WORKING DISTANCE	100 mm
BINOCULAR TUBE	Inclined 45°
INTERPUPILLARY DISTANCE	Adjustable 47-70 mm
DIOPTER ADJUSTMENT	±5 Diopter (both eyepiece tubes)
MICROSCOPE BODY	Rotatable 360°
AUXILIARY LENSES	
Total magnification	2× - 225×
Biggest Field of View	Ø 110 mm
Working Distance	26-287 mm
SHIPPING WEIGHT	23 lh

Inverted Trinocular Microscope



Accessories Included

Green and neutral filters, dust cover.

Options

Multiple photo options

Weight

24 lb (10.9 kg)

INV-101	VV-101 Trinocular Inverted Microscope				
503510	30 mm 10x Eyepiece with 100/10 reticle				
503520	Replacement lamp				
503512	Deluxe Optical Cleaning Kit				

All necessary tools are included for routine adjustments, alignments, and assembly: 8 oz. Air Duster, 280 sheets lint-free lens paper, 1 oz. "no-residue" lens cleaning fluid, 100 cotton-tipped applicators, 9x9 microfiber lens cloth, soft-bristled dusting brush, micro-glide gear lubricant, allen wrenches, double-sided friction collar wrench, precision screwdriver set, nylon carry case

3000-LED **Series**

Superb Optics, Durable and High Performance **Microscopes**



3000-LED SERIES SPECIFICATIONS

OPTICAL SYSTEM Infinity Optical System, f=180mm, Anti-Mold

VIEWING HEAD Siedentopf type, inclined 30°; interpupillary distance adjustment

48-75mm

EYEPIECES HWF Plan 10× eyepiece, 20mm field of view with built-in diopter

adjustment; a pointer is standard in one eyepiece

NOSEPIECE Rear facing quadruple

Infinity Plan achromat 4× (N.A. 0.10), 10× (N.A. 0.25), 40×R* (N.A. 0.65), **OBJECTIVES**

100×R oil* (N.A. 1.25) are standard

STAGE Stage size 140mm × 132mm with X-Y movement range of

76mm × 50mm

PHASE CONTRAST (OPTION) Phase sliders for 10×/BF/40× CONDENSER N.A. 1.25 Abbe condenser

ILLUMINATION 3 watt LED with variable intensity control

ACCESSORIES Dust cover, immersion oil and instruction manual, universal power

supply 110v-240v

STAND Cast alloy aluminum; coaxial coarse and fine focusing controls.

DIMENSIONS $15.25 \times 7.75 \times 15.4$ in. (387 × 196 × 391mm)

WEIGHT 16 lb (7.26 kg)

> The ACCU-SCOPE® 3000-LED Microscope Series delivers outstanding optical performance, value and resolution to meet the exacting standards of life science professionals and students. With a newly designed infinity plan optical system, a best-in-class 20mm field of view and a super-bright 3-watt LED illuminator, the 3000-LED Series provides high contrast images with outstanding resolution, precision design and enhanced illumination.

504221	Binocular, Infinity Plan achromat 4×, 10×, 40×R and 100×R oil objectives
504443	Trinocular, Infinity Plan achromat 4×, 10×, 40×R and 100×R oil objectives
504445	Binocular, Plan phase contrast 10× and 40×R objectives
504444	Trinocular, Plan phase contrast 10× and 40×R objectives
504416	0.50× c-mount adapter for 1/2" sensors, adjustable focus
504417	0.35× c-mount adapter for 1/3" sensors, adjustable focus
500828	Stage Micrometer, 1mm scale, 200 div. at 10 µm
504606	Stage Micrometer, 50-0.5mm scale and 10um scale

W30 Professional-Grade Microscope



The W30 professionalgrade microscope is a best-seller in universities, medical schools, and reseach laboratories. Equipped for performance, its features include titanium-finished DIN or Semi-Plan optics and a 30-year anti-fungal coating. The W30 is the choice for superior performance at a great price.

W30S SPECIFICATIONS

HEAD Binocular (Seidentopf)

True Trinocular

Inclined 30°, rotates 360° Dual diopter adjustment, Interpupillary

range 55-75 mm

10×/18 wide field eyepieces

NOSEPIECE Quadruple forward-facing nosepiece

OBJECTIVES DIN Plan, anti-fungal

4×, 10×, 40×, 100×R (oil) Parfocal, parcentric, color-coded

STAGE Mechanical stage (140 mm × 140 mm)

Coaxial drive controls

XY Movement: 73 mm × 43 mm

FOCUS Coarse adjustment: range of 30 mm

Fine adjustment: graduation of 2 µm

Tension control knob

ILLUMINATION Moveable Abbe condenser, NA 1.25,

Iris diaphragm

Variable LED light source (3W bulb) 110V/220V switchable electronics

ACCESSORIES INCLUDED

attachment cover.

distance

Replacement 0.5 amp fuses, mirror

(for field use), blue and green filters, dust

DIMENSIONS AND WEIGHT 15" (38 cm) × 9" (23 cm) × 7" (17.8 cm)

14 lb (6.4 kg)

W30S-LED	Binocular Microscope
W30ST-LED	Trinocular Microscope
503513	21 mm 10X Eyepiece with 100/10 reticle

Precision SurgioScope

Ideal for small animal surgery

WPI's improved precision SurgioScope (now with five magnification steps) is a portable high quality surgical microscope offering outstanding image quality and value. Incorporating an agile extension arm and excellent working distance objectives, the SurgioScope provides convenient movement and maneuverability necessary for accurate positioning. These important features, together with a high quality optical system, provide sharp image contrast and enhanced large field of vision. The SurgioScope comes fully equipped with a foot-controlled motorized focusing system, normally only found in more expensive surgical microscopes. A unique dual lamp housing enables safe and rapid changing of the lamp during an operation, without the need to power down. The optional video port on the "T" version permits operational procedures to be monitored or recorded simultaneously using a video recorder and a COLCAM video camera or digital stills with USBCAM50.

PSMB5N	Binocular SurgioScope, F200 objective (Specify post height.)				
PSMT5N	Trinncular SurginScone heam solitter std video adante				

F200 objective (Specify post height.)

Specify 89 cm or 103 cm post Specify line voltage

OPTIONS AND ACCESSORIES

501636	¹ / ₂ " CS-mount Adaptor (requires			
	Beam Splitter 501637)			
501637	Beam splitter			
504284	F100 Objective			
504285	F250 Objective			
504286	F300 Objective			
504287	F350 Objective			
500162	Replacement lamp, 12V, 100 W			

RANGE OF MOTION

=
Maximum Stretch Radius of Arm870 mm
Vertical Movement Range700-1100 mm
ILLUMINATION
Spot = 42 mm
Dual lamp housing with quick-change spare and internal
coaxial fiber optic cable.
HALOGEN-TUNGSTEN LAMP12V, 100W, with cold reflection
OPTIONAL CAMERACOLCAM, USBCAM50 (½" CCD)
USBCAM33 (1/3" CCD)
POWER110V, 50-60 Hz, or 220V, 50-60 Hz
CHIDDING WEIGHT 04 lb (42 kg)



SURGIOSCOPE Specifications

TOTAL MAGNIFICATION (F200)	3.2× — 25×
ADJUSTABLE DIOPTER	± 6 Diopter
ADJUSTABLE INTERPUPILLARY DISTANCE	min. 50 mm — max. 70 mm
EYEPIECE	12.5x

Objective	Working Distance	Magnification step	Visual Field of view (mm)	Camera field 1/2" CCD (mm)	Camera field 1/3" CCD (mm)
F100 #504284	90mm	6.4, 10, 16, 26, 40x	25, 15.5, 10, 6, 4	25, 15.5, 10, 6, 4.5	17.5, 11.5, 7, 4.6, 2.8
F200 (included)	190mm	3.2, 5, 8, 13, 20x	50, 31, 20, 12, 8	50, 31, 20, 12, 8	35, 23, 14, 9, 5.5
F250 #504285	240mm	2.6, 4, 6.4, 10.4, 16x	65, 40, 25, 16, 10	63, 40, 25, 16, 10	45, 28, 18, 11, 7
F300 #504286	290mm	2.1, 3.3, 5.3, 8.7, 13x	75, 46.5, 30, 18, 12	75, 46.5, 30, 18, 12	52.5, 34.5, 21, 13.5, 8.3
F350 #504287	340mm	1.8, 2.9, 4.6, 7.4, 11x	91, 57, 36, 22, 14	88, 55, 35, 21, 13	60, 38, 24, 15, 9.5

For additional objectives and specifications, please go to WPI website—www.wpiinc.com

FINE FOCUS ADJUSTMENT RANGE......30 mm

WORKING HEIGHT (Arm Movement Range Above Floor)

89 cm Post...... Focus on specimens 34.5" (88 cm) to 51" (130 cm) above floor * 103 cm Post Focus on specimens 40.5" (103 cm) to 57" (146 cm) above floor *

* Subtract Working Distance for height above specimen, 103 cm post recommended for F350

High intensity source for fiber optic illumination

The **Z-LITE** Fiber Optic Illuminator provides reliable, uninterrupted high-intensity light for microscopes. **Z-LITE** allows a continuous range of subdued or concentrated lighting controlled by a rotary dimmer on the front panel. **Z-LITE** may be used with a ring light and single or bifurcated flexible fiber bundles, enabling the light beam to be placed exactly where needed. Forced air cooling prolongs lamp life. Lamp color temperature is 3350°K. An interlock switch automatically cuts off power when front panel is opened to replace bulb.



Z-LITE-186	Z-Lite & Bifurcated Light Guide (115 V, 60 Hz, beige case)
Z-LITE-Z186	Z-Lite & Bifurcated Light Guide (230 V, 50 Hz, black case)
Z-LITE	Z-Lite Fiber Optic Illuminator (115 V, 60 Hz, beige case)
Z-LITE-Z	Z-Lite Fiber Optic Illuminator (230 V, 50 Hz, black case)


LED Ringlight for PZMIII and PZMIV stereo microscopes • "White" light illumination — 72 LED bulbs Maximum opening 61 mm Ring light is divided into four areas and each area is turned on and off separately Brightness adjustable ESD safe Power supply AC 90-264V, 50/60 Hz, US plug only 504134 LED Ringlight

LIGHT GUIDES AND ACCESSORIES

500186	Bifurcated Light Guide (with lenses)
R-8-8-WPI01	Ring Light Guide for PZM and PZMIII Series*
13338	Ring Light Adapter (48 mm Ø) for PZM, PZMII, PZMIII
502015	Ring Light Adapter for PZMIV
5475	Adapter for SMA-terminated Fiber Optic Cables
EJA	Replacement Halogen Lamp, 150W, 3350°K, 40-hour
EKE	Replacement Lamp, 150W, 3250°K, 200-hour

*Ring Light Guide requires adapter #13338 for use with PZM, PZMII and PZMIII, included with each PZMIII and PZMIV microscope.





Environmental Control Systems for Live Cell Microscopy

This Live Cell Microscope Incubator was extensively tested in laboratories. When compared with other systems, it offers dramatic advantages. For example, other incubators for live cell microscopy rely on passive, random diffusion of heated air from a single source to maintain the desired temperature setpoint. With no hot air return vent, the heated air

escapes from the system through cracks at the microscope/incubator junction in an uncontrolled, random fashion. These systems offer no temperature uniformity, suffer from focus drift and often experience electrical and vibrational interference from the heater. You will also notice dramatic temperature drifts when the imaging environment is disturbed.

Features

- Unique, diffusion grid, combined with air input and return vents provide an air flow pattern for consistent, even heating, with no hot or cold spots in the chamber
- External heater that can be placed far enough from the system to eliminate electrical and vibrational interference from the heater
- High degree of temperature precision and stability
- Minimal focal drift after equilibrium is achieved Accuracy ±0.1°C at the sample itself, and 0.2°C across the microscope stage (allowing for uniform heating of multiwell dishes)
- Airflow pattern and temperature uniformity eliminate dramatic changes in environmental temperature when the incubator door opens
- Ergonomic design for ease of use– The focus and x/y stage controls are outside of the incubator itself. Large doors allow easy access to the specimen and small ones for cords. tubing, etc.
- Precision, shielded temperature probe
- Simple, one person setup of the system

Live Cell Microscopy System Components

This unique, acrylic Live Cell Imaging chamber, combined with an Environmental Control Unit (ECU) and an AirTherm controller, ensured precision control of your incubator environment.

The ECU comes in four varieties so you have all the control you require. With the ECU-H5, you can control air flow and heating. In addition to air flow and temperature control, the ECU-HC and ECU-HCP let you control the CO₂ level. One has an internal sensor, and the other has a remote sensor that can be positioned inside the microscope chamber. Finally, the ECU-HOC adds control of the O₂ level, which is accomplished by displacing the oxygen with nitrogen. The first three ECU units are capable of controlling a simple, external heater, like the AirTherm Satelite (AirTherm-SAT) or a microscope lens warmer. The AirTherm-SMT can monitor and control both temperature and humidity inside the microscope chamber.

Consistent Air Flow

Air flow affects the temperature uniformity of incubators. The red arrows on Fig. 1 and Fig. 3 indicate air flow. The Live Cell Microscope Incubator uses a diffuser grid and proper venting to insure consistent air flow. Traditional incubators with poor air flow suffer with hot and cold spots in the incubator, as seen in thermal images (Fig. 2 and Fig. 4). Warmer temperatures are indicated by red and cooler temperatures by blue.

Ordering an Environmental Chamber

Acrylic enclosures are essentially custom-built. When ordering a system, you will need to provide the following information:

Microscope Perfect Focus Right Port Tirf Tilting Head Emission Cells

Stage Camera Analyzer White Light Tirf Filter Wheels: Dual Lamphouse 35mm/60mm Coverslips

Stage-Up Left Port Fluor Attachment Binocular D Head Excitation Transmitted Light Shutter Wells



Fig. 1-Single air input and no venting causes random air flow in a traditional incubator.



Fig. 2-Hot and cold spots result from inconsistent flow.



Fig. 3-A diffusion grid with air input and exhaust vents yields consistent air flow.



Fig. 4–Consistent air flow means uniform heating.


IN VIVO FULL CHAMBER SYSTEM WITH ECU CONTROLLER & NEW AIRTHERM-SAT/SMT

	All Systems Include: Proprietary Humidification Module, Stage Adaptor, Stage Dish with Optical Grade Glass.
IV-100SMT	Environmental Chamber, Heat controller. Requires pre-mixed 5% CO ₂ gas supply
IV-200ECU	Environmental Chamber, Heat & CO ₂ controller. Requires 100% CO ₂ & ambient air supplies
IV-200-OX	Microscope Environmental Chamber, Heat Controller, CO ₂ and O ₂ Controller.
IV-300ECU	Environmental Chamber, Heat & CO_2 controller. Requires 100% CO_2 & ambient air supplies. Features a CO_2 probe to detect concentration at the sample.
IV-CUST	Custom Design includes but not limited to Tilting Heads, Black Chamber, Confocals, Camera on Right, HMX Lamphouse, and Manipulators. Fee added for design not included in standard system

Stagetop **Environmental** Control

Control temperature and CO2 in a microscope stagetop environment

- Four programmable digital control loops
 - Independent incubator base temperature PID control with ±0.1°C precision
 - Independent incubator lid temperature PID control with ±0.1°C precision
 - CO₂ digital PID control with ±0.1% precision Airflow digital PID control from 0-900 SCCM
- USB-based remote control and data logging
- Electronic flow meter
- Programmable alarm for out of tolerance condition on all four channels
- Compact and lightweight



For short term or long term studies of living cell cultures under a microscope or for time lapse video research, a microscope stagetop incubator is essential.

Perfect for Live Cell Imaging, STEV (the stagetop environmental control system) is a compact environmental case that houses your culture wells and fits on a microcope stage inside the live cell microscope incubator.

This system offers precision control of both temperature and carbon dioxide, as well as remote control and data logging via a USB connection. The system is flexible and easy to configure for

a variety of experimental conditions.



The system includes the **Environmental Control Unit** electronics which use four programmable loops to control the temperatures of the case and the lid, CO₂ within the environmental case and airflow within the incubator.

STAGETOP ENVIRONMENTAL CONTROL SYSTEMS

IV-ECU-H5	Chamber, Controller, Heat, Digital flow control for bottle gas
IV-ECU-HC	Chamber, Controller with CO2 and Heat (using CO2 internal sensor)
IV-ECU-HCP	Chamber, Controller with CO2 and Heat (CO2 probe sensor)
IV-ECU-HOC	Chamber, Controller with CO2, O2 and Heat (using CO2 internal sensor)

Precision Heat Control

Smart, electrically quiet air heater for live cell imaging systems and custom incubators



The **AirTherm™ SMT** is a new generation of heat control system from WPI designed to be used in Live Cell Imaging applications with microscopes fitted with a full microscope environmental chamber enclosure. The standard **AirTherm™ SMT** controls temperature and, as an option,

The **AirTherm™ SMT** uses a PID control algorithm to maintain tightly controlled loops of heat and humidity environment control.

With AirTherm™ SMT, the temperature of the sample and microscope optics can be controlled to within 0.2°C. During operation, air is drawn out of the chamber through a flexible hose, heated by the **AirTherm™ SMT** heater and re-circulated to the chamber by the return hose.

- Precision heat controller for use in live cell imaging and custom incubators
- Control heat and humidity (optional) with a single controller
- Electrically and acoustically quiet
- Quick, precise response to thermal change

The system is typically used in a closed loop configuration.

The AirTherm™ SMT system includes:

98727


98728

- Two coil-reinforced heater hose pieces and hose clamps.
- Temperature sensor for remote placement in environmental

Humidity probe for monitoring chamber humidity available as an option.

riairiiaity pro	DC TOT THOTHE	orning charmoer marmately available as air option.
AIRTHERM-	SMT-1W	Airtherm™ SMT Heater, 110V
AIRTHERM-	SMT-2W	Airtherm™ SMT Heater, 230V
AIRTHERM-	SAT-1W	AirTherm Satellite Heater, 110V
AIRTHERM-	SAT-2W	AirTherm Satellite Heater, 230V
OPTIONAL ACCESSORIES		
15590	Replacemer	nt Hoses, 2.5" diam., 4.5 ft
300276	Replacemer	nt Platinum Temperature Probe

A typical AirTherm™ SMT installation places the heated air inflow at the bottom and the cold air return at the top of the microscope chamber.



Air-Therm SMT

AIRTHERM SPECIFICATIONS

0.2°C

AIR FLOW RATE	20-50 CFM (0.55-1.4 m³/minute)

CONTROL TEMPERATURE RANGE Ambient to 60°C TEMPERATURE RESOLUTION 0.1°C

HUMIDITY CONTROL RANGE Ambient-90%

HUMIDITY RESOLUTION 0.1% **HUMIDITY ACCURACY** 5% **HUMIDITY REPEATABILITY** 0.5%

TEMPERATURE ACCURACY

ANALOG OUTPUT 0.5°C resolution; FOR CHART RE-CORDER 0-10V represents 0-100°C

HEATING VOLUME Less than 50 CF (1400L), re-circulating

TEMPERATURE SENSOR TYPE Platinum RTD 1000 W

HUMIDIFIER TYPE Ultrasonic **HUMIDIFIER TANK CAPACITY** 0.5 gallons **HUMIDIFIER DAILY OUTPUT** 2 gallons

> For 120 VAC, 8A 250 V 5x20 mm metric For 230 VAC, 4A 250 V 5x20 mm metric 450 W. 95-135 V or 220-240 V. 50/60 Hz

DIMENSIONS 6½ x 8 x 7½ in. (15.5 x 21 x 19 cm)

Replacement Temperature Probe

Humidity Probe

FUSE

POWER

USBCAM152

1/2" Interline SXGA color progressive CCD 4.65 (H) x 4.65 (V) µm, 8mm diagonal



- Ultra-Compact USB cameras with color CCD
- Hardware & Software Trigger

USBCAM133 with a 1/3-inch color CCD, 6mm diagonal. USBCAM152 with a 1/2-inch color CCD, 8mm diagonal. USBCAM202 with a 1/1.8-inch color CCD, 7mm diagonal.

All three cameras are cased models with a rear mini-USB connector; hardware and software triggering image capture, digital zoom and a feature-rich user based menu setup and control. Sentech USB cameras include a SDK, DirectX, Twain and Linux driver, as well as the Sentech Viewing Software.

USBCAM133	Digital Microscope Camera
USBCAM152	Digital Microscope Camera
USBCAM202	Digital Microscope Camera
504570	Replacement USB Cable, 3m (10 ft)
504574	CONN 6-PIN (external camera triggering control)

USBCAM202		
	USBCAM133	
Image Sensor	1/3" Interline SXVGA color progressive CCD	
Cell Size	3.75 (H) x 3.75 (V) µm, 6mm diagonal	
Scanning System	Progressive	
Resolution	1280 (H) x 960 (V)	
Min. Scene Illumination	11 Lux at F1.2	
Speed	22.4 Frames per Second	
Electronic Shutter	Auto / Manual (software selectable)	
Gain	Auto / Manual (software selectable)	
Gamma	Manual (software selectable)	
White Balance	Auto / Manual / One shot (software selectable)	
Input / Output	USB 2.0 High Speed	
Power	+5 Vdc through USB connector, < 300 mA	
Dimensions	28 (W) x 28 (H) x 37 (D) mm (excluding connecto	
Lens Mount	CS mount	
Weight	Approximately 45g	

Progressive 1360 (H) x 1024 (V) 18 Lux at F1.2

19.26 Frames per Second Auto / Manual (software selectable) Auto / Manual (software selectable) Manual (software selectable) Auto / Manual / One shot (software selectable) USB 2.0 High Speed +5 Vdc through USB connector, < 420 mA 28 (W) x 28 (H) x 42 (D) mm (excluding connector) Approximately 45g USB: mini-B USB connector IO signal: 6pin connector (HR10A-USB: mini-B USB connector IO signal: 6pin connector (HR10A-7R-7R-6PB or equivalent) 6PB or equivalent) RoHS Compliant RoHS Compliant

USBCAM202 1/1.8" Interline UXGA color progressive CCD 4.40 (H) x 4.40 (V) µm, 7mm diagonal Progressive 1600 (H) x 1200 (V) 7.7 Lux at F1.2 15.3 Frames per Second Auto / Manual (software selectable) Auto / Manual (software selectable) Manual (software selectable) Auto / Manual / One shot (software selectable) USB 2.0 High Speed +5 Vdc through USB connector, < 450 mA 28 (W) x 28 (H) x 42 (D) mm (excluding connector) Approximately 45g USB: mini-B USB connector IO signal: 6pin connector (HR10A-7R-6PB or equivalent) RoHS Compliant

USBCAM33 / USBCAM50

Record images directly to your computer. These digital microscope cameras offer flexibility, with a range of configurations for image capture, a choice of mount option (C or CS) and file output alternatives. Since both cameras connect via the USB port, installing the image capture

Interface Connector

software is simple. Either camera can be used on WPI's stereo microscopes PZMTIV, PZMTIII, compound microscopes

W30ST and GPL-T and also the PSMT5 Microscope. Choose from the one third-inch CCD with 1024×768 resolution and 30 frames per second (USBCAM33) or one half-inch CCD

USBCAM33 Digital Microscope Camera, 1/3-in. CCD **USBCAM50** Digital Microscope Camera, 1/2-in. CCD Cable, USB Extension (male-female) 503536

with 1280×960 resolution and 15 frames per second (USBCAM50).

These cameras include **IC Imaging Control** 3.0 software that features:

- Real-time video preview
- Text and graphics can be drawn on a live video stream
- Scroll and Zoom
- Acquisition of single frames
- Capture pause, for intermittent image capture
- Timestamps

	USBCAM33	USBCAM50
IMAGE SENSOR	1/3" Sony CCD, progressive scan	1/2" Sony CCD, progressive scan
MAX RESOLUTION	1024 x 768	1280 x 960
SIZE	4.65 μm x 4.65 μm, 6mm diam.	4.65 μm x 4.65 μm, 8mm diam.
SPEED (PC DEPENDENT)	30fps, 15fps, 7.5fps or 3.75fps	15fps, 7.5fps or 3.75fps
SENSITIVITY	0.5 lux @ 1/15 s	0.5 lux @ 1/7.5 s
EXPOSURE, SHUTTER CONTROL, WHITE BALANCE	Automatic/Manual	Automatic/Manual
INTERFACE	USB 2.0 cable	USB 2.0 cable
SYSTEM REQUIREMENT	Windows Vista (32 & 64 bit) or Windows 7 (32 & 64 bit)	Windows Vista (32 & 64 bit) or Windows 7 (32 & 64 bit)
SOFTWARE	IC Imaging Control Software	IC Imaging Control Software
LENS MOUNT	C/CS-Mount	C/CS-Mount
CAMERA BODY	50.6 x 50.6 x 50 mm	50.6 x 50.6 x 50 mm
WEIGHT	265 g (9.5 oz)	265 g (9.5 oz)

Color Video Cameras for Microscopy

IMAGER

ROHS

16:9 Aspect Ratio, 1:1, no scaling



DVI signal output via HDMI cable

> Improved design of COLCAM-HD with better low light sensitivity

 ϵ

COLCAM-HD and **COLCAM-HD1080P**

The **COLCAM-HD** is a CCD-based camera that outputs a true HD 720P at 60fps in the 16x9 format. The COLCAM-HD1080P is a CMOS-based camera that outputs a true HD 1080P or 720P image at 60fps in the 16x9 format. Both models feature the capability to program individual DSP profiles accessed via remote hand held controller.

COLCAM-HD	Color Video Camera
-----------	--------------------

Includes 3m DV cable, power supply, C/CS mount.

COLCAM-HD1080P Color Video Camera

Includes 3m DV cable, power supply, C/CS mount.

3 meter HDMI cable 504136

Hand held control pad for COLCAM-HD* 504138

504137 C/CS Adaptor

*Controller for COLCAM-HD1080P available soon.

Low cost alternative to existing HD and 3 CCD cameras, with outstanding image quality

COLCAM SPECIFICATIONS

COLCAM-HD	COLCAM-HD1080P
1/3" Interline SXGA CCD:	1/2.8" 230 Mega pixel CMOS

ICX445AQA (SONY: IMX136), Rolling Shutter HD ACTIVE PICTURE 1280 (H) x 720 (V) 1920 (H) x 1080 (V) **ELEMENT**

CHIP SIZE 6mm diagonal 6.4 mm diagonal MINIMUM SCENE

RGB, 1280 (H) x 720 (V); 60 Hz /

4 Lux at F1.2 ILLUMINATION SYNC SYSTEM Internal Internal

VIDEO OUTPUT DVI 1.0 compliant; 720P DVI 1.0 conformity 1080P

> 1080P30,1080P29.97, 1080P25, 720P60, 59.94 Hz / 50 Hz 720P59.94, 720P50

RGB, 1080P60,1080P59.94, 1080P50,

40 (W) x 40 (H) x 45.8 (D) mm

GAIN AGC * AGC * SHUTTER SPEED Auto * Auto *

GAIN AGC or Fixed gain * AGC or Fixed gain * GAMMA presets or manual gamma * presets or manual gamma * WHITE BALANCE Auto / Manual / Push-to-set * Auto / Manual / Push-to-set *

POWFR 12V power jack; 5.5x2.1mm **DIMENSIONS** 40 (W) x 40 (H) x 45.8 (D) mm

OPTICAL FILTER IR cut filter included IR cut filter included LENS MOUNT C/CS mount C mount VIDEO OUTPUT HDMI connector HDMI connector

WEIGHT Approximately 120g TBD

> RoHS compliant RoHS compliant * Selectable via the UART communication

C-Mount Eyepiece Adapters

For 1/4-inch (4mm diagonal), 1/3-inch (6 mm diagonal) and 1/2-inch (8 mm diagonal) video cameras and eyepiece camera conversion, this lens and its accessories make it possible to connect a typical video or C-mount camera to almost any microscope on the market. The lens fits right into the ocular socket of standard 23.2 mm microscopes and the 30 mm adapter allows for use on the typical stereo zoom microscope. If you already have a trinocular microscope you can add the included C-Adapter to the top of an existing 1X C-mount (no lens) adapter. For a list of 1X C-mounts for popular microscopes, visit WPI's website and search for "C-mounts".

503097	Adapter, 0.45X for 1/3-in. and 1/2-in. video cameras,	
	30 mm Stereo adapter, 1X C-mount adapter	
503098 *	Adapter, 0.28X for 1/4-in. video cameras, 30 mm	

Stereo adapter, 1X C-mount adapter

^{*} For adapter 503098 to mount properly, the 1/4-inch CCD camera must have a minimum of 13 mm clearance, measured from the front of the C-mount ring to the CCD surface.



SLR Digital Camera-to-Microscope **Eyepiece Adapter**

This adapter connects T-mount SLR digital cameras to almost any microscope on the market. The adapter is built to 23.2 mm ocular tubes that are found on most high magnification (upright, inverted, standard) microscopes. The 30 mm adapter allows mounting on most Stereo zoom microscopes that use 30 mm oculars. If you already have a trincocular microscope you can add this adapter to the top of an existing 1X C-mount adapter.

The 2X magnification of this microscope adapter yields an approximate 65% field of view from the visual field as measured on a Canon 10D Digital camera. (CCD Sensor size = 22.7 x 15.1 mm). 35 mm film reference size is 24 x 36 mm.

Please contact your camera dealer for a suitable T-mount to bayonet adapter for your camera.

503099 Adapter, 2.0X for SLR digital Cameras (includes 30 mm Stereo adapter, 1X C-mount adapter)



FluoroDish"

Cover-glass bottom for observing and growing cells in imaging related research

- Optical quality glass bottom for better imaging quality (RI=1.525)
- Low sample volume for expensive chemicals
- Lowest access angle for micropipette
- Low toxicity adhesive for embryo research

WPI's FluoroDish™ tissue culture dishes are now available in a larger range of sizes and coatings. These dishes provide exceptional imaging quality for many applications requiring the use of inverted microscopes such as high resolution image analysis, microinjection and electrophysical recording of fluorescent-tagged cells. Taking advantage of WPI's extensive experience with low toxicity adhesives, FluoroDish uses a specially formulated adhesive that is optically clear, durable and with extremely low toxicity. Tests by an independent laboratory have shown that the 96-hour surviving rate of embryos is 100% when kept in FluoroDish, substantially better than some other brands. The bottom glass has superior UV transmission (30% transmission at 300 nm, compared to less than 7% for the most popular German glass). Stringent quality control ensures that glass thickness stays within the 0.17 ±0.01 mm range.

Conventional plastic dishes and chambers limit the utility of the inverted scope for many applications because the thick plastic bottom requires a long working distance objective available only in lower magnifications. Each WPI dish has a flat (0.17 mm thick) optical quality glass bottom, allowing the use of a much shorter objective working distance, larger numerical aperture (NA), and a higher magnification (up to 100x). The larger NA and

higher magnification provide superior quality imaging for both classical and fluorescence microscopy. Higher effective NA yields brighter images for fluorescence and higher resolution

in Image Analysis. The glass bottom permits the use of immersion objectives with medium such as water, glycerin or oil for the highest magnification possible. To optimize heatexchange, WPI's glass-bottom dish is designed to be flush (flat) with the microscope stage or heating unit, therefore eliminating the air gap that exists with modified plastic dishes in which a glass cover slip has been inserted.


Three different sizes of FluoroDish are offered, one type of 50 mm diameter dish and two types of 35 mm diameter dishes. An inner well is created within the dish by the glass bottom and the tissue culture grade polystyrene which forms the sides of the dish. All WPI dishes have the advantages of low toxicity and good UV transmission bottom glass. They are individually packed and gamma sterilized.

The 35 mm dish has outside dimensions similar to that of a Corning 35 mm dish and has ø23.5 mm glass window (FD35) or Ø10 mm glass window (FD3510). Most heaters and perfusion adapters designed for the Corning 35 mm dish will also fit this dish. The 23.5 mm glass window dish is available uncoated or poly-Dlysine-coated. Certain types of cell lines (e.g., PC3 and HEK) adhere well to the uncoated glass bottom dish. The poly-D-lysine coating has been reported to improve the adhesion of neuron cells. The users can also apply to the




The 10 mm glass window dish (FD3510) has low sidewall for easy microelectrode access and low solution volume. The low microelectrode access angle is the lowest among all of 35 mm dishes on the market (very close that of a 50 mm dish). The dish needs only about 115 µL to cover the bottom well, an important feature when using expensive drugs and chemicals.

The 50 mm dish (FD50) has a large growth area (35 mm well diameter), a low access angle for microelectrodes, and grips for easy handling.



Standard Fluorodish geometry. See the table below.



FD3510

Part Number	ID (mm	OD (mm)	Glass Ø (mm)	height (inside)	Height (outside)	Access Angle
FD35	33	35.5	23.5	7.8	9	29°
FD3510	10	35.5	10	1.5	4.65	17°
FD5040	47.5	49.82	35	7.25	7.4	17°

Choose from poly-D-lysine coated or uncoated

FD35-100	FluoroDish Sterile Culture Dish, clear wall, 35 mm, 23 mm well, box of 100
FD35PDL-100	FluoroDish Sterile Culture Dish, Poly-D-Lysine Coated, clear wall, 35 mm, 23 mm well, box of 100
FD3510-100	FluoroDish Sterile Culture Dish, clear wall, 35 mm, 10 mm well, low sidewall, box of 100
FD5040-100	FluoroDish Sterile Culture Dish, clear wall, 50 mm, 35 mm well, box of 100



Cover Slips

These cover slips made of German glass can be used for growing and culturing cells that normally have poor adhesion to plastic surfaces. They are small enough to be placed in the micro plate or other cell culture devices. The 5 mm size will fit inside the 96-well culture plate and leave enough room to pick it up from the bottom of the well with forceps. The 8 mm size fits inside the 24-well plates.

Order Number	Diam.	Thickness	Quantity
502040	5 mm	#1.5 (0.16 - 0.19 mm)	100
502041	8 mm	#1.5 (0.16 - 0.19 mm)	100
503508	25 mm	#1.5 (0.16 - 0.19 mm)	100





Anatomy Kit



Entomology Kit

Student Entomology Kit

· 2 Teaser Needles Straight

· 2 Teaser Needles Angled

• 1 Adson Micro Forceps, 12cm

· 2 Entomology Forceps, 10cm

· 1 Stainless Steel Scale, 6"

· 2 Micro Dissection Forceps, 12cm Str

· 2 Micro Dissection Forceps, 12cm Cvd

• 1 Tweezers with Magnifying Mirror

• 1 Iris Scissors, 11.5cm • 1 Scalpel Handle #3

504166

Kit Contains:

1 Dropper

· 10 Blades #10

· 1 Leather Pouch

504164 Student Anatomy kit

Kit Contains:

- · 1 Plastic Ruler, 6"
- · 1 Mall Probe
- 1 Teaser Needle Angled
- · 1 Teaser Needle Straight
- 1 Dressing Forceps, 16cm
- · 1 Mosquito Forceps, 12.5cm
- 1 Scalpel Handle #3
- · 4 Blades for Scalpel Handle #3
- 1 Scalpel Handle #4
- · 4 Blades for Scalpel Handle #4
- 1 Dissecting Scissors, 14cm Sharp/Blunt
- · 1 Hooks with Chain
- · 1 Leather Pouch

Botanical Kit

504165 Student Botanical Kit

Kit Contains:

- 1 Scissors, 14cm
 - 1 Dropper
 - 1 Scalpel Handle #3
 - 4 Scalpel Blades
 - 1 Teaser Needle Angled
 - 1 Teaser Needle Straight
 - · 1 Mariam Tweezers Angled
 - 1 Tweezers Straight
 - 1 Leather Pouch

503507	503506
	503507


503505	Plain Glass Microscope Slides, Box of 144
503506	Frosted Glass Microscope Slides, Box of 144
503507	Red frosted Glass Microscope Slides, Box of 144

Slides

These clean glass microscope slides are 25 x 75 mm, 1.0~1.2 mm thick with 90° grounded edges. and are available plain, frosted, and red ended. The frosted end slides feature a fine 20 mm frosted area on both sides of one end for easy marking. The red frosted slides feature a 20 mm colored end useful for identifying hazardous materials.

Pumps & Fluid Handling

	Fluid Range	Channels	Special features	Pa
Perista	ltic Pumps			
MINISTAR	0.006- 37 mL/min	1	Compact design, remote control	15
Peripro-2HS	0.8 - 300 mL/min	2	Calibrated output, replaceable tubing cartridges	1!
Peripro-4HS	0.8 - 300 mL/min	4	Calibrated output, replaceable tubing cartridges	1
Peripro-4LS	0.01-80 mL/min	4	Calibrated output, replaceable tubing cartridges	1
Peripro-8LS	0.01-80 mL/min	8	Calibrated output, replaceable tubing cartridges	1
<u>-</u>				
			Solution Exchange	
MPS2	Up to 250 μL/min	8	Programmable control; low dead volume	1
	tory Syringe Pum	ps		
AL-1000	0.000073-1699 mL/hr	1	Push/pull	1
AL-2000	0.000073-1699 mL/hr	2	Push/pull (2 networked pumps)	1
SPLG100	1.26 pL/min to 88.32 mL/min	1	Infuse only	1
SPLG101	1.26 pL/min to 25.99 mL/min	2	Infuse only	1
SPLG110	1.26 pL/min to 88.28 mL/min	1	Infuse/Withdraw	1
SPLG210	0.5 pL/min to 220.97 mL/min	2	Infuse/withdraw	1
SPLG212	0.5 pL/min to 220.97 mL/min	2	Infuse/withdraw programmable	1
SPLG270	0.5 pL/min to 220.97 mL/min	2	Push-pull	1
SPLG272	0.5 pL/min to 220.97 mL/min	2	Push-pull programmable	1
SPLG200	0.5 pL/min to 220.97 mL/min	2	Infuse only	1
SP100i	0.0001-519 mL/hr	1	Basic single channel	1
SP101i	0.001 μL/hr – 35mL/min	2	Micro dialysis application	1
SP120p	0.1 μL/hr – 127 mL/hr	1+1	Push pull, single cycle	1
SP200i	0.001 µL/hr - 145 mL/min	2	RS232 TTL/Footswitch	1
SP210c	0.001 μL/hr - 86 mL/min	2+2	RS232 push pull, continuous	1
SP210iw	0.001 μL/hr - 145 mL/min	2	RS232 Infuse/Withdraw	1
SP220i	0.001 μL/hr - 21 mL/min	10	RS232 Infuse Only	1
SP230iw	0.001 μL/hr - 21 mL/min	10	RS232 Infuse/Withdraw	1
SP250i	0.001 μL/hr - 21 mL/min	4	RS232 Infuse Only	1
SP260p	0.001 μL/hr - 86 mL/min	2+2	RS232 push pull, single cycle	1
licro Syrin	ige Pump / Stereo	taxic lı	njection	
UMP3	0.03nL/min - 10 μL/sec	1	Ultra micro infuse/withdraw RS232	1
MMP	Manual 100 μL-1mL syringe	1	Manual	1
DMP	Manual 100 μL-1mL syringe	1	Digital readout micrometer	1
licroinject	tion			
PV820	Injected volumes from picoliters	1	Injection pressure and holding pressure	1
1 1020	to nanoliters Injected volumes from picoliters	'	3 1	
PV830	to nanoliters	1	Injection pressure and holding pressure and vacuum	1
ANOLITER2010	Bolus, 2.3-69 mL/Injection	1	Oocyte injector, infuse only	1
_	on Systems (Zebrafish, C. Elegans,	Drosophila,)	Xenopus oocytes)	1
licrofluidi	CS			
ExiGo	50 nL/min - 10 mL/min	1	Infuse only, feedback via integrated flow sensor, includes iPad	1
Mirus	100nL/min – 10mL/min ±1%	8	mini which can control up to four pumps Infuse only, reversible flow, ~600μL dead volume, PC control	1
Kimo	15 – 35 mL/hr ±4%	1	Infuse only, recirculating pump controlled by iPod Touch, Wi-Fi	1
		'	communication, <300μL dead volume	
upplies &				
	tandard and Autoclavable			1
	icrovolume, for UltraMicroPump III			1
Miniature Vac				1
	roinjection and Transplantation Mo	lds		1
	ecialty Microsyringe			1
NanoFil™ App	olication Kits			1
Pressure Man				1
	World Precision	an Inctrum <i>o</i>	ents www.wniinc.com	



Mini★Star™

Miniature DC Peristaltic Pump

This compact and lightweight peristaltic pump fits just about anywhere. It can be mounted directly on the bench, in a regular rack or to a post. The MiniStar's speed can be adjusted from 1rpm to 50 rpm. With recommended silicone tubing, the volume can be set from 0.06 mL/min to 14.0 mL/min. The MiniStar also features a hand held remote control that allows users to start and stop the pump, purge or adjust its speed and direction.

MINISTAR	Miniature Peristaltic Pump, 1-channel
504011	MiniStar and Stand (as pictured above)
503120	TTL Control Module
503121	Silicone Tubing w stops, 2.4mm ID x 0.8mm wall x 1 m (5-pk)
503122	Silicone Tubing w stops, 1mm ID x 1mm wall x 1 m (5-pk)

MINISTAR SPECIFICATIONS

CHANNEL SPEED 1-50.0 rpm, forward/reverse FLOW RANGE 0.06~14.0 mL/min RESOLUTION 1 rpm (0.1 rpm computer control)

SPEED CONTROL Remote control

DISPLAY Indicators for status and speed POWER 12 V DC (110/220 VAC adapter incl.) WORKING CONDITION Temperature 0-40°C, humidity <

US: Tel: 941-371-1003 • sales@wpiinc.com

TUBING Two-stop Silicone 0.8~1.0 mm Wall Thickness Outer Diameter ≤ 4.8mm

DIMENSION OF DRIVER 135×72×72 mm (L×W×H) DIMENSION OF REMOTE CONTROL 105×50×16 mm (L×W×H)

WEIGHT OF DRIVER 0.5 Kg

80%


Peri-Star™ Pro

High performance digital peristaltic pump at an affordable price!

- Display either rotation speed (RPM) or flow rate (mL/min)
- Wide flow range: 0.01 -280 mL/min
- Accuracy of flow rate: 0.5% using self calibration function
- Accuracy of speed: 0.1 rpm
- Large backlit digital LCD display
- Programmable for all tubing sizes between 0.8 mm and 6.4 mm ID
- Easy and fast tubing replacement using snap-on cartridges
- Membrane keypad allows easy programming while protecting controls from fluid entry
- Actively driven rollers by planetary gears for long lasting tubing life

Peri-Star™ Pro peristaltic pumps provide accurate and precise pumping with convenience and versatility. Peri-Star Pro can be run in either flow rate mode (mL/ min) or rotation speed mode (rpm). For good laboratory practice, pumps must be calibrated after changing the tubing and solution. Users can easily calibrate Peri-Star Pro to deliver flow as accurate as 0.5% in a wide flow range from 0.01 mL/min to 280 mL/min. Under rotation speed mode, the digitally controlled stepping motor provides accurate and reproducible operation with 0.1% rpm both forward and in

Large backlit digital LCD display provides readouts of rotation direction, flow rate or rotation speed, tubing ID, drive status and



remote control mode simultaneously. Water resistant membrane keypad allows easy programming while protecting LCD display and controls from fluid entry.


Built-in Human Machine Interface (HMI) with screen instructions in plain English steps users through initial setup, calibration and operating procedures. The user-friendly interface reduces the need to frequently check the printed manual for instruction and reference.


Peri-Star Pro is available in two versions: a

4-roller version for high flow and an 8-roller version for lower volumes which provides high pressure with minimal pulsations.

A unique planetary gear design with eight actively driven rollers (four rollers for higher flow rate model), together with independent tubing compression fine adjustment, greatly increases flow accuracy and prolongs tubing life. Snap-on cartridges allow tubing to be changed quickly without cross contamination of solutions.

Available in 2-, 4- and 8-channel versions









D	EDI	_CTAD	DDO	CDECIEL	CATIONS

	Peri-Star Pro 2H / 4H (High Rate)	Peri-Star Pro 4L / 8L (Low Rate)
NUMBER OF ROLLERS	4	8
NUMBER OF CHANNELS	2–4	4-8
ROTOR SPEED RANGE	1–100 rpm	1–100 rpm
FLUID FLOW RANGE	0.8–280 mL/min #17 Tubing: 3.5–280 mL/min	0.01-80 mL/min #14 Tubing: 0.2-18 mL/min
TUBING RANGE	3.1-6.4 mm ID	0.5-2.4 mm ID
SELF-CALIBRATION	Yes	Yes
WORKING ENVIRONMENT	0-45°C, Humidity < 80%	0-45°C, Humidity < 80%
POWER	110 V or 220 V AC, 50 - 60 Hz	110 V or 220 V AC, 50 - 60 Hz
DIMENSIONS	190 x 162 x 275 mm	190 x 162 x 275 mm
SHIPPING WEIGHT	11 lb / 5 kg	11 lb / 5 kg

NEW PERI-STAR PRO PUMPS

INCHES I ENTER ST	AN THO TOWN 5	
PERIPRO-2HS	Peri-Star™ Pro, 2-channel, High Rate, Large Tubing (110-220V)	
PERIPRO-4HS	Peri-Star™ Pro, 4-channel, High Rate, Large Tubing (110-220V)	
PERIPRO-4LS	Peri-Star™ Pro, 4-channel, Low Rate, Small Tubing (110-220V)	
PERIPRO-8LS	Peri-Star™ Pro, 8-channel, Low Rate, Small Tubing (110-220V)	
OPTIONAL A	CCESSORIES	
503049	Replacement Tubing Cartridge, Large	
503050	Replacement Tubing Cartridge, Small	
503022	Replacement Silicone Tubing, 1m, 1.6 mm I.D., #14, with stops	
503023	Replacement Silicone Tubing, 1m, 6.4 mm I.D., #17	
503120	TTI Control Module	


US: Tel: 941-371-1003 • sales@wpiinc.com

159





- Manual or Programmable PC control with user-friendly GUI interface
- Fast LAFF solenoid valve
- Color-coded polyurethane tubing for easy identification
- Super low dead volume (<100 nL) micromanifold
- Economically priced



Micromanifold closeup: Fluid-filled passages are shown in magenta.

MPS-2 is a programmable 8-channel perfusion system designed for single channel and whole-cell patch preparations. Offering the best combination of performance and value, the MPS-2 incorporates the same high quality solenoid valves found on similar but much more expensive systems. Unlike other perfusion systems on the market, which often compromise performance to fit every possible application, the MPS-2 is the only perfusion system designed and optimized specifically for single-channel and whole-cell patch perfusion applications.

The system can be controlled manually via membrane switches on the front panel or through a PC. Two different manual control modes are offered. One controls each channel independently and the other mode allows the user to assign a master channel that will keep the system flow when all other channels are switched off. User-friendly graphic timing software is included, and the programmed perfusion sequence can be started by computer, a patch clamp amplifier or other external trigger, or manually by the user

The perfusion fluid flows through specially designed color-coded

polyurethane ribbon style tubing. The color-coding allows the user to easily trace each channel for diagnostic checks or set up and the ribbon style of tubing keeps the system very neat and organized. Unlike PVC based tubing, polyurethane tubing contains no plasticizer, which can cause contamination. The tubing ribbon is designed as an economical disposable item, which is often critical when cleanness is needed.

The most unique feature of the MPS-2 is its perfusion micromanifold. Using the latest microfluidic techniques, the injection molded micromanifold provides the least flow resistance and dead volume of any product on the market. The flow channel inner diameter is approximately 1.0 mm, except for the last 5 mm before the junction point. This design allows a fast flow rate without using a pressured system. The maximum flow rates are 1 and 16 microliter per second for the 50 mm long 100 μ m and 250 μ m ID tips, respectively. Small channels and a unique design at the merging point further reduce the chance of cross contamination. Dead volume is less than 100 nL.

MPS-2 Multichannel Perfusion System & Control Software

REPLACEMENT PARTS

502109-15 Color-coded Polyurethane Tubing, 1/16" ID x 8 Channels, 15 ft

502110 Micromanifold, 100 µm ID tip, 2 pcs/pk 502125 Micromanifold, 200 µm ID tip, 2 pcs/pk

Specify line voltage and Micromanifold tip OD when ordering.

MPS-2 SPECIFICATIONS

CHANNELS 8
VALVE RESPONSE TIME 2

VALVE CONTROL USB, TTL, external start via software

SYRINGE RESERVOIR VOLUME 10 mL

MANIFOLD 8 to 1

TIP ID 250 micron and 100 micron.

MAXIMUM FLOW RATES (gravity fed) 100 µm ID tip, 8 µL/min. at 50 cm

 $250 \ \mu m \ ID \ tip, \ 250-500 \ \mu L/min. \ at \ 50 \ cm$ DEAD VOLUME < 100 nL excluding the single outlet tubing

161



Programmable Syringe Pump



More Features

Dual Pumping Action

continuous pumping with optional check valve set.

ALADDIN SPECIFICATIONS

ALADDIN SPECIFICATIONS				
	AL-1000	AL-1000HP		
SYRINGE SIZES	Plastic syringes up to 60 mL and selected glass micro syringes from 0.5 to 500 μ L.	Plastic syringes up to 60 mL and selected glass micro syringes from 0.5 to 500 $\mu L. $		
NUMBER OF SYRINGES	1	1		
MOTOR TYPE	Step Motor, 1/8 to 1/2 step modes			
STEPS PER REVOLUTIONS	400	200		
STEPPING (max./min.)	0.21 μm to 0.850 μm			
MOTOR TO DRIVE SCREW RATIO	15/28	15/28		
SPEED (max./min.)	5.1 cm/min / 0.0042 cm/hr	18.36964 cm/min / 0.008409 cm/hr		
PUMPING RATES	1699 mL/hr with 60 mL syringe, to 0.73 μL/hr with 1 mL syringe	6120 mL/hr with 60 mL syringe, to 1.459 mL/hr with 1 mL syringe		
MAXIMUM FORCE	35 lb at min. speed, 18 lb at maximum speed	100 lb at minimum speed, 18 lb at maximum speed		
NUMBER OF PROGRAM PHASES	41	41		
RS-232 PUMP NETWORK	100 pumps maximum	100 pumps maximum		
POWER SUPPLY	Wall adapter 12V DC @ 850 mA	Wall adapter 12V DC @ 1000 mA		
DIMENSIONS	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)		

1.6 kg (3.6 lb)

WEIGHT

Need a pump for two syringes? Two Aladdin pumps when daisychained are more efficient and affordable than any competitor's dual syringe models. Two Aladdins (AL-2000) will perform as a dual infusion/withdrawal pump, a double pump for infusing at different rates, a push/pull pump with one infusing and one withdrawing at the same or different rates, two independent pumps, or a master/slave pump. One Aladdin can even control the second for

The Aladdin pump series will accept syringes from Becton Dickinson, Monoject, Terumo, and Air-Tite.

AL-1000	Programmable Syringe Pump
AL-1000HP	Programmable Syringe Pump, High Pressure
AL-2000	Two AL1000 Syringe Pumps

Includes GN-TTL Interconnecting Cable for push/pull or continuous pumping. Valves not included.

Specify line voltage

When ordering 220V models, specify UK, Euro or Australian line cord.

OPTIONAL ACCESSORIES

OI HOUTE A	CESSORIES
GN-PC7	PC to pump cable, 7 ft
GN-PC25	PC to pump cable, 25 ft
GN-NET7	Pump-to-pump network cable, 7 ft
GN-NET25	Pump-to-pump network cable, 25 ft
GN-TTL	Pump-to-pump reciprocating cable
ADPT2	Footswitch

World Precision Instruments

1.6 kg (3.6 lb)

Legato Syringe Pumps

The large touch screen color display lets you see all of the pump's operating parameters to ensure proper operation during the experiments. Syringe size and flow rate are easily displayed, as well as the volume delivered and elapsed time. Set up is easy using the icon-driven software. An LED on the front panel makes it easy to see if the pump is running. Advanced micro stepping techniques are employed to further reduce the step angle to eliminate flow pulsation. Accuracy is

±0.5%. A wide dynamic flow range from picoliters per minute to millimeters per minute can be programmed into the pump. These versatile pumps can be connected through an RS485 interface. Add the new *Adagio* software to maximize the use of the pump's functions and features. Adagio allows you to configure the pump through the software, as well as operate one or multiple pumps. LabVIEW drivers are available on the National Instruments website.



· Better flow performance

Features

- Accuracy ±0.35%
- Holds one or two syringes from 0.5 μ L to 140 mL
- High resolution color touch screen
- · Real time clock
- · Unparalleled ease of use
- · Touch pad "lock" feature
- · LED light on front panel
- · Full metal chassis

Benefits

- · Automatic dispensing of small volumes
- · Constant delivery of fluids
- · Hands free operation
 - · Built in syringe table
 - Up to 75 lb linear force
 - Advanced microstepping techniques
 - Built in RS-485 interface to link multiple pumps
 - · USB port & RS232 Interface
 - I/O & TTL interface
 - · Continuous mode of operation
 - Protection with a spill dam
 - · Analog control option
 - CE, UL, CSA, CB Scheme, EU RoHS

including gl	ass, plastic or stainless steel, are held securely in plac	e.
SPLG100	Legato 100 Syringe Pump, Infuse-Only	
504578	Software Adagio/USB Key	

The SPLG100 is the world's first single-syringe infusion-only pump

with a touchscreen interface. The *SPLG100* has a wide flow rate range from 1.26 pL/min to 88.32 mL/min, depending on syringe size. It

accommodates a single syringe from 0.5 µL to 60 mL. Any type of syringe,

SPLG100 SPECIFICATIONS

SYRINGE SIZE
POWER
MOTOR DRIVE CONTROL
LINEAR FORCE (MAXIMUM)
NUMBER OF MICROSTEPS PER ONE
REVOLUTION OF LEAD SCREW
STEP RATE (MINIMUM)
STEP RATE (MAXIMUM)
DRIVE MOTOR

PUSHER TRAVEL RATE (MINIMUM)
PUSHER TRAVEL RATE (MAXIMUM)

FLOW RATE (MINIMUM)
FLOW RATE (MAXIMUM)
DIMENSIONS
WEIGHT

CONNECTORS

0.5 µL to 60 mL

100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse Microprocessor with 1/16 microstepping 13.6 kg (30 lb) @ 100% Force Selection

15,360

27.5 sec/µstep 26 µsec/µstep

0.9 degree Stepper Motor

0.15 µm/min 159 mm/min

1.26 pL/ min (0.5 µL syringe) 88.32 mL/min (60 mL syringe) 22.6 x 19.05 x 15 cm (9 x 7.5 x 5 in) 2.66 kg (5.9 lb)

RS485 - IEEE-1394 6 pos, USB Type B

FLOW RATES Diameter Minimum Syringe Maximum 0.5 μL 0.103 mm 1.260 pL/min 1.325 µL/min 1 μL 0.146 mm 2.520 pL/min 2.651 µL/min 2 μL 0.206 mm 5.100 pL/min 5.299 µL/min 14.690 µL/min 5 μL 0.343 mm 14.100 pL/min 10 μL 0.485 mm 28.260 pL/min 29.380 µL/min 25 μL 0.729 mm 63.900 pL/min 66.370 µL/min 50 µL 1.03 mm 132.500 µL/min 127.600 pL/min 100 µL 1.457 mm 255.20 pL/min 265.100 µL/min 250 µL 2.304 mm 638.300 nL/ min 662.900 µL/min 500 µL 3.256 mm 1.275 nL/min 1.324 mL/min 1000 µL 4.608 mm 2.553 nL/min 2.652 mL/min 4.699 mm 2.655 nL/min 2.757 mL/min 1 mL 3 mL 8.585 mm 8.863 nL/min 9.204 mL/min 11.989 mm 5 mL 17.290 nL/min 17.950 mL/min 10 ml 14.427 mm 25.030 nL/min 25.990 mL/min 20 mL 19.05 mm 45.320 mL/min 43.640 nL/min 30 mL 21.59 mm 56.050 nL/min 58.210 mL/min 60 mL 26.594 mm 85.050 nL/ min 88.320 mL/min



SPLG101

Dual Infuse-Only Syringe Pump

The SPLG101 is ideal for applications where dual syringes are required with small volumes under 10mL. It accommodates two syringes from 0.5 µL to 10 mL. The SPLG101 has a wide flow rate range from 1.26 pL/ min to 25.99 mL/min, depending on syringe size.

SPLG101 Legato 101 Syringe Pump, Dual Infuse-Only

504578 Software Adagio/USB Key

SPLG101 SPECIFICATIONS

SYRINGE SIZE 0.5 µL to 10 mL POWER 100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse MOTOR DRIVE CONTROL Microprocessor with 1/16 microstepping LINEAR FORCE (MAXIMUM) 13.6 kg (30 lbs) @ 100% Force Selection NUMBER OF MICROSTEPS PER ONE 15.360

STEP RATE (MINIMUM) 27.5 sec/ustep STEP RATE (MAXIMUM) 26 µsec/µstep

REVOLUTION OF LEAD SCREW

DRIVE MOTOR 0.9 degree Stepper Motor

PUSHER TRAVEL RATE (MINIMUM) 0.15 um/min PUSHER TRAVEL RATE (MAXIMUM) 159 mm/min

FLOW RATE (MINIMUM) 1.26 pL/ min (0.5 µL syringe) 25.99 mL/min (10 mL syringe) FLOW RATE (MAXIMUM **DIMENSIONS** 22.6 x 19.05 x 15 cm (9 x 7.5 x 5 in)

WEIGHT 2.66 kg (5.9 lbs)

CONNECTORS RS-232 - 9 Pin D-Sub Connector,

RS485 - IEEE-1394 6 pos,

USB - Type B

CE

SPLG110

Infuse/Withdraw Syringe Pump

The SPLG110 offers infuse/withdraw flow control and programmability for up to two multi-step programs of 50 steps each. The SPLG110 has a wide flow rate range from 1.26 pL/min to 88.28 mL/min, depending on syringe size. The SPLG110 accommodates a single syringe from 0.5 µL to 60 mL. Any type of syringe can be used in the unit including glass, plastic or stainless steel. The pump is ideal for more complex multi-step dosing and has all multi-mode operation including infusion only, withdrawal only, infusion and withdrawal and withdrawal/infusion modes.

SPLG110 Legato 110 Syringe Pump, Infuse/Withdraw

504578 Software Adagio/USB Key

SPLG110 SPECIFICATIONS

SYRINGE SIZE 0.5 µL to 60 mL **POWER** 100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse MOTOR DRIVE CONTROL Microprocessor with 1/16 microstepping

LINEAR FORCE (MAXIMUM) 13.6 kg (30 lbs) @ 100% Force Selection NUMBER OF MICROSTEPS PER 15,360

ONE REVOLUTION OF LEAD SCREW

STEP RATE (MINIMUM) 27.5 sec/ustep

STEP RATE (MAXIMUM) 26 µsec/µstep

0.9 degree Stepper Motor **DRIVE MOTOR**

PUSHER TRAVEL RATE (MINIMUM) 0.15 µm/min PUSHER TRAVEL RATE (MAXIMUM) 159 mm/min

FLOW RATE (MINIMUM) 1.26 pL/ min (0.5 µl syringe) FLOW RATE (MAXIMUM) 88.28 mL/min (60 ml syringe) DIMENSIONS 22.6 x 19.05 x 15 cm (9 x 7.5 x 5 in)

WEIGHT 2.66 kg (5.9 lbs)

CONNECTORS RS-232 - 9 Pin D-Sub Connector,

RS485 - IEEE-1394 6 pos,

USB - Type B

36.003 ml/min E LEVEL PLES SERECTION: 00:00:09 THE ELAPSED 00:00:08 THE REHADING: 4.741 ml

Intuitive Run Screen — Combining multiple parameters simultaneously with internationally recognizable icons allow the Legato™ Series to provide a new level of intuitive syringe pump operation.

Adagio Software

 Low Cost Simple Installation
 Flow Evolution Graph
 Import & Export Programs
 Quick & Easy Manual Pump Control ● Monitor One or More Pumps ● Program Data Logging

The manual pump control tool allows easy direct control of the pump. Pump commands can be entered directly into the log. Multiple programs can be opened at the same time. The program's progression is tracked, and can be stored in a file for later access.

System Requirements:

• 1 GHz Pentium processor or higher • 512 MB of RAM (1 GB recommended) ● Windows XP SP3 or Vista (XP recommended) ● Free RS232 or USB 2.0 ports (depending on the quantity, model and connectivity of the controlled pumps; daisy chained pumps require a single port; direct PC-to-Pump connections require one free port per pump)

Microsoft Excel 97 or higher





Infuse/Withdraw Continuous

The SPLG270 is a Push-Pull syringe pump. It accommodates two syringes from 0.5 μL to 140 mL for infusion and two syringes for withdrawal. This model supports infusion and withdrawal simultaneously at user-defined flow rates and with selectable target volumes to control the total volume pumped. It also supports infuse only, withdraw only, infuse/withdraw, withdraw/infuse and continuous mode. The touch screen interface lets you quickly create configurations and recall them for easy use. The 4.3-inch TFT color display with touch pad interface presents all the pump operating parameters on one easy-to-view run screen. Protective cover over the display prevents leakage into the display.

SPLG270	SPL Syringe Pump, Push-Pull
SPLG272	SPL Syringe Pump, Push-Pull Programmable
504576	Small Syringe Multi Rack
	(for six 30-60mL syringes or ten 0.5µL-20mL syringes)
504577	Large Syringe Multi Rack
	(for up to four 60-140 mL plastic syringes)
504578	Software Adagio/USB Key

SPLG270 SPECIFICATIONS

SYRINGE SIZE 0.5 uL to 140 mL 100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse **POWER** MOTOR DRIVE CONTROL Microprocessor with 1/16 microstepping LINEAR FORCE (MAXIMUM) 34 kg (75 lbs) @ 100% force selection NUMBER OF MICROSTEPS PER ONE 6400 REVOLUTION OF LEAD SCREW

STEP RATE (MINIMUM) 27.5 sec/ustep STEP RATE (MAXIMUM) 26 µsec/µstep

DRIVE MOTOR 1.8 degree Stepper Motor

PUSHER TRAVEL RATE (MINIMUM) 0.36 µm/min PUSHER TRAVEL RATE (MAXIMUM) 190.80 mm/min

FLOW RATE (MINIMUM) 5 pL/min (0.5 µL syringe)

FLOW RATE (MAXIMUM) 215.803 mL/min (140 mL syringe) DIMENSIONS 8.89 x 25.4 x 27.94 cm (3.5 x 10 x 11 in)

WEIGHT 4.9 kg (10.75 lb)

CONNECTORS RS-232 - 9 Pin D-Sub Connector, RS-485 - IEEE-1394 6 pos,

USB - Type B,

I/O & TTL - 15 Pin D-Sub Connector



SPLG210

Infuse/Withdraw Syringe Pump

The SPLG210 Infuse/Withdraw syringe pump offers unparalleled ease of use through the high resolution touch screen. The basic model works with one syringe or two (from 0.5µL to 140mL) and can be reconfigured in the field to be used with multiple syringes. rotective cover over the display prevents leakage into the display. To optimize your bench space the SPLG210 can be placed on its side to reduce the footprint to only 3.5 in. x 9.75 in. The display also tilts with the change to allow the user to operate the pump vertically.

The programmable model offers maximum flexibility for configuring and running different programs. Up to 40 programs of 20 steps each can be configured and stored in the unit and recalled quickly with the touch of a button.

SPLG210	SPL Syringe Pump, Infuse/Withdraw
SPLG212	SPL Syringe Pump, Infuse/Withdraw Programmable
504576	Small Syringe Multi Rack
	(for six 30-60mL syringes or ten 0.5µL-20mL syringes)
504577	Large Syringe Multi Rack
	(for up to four 60-140 mL plastic syringes)
504578	Software Adagio/USB Kev

SPLG210 SPECIFICATIONS

SYRINGE SIZE	0.5 μL to 140 mL
POWER	100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse
MOTOR DRIVE CONTROL	Microprocessor with 1/16 microstepping
LINEAR FORCE (MAXIMUM)	34 kg (75 lbs) @ 100% force selection
NUMBER OF MICROSTEPS PER ONE	6400
REVOLUTION OF LEAD SCREW	
STEP RATE (MINIMUM	27.5 sec/µstep
STEP RATE (MAXIMUM)	26 µsec/µstep
DRIVE MOTOR	1.8 degree Stepper Motor
PUSHER TRAVEL RATE (MINIMUM)	0.36 µm/min
PUSHER TRAVEL RATE (MAXIMUM	190.80 mm/min
FLOW RATE (MINIMUM)	5 pL/min (0.5 μL syringe)
FLOW RATE (MAXIMUM)	215.803 mL/min (140 mL syringe)
DIMENSIONS	8.89 x 25.4 x 27.94 cm (3.5 x 10 x 11 in))
WEIGHT	4.9 kg (10.75 lb)
CONNECTORS	RS-232 - 9 Pin D-Sub Connector, RS-485 - IEEE-1394 6 pos

I/O & TTL - 15 Pin D-Sub Connector

USB - Type B,



SPLG200	SPL Syringe Pump, Infuse Only
504576	Small Syringe Multi Rack
	(for six 30-60mL syringes or ten 0.5µL-20mL syringes)
504577	Large Syringe Multi Rack
	(for up to four 60-140 mL plastic syringes)
504578	Software Adagio/USB Key

configured and stored in the unit and recalled quickly with the touch of a

button.

SPLG200 S	PECIFICTIONS
SYRINGE SIZE	0.5 μL to 140 mL
POWER	100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse
MOTOR DRIVE CONTROL	Microprocessor with 1/16 microstepping
LINEAR FORCE (MAXIMUM)	34 kg (75 lbs) @ 100% force selection
NUMBER OF MICROSTEPS PER ONE REVOLUTION OF LEAD SCREW	6400
STEP RATE (MINIMUM)	27.5 sec/µstep
STEP RATE (MAXIMUM)	26 µsec/µstep
DRIVE MOTOR	1.8 degree Stepper Motor
PUSHER TRAVEL RATE (MINIMUM)	0.36 μm/min
PUSHER TRAVEL RATE (MAXIMUM)	190.80 mm/min
FLOW RATE (MINIMUM)	5 pL/min (0.5 μL syringe)
FLOW RATE (MAXIMUM)	215.803 mL/min (140 mL syringe)
DIMENSIONS	8.89 x 25.4 x 27.94 cm (3.5 x 10 x 11 in.)
WEIGHT	4.9 kg (10.75 lb)
CONNECTORS	RS-232 - 9 Pin D-Sub Connector,
	S485 - IEEE-1394 6 pos,
	USB - Type B,
	I/O & TTL - 15 Pin D-Sub Connector

Syringe Diameter Minimum Maximum 0.5 μL 0.103 mm 3.12 pL/min 1.589 μL/min 1 μL 0.146 mm 6.18 pL/min 3.180 μL/min 2 μL 0.206 mm 12.301 pL/min 6.358 μL/min 5 μL 0.343 mm 33.96 pL/min 17.630 μL/min 10 μL 0.485 mm 67.72 pL/min 35.249 μL/min 25 μL 0.729 mm 153.42 pL/min 79.640 μL/min 50 μL 1.03 mm 306.24 pL/min 158.984 μL/min 100 μL 1.457 mm 612.72 pL/min 318.126 μL/min 250 μL 2.304 mm 1.533 nL/ min 795.51 μL/min 500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min	SPLG200 FLOW RATES			
1 μL 0.146 mm 6.18 pL/min 3.180 μL/min 2 μL 0.206 mm 12.301 pL/min 6.358 μL/min 5 μL 0.343 mm 33.96 pL/min 17.630 μL/min 10 μL 0.485 mm 67.72 pL/min 35.249 μL/min 25 μL 0.729 mm 153.42 pL/min 79.640 μL/min 50 μL 1.03 mm 306.24 pL/min 158.984 μL/min 100 μL 1.457 mm 612.72 pL/min 318.126 μL/min 250 μL 2.304 mm 1.533 nL/ min 795.51 μL/min 500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min				
2 μL 0.206 mm 12.301 pL/min 6.358 μL/min 5 μL 0.343 mm 33.96 pL/min 17.630 μL/min 10 μL 0.485 mm 67.72 pL/min 35.249 μL/min 25 μL 0.729 mm 153.42 pL/min 79.640 μL/min 50 μL 1.03 mm 306.24 pL/min 158.984 μL/min 100 μL 1.457 mm 612.72 pL/min 318.126 μL/min 250 μL 2.304 mm 1.533 nL/ min 795.51 μL/min 500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min				
5 μL 0.343 mm 33.96 pL/min 17.630 μL/min 10 μL 0.485 mm 67.72 pL/min 35.249 μL/min 25 μL 0.729 mm 153.42 pL/min 79.640 μL/min 50 μL 1.03 mm 306.24 pL/min 158.984 μL/min 100 μL 1.457 mm 612.72 pL/min 318.126 μL/min 250 μL 2.304 mm 1.533 nL/ min 795.51 μL/min 500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min				
10 μL 0.485 mm 67.72 pL/min 35.249 μL/min 25 μL 0.729 mm 153.42 pL/min 79.640 μL/min 50 μL 1.03 mm 306.24 pL/min 158.984 μL/min 100 μL 1.457 mm 612.72 pL/min 318.126 μL/min 250 μL 2.304 mm 1.533 nL/ min 795.51 μL/min 500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min				
25 μL 0.729 mm 153.42 pL/min 79.640 μL/min 50 μL 1.03 mm 306.24 pL/min 158.984 μL/min 100 μL 1.457 mm 612.72 pL/min 318.126 μL/min 250 μL 2.304 mm 1.533 nL/ min 795.51 μL/min 500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min				
50 μL 1.03 mm 306.24 pL/min 158.984 μL/min 100 μL 1.457 mm 612.72 pL/min 318.126 μL/min 250 μL 2.304 mm 1.533 nL/ min 795.51 μL/min 500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min				
100 μL 1.457 mm 612.72 pL/min 318.126 μL/min 250 μL 2.304 mm 1.533 nL/ min 795.51 μL/min 500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min				
250 μL 2.304 mm 1.533 nL/ min 795.51 μL/min 500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min	l			
500 μL 3.256 mm 3.06 nL/min 1.588 mL/min 1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min	1			
1000 μL 4.608 mm 6.129 nL/min 3.181 mL/min				
1 ml 4 600 mm 6 272 nl /min 2 200 ml /min				
1 mL 4.699 mm 6.373 nL/min 3.308 mL/min				
3 mL 8.585 mm 21.272 nL/min 11.044 mL/min				
5 mL 11.989 mm 41.485 nL/min 21.539 mL/min				
10 mL 14.427 mm 60.073 nL/min 31.19 mL/min				
20 mL 19.05 mm 104.74 nL/min 54.383 mL/min				
30 mL 21.59 mm 134.533 nL/min 69.852 mL/min				
50 mL 26.594 mm 204.122 nL/ min 105.985 mL/min	n			
100 mL 35.7 mm 367.839 nL/min 190.992 mL/min	n			
140 mL 38.4 mm 415.623 nL/min 215.803 mL/min	n			

SP Series Syringe Pumps

Syringe pumps for high metering precision at low, pulse-free rates

SP Pumps are sturdy and reliable, extremely simple to set up and use and surprisingly affordable. Liquid crystal displays (LCDs) prompt you through setup:

- 1. Select syringe from table stored in the pump's memory and displayed on the LCD.
- 2. Enter the volume to be dispensed.
- 3. Enter the flow rate and press Start.

It's fast and simple. Your settings are permanently stored in memory there's no need to re-enter them each day. SP pumps feature preset rate and volume control. Just set the volume you want dispensed. Volume is tracked continuously on the LCD display. Then, when the preset volume has been dispensed, the pump shuts off automatically. The easy-to-read digital display provides real-time readings using both parameters and values for clearer, mistake-free readings. The SP200 Series pumps offer TTL and RS-232C interfaces and automatic shutoff under stall conditions..

Two-Syringe Microdialysis Pump

See www.wpiinc.com/sppumps for specifications on the SP Series Syringe Pumps.

Single-Syringe Infusion Pump


SP100i

This inexpensive single-syringe infusion pump combines precision and simplicity with outstanding ease of use and durability.

- Holds any size syringe, 10 μL to 50 mL
- Automatic volume control and
- Simple menu-driven setup: dispense volume, dispense flow rate, syringe diameter
- Last settings stored in permanent memory



microdialysis experiments, this pump produces very fine syringe movement. The modified SP100i gearing features a sixfold gear reduction compared to standard models, allowing pumping at much smaller flow rates. (See the Micro-C for detection of dialysates.)



• Holds two syringes, 10 μL to 10 mL



Continuous dispense volume display Two-Syringe Push-Pull Pump

SP120p



A second syringe mount has been added to the basic SP100i, with both syringes activated by a single pusher block for simultaneous infusion and withdrawal.

- All the features of SP100i
- Holds two syringes, from 10 μL to 10 mL.

Two-Syringe Infusion Pump

SP200i

This feature-laden, two-syringe, infusion pump combines a broad speed range and holds a wide range of syringe sizes to meet the requirements of virtually any laboratory application.

- Holds two syringes, 10 to 140 mL
- Knob locks/unlocks drive block for effortless, drag-free adjustment
- Simple menu-driven setup: Syringe diameter, Dispense volume, Dispense flow
- Continuous dispense volume display
- Preset volume control and automatic shutoff
- Review or change settings during operation
- Optical encoder stall detection
- Choice of unit selection
- Last settings stored in permanent memory
- Built-in RS-232C interface for computer linking or "daisy chaining" up to 100 pumps.



Continuous Cycle Syringe Pump SP210c

The SP210c holds up to four syringes and can cycle continuously back and forth in a push-pull action.

As two syringes are infusing, two other syringes are withdrawing at the same rate. At the end of the set volume the direction is automatically reversed and the next cycle begins. With the use of 2-way valves, the pump can empty and refill syringes for continuous dispensing.

 Holds four syringes, 10 ml to 60 ml each



Two-Syringe Infusion/ Withdrawal Pump

SP210iw

The SP210iw offers more advanced features than any other infusion/withdrawal pump in its price rangeincluding five operating modes plus independent rate and volume settings for both infusion and withdrawal.

- All features of SP200i
- Independent rate and volume settings for infusion and withdrawal
- Multiple mode selection: infusion, withdrawal, infusion then withdrawal, withdrawal then infusion or continuous cycle.



Multi-Syringe Infusion/

Multi-Syringe Infusion Pump

Ideal for applications requiring multiple syringes, the SP220i is an adaptation of the SP200i and has been modified to hold up to 10 syringes.

All features of SP200i

SP220i

 Accommodates 10 syringes up to 10 mL, or six syringes up to 50 mL, or four syringes up to 140 mL.



Four-Syringe Nanoliter Infusion Pump SP250i

Each syringe can be sized differently and is clamped independently.

- Holds four syringes, up to 10 mL each
- Separate clamping accommodates different sizes
 - Syringes may be positioned independently for sequential dispensing by the pusher block



Withdrawal Pump SP230iw

Ideal for applications requiring multiple syringes, the SP230iw is an adaptation of the SP210iw, and has been modified to hold up to 10 syringes.

- All features of SP200iw
- Multiple syringe holder accommodates 10 syringes up to 10 mL, or six syringes up to 50 mL, or four syringes up to 140 mL.

Four-Syringe Push-Pull Pump

CE

SP260p single cycle

The SP260p can hold up to four syringes. As two syringes are infusing, two other syringes are withdrawing at the same rate. The SP260p is used for single-cycle applications only.

- All the features of SP200i
- Holds up to 4 syringes



SP100i	Syringe Pump, Infusion (single), 95-135 V
SP100iZ	Syringe Pump, Infusion (single), 220-240 V
SP101i	Syringe Pump, Microdialysis (double, slow speed), 95-135 V
SP101iZ	Syringe Pump, Microdialysis (double, slow speed), 220-240 V
SP120p	Syringe Pump, Infusion-Withdrawal (double), 95-135 V
SP120pZ	Syringe Pump, Infusion-Withdrawal (double), 220-240 V
SP200i	Syringe Pump, Infusion (double), 95-135 V
SP200iZ	Syringe Pump, Infusion (double), 220-240 V
SP210c	Syringe Pump, Infusion & Withdrawl (Continuous Action), 95-135 V
SP210cZ	Syringe Pump, Infusion & Withdrawl (Continuous Action), 220-240 V
SP210iw	Syringe Pump, Infusion & Withdrawal (double), 95-135 V
SP210iwZ	Syringe Pump, Infusion & Withdrawal (double), 220-240 V
SP220i	Syringe Pump, Infusion (multiple), 95-135 V
SP220iZ	Syringe Pump, Infusion (multiple), 220-240 V
SP230iw	Syringe Pump, Infusion & Withdrawal (multiple), 95-135 V
SP230iwZ	Syringe Pump, Infusion & Withdrawal (multiple), 220-240 V
SP250i	Syringe Pump, Infusion (multiple, mixed volumes), 95-135 V
SP250iZ	Syringe Pump, Infusion (multiple, mixed volumes), 220-240 V
SP260p	Syringe Pump, Infusion-Withdrawal (double) Single Cycle Action, 95-135 V
SP260pZ	Syringe Pump, Infusion-Withdrawal (double) Single Cycle Action, 220-240 V
	All 240-volt pumps are CE-approved.
####-A	Audible Alarm (add "A" to pump part number when ordering)
####-P	Programmable Ramp Option (SP200 Series)
OPTIONA	L CABLES
15623	Serial cable, SP Pump-to-IBM 9-pin "D" connector
13685	SP Pump-to-Pump "Daisy-Chain" linking cable, 7 ft
13962	Footswitch for SP200 Series Pumps

World Precision Instruments



- Lightweight and conformable ergonomic design
- Easy calibration using provided tool
- Made from biologically inactive and chemical inert polymers
- Easy for cleaning and parts replacement
- CE and ISO13485 Certified

Model	Volume Range µL	Increment µL	Nominal Volume µL	Tolerance %	Repeatability %
REG2	0.2 ~ 2	0.01	0.2 0.5 2	±12.0 ±5.0 ±2.0	≤6.00 ≤2.50 ≤0.70
REG10	1 ~ 10	0.1	1 5 10	±3.0 ±1.5 ±1.0	≤1.50 ≤0.60 ≤0.40
REG20	2 ~ 20	0.1	2 10 20	±3.0 ±1.0 ±1.0	≤1.50 ≤0.50 ≤0.30
REG50	5 ~ 50	0.5	5 20 50	±2.0 ±1.2 ±1.0	≤1.50 ≤0.40 ≤0.20
REG100	10 ~ 100	1	10 50 100	±2.0 ±0.8 ±0.8	≤0.50 ≤0.30 ≤0.15
REG200	20 ~ 200	1	20 100 200	±2.0 ±0.8 ±0.8	≤0.50 ≤0.30 ≤0.15
REG1000	100 ~ 1000	5	100 500 1000	±1.5 ±0.8 ±0.8	≤0.30 ≤0.30 ≤0.15
REG5K	1000 ~ 5000	50	1000 2000 5000	±1.0 ±0.7 ±0.7	≤0.50 ≤0.25 ≤0.15
REG10K	1000 ~ 10000	100	1 mL 5 mL 10 mL	±3.0 ±0.7 ±0.7	≤0.30 ≤0.20 ≤0.15



GERS5	Regal Pipetters (set of any 5) & stand
GERS6	Regal Pipetters (set of any 6) & stand
GERS7	Regal Pipetters (set of any 7) & stand
504591	Stand for Regal Pipetters (holds 8)



Universal Filter Tips (sterile)						
Tip Volume	For Pipetter	Rack	Part No.			
0.1 - 10 μL	REG2 REG10 REG20	960 (10 racks of 96)	500199			
REG20 10 - 200 µL REG50 REG100 REG200		960 (10 racks of 96)	500200			
1000 - 10,000 μL	REG10K	250 (10 racks of 25)	504590			

Universal Tips						
Tip Volume	For Pipetter	Bulk	Part No.	Rack	Part No.	
0.1 - 10 μL	REG2 REG10 REG20	Bag of 1000	500191	960 (10 racks of 96)	500192	
5 - 200 μL	REG20 REG50 REG100 REG200	Bag of 1000	500193	960 (10 racks of 96)	500194	
100-1000 μL	REG1K	Bag of 1000	500195	1000 (10 racks of 100)	500196	
500 - 5000 μL	REG5K	Bag of 250	500197 *	500 (10 racks of 50)	500198 *	
1000 - 10,000 µL	REG10K	Bag of 1000	504588	250 (10 racks of 25)	504589	

WPI's Universal Pipette Tips are for use with Regal and most other pipetters, including Gilson, Oxford Benchmate, Socorex, and

* Tips 500197 and 500198 fit Regal, Eppendorf, and BioHit pipetters.

same as leading brands the price!



UltraMicroPump III



This versatile injector uses microsyringes to deliver picoliter volumes

Perfect for a wide range of applications including intracellular injection, micro delivery of biochemical agents or dyes, cell separation, and in vitro fertilization.

See "Reproducible and Efficient Murine CNS Gene Delivery Using a Microprocessor Controlled Injector," A.I. Brooks et al., Journal of Neuroscience Methods, 80 (1998) 137-147.

ULTRAMICROPUMP SPECIFICATIONS

(based on 10 μL syringe)

NORMAL MODE

TRAVEL

MINIMUM DISPENSING VOLUME 0.58

LINEAR MOTION PER STEP

WEIGHT

MOUNTING ROD DIAMETERS
MAINS POWER SUPPLY

DIMENSIONS

62 mm

 $0.58~\text{nL}\,/\,\text{step}$ (10 μL syringe)

3.175 microns 325 g (11.5 oz)

7.9 mm (0.31 in.) 90-264VAC @ 47-63Hz

Ø 32 mm x 190 mm (Ø 1.3 in. x 7.5 in.)

MICROSTEPPING MODE

Precision is increased eight-fold

The world's leading ultramicropump!

With its digital controller, **UltraMicroPump III** can dispense as little as 600 picoliters per incremental advance of the syringe piston (using a 5 μ L syringe). Syringes may be filled externally and then inserted into the pump or filled while mounted in the pump. Fluids injected or withdrawn are held entirely within the micro syringe to maintain a low fluid dead volume.



For positioning, the **UltraMicroPump III** may be attached to any of several WPI micropositioners such as the **M3301** (manual), **DC3001** (motorized), or any manual stereotaxic manipulator.

UMPIII shown mounted to stereotaxic frame (not included).

Smart Controller

An Integral component in the **UMPIII** system is a microprocessor-based controller, **SYS-Micro4**, which provides an "intelligent" and easy-to-use interface to up to four syringe pumps. Operating parameters are set with the membrane keypad and LCD display. From the keypad the user can select the following functions: set pump to infusion or withdrawal mode, enter the volume to be infused or withdrawn, rate of delivery, and syringe type as well as synchronize the starting and stopping of any combination of syringe pumps.

User parameters can be stored in the device's "non-volatile" memory for instant recall when the unit is powered on.

An optional footswitch can be plugged into a connector on the rear of the controller for "hands free" start-/-stop operation.

Computer Control—An **RS-232 port** on the rear of the controller can be used to connect it to a computer for use with computer control programs.

UMPIII ACCEPTS: glass syringes with barrel diameters from 5.5 to 9 mm.				
UMP3-1	UltraMicroPump III (one) and Micro4 Controller			
UMP3-2	UltraMicroPump III (two) and Micro4 Controller			
UMP3-3	UltraMicroPump III (three) and Micro4 Controller			
UMP3-4	UltraMicroPump III (four) and Micro4 Controller			
UMP3	UltraMicroPump III (without controller)			
SYS-MICRO4	Micro4 Controller, Four-Channel			
OBTIONS AND ACCESSORIES				

OPTIONS AND ACCESSORIES			
15867	Footswitch for Micro4		
40500	RS-232 Cable, 9-pin "D" connector		
502201	V-clamp for Stereotaxic Frame		
503301	Extension Cable, miniDIN (male-female) 10 ft		
503207	Small Base Stand and Clamps		

Microvolume Syringes

Syringes with Luer Fitting (no needle)

Order No.	Volume	Description	O.D.	SCALE LENGTH	UMP3	UMP2
ILS005LT	5 μL	ILS 5 μL Gas-tight Luer tip	6.5 mm	54.1 mm	Υ	Υ
ILS010LT	10 μL	ILS 10 μL Gas-tight Luer tip	6.5 mm	54.1 mm	Υ	Υ
ILS025LT	25 μL	ILS 25 μL Gas-tight Luer tip	8.0 mm	60 mm	Υ	Υ
SGE050TLL	50 μL	SGE 50 µL Gas-tight Teflon Luer Lock	8.0 mm	60 mm	Υ	Υ
SGE100TLL	100 μL	SGE 100 µL Gas-tight Teflon Luer Lock	8.0 mm	60 mm	Υ	Υ
SGE250TLL	250 μL	SGE 250 µL Gas-tight Teflon Luer Lock	8.0 mm	60 mm	Υ	N

Syringes with Replaceable Beveled Needles

Order No.	Volume	Description	O.D.	SCALE LENGTH	UMP3	UMP2
SGE0005RN*	0.5 μL	SGE 0.5 µL 23 ga (0.63 mm), 70 mm long needle	8.0 mm	54.1 mm	Υ	Υ
SGE001RN*	1.0 µL	SGE 1.0 µL 26 ga (0.47 mm), 70 mm long needle	8.0 mm	54.1 mm	Υ	Υ
SGE005RN	5 μL	SGE 5 µL 23 ga (0.63 mm), 50 mm long needle	8.0 mm	54.1 mm	Υ	Υ
SGE010RNS	10 μL	SGE 10 µL 26 ga (0.47 mm), 50 mm long needle	8.0 mm	54.1 mm	Υ	Υ
SGE025RN	25 µL	SGE 25 µL 25 ga (0.50 mm), 50 mm long needle	8.0 mm	60 mm	Υ	Υ
SGE050RN	50 μL	SGE 50 µL 25 ga (0.50 mm), 50 mm long needle	8.0 mm	60 mm	Υ	Υ
SGE100RN	100 µL	SGE 100 µL 25 ga (0.50 mm), 50 mm long needle	8.0 mm	60 mm	Υ	Υ

^{*} The capacity of this syringe is so small that the entire sample is contained within the needle. The plunger extends to the tip of the needle, displacing the full sample during injection — which gives the syringe zero dead volume.

SGE and ILS are respective trademarks of Scientific Glass Engineering and Innovative Labor Systeme.

Replacement Needles

RN0005	For syringe SGE0005RN, 23 ga (0.63 mm) 70 mm long
RN001	For syringe SGE001RN, 26 ga (0.47 mm) 70 mm long
RN005	For syringe SGE005RN, 23 ga (0.63 mm) 50 mm long
RN010	For syringe SGE010RN(S), 26 ga (0.47 mm) 50 mm long, 5-pack
RN025	For syringes SGE025RN, SGE050RN, SGE0100RN, 26 ga (0.47 mm) 50 mm long, 5-pack



Pneumatic PicoPumps

Repeatable microinjection in volumes ranging from picoliters to nanoliters

Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's PicoPumps use carefully regulated air pressures for securing cells and injecting them with fluid. Injected volumes range from picoliters to nanoliters. Separate ports supply positive and negative pressure—positive pressure for high-pressure ejection, and suction for supporting the cell or for filling the pipette from the tip. A second pressure port maintains a low positive "holding" pressure to the injecting pipette between injection pulses, to prevent fluid uptake through capillary action or diffusion. Timing, ejection pressure, holding

pressure, and suction are adjusted independently by control knobs and indicator gauges on the front panel. Injection pressure is controlled by a 20-turn regulator on the front panel. A built-in timing circuit allows precise control of the amount of time that the injection pressure is applied to the output port. Time intervals can range from 10 seconds down to 10 ms or less, depending on the eject pressure setting. The injection pressure interval can be triggered manually on the front panel, by footswitch, or by computer controlled TTL pulse. A 5-volt monitor output provides a logiclevel pulse for your computer or other monitoring device.

The most recognized picopump in the world!



For a complete list of pre-pulled micropipettes, see µTips™, or call us with your special requirements.

PV830 — Eject pressure, Hold pressure, and Vacuum are all available, controlled by separate regulators on the front panel. Eject pressure supplies a high-pressure pulse for injecting fluid. Hold pressure, which is not sufficient to cause fluid ejection, is used to prevent back filling of the pipette by capillary action or diffusion when the solenoid is inactive. Pressure in the injection pipette is automatically switched between

Eject and Hold pressure by a precision timing circuit that controls a solenoid valve. Vacuum is used to fill pipettes from the tip or to secure a floating cell during microinjection. Vacuum is regulated the same way, by a 20-turn knob on the front panel. Vacuum may be switched from regulated vacuum to atmosphere by using the pneumatic toggle on the front panel. Vacuum can also be routed to the eject port.



New PicoNozzle Kit 5430-ALL (included) allows micropipettes to be securely mounted in micropositioners for stable axial air delivery. Because air enters the pipette axially, lateral whipping during injection is eliminated.

Each PicoPump is supplied with a 5430-ALL kit that includes two PicoNozzles and tubing

to connect the holders to the pressure and vacuum ports. **SYS-PV820** PicoPump w/ hold pressure

Specify line voltage All PicoPumps require external vacuum source — see below.

OPTIONAL	ACCESSORIES
3260	Foot Switch
2932	Rack Mount Kit, 3½-in. high (PV820)
2933	Rack Mount Kit, 5¼-in. high (PV830)
5430-10	PicoNozzle Kit (MPH6S for 1.0 mm pipette & 5-ft tubing assembly)
5430-12	PicoNozzle Kit (MPH6S for 1.2 mm pipette & 5-ft tubing assembly)
5430-15	PicoNozzle Kit (MPH6S for 1.5 mm pipette & 5-ft tubing assembly)
5430-20	PicoNozzle Kit (MPH6S for 2.0 mm pipette & 5-ft tubing assembly)
5430-ALL	PicoNozzle Kit (for 1.0, 1.2, 1.5, and 1.65 mm pipettes & 5-ft tubing assembly)
75122-110	Replacement Gaskets for PicoNozzle — 1.0 mm, green, package of 10
75122-210	Replacement Gaskets for PicoNozzle — 1.2 mm, black, package of 10
75122-310	Replacement Gaskets for PicoNozzle — 1.5 mm, blue, package of 10
75122-410	Replacement Gaskets for PicoNozzle — 1.65 mm, red, package of 10
MPH6S	Micropipette Holder (specify 1.0, 1.2, 1.5 or 2.0 mm)
MPH6R	Micropipette Holder (specify 1.0, 1.2, 1.5 or 2.0 mm)
3316	Replacement Input Kit

SYS-PV830 PicoPump w/ hold pressure and vacuum



PV820 offers separate regulated Hold and Ejection pressure, used to maintain a low pressure in the pipette between injections to prevent unwanted fluid uptake by capillary action or diffusion. A precision timing circuit switches from Eject pressure to Hold pressure automatically, once

timing has been set. Although regulated vacuum is not provided in this model, suction can be provided by connecting a vacuum source to the vacuum port on the rear panel. Suction is then available through the pressure ports.

PICOPUMP SPECIFICATIONS

PRESSURE

PRESSURE INPUT
PRESSURE OUTPUT
PULSE WIDTH (10-turn dial)
REGULATOR ACCURACY
REGULATOR REPEATABILITY
GAUGE ACCURACY
INPUT CONNECTOR
OUTPUT CONNECTOR
CONTROL

VACUUM

VACUUM INPUT
VACUUM OUTPUT
LOWEST REGULATED VACUUM
REGULATOR ACCURACY
REGULATOR REPEATABILITY
GAUGE ACCURACY
INPUT CONNECTOR
OUTPUT CONNECTOR
CONTROL
VENT

CONNECTIONS INCLUDED

INPUT KIT OUTPUT KIT

PHYSICAL SPECIFICATIONS

POWER
DIMENSIONS
SHIPPING WEIGHT

PV820

0 to 150 psi 0.3 to 90 psi * 10 ms to 10 s in Timed Mode 0.1% (20-turn dial) * 0.05 psi * 3% at full scale * Quick Connect (¼ in. OD Tubing) Barbed (¼-in. ID Tubing) Solenoid * Both Hold and Eject Pressures

0 to 30.0 in. Hg Unregulated Unregulated Unregulated Unregulated None

Quick Connect (¼ in. OD Tubing) Barbed (¼ in. ID Tubing) Manual

Atmosphere

PV830

0 to 150 psi 0.3 to 90 psi 10 ms to 10 s in Timed Mode 0.1% (20-turn dial) * 0.05 psi * 3% at full scale * Quick Connect (¼ in. OD Tubing) Barbed (¼6-in. ID Tubing) Solenoid * Both Hold and Eject Pressures

0 to 30.0 in. Hg
0.2 to 29.9 in. Hg
3 in. water
0.1% (20-turn dial)
0.03 in. Hg
3% at full scale
Quick Connect (¼ in. OD Tubing)
Barbed (¾ in. ID Tubing)
Manual

Atmosphere

10-ft nylon tubing (0.25-in. OD, 1000 psi), one ½-inch female NPT adapter Two PicoNozzle assemblies, each consisting of one MPH6S pipette holder, 60-in. of PVC tubing

95-135 V or 220-240 V, 50/60 Hz 17 x 3.5 x 9.5 in. (43 x 9 x 24 cm)

11 lb (5 kg)

95-135 V or 220-240 V, 50/60 Hz 17 x 5.25 x 9.5 in. (43 x 13 x 24 cm) 14 lb (6.3 kg)

(200 psi), and a luer-fitted aluminum handle.

Miniature Vacuum Pump

- Oil free
- Maintenance free
- Minimal vibration
- Low noise
- Extremely long life time
- Compact (18 x 7 x 7cm)
- Durable aluminum exterior



This miniature vacuum pump is durable and accurate, ideal for any application requiring a small, reliable pump that provides vacuum pressure up to 250 mbar. The industrial-strength aluminum exterior, neoprene diaphragm and neoprene/ silicone valves ensure this pump will stand up to daily use.

MINI VAC SPECIFICATIONS

230 (50 Hz) 120 (60 Hz) **Power Source** Motor Type Vibrating Power 4.0 W 4.0 W Free Flow 4.0 L/min. 3.0 L/min. At -100 mbar 2.0 L/min. 1.5 L/min. **Maximum Pressure** Maximum Vacuum -250 mbar -250 mbar Pump Head Construction Aluminum CR-neoprene Diaphragm Valves CR-neoprene/FPM (Viton)/Silicone 185 x 72 x 72 mm Dimensions Weight 850 g

801566	Mini Vacuum Pump (110V)
801963	Mini Vacuum Pump (220V)

Microinjection System

Zebrafish, C. Elegans, Drosophila, Xenopus Oocytes



ebrafish (*Danio Rerio*) are rapidly gaining in popularity as bio-medical research subjects because of the ability to generate high resolution, *in vivo* images of the embryos. Zebrafish are easy to maintain and produce a large number of offspring. Additionally, the embryos have a nearly transparent skin, making their development easily visible. These fish are used for a variety of disciplines, including neuroscience, genetics and aging studies.

Serving scientist for over 45 years, WPI offers a variety of instruments for microinjection including pumps, pipetters, microscopes and more. One of our most popular pumps for microinjection is the PV820 Pneumatic PicoPump.

The *PV820* and *PV830* Pneumatic PicoPumps were designed to simplify intracellular injection. You get repeatable microinjection in volumes ranging from picoliters to nanoliters. *PV820* offers eject and hold pressure. The hold pressure prevents backfilling of the pipette by capillary action. In addition, the *PV830* also has vacuum pressure which allows you to securely hold a cell with one pipette while you inject it with another. The volume injected is controlled by the inside diameter of the glass tip, the pressure and the time.

Recently, WPI introduced its customizable Microinjection System with everything you need to get started. The basic system is shown here. Starting below, you will find many options and accessories to customize your system.

Options for Customizing Your System

INJECTOR

- PV820 Pneumatic PicoPump with Hold Pressure
- PV830 Pneumatic PicoPump with Hold Pressure and Vacuum
- UMPIII UltraMicroPump
- Nanoliter2010



Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's **PicoPumps** use carefully regulated air pressures for securing cells and injecting them with fluid. Injected volumes range from picoliters to paneliters.



The versatile **UMP3** injector uses microsyringes to deliver picoliter volumes.



Microprocessorcontrolled **Nanoliter 2010** uses direct piston displacement.

Complete system for microinjection

MICROSCOPE

- **PZMIII** Precision Stereo Zoom Microscope on Track Stand
- PZMIV Precision Stereo Zoom on Track Stand
- 504928 LED Lighted Microscope Stand, 12.5"
- 504929 LED Lighted Microscope Stand, 10.5"
- 504596 Halogen Lighted Microscope Stand



PZMIII Precision Stereo Zoom Microscope on Track Stand



504928 LED Lighted Microscope Stand



504596 Halogen Lighted Microscope Stand

PULLERS

- PUL1000 Microprocessor-controlled 4-Step micropipette puller
- Nanofil Microliter syringes
- *MicroFil* for backfilling glass needles
- Glass capillaries
- **Pipetters**
- *MicroTip* pre-pulled pipettes
- E2XX Micropipette Storage Jar

Use the **E2XX** jars to store up to 30 micropipettes, filled or unfilled, up to three inches in length.



PUL-1000 is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, microperfusion or microinjection. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. This allows graduated cycles for applications like patch clamp recording.



MANIPULATOR

- M3301 Manual Micromanipulator
- KITE Manual Micromanipulator
- MMI Joystick-Controlled Micromanipulator
- DC3001 Motorized Micromanipulator
- SN-PCZ-50 Miniature Piezo Micromanipulator with controller



Weighing just 550 grams, the **M3301** is a well-built micromanipulator that outsells all others worldwide for high precision experiments where magnification is in the range of up to 250×.

> The MMJ is specially adapted for use with the *Nanoliter Injector* for oocyte injection and similar applications, this joystick-controlled micromanipulator allows an easy "steering" motion that translates normal hand movement into smooth submillimeter shifts.



LIGHTS

- Z-LITE-186 Fiber Optic Illuminator with (500186) Bifurcated Light Guides
- **LED-Lite** Modular LED Light Source with Exchangeable LEDS



LED-lite is a power supply for WPI's ELS LED modules for monochromatic light excitation.



The **Z-LITE** Fiber Optic Illuminator provides reliable highintensity light for microscopes.

ACCESSORIES

- 801566/801963 Vacuum Pump for use with the PV830
- Fluorodish Optical glass bottom dishes
- 504134 LED Ring Light
- M10 or M-3 manipulator base
- **Z-MOLDS** Microinjection and Transplantation Molds
- Many surgical instruments, including: 77020 Glass tweezers

500342 Dumont Tweezers #5 **504507** Fine pointed Swiss tweezers **14003** Vannas spring scissors



Z-MOLDS Microinjection and Transplantation Molds create grooves in agarose gel. Pipette the embryos into the grooves. The embryos self-align.

Microinjection and Transplantation Molds

Are you looking for a simplier method for handling zebrafish embryos?


To speed up your process and accuracy, try Z-MOLDS. Mold your agarose, pipette in your embryos, and watch them auto aligning in the grooves.

Z-MOLDS Microinjection and Transplantation Molds (4 per kit) are designed for zebrafish research. Using the molds is easy.

Turn the molds up-side down (ridged side down) and place them in liquid agarose gel.



Allow the agarose to solidify. Then, remove the



Pipette your embryo eggs into the grooves formed by the **Z-MOLDS**. The embryos self-



It's truly that easy. Now you are ready for microinjection.

See the video: www.wpiinc.com/z-molds



Proteomics and Large Screening

This mold is designed for injecting many embryos-up to 1000. The grooves made by the mold in the agarose gel enable the embryos to self align.







Xenograft and Larval Injection

This mold is designed for larval injections. The sloped ridges make perfect angles in the agarose gel, which then makes it easier to do microinjections in the larvae.







Transplantation

This mold is designed for increasing the speed of doing microinjections. Simply turn the petri dish as you are injecting.







Standard Microinjection

This mold is designed for blastomere transplantation.









Z-MOLDS Microinjection & Transplantation Molds



Nanoliter Injector For oocyte injection and applications in the 2 to 70 nanoliter range



Micromanipulator not included.

IN RE

GL ST IN

FIL

FΝ

SMALLEST VOLUME

SHIPPING WEIGHT

TO CHANGE VOLUME

NANOLITER 2010

INJECTIONS PER FILLING, MAX. 100 injections

WPI's microprocessor-controlled Nanoliter 2010 uses direct piston displacement. By either pushing the injection button on the control box or pressing on the optional footswitch, a calibrated volume will be smoothly injected. The process is quiet and vibration free. Capillary filling and injection speeds are 23 nL/sec and 46 nL/sec (emptying speeds are 92 nL/sec and 230 nL/ sec). Maximum fluid ejection is 5 μL. Each unit comes with sufficient glass to pull at least 300 tips. Glass is 1.14 mm O.D. (nominal) and 0.5 mm I.D.

By setting the DIP switch, the injection volume can be changed from 2.3 to 69.0 nL in 16 steps. Up to 100 injections may be triggered per filling. Since the volume of a normal Xenopus oocyte is about 500 nL, the instrument has the capability to inject from less than 1% to over 10% of the total volume of the oocyte in one preset step increment.

Included: 1 vial 3.5 in. capillaries (300); replacement "O" rings; Allen wrench; MicroFil™ MF34G backfilling needle; and two sample µTip™ prepulled micropipettes.

Smart Controller: Micro4, an optional microprocessor-based controller, can provide an "intelligent" and easy-to-use interface to up to four Nanoliter Injectors. Operating parameters are set with the membrane keypad and LCD display. From the keypad the user can set pump to infu-

sion or withdrawal mode, enter the volume to be infused or withdrawn, and rate of delivery, as well as synchronize the starting and stopping of a combination of Injectors. User parameters can be stored in the device's "non-volatile" memory for instant recall when the unit is powered on. An optional footswitch can be plugged into a connector on the rear of the controller for "hands free" start/stop



Nanoliter 2010 operation. An RS-232 port on the rear of the controller can be used to connect it to a computer for use with computer control programs.

US: Tel: 941-371-1003 • sales@wpiinc.com

WPI recommends the following equipment for a complete system,: PZMIII stereo microscope · MMJ joystick micromanipulator · M10 magnetic base · μTip™ micropipettes.

IVAIVOLITE			
SPECIFICATIONS		NANOLITER2010	Nanoliter 2010
		NL2010MC4	Nanoliter Injector & Micro4 Controller (small controller not included)
NJECTION VOLUME	Variable	OPTIONAL ACCES	CODIEC
REMOTE CONTROL	Yes		
GLASS OD	1.14 mm	13142	Footswitch for Nanoliter 2010
GLASS ID	0.5 mm	15867	Footswitch for Micro4 Controller
STEP	12.7 µ/step	4878	Replacement 3.5-in. glass capillaries (300)
NIECTION SPEED	p. step	4879	Replacement 7-in. glass capillaries (300)
Slow	23 nL/sec	TIP10XV119	Micropipettes for Nanoliter Injector (10)
Fast	46 nL/sec	SYS-MICRO4	Micro4 Controller, Four-Channel
ILL SPEED		300033	Adapter for Micro4
Slow	23 nL/sec	300521	Replacement O-rings (five)
Fast	46 nL/sec	500778	Replacement Nanoliter Injector Universal Adapter
MPTY SPEED	92 nL/sec	500299	Pistons, 5-pack
/ARIABLE VOLUME RANGE	2.3 - 69.0 nL		

2.3 nl

Set switch

3 lb. (1.1 kg)



Sub-microliter injection system for small animal research



- The world's smallest dead volume injection syringe
- Comes with various needle sizes from 26 ga. to 36 ga.
 - Versatile research applications RPE and IO Kits
- Custom needle shapes available blunt, sharp, beveled
- Compatible with WPI's UMP3 and PV800 series microinjection systems

NanoFil is a specially designed 10 microliter syringe developed in response to customer requests for improved microinjection in mice and other small animals. It makes quantitative nanoliter injection much easier and more accurate than any other method currently in use.

NanoFil's low dead volume eliminates the need for oil backfilling, a messy process which risks contamination of the injected sample.Injection is now simpler, and less messy, and there is no possibility of oil contamination in critical applications such as ophthalmology research (see the Retinal Pigment Epithelial (RPE) and Intra Ocular (IO) injection kits listed below).

When the inner tip diameter of a conventional syringe is reduced to less than 100 micron, it is very difficult to backfill the solution at a reasonable speed. NanoFil solves this problem by using a tip coupling mechanism that makes it possible to *change the syringe tip during the experiment*. Simply load the sample using a larger tip, such as the 26 gauge needle

provided with the syringe, and then replace it with a micro tip for sample injection. On a conventional 10 microliter syringe, a solid ring or bushing is permanently bonded to the tubing. Replacing the tip in middle of the experiment is not practical. With the NanoFil, tips can be exchanged by a simple twist of the brass lock, gently pulling out the tip, and replacing with the desired new tip. To secure the tip, NanoFil uses an olive shaped silicon gasket that is similar to, but much sturdier than, some of the microelectrode holders used for electro physiology recording. The silicone gasket makes it possible to hold not only metal tips but also glass and quartz tubing. Many types of tubing

can be easily connected to the syringe as long as the outer diameter (OD) is close to, but not more than, the inner diameter (ID) of the inside barrel. Flexible quartz capillaries used in Gas Chromotography (GC) and Capillary Electrophoresis (CE) can also be easily coupled to the syringe.

Specially designed tips as small as 36 gauge (110 micron OD) are offered in both blunt and beveled styles. Our studies have shown that these tips will cause less trauma to the tissue than any other form of micro syringe currently in use. NanoFil has a unique coupling mechanism that allows many different forms of small tubing and tips to be coupled with the syringe barrel.

NANOFIL NanoFil Syringe, 10 microliter
NANOFIL-100 NanoFil Syringe, 100 microliter

NanoFil syringe does not contain any injection tips, those must be purchased separately. It does include a 26 gauge beveled needle for backfilling.

REPLACEMENT BACKFILL NEEDLES

NF26BV-2 26G Beveled Needle, 460 μm nominal diameter (package of 2)

Using NanoFil in different configurations

Direct injection by hand: This is the simplest and most economical way to inject. Any of our tips can be inserted directly into the NanoFil syringe. Even the SilFlex tubing can be inserted to switch from hand injection to the other methods listed below. The limitation of this method is the difficulty achieving sub microliter resolution.

Installed on WPI's UMP-III microsyringe pump: This will allow the user to achieve nanoliter resolution and reproducibility. For neural system injection, mount the UMPII on a stereotaxic frame.

SilFlex tubing and holder: The needle is mounted on a small plastic holder that is connected to the NanoFil by a 35 cm length of flexible tubing. The NanoFil is mounted on the UMP II pump. This configuration allows the user to hold the animal in one hand and insert the needle with the other. When the needle reaches the desired location, activate the pump using the footswitch and the pre-programmed injection volume will be delivered. This configuration gives a nanoliter level of accuracy and reproducibility. It is best suited for applications such as the RPE and IO injection.


Selecting the correct tip for your application

The replaceable needles used with the NanoFil are available with either blunt or beveled tips. The blunt tip is used for injection into soft tissue and when a uniform solution distribution is needed. The beveled style is used for applications that involve the penetration of a tough tissue.

One of the main factors that limit the resolution and accuracy of conventional micro syringes to the upper tens of nanoliters range is diffusion in the large tip ID. When the tip ID is equal or larger than 100 micron, the error caused by tip diffusion is in the nanoliter range level $[(100 \text{ micron})^3 = 1 \text{ nanoliter}]$. With a 36 gauge needle installed on the NanoFil, the error caused by diffusion will be reduced to the subnanoliter level, making accurate injection of a nanoliter possible.

All of WPI's beveled tips have a unique 25 degree tri-surface bevel that is optimized for microinjection. A 10 degree single-surface beveled tip penetrates better than one with a 25 degree angle, however the distance between the upper opening to the tip (the dimension "F" in the drawing above) is longer. As a result, it requires a deeper penetration of the tip to achieve the same level of liquid delivery. Deeper penetration means more tissue damage. WPI's unique 25 degree beveled tip solves





this problem with two extra beveled surfaces. The tip of a single surface beveled tip is actually a blade instead of a point. It dulls very quickly. In contrast, the tri-surfaced tip has a real point. It not only penetrates much better but is also much more durable. Our tests show that our 33 gauge, 25 degree beveled tip penetrates easier and lasts longer than other manufacturers' 33 gauge, 10 degree single beveled tips. With a 35 gauge tri-surface beveled tip, the resistance to the penetration becomes even less. Each of our tips undergo a penetration test before leaving the factory to guarantee the best results for our customers.

		 •
Avai	Ish	Inc
Avai	Iav	IDS

33 gauge: This tip is similar to Hamilton's 7762 and 7803 series removable needles in both tip length and outer diameter. However, our beveled tip version is shorter, more durable, and penetrates better due to the special tri-surface grinding technique. In the past, 33 gauge tips were the smallest size sold by other manufacturers and were frequently cited in literature. However, our new 35 gauge tip is much better for injections involving small animals, especially mice. Compared with Hamilton's 33 gauge, 10 degree beveled tip, our 35 gauge 25 degree beveled tip can reduce the depth of penetration by almost 80%. The distance between the tip and the upper rim of the opening (dimension F on the drawing) is 348 microns for the 33 gauge tip. The distance for our 35 gauge tip is only 230 microns. In addition, the smaller tip size significantly reduces the required penetration force. In nearly all applications, a 33 gauge tip can be replaced with our 35 gauge tip and produce better results.

34 gauge: This is a transitional size between the 33 gauge and 35 gauge. If the 35 gauge is too weak and the 33 gauge is too large, this makes a good alternative.

35 gauge: This was the most popular and preferred tip of most scientists during our field trial. The combination of its strength, length, durability, and clogging resistance creates a balance with very little compromising of the individual properties. It is much smaller than the 33 gauge tip offered by other

Tip Order Number	Tip O.D. "A"	Tip I.D. "B"	Tip Length "C"	Total Length "D"	Shank O.D. "E"	Bevel Length "F"	Total Dead Volume	Tip Material
NF33BV-2	210 µm	115 µm	10 mm	40 mm	460 µm	≈348 µm	0.416 μL	Stainless Steel
NF34BV-2	185 µm	85 µm	5 mm	35 mm	460 µm	≈290 µm	0.199 µL	Stainless Steel
NF35BV-2	135 µm	55 µm	5 mm	35 mm	460 µm	≈204 µm	0.435 μL	Stainless Steel
NF36BV-2	110 µm	35 µm	3 mm	33 mm	460 µm	≈156 µm	0.340 µL	Stainless Steel
NFQ34-5	160 µm	100 µm	55 mm	75 mm	460 µm	n/a	0.589 µL	Quartz
NF33BL-2	210 µm	115 µm	10 mm	34 mm	460 µm	≈348 µm	0.416 µL	Stainless Steel
NF34BL-2	185 µm	85 µm	5 mm	29 mm	460 µm	≈290 µm	0.199 μL	Stainless Steel
NF35BL-2	135 µm	55 µm	5 mm	29 mm	460 µm	≈204 µm	0.435 µL	Stainless Steel
NF36BL-2	110 µm	35 µm	3 mm	27 mm	460 µm	≈156 µm	0.340 μL	Stainless Steel
Silflex		100 µm		35 cm			2.749 µL	
NF26BV-2	460 µm	110 µm	3 mm	40 mm	460 µm		0.380 µL	

manufacturers. It is only slightly larger than the 36 gauge tip but is much stronger and less likely to be clogged. Samples can be directly loaded with this tip. Its 5 mm length is sufficient enough for almost all injection applications in mice.

36 gauge: This is the smallest tip that is commercially available. The tip is so small that it can be inserted into the opening of the 33 gauge needle tip. Because this is pushing the limits of what current technology can produce, there are some limitations to consider before using. Its thin diameter makes it necessary to limit its length to 2.5 to 3 mm in order to maintain a usable strength. Since the tip ID is in the 25 to 50 micron range, it is very easily clogged. Therefore, only well filtered solutions can be used. Depending on the viscosity of the

sample, the user might also need to pre-load the syringe with a regular tip before switching to this tip for injection. We recommend using the 35 gauge tip instead of the 36 gauge unless it is absolutely necessary.

Flexible Quartz Tubing: The flexible quartz tubing tip is made of 160 micron OD polyimide coated quartz tubing with a special adapter sleeve mounted at the end. It is designed for filling glass capillary electrodes or pipettes, just like WPI's traditional MF34G Microfil. However, unlike the traditional MicroFil, which has about 50 microliters of dead volume in its luer hub, the dead volume of this tip is less than 0.4 microliters. It is useful for loading electrodes with solutions that have a limited volume or are too expensive to waste.

US: Tel: 941-371-1003 • sales@wpiinc.com

NANOFIL NEEDLES

NF33BL-2	33 G blunt NanoFil needle (pkg of 2)
NF34BL-2	34 G blunt NanoFil needle (pkg of 2)
NF35BL-2	35 G blunt NanoFil needle (pkg of 2)
NF36BL-2	36 G blunt NanoFil needle (pkg of 2)
NF33BV-2	33 G beveled NanoFil needle (pkg of 2)
NF34BV-2	34 G beveled NanoFil needle (pkg of 2)
NF35BV-2	35 G beveled NanoFil needle (pkg of 2)
NF36BV-2	36 G beveled NanoFil needle (pkg of 2)
NF33-36BL	Assortment of 4 blunt NanoFil needles
NF33-36BV	Assortment of 4 beveled NanoFil needles

REPLACEMENT PARTS & ACCESSORIES

NFINHLD	NanoFil Injection Holder
SILFLEX-2	SilFlex tubing 35 cm long (pkg of 2) (dead volume = 2.74 μL)
NFGSK-5	Spare Silicone Gasket for NanoFil & Holder (pkg of 5)
NFQ34-5	34 Gauge Flexible Quartz Tubing for filling (pkg 5)





Holder assembly for RPE-KIT.

The optional UMP3 stand in the photo below includes the small base (503084), openside clamp (14073-4) and 25cm rod (503070).

These kits are specially designed for eye research for injecting retinal pigment epithelium (RPE) and intraocular (IO) in addition to brain injection in mice. They need to be used with a NanoFil syringe and UMP3 to achieve accurate, repetitive, and oil free injection in the submicroliter range. Each kit includes two pieces of Silflex tubing (one for a spare), a holder assembly, spare gaskets, and an assortment of four tips — blunt for the RPE kit and beveled tips for the IO kit. Each kit comes with one each of 33, 34, 35 and 36 gauge tips so that first time users can find the best size for their application.

The Silflex tubing is the most critical component of the kit. This 35 cm long, flexible tubing has a very precise outer diameter for airtight fitting with the syringe. It also has a small inner diameter for minimum dead volume, and is very durable when handled correctly. The SilFlex

is coupled to the injection tip with a mechanism similar to that of the NanoFil. The dead volume of the entire kit (including the tubing) is less than 3 microliters. All of the components in the kit are constructed of inert, solvent resistant materials for easy cleaning after viral injection.

RECOMMENDED ACCESSORIES


RPE-KIT	Retinal Pigment Epithelium (RPE) injection kit (SilFlex
	tubing, gasket, holder, and blunt tip mix)
IO-KIT	Intraocular (IO) injection kit (SilFlex tubing, holder,
	gasket, and beveled tipmix)
503207	Stand & Clamps



Manual Microsyringe Pump

The MMP and DMP are convenient tools for precise manual injection of fluid using glass pipettes or similar injection devices. The design allows visual feedback of flow at the pipette tip. They can also be used as a manual micro syringe pump for perfusion or withdrawal of liquids. The resolution of the injection volume can be continuously varied from 10 nanoliters to the microliter range, depending on the syringe used. Either oil or air can be used as the transfer media to assist the injection of fluid. The DMP comes with an exclusive digital micrometer that will allow the reading of piston advancement easily with a 0.001-millimeter resolution. The optional software and cable kit can transmit advancement data directly into computer. Model MMP has the

traditional mechanical micrometer head with a resolution of 10 microns per division and advances 500 micrometers per revolution. The entire frame body of the injector is constructed with polished stainless steel for excellent stability



DMP & MMP SPECIFICATIONS

Travel Distance 25 mm

Advances Resolution 0.001 mm for DMP and 0.01 mm for MMP 10 µL to 1 mL gas tight luer tip syringe Syringe Size 1.5 m of PTFE tubing with 0.5 mm ID Tubing

0.24" x 5.2" Pipette Holder

Pipette Holder Fits 1.0 to 1.5 mm OD pipette and durability. The piston of the micrometer can be slid across the rail to the syringe's plunger position. Small diameter PTFE tubing is used to improve the accuracy and solution compatibility. The unique design of the pipette holder can securely hold any pipette with an outer diameter of between 1.0 mm and 1.5 mm. All necessary accessories for removing air and filling the syringe and tubing with liquid are included. The system comes complete with a 100 µL gas tight syringe and other syringe sizes can be purchased.

MMP	Manual Microsyringe Pump
DMP	Manual Microsyringe Pumn with Digital Dist

ACCESSORIES

MMP-KIT **Injection Assembly Parts Kit**

Not including valve—see #14057-10, page 97)

Pressure Manometer

For measuring hydrostatic pressures



Hand-held and battery operated, PM Series pressure manometers monitor vacuum and pressure in nonaqueous fluids. An integral transducer and digital display allow easy and accurate pressure readings.

Three versions measure pressures in the range of ±1 psi, ±15 psi or ±100 psi.

allows measurement in units of psi or kPa for the 100 psi

version, and psi or mmHg for the 15 psi version. Pressure can be read on the built-in LCD display or relayed to a chart recorder, oscilloscope, or

PM Series pressure manometers come with 4 feet of 1/8-inch ID soft vinyl tubing. A mini-phone-to-BNC cable for the recorder output is also available (Part #CBL102). Standard versions are equipped with a nine-volt alkaline battery.

PRESSURE RANGES ±1 psi (±52 mm Hg) MAX. PRESSURE 20 psi (1035 mm Hg) RESOLUTION 0.001 psi (0.1 mm Hg) OUTPUT 1 V/psi **OUTPUT RANGE** ±1.0 V LINEARITY 0.5% full-scale

TEMPERATURE EFFECT **ZERO RESPONSE TIME POWER BATTERY LIFE** RECORDER OUTPUT **OUTPUT IMPEDANCE**

PNEUMATIC CONNECTORS **DIMENSIONS** SHIPPING WEIGHT

SYS-PM01D Pressure Manometer (1 psi) SYS-PM01R Pressure Manometer (1 psi), Rechargeable* Pressure Manometer (15 psi) SYS-PM015D SYS-PM015R Pressure Manometer (15 psi), Rechargeable* SYS-PM100D Pressure Manometer (100 psi) SYS-PM100R Pressure Manometer (100 psi), Rechargeable* **CBL102** Mini-Phone-to-BNC Cable

Specify line voltage

*Rechargeable versions come with nickel/cadmium battery and charger

PRESSURE MANOMETER SPECIFICATIONS

PM100

±100 psi (±690 kPa) ±15 psi (±775 mm Hg) 30 psi (1550 mm Hg) 150 psi (1035 kPa) 0.01 psi (1 mm Hg) 0.1 psi (1 kPa) 100 mV/psi 10 mV/psi ±1.5 V ±1.0 V

1.0% full-scale (0-70°C) Screwdriver-adjust 30 ms

Nine-volt battery

Alkaline, 200 hours; rechargeable, 25 hours Mini-phone jack, 0.141 inch (3.5 mm)

Barbed, for 1/8-inch or 3/16-inch ID soft tubing

US: Tel: 941-371-1003 • sales@wpiinc.com

 $3 \times 6 \times 1$ inches $(8 \times 15 \times 4 \text{ cm})$

3 lb (1.4 kg)

Microfluidic Pumps

ExiGo is a precision syringe pump based on the 5-phase stepper motor drive that has more microsteps per revolution of the lead screw vs. standard syringe pumps on the market. ExiGo has 250,000 microsteps/ revolution and even at low rotational speed / low flow rates, it has a very low pulsation and high accuracy. When coupled with the flow sensor and active PID feedback; this results in very fast response times for changing flow rates. A standard syringe pump typically has a smaller number of microsteps and so usually the only way a standard "microfluidic" pump can achieve pulse-free flow control is to use small syringes; e.g. 0.5µL; 1µL; 5µL etc. to achieve non-pulsatile stable flow rates in the nanolitre/ minute range. By comparison, the ExiGo pump with the flow sensor can use a standard 250µL glass syringe to produce stable non-pulsatile flow rates of 10nL/min – 1ml/min; or a 5mL plastic syringe to produce stable non-pulsatile flow rates of 100nL/min - 20mL/min. ExiGo can be used in conjunction with expandable (flexible) element and fluidic resistance in order to dampen any pulsation occurring during the stepper motor operation. As it employs active feedback, the response time of the pumps still remains fast.

Mirus is a precision syringe pump, which uses the combination of the expandable (flexible) element and fluidic resistance in order to dampen the pulsation of the syringe pump stepping. The main difference is that the ExiGo has feedback control; the Mirus does not. Because the Mirus does not have feedback control, it also does not have a fast reaction time and so takes a long time to change flow rates. Conversely due to the expandable element, the Mirus has a very stable flow profile.

Additionally, Mirus is equipped with 3-way valve allowing automatic recharging and washout of syringe. Mirus is also provided with an 8-way flow splitter, allowing multiple executions(8 parallel experiments) simultaneously.

Kima pump is a microfluidic pump designed to aid cell culture (e.g. endothelial cells) under physiological conditions (shear flow) in various biochips and flow chambers, including Cellix's Vena8 Endothelial+ biochips where it is possible to culture 8 cell monolayers simultaneously over 24–48 hours. The pump is a solenoid pump and so it is not pulse-



CX-MIRUS-P	Microfluidic syringe pump for shear flow studies. Includes: tubing kit, VenaFlux Assay Software PC control and MultiFlow8
CX-M1000	Microfluidic syringe pump for shear flow studies. Includes: tubing kit, VenaFlux Assay Software PC control and MultiFlow8. Not included: MultiFlow8.
CX-M8000	MultiFlow8 attachment for Mirus Evo Nanopump: for higher throughput enabling 8 assays in parallel in the microfluidic biochips.



CX-EXIGO-PFS	Microfluidic syringe pump with precise flow control with active feedback via integrated flow sensor.
	Flow rate: 50 nL/min - 10 mL/min ±0.5%
	Includes: iPad mini which can control/program up to 4 pump modules independently
CX-E1000	Microfluidic syringe pump with precise flow control. For customers who have bought ExiGo Pro 1.0 and already have iPad mini.
	DOES NOT INCLUDE IPAD MINI.
CX-EXIGO-P	Microfluidic syringe pump with precise flow control. Includes: iPad mini which can control/ program up to 4 pump modules independently
CX-E2000	Microfluidic syringe pump with precise flow control with active feedback via integrated flow sensor.
	Pump only — for customers who have bought ExiGo Pro 1.0 and already have iPad mini.



CX-KIMA-P	Microfluidic pump for continuous cell culture under shear flow mimicking physiological flow in the human vasculature. Includes: iPod Touch and controller.
CX-K1000	Microfluidic pump for continuous cell culture under shear flow mimicking physiological flow in the human vasculature. KIMA pump only.

US: Tel: 941-371-1003 • sales@wpiinc.com

			Biochips				
	Vena8 Fluoro+	Vena8 Endothelial+	VenaT4	Vena8 Glass Coverslip		Vena Delta	
				Low Flow Rates	High Flow Rates	Y1	Y2
Dimensions			Hillian Comment				
Channel width, b (cm):	0.04	0.08	0.08	0.16	0.08	0.008	0.008
Channel height, h (cm):	0.01	0.012	0.01	0.016	0.008	0.012	0.012
Channel length, L (cm):	2.8	2.8	2.8	2.8	2.8	2.8	2.8
Channel volume (cm³):	0.00112	0.00269	0.00224	0.00717	0.00179	0.00269	0.00269
Channel volume (µL):	1.12	2.69	2.24	7.17	1.79	2.69	2.69
Thickness of bottom substrate (mm)	0.17	0.5		0.17	0.17	0.5	0.5
# of channels / chip	8	8	4	8	8	4	4
# of assays / pack	40/80	40/80	40	80	80	40	40
Pack of 5	CX-002 \$768	CX-004 \$768	10	30	30	.0	
Pack of 10	CX-001 \$1532	CX-003 \$1532	CX-005 §1532	CX-009 \$1532	CX-010 \$1532	CX-007 \$1532	CX-008 \$1532
Specifications							
Brightfield, phase contrast & immunostaining	1	✓	✓	✓	✓	✓	✓
Confocal microscopy	1			1	1		
Cell types: whole blood (human, animal); PRP; platelets; PBMC; monocytes; T-cells (primary 8 cell lines); eosinophils; neutrophils etc.	1	✓	1	1	1	✓	✓
Cell types: adherent cells e.g. endothelial – HUVECs; HMVECs etc.; HepG2; stem cells; muscle cells etc.	1	1	1	1	1	1	✓
Protein coatings: collagen, fibronectin, fibrinogen, vWF, VCAM, ICAM, selectins, MadCAM etc.	✓	√	✓	✓	✓	✓	✓
Applications							
Platelet adhesion, aggregation 8 thrombi formation; Leukocyte rolling, adhesion 8 migration; Thrombosis; Immunology (Inflammation); Infectious diseases (e.g. Malaria); Sickle cell disease; Respiratory (Asthma & COPD)	✓			√	√		
Cell adhesion and culture under perfusion / shear flow; Leukocyte cell-cell rolling, adhesion & migration; Oncology (Melanoma, Breast Cancer etc.); Cardiovascular (Atherosclerosis, Drug eluting stents); Immunology (Inflammation); Respiratory (Asthma & COPD)		✓					
Biofilm assays, microbe seeding & culturing; Biochips with glass coverslips (attached / not attached; treated / non-treated); Biochips for the attachment of coupons for biofilm studies				1	1		
Chemotaxis, transmigration and invasion assays; 2D and 3D cell culture; Mimicking tumour microenvironment with gels (ECM gel, hydrogel, matrigel, collagen gel)			1				

All biochips are:

- Disposable plastic; some with glass coverslips.
- \bullet Require no assembly; unlike many standard perfusion chambers / flow chambers.
- Require no luer lock connections which increase dead volume. Cellix's biochips have a unique plug & play connection with tubing connections which are autoclaveable and reuseable.

US: Tel: 941-371-1003 • sales@wpiinc.com

Optical Detection Systems




Bio Photometric Detection System

Choose the wavelength data you want to see. In many nutrient, water purity and process applications, full spectrum analysis is not required. With LEDspec you can eliminate the extraneous data and focus on those wavelengths you need to see. You can conduct flow analysis and single-scan applications with high precision over a large dynamic range.

- Affordable spectroscopy
- Dynamic range and baseline noise outperform CCD and photodiode array based spectrometers
- LED light source eliminates costly replacement lamps
- Integrated reference channel eliminates signal drift
- Full computer control
- Integrated math functions allow for baseline correction at a second wavelength, signal ratio and
- Simplified display with just the data you want to view
- Your choice of three wavelengths included

Two models are available: 2-channel and 4-channel. LEDspec-2 (2-channel) comes with your choice of three LED modules (wavelengths). LEDspec-4 (4-channel) also includes your choice of three LED modules, however, up to four additional wavelengths can be added, if desired.



DNA Calibration Curve using WPI's **V-Vette** combined with a **LEDspecUV** and pharmaceutical compliant spectrophotometer.

spec

- Measures visible wavelengths
- Sample cells: LWCC, Fiber Optic Cuvette Holders, V-Vette
- Wavelength range (nm): 400, 450, 540, 560, 600, 650, 700, custom
- Applications include:

Environmental/Oceanography

Ammonia at 650nm Iron at 560nm Nitrite/Nitrate at 540nm Phosphate at 700nm

Pharmaceutical

Process Control

Semiconductors

Water purity, trace metal analysis (Fe, Pd, Cu, U)

Data You Want to See

Many **biochemistry** applications require information at specific, important wavelengths, instead of a full spectrum analysis. For example, the Bradford, BCA and Lowry assays for protein analysis rely on specific wavelengths.

LEDspec is ideally suited for **oceanographic**

applications such as detecting nM concentrations of nitrite/ nitrate, phosphate and iron using WPI's LWCC sample cells. Two or four independent channel FIA detection systems can be assembled using a **LEDspec-2** or **LEDspec-4**, respectively.

LEDspec is a stand-alone LED-based bio-photometric detection system designed to give you the information you want to see. Now you can conduct flow analysis and single-scan applications with **high precision** and a large **dynamic range**.

LEDspec can be equipped with **up to 7 LEDs** of different wavelengths.

Its noise (< 0.1 mAU peak to peak) and drift performance (<0.5 mAU/h) exceeds that of a CCD or photodiode array detection system **at a fraction of the cost**.

LEDspec uses dual-beams to **reduce light source drift**. Conventional single beam spectrometers notice baseline drift caused by warm up, temperature stability and bulb aging. An internal reference channel in the LEDspec corrects for baseline while you make sample measurements.

Data Collection and Analysis

Now, you can analyze output data with **LEDspec**'s easy-to-use software and export chromatographs directly to your PC (via USB) in Microsoft® Excel format. Software provides:

- Full computer control of **LEDspec**
- Continuous flow or single-shot analysis of up to four independent channels simultaneously or sequentially.
- Immediate calibration and analysis (mean and standard deviation) of up to four channels



LEDSPEC SPECIFICATIONS

OPTICAL BASICS

LED-based multiple wavelength detector with build-in reference channel

CHANNELS 2 or 4
DETECTOR Photodiode

SPECTRAL BANDWIDTH (FWHM) 10 nm (LEDs >400nm) 4 nm (260, 280, 340nm LEDs)

DYNAMIC RANGE 0-3 AU

DETECTOR RESOLUTION 24 Bit

NOISE (PEAK TO PEAK) < 0.1 mAU

WARMUP TIME Instant

FIBER OPTIC INPUT 600 µm

DRIFT < 0.5 mAU/h

DIGITAL INPUTS AND OUTPUTS 8/8

ANALOG OUTPUT +/- 10 V, scaleable output

DIMENSIONS (W*H*D) 290 x 80 x 250 mm
(11.4" x 3.2" x 9.9")

WEIGHT 2 kg (2.2 lbs)
INTERFACE USB 2.0

MAINS 100 – 240 V / 50 - 60 Hz



INTRODUCTORY PRICE

LEDSpec-2	LEDspec biophotometric detection system, 2 channel, 3 VIS LED modules (choose when ordering)
LEDSpec-4	LEDspec biophotometric detection system, 4 channel, 3 VIS LED modules (choose when ordering)
89273	UV LED module, 260 nm
89272	UV LED module, 280 nm
89274	UV LED module, 340 nm
89245	VIS LED module, 400 nm
89246	VIS LED module, 450 nm
89247	VIS LED module, 540 nm
89248	VIS LED module, 560 nm
89275	VIS LED module, 600 nm
89276	VIS LED module, 650 nm
89249	VIS LED module, 700 nm
PERIPRO-4LS	Peri-Star Pro, 4-channel, low rate, small tubing (see page 162)
MINISTAR	Miniature Peristaltic Pump, 1-channel (see page 164)

Optical Detection Systems UltraPath A unique multiple long pathlength sample cell for absorbance spectroscopy 18.14

Process Control & OceanographyRugged system for laboratory and onboard measuring • Portable & easy to use • User-selected optical path lengths: 2, 10, 50 & 200 cm • Highly sensitive and stable

UltraPath™ is a unique high-performance spectrophotometer system offering user-selectable optical path lengths of 2, 10, 50 and 200 cm. The instrument operates in the wavelength range of 250 to 730 (UPUV) or 380 to 730 nm (UPVIS) and has an exceptional dynamic range. Designed for the detection of low absorbing species in aqueous solutions, UltraPath is an ideal tool for any study requiring precise and highly sensitive spectroscopic determination of analytes, either in the lab or in the field.

Appearance of instruments presently in these systems differ from those pictured here

Background

UltraPath was developed by WPI under a collaborative agreement with NASA (Stennis Space Center) for the spectroscopic determination of colored dissolved organic matter (CDOM) in seawater and fresh water environments. It can be used in the laboratory and in the field (i.e., at sea). CDOM concentrations vary significantly between open ocean samples with low CDOM (e.g., 0.007 m⁻¹ at 380 nm), and high CDOM freshwater environments (e.g., 10-20 m⁻¹ at 380 nm). To address these problems the design requirements of UltraPath mandated the development of a rugged portable system capable of high sensitivity measurements across a wide dynamic range. The UltraPath system meets these stringent design criteria and enables reliable measurement of CDOM in the range of 0.002 m⁻¹ to 200 m⁻¹ (250 to 730 nm).

Design

UltraPath has four optical pathlengths contained within a single sample cell (i.e., 2 cm, 10 cm, 50 cm and 200 cm). The pathlengths are userselectable, offering a very high sensitivity and an extended dynamic range for UV and VIS absorbance measurements. The fluid path of the sample cell is optimized to produce a laminar flow that is virtually free of interference from trapped air bubbles and adherence of dissolved

substances to the cell wall. In particular, the design greatly minimizes the problems commonly found with flow cells of long optical pathlengths: the risk of trapping dust particles, fibers or particulate matter inside the cell. The UltraPath system includes a low noise photodiode array-based spectrometer module (TIDAS I: FWHM = 5 nm, noise <0.2 mAU) and a light source (D4H with UPUV; FO6000 with UPVIS) to measure sample absorption. Light is coupled from the light source to the sample cell and from the sample cell to the detector via two fused silica fibers. A peristaltic pump (PeriPro-4LS) is utilized to draw the sample into the UltraPath sample cell.




Fig. 1 — Two typical absorption spectra measured using UltraPath. The sample labeled "Sarasota Bay" is a CDOM sample with 34 PSU salinity collected from Sarasota Bay (Nov. 2007), and the sample labeled "Pond" is a highly concentrated CDOM sample collected from a local pond in Sarasota, Florida (Nov. 2007).




Fig. 2 — CDOM Sample "Mayagüez Bay" was collected from the high salinity oligotrophic waters of Mayagüez Bay on the west coast of Puerto Rico (2001). Data courtesy of NASA Stennis Space Center.

A standard PC or laptop (not included) is connected to the TIDAS E via a RJ-45 Ethernet interface. For spectrometer requirements and software options, see **TIDAS-E**.

Mobility

The system is designed for mobility. The components of the UltraPath system are designed to function over a broad range of laboratory and field environments.

Samples

Two typical absorption spectra recorded with an UltraPath (UPUV) of a seawater and a fresh water sample collected in November 2007 are shown in Fig. 1. Due to their high absorbance, both samples were analyzed in the 10 cm pathlength. The CDOM sample labeled Mayagüez Bay in Fig. 2 is from oligotrophic, low productive waters with high salinity collected off the west coast of Puerto Rico in the Mayagüez Bay. Special attention should be drawn to the exceptional sensitivity of UltraPath enabling detection of CDOM absorption below 0.03 m-1. To exemplify the performance of the UltraPath in Laboratory Chemistry and Process Control, Ponceau S absorbance was measured with the 200 cm pathlength of an UltraPath. Normalizing the Ponceau absorbance graph to AU/cm, the range of this measurement is 150 µAU with a noise level below 2 µAU peak to peak. Sub-nanomolar concentration of this dye can clearly and reliably be detected, which is a novelty in absorbance based spectroscopy.

Particulate Absorption

Particulate absorption can be measured by the well established Quantitative Filter Technique (QFT). WPI now offers a fiber optic filter holder for Glass Fiber Filters (QFT1, page 206) which can be used with the spectrometer (TIDAS E) and light source (D4H or FO6000) supplied with the UltraPath. With this accessory, particulate absorption can be measured on site, avoiding loss of spectral information due to freezing and shipping particulate samples to a laboratory.

Reference

N. B. Nelson, D. A. Siegel, C. A. Carlson, C. Swan, W. M. Smethie Jr. and S. Khatiwala. 2007. Hydrography of chromophoric dissolved organic matter in the North Atlantic. Deep-Sea Res. I. 54: 710 – 731.

V. Kitidis, A. P. Stubbins, G. Uher, R. C. Upstill Goddard, C. S. Law, E. M. S. Woodward, "Variability of chromophoric organic matter in surface waters of the Atlantic Ocean", Deep Sea Research Part II: Topical Studies, Vol. 53, Issue 14-16, 2006, p. 1666-1684.

R. L. Miller, M. Belz, C. Del Castillo, R. Trzaska,

"Determining CDOM Absorption Spectra in Diverse Coastal Environments Using a Multiple Pathlength, Liquid Core Waveguide System", Continental Shelf Research, July 2002, 22:9, p 1301-1310.

"System Analyzes Water Samples at Sea", NASA Aerospace Technology Innovation, 2001, 9 (5). http://nctn.hq.nasa.gov/innovation/innovation95/3-techtrans2.html

R. L. Miller and E. D'Sa. "Evaluating the influence of CDOM on the remote sensing signal in the Mississippi River Bight". In Eos Transactions AGU Ocean Sciences, 2002. Honolulu, HI, p. 171.

E. D'Sa, R.L. Miller and R. Trzaska. "Aparent Optical Properties in Waters Influenced by the Mississippi River", Proceedings of the Seventh Thematic Conference, Remote Sensing for Marine and Coastal Environments, 2002, 6 pg, Miami, FL.

R. L. Miller, C. Hall, C. Del Castillo, B. McKee and M. Dagg. "Bio-optical Properties of the Mississippi River Plume and Adjacent Shelf." ASLO Aquatic Sciences, Albuquerque, NM, 2001.

R. L. Miller, M. Belz and S. Y. Liu, "Measuring the absorption of CDOM in the field using a multiple pathlength liquid waveguide system", Ocean Optics XV, paper 1308, Monaco, October 2000.




Fig. 3 – Ponceau S absorption measured with UltraPath (200 cm cell). Ponceau S was dissolved in Millipore water.

UltraPath Specifications

DYNAMIC RANGE 5 µAU/cm to 1 AU/cm 0.002 m⁻¹ to 200 m⁻¹

WAVELENGTH RANGE 250 nm - 730 nm (UPUV) 380 nm - 730 nm (UPVIS)

WAVELENGTH RESOLUTION (FWHM) 5 nm

 NOISE (PEAK TO PEAK)
 < 0.2 mAU</td>

 DRIFT
 < 1 mAU/h</td>

 OPTICAL PATHLENGTH
 2, 10, 50 & 200 cm (user selectable)

SAMPLE CELL INNER DIAMETER 2 mm

AMPLE CELL INNER DIAMETER 2 MM

CELL VOLUME 10 mL (at 200 cm pathlength)

SAMPLE INLET / OUTLET 1/8"
FIBER INPUT/OUTPUT 600 µm

SOLVENT RESISTANCE Most organic and inorganic solvents

SHIPPING WEIGHT UPUV: 44 lb (20 kg)
LIPVIS: 33 lb (15 kg)

UPVIS: 33 lb (15 kg)

UPVIS Ultrapath System, Visible Light UPUV Ultrapath System, Ultraviolet & Visible Light

The UltraPath system includes: Multiple pathlength cell, Tidas E with TidasDAQ/SpectraView software, FO-6000 light source (UPVIS) or D2H light source (UPUV), two FO-600-SMA1M optical fibers, PeriStar Pro peristaltic pump, silicone tubing, sample injector and Waveguide Cleaning Kit.

Specify line voltage

501609	Waveguide Cleaning Kit
KIT-UPVIS-STARTUP	FO-600-SMA1M, 501609, 72100, 800120, 15807
KIT-UPUV-STARTUP	FO-600-SMA1M, 501609, 72100, D2H-DB, D2H-HB, 15807
89575	QFT1, Fiber Optic Holder for Glass Fiber Filters

Optical Detectors



- Photodiode array spectrometer module
- Low noise detection (<0.1 mAU peak to peak)
- Wavelength range 190 nm to 720 nm
- Fiber optic design



Tidas E

High performance fiber optic spectrometer systems

WPI's **Tidas E** is a high end photodiode array-based spectrophotometer including a light source, a cuvette holder and fiber optic coupling. Its companion is the **TIDAS E Base**, a fiber optic spectrometer module with a fiber optic connector for modular spectrometer systems. The Tidas family of spectrophotometers and spectrometer modules outperforms conventional bench-based spectrophotometers and CCDbased spectrometer modules, when it comes to high precision fiber

optic sampling. It relies on a monolithic optical bench made by Zeiss, which is optimized for fiber optic applications. Most cuvette-based standard spectrometers lose more than 90% of light through expensive prism decoupling. The Tidas E is designed for fiber optic sampling cells. Using suitable light sources and sample cells, spectral detection in the wavelength range of 190 to 720 nm can be performed at noise levels <0.04 mAU peak to peak.

Applications

The Tidas E is ideally suited for WPI's fiber optic sampling equipment. High sensitivity detection systems for flow analysis can be assembled using WPI's Liquid Waveguide Capillary Cells (LWCC) with effective pathlengths ranging from 50 to 500 cm. These setups are frequently used in fluid injection analysis systems for nutrient analysis (nitrite, nitrate, phosphate, iron) in oceanographic applications. Microliter sampling systems for UV/VIS applications can be assembled using WPI's V-Vette or DipTip™ dipping probes.

Software

TidasDAQ 3 software is included with each instrument for data collection and data analysis. TidasDAQ is used to run the spectrometer module, collect spectra in either single or continuous mode, control the digital I/Os, save the experimental data to disk, and analyze the data. Further, TidasDAQ can export data directly into GRAMS/AI, a feature very useful for advanced data analysis for pharmaceutical applications and

Wavelength Range	Resolution	Wavelength Accuracy	TIDAS E BASE	Light Source	TIDAS E	Light Source
UV 190-390 nm	< 3 nm	±1 nm	504717	no	504720	UV
UV/VIS 190-720 nm	< 7 nm	±1 nm	504718	no	504721	UV / VIS
VIS/NIR 300-1100 nm	< 10 nm	< 3 nm	504719	no	504722	VIS
Integrated Cuvette Holder			no			yes
Number of Pixels	256					
Dimensions	260 x 150 x 1	40 mm				
Weight	2.5 kg					
Power Supply	100 - 240 VA	C / 47 - 63 Hz				
Digital I/O	2 x IN / 2 x O	UT				
Optical Fiber Connection	SMA 905					

Includes power supply, TIDAS DAQ software, RJ-45 cable, and manual




Figure 1: TIDASDAQ acquisition window, showing an absorbance baseline.




Figure 2: Spectra can be displayed and analyzed in 2D and 3D format. This allows the user to conveniently interpret "time acquisition" data typically done with a TIDAS-E-BASE-LWCC flow system.

TIDAS E SPECIFICATIONS

OPTICAL BASICS Monolithic Spectrometer Module; Concave Aberration Corrected

Holographic Grating; Fiber optic cross section converter for increased light throughput; 2nd order multilayer filter

DETECTOR ARRAY Hamamatsu photodiode array, 256 pixel

DETECTOR RESOLUTION 16 Rit

NOISE (PEAK TO PEAK)* < 0.04 mAU @ 254 nm

WAVELENGTH ACCURACY <1 nm WAVELENGTH REPRODUCIBILITY < 0.1 nm FIBER OPTIC INPUT 600 µm SYSTEM REQUIREMENTS Windows XP, 7, 8 SOFTWARE (INCLUDED) TIDASDAQ

DIMENSIONS (WxHxD) $260 \text{ mm} \times 150 \text{ mm} \times 140 \text{ mm}$ $(10.25" \times$

 $5.9" \times 5.5"$)

WEIGHT 2.5 kg (5.5 lb) INTERFACE External (RJ-45)

POWER 100 - 240 V / 50 - 60 Hz


TidasDAQ: Data Collection & Instrument Control

With TidasDAQ, high precision intensity, absorbance, transmittance or normalized spectra can be obtained in less than a second. Only a few parameters need to be adjusted to obtain spectral data. Sampling of single scans, continuous full spectra scans or triggered scans is possible. Chromatograms can be displayed and logged to disk at up to four wavelengths. Data Export of 2D and 3D Spectrograms, as well as Chromatograms is supported in ASCII, Spectralys/SpectraView, Excel and Grams/Al formats. Light sources and other sampling instrumentation can be controlled via the TTL level digital outputs, as well as data collection can be triggered by TTL leveled external inputs of the TIDAS E.


Spectra can be recorded in 2D and 3D view. Mathematical computation, Derivation, Smoothing, Quantification and other functions are available to work with your data. The Quantification module allows single point and multiple point analysis, multiple linear regression, partial least square and principle component analysis. Data can be exported out of a 3D analysis file into separate scans. Further, chromatograms as well as spectrograms can be copied directly into Excel for further data analysis.

The TIDAS E Optical Path


With help of the patented cuvette holder and the mirror cuvette, various measurement setups are possible in conjunction with either internal or external light source. The patented cuvette holder and the external measurement setup can be simultaneously used. The flexibility of the design of TIDAS E is manifested below in various possible options of measurements using different optical paths formed by combinations of cuvette holder, mirror cuvette, optical fiber and light sources.




The light radiates from the internal light source, and is conducted from the redirected mirror, through the cuvette holder, into the correct detector. The measuring sample can be located in the second cuvette chamber.



The light from the external light source is coupled over the cuvette holder and redirected over the mirror cuvette to the detector; thereby the sample in the standard cuvette chamber can be measured.



From the internal light source, radiated light passes through the cuvette holder and illuminates through the external measuring cell. The light is redirected through the cuvette holder and channelled through the detector.



The light from the external light source is coupled over the cuvette holder and redirected over the mirror cuvette to the detector.

Biofluorometer



The new SI-BF100 is an LED-based fluorometer for life science applications. It is ideally suited for ratiometric calcium detection (FURA-2) and ATPase detection (via NADH fluorescence). With up to seven LED modules (wavelengths), the SI-BF100 covers many fluorometric applications in neuroscience and cell biology. Recent advancements in optics and LED technology simplify ratiometric calcium imaging, making this equipment more affordable. A breakthrough in WPI patented technology allows the SI-BF100 to use wavelengths below 380nm and produce more light in those spectra. This technology significantly cuts the cost of photometric calcium imaging without sacrificing resolution or quality.

LED light sources require less power, give off less heat and are more compact and affordable

Sampling rates up to 1kHz (1,000 ratios/ second maximum). At lower speeds, signal averaging is used for noise reduction.

Two auto ranging photomultiplier inputs allow you to monitor multiple wavelengths from a single emission output that can be comprised of any wavelength of light for which an LED module is available

Using a separate reference channel, ultrastable, continuous ratio calculations automatically compensate for LED intensity drift. This ensures less noise and produces more accurate measurements.

Application-specific probes are available for existing tissue baths and cuvette systems.

Ratio noise is <0.05 peak to peak, drift is less than 0.1 unit/hour

The warm up time of less than one minute is a dramatic improvement over the common 20-60 minutes required by xenon or mercury light sources

Replace the emission filter easily or change the LED modules to transform the SI-BF100 into a general purpose fluorometer for many other applications

How it Works

Up till now, calcium imaging systems have been required to compensate for errors and noise introduced by the complexity of their design. The systems require mechanical filters and use expensive xenon or mercury light sources. The beauty of the SI-BF100 is its simplicity. The elegance of its design reduces the noise introduced into the system and the errors inherent in traditional systems.

Monochromatic LED light sources using WPI patented technology eliminate the need for complex and expensive white light

Because the LED modules can be pulsed, sampling frequencies up to 1,000 cycles per second are possible.

sources and filter wheels

The LED light source emits specific excitation frequencies which travel through the probe. The excitation light can be comprised of any wavelength of light for which an LED module is available. The probe returns a single emission

output to one or two photomultiplier inputs on the front of the SI-BF100, which are independently filtered for specific wavelengths. This design allows you to monitor multiple wavelengths from a single emission output.

The LED light source in the WPI design makes this ratiometric fluorometer more compact, energy efficient and affordable. As added benefits, the low-power light source produces much less heat, and it warms up in less than one minute!

This incredible design is not limited to calcium imaging either. By simply replacing the 510nm emission filters Unique fiber optic coupling probe allows in front of the

photomultipliers with the desired wavelength filters,

your SI-BF100 becomes

a general purpose fluorometer for any application you can imagine. Changing a filter involves removing the two screws that hold the filter carriage on the face of the SI-BF100, swapping the filter and reinstalling the integral SMA/filter carriage.

SI-BF-100	Biofluorometer
OPTIONAL COM	PONENTS
M3301	Manual Manipulator for securing the probe
M10	Magnetic Base

for highly efficient transfer of light and

ease of placement. Custom probes fit

your existing systems.

SI-BF100 Specifications

CALL FOR APPLICATION

FIBER OPTIC LIGHT INPUT/OUTPUT SMA terminated RATIO NOISE < 0.1 peak to peak

ANALOG OUTPUT RANGE0–10V (continuous, equivalent to a ratio 0–10)

Excitation and Emission fibers in a single probe

Single and dual emission configurations available for both the small probe and the tissue probe



All probes use fibers with a 300µm core diameter. Excitation fibers have 1000 micrometer SMA connectors for Excitation and Double Emission Probes and 1500 micrometer SMA connectors for Single Emission Probes.

A smaller probe, optimized for use with small tissues like trabecula or in multi-well plates, is in development.

94650	Single Emission, Small Tissue Probe
94689	Dual Emission, Small Tissue Probe

World's Smallest Fiber Optic Dipping Probe

for UV/Vis Spectroscopy




Mini DipTip™ is a miniature transmission probe for microliter spectroscopic sampling. Mini DipTip's tip diameter is only 1.5 mm the size of a 17-gauge needle. It will fit into all micro centrifuge tubes on the market and is a useful tool for measuring protein and DNA samples. It can also be used for a dissolution system.

Microliter samples can be analyzed cost effectively when you combine the Mini DipTip with one of the following:

- The Fiber optic-based spectrometer (Tidas I) and a light source (D4H and FO-6000)
- WPI's biophotometric detection system (LEDspec)

The Mini DipTip is ideal for multi-channel applications.

Compatible with many standard spectrophotometers (600µm fiber optic coupler connections)



DIPTIP Specifications


	DIP-UV-MINI
TIP DIAMETER	1.5 mm
LIGHT PATHLENGTH	2, 5, 10mm
WAVELENGTH RANGE (nm)	200-1000
SAMPLE VOLUME REQUIRED	20-50 μL
DISTANCE FROM TIP TO UPPER EDGE OF SAMPLE WINDOW	7 mm
FIBER LENGTH	1.0 m
FIBER OPTIC CONNECTION	SMA 905
LAUNCH FIBER BUNDLE (7 \times 200 μ m)	680 µm*
RETURN FIBER BUNDLE (7 x 200μm)	680 µm*

*Circular packaging of the fiber bundle results in an active area equivalent to a fiber with a core diameter of 680 μm . Using a 600 μm connection is recommended and will result in negligible light loss.

DIP-UV-MINI-2	Mini DipTip™ for UV/VIS/NIR (2mm path)
DIP-UV-MINI-5	Mini DipTip™ for UV/VIS/NIR (5mm path)
DIP-UV-MINI-10	Mini DipTip™ for UV/VIS/NIR (10mm path)

Light Sources

Deuterium halogen light source with integrated TTL shutter






Replacement Deuterium . Lamp 503847



Replacement Halogen Lamp 503848

The **D4H** is a combined deuterium and halogen light source for UV/ VIS and NIR applications. This light source is ideally suited to work with WPI's spectrometer modules and sample cells. It supplies a continuous spectrum in the UV, VIS and NIR range from 200 nm to 1100 nm. The D4H is equipped with an integrated electrical shutter, which can be controlled by a switch or a TTL signal.



D4H	Deuterium Halogen Light Source (200 nm-1100 nm)
503848	Halogen Replacement Lamp for D4H
503847	Deuterium Replacement Lamp for D4H (> 215 nm)

LIGHT SOURCE SPECIFICATIONS

	D4H	FO-6000
APPLICATION	UV/VIS/NIR	VIS/NIR
SPECTRAL RANGE	200—1100 nm	380—1700 nm
DEUTERIUM LAMP LIFE	2000 hr	NA
TUNGSTEN/HALOGEN LAMP LIFE	2000 hr	3000* hr
STABILITY	1-2 mAU/h	<0.5 mAU/h
POWER CONSUMPTION	140 W	6 W
POWER REQUIREMENTS 1	10/240V, 50-60Hz, 1A	12VDC/1A
SHUTTER/TTL TRIGGER	Yes	Yes
MAX. FIBER OUTPUT	1000 μm	1000 µm
CONNECTIONS	SMA	SMA
SHIPPING WEIGHT	13.2 lb (6 kg)	1.3 lb (0.6 kg)
DIMENSIONS (W/H/L)	7 x 6.2 x 9.8 in.	4.8 x 2.8 x 7.5 in.
((17.8 x 15.7 x 25 cm)	(12 x 7 x 19 cm)
*I amp life is depende	ant upon internal nower	cattings

*Lamp life is dependent upon internal power settings

Replacement Lamps for Tidas II

TIDAS-D2	Replacment Deuterium Lamp	
TIDAS-H	Replacment Halogen Lamp (Type 1)	
TIDAS-H2	Replacment Halogen Lamp (Type 2)	
Replacement Lamps for D2H		


D2H-DB	Replacment Deuterium Lamp
D2H-HB	Replacment Halogen Lamp
D2H-HBER	Replacment Deuterium Lamp, Extended Range

FO-6000

High color temperature tungsten light source



This new miniature tungsten light source has been developed for high precision portable and low-power spectroscopy applications. Its advantage lies in its high light power output, its effective color temperature of 6000K and its exceptionally low drift below 0.5 mAU/h. The FO-6000 is a light source for the extended visible



part of the light spectrum (380 nm - 1700 nm). It has a SMA type output connector. Both shutter and lamp can be operated via TTL external triggering. This light source offers a wide assortment of field applications in analytical chemistry as well as environmental and life science.

A significant problem with tungsten light bulbs is their inherent low light output at wavelengths below 430 nm. The FO-6000 was developed to overcome this limitation. The light intensity of a tungsten light bulb (2760K) drops below 10% at 420 nm wavelength. However, using FO-6000, the light intensity drops below 10% at 370 nm, where the intensity of the conventional tungsten light bulb is at approximately 2% relative light output. The inherent low noise and low drift of the FO-6000

makes it particularly suitable for low-noise detection systems.



FO-6000-FILT

The FO-6000-FILT inline filter holder directly attaches to the FO-6000 light source. This allows a virtual light loss free insertion of optical filters with outer diameters from 8 to 25.4 mm and thickness ranging from 2 to 10 mm into the light path of the FO-6000. With this filter holder and an optical filter, a highly stable monochromatic light source can be assembled.

FO-6000	Fiber Optic Light Source
FO-6000FILT	Inline Filter Holder Adapter for FO-6000
800120	Replacement Lamp for FO-6000

LED-Lite[™]

Modular LED Light source with exchangeable LEDs

The **LED-lite** is a power supply for WPI's **ELS** LED modules for monochromatic light excitation. Each ELS module has an SMA bulk head fitting and allows direct attachment of SMA terminated fibers.



ELS Spectral Distributions


Color	λ_{max}	Spectral Line Half Width	Estimated Output
UV	370 nm	12 nm	85 µW
Blue	430 nm	65 nm	15 μW
Blue	450 nm	70 nm	119 µW
Blue	470 nm	20 nm	140 µW
Blue-Green	495 nm	35 nm	227 µW
Green	525 nm	40 nm	80 µW
Yellow	590 nm	13 nm	60 µW
Orange	623 nm	15 nm	114 µW
Red	660 nm	35 nm	275 μW

Estimated output is after light has passed through a 1 mm fiber.

LED-LITE ELS Power Supply (requires ELS module) Includes transformer and AC adapter. Specify line voltage

ELS-xxx	External Light Source Module (specify wavelength)
ELS-370	ELS Module (370 nm)
300051	Fiber Optic Collimator (SMA)
300052	Fiber Optic Collimator (ST)

To order ELS, use wavelength as suffix to part number (e.g. ELS-430).



Flow Cells



Liquid Waveguide Capillary Cell

Liquid Waveguide Capillary Cells (LWCC) are fiber optic cells that combine an increased optical pathlength (50–500 cm) with small sample volumes (125–1250 µL). They can be connected via optical fibers to a spectrophotometer with fiber optic capabilities. Ultra-sensitive absorbance measurements can be performed in the ultraviolet (UV), visible (VIS) and near-infrared (NIR) to detect low sample concentrations in a laboratory or process control environment. According to Beer's Law the absorbance signal is proportional to chemical concentration and light path length. Compared with a standard 1 cm cell, a 1 mAU signal is enhanced fifty-fold with a


50 cm cell to 50 mAU, using WPI's patented aqueous waveguide technology*. The LWCC can be connected directly to a pump or can

even be filled using a syringe. Based on fiber optics, the LWCC is designed for use with WPI's **LEDspec** photometric detector (see page 182).

LWCC Key Features

Pathlength, internal volume, and wavelength range (measured with ultrapure water and a Tidas II spectrophotometer

	Pathlength [cm]	Internal Volume [μL]	Wavelength Range [nm] measured with Tidas II
LWCC-3050	50	125	230-800
LWCC-3100	100	250	230-730
LWCC-3250	250	625	250-730
LWCC-3500	500	1250	280-730



These spectra show the optimal detection limits for LWCCs of varying pathlength.

LWCC specifications

WAVEGUIDE MATERIAL Fused silica tubing coated with a low refractive index polymer

OPTICAL PATHLENGTH 50-500 cm INNER DIAMETER 550 µm

INTERNAL VOLUME $\approx 125 - 1250 \,\mu L$

SAMPLE INLET/OUTLET

COMPRESSION FITTING 1/16", 1/32" SMA, ID = $400 \mu m$ FIRER INPLIT

MINIMUM PRESSURE* 1.5 - 3 PSI

SOLVENT RESISTANCE Most organic & inorganic solvents

SHIPPING WEIGHT 1.4 kg (3 lb)

*A one-meter Type II waveguide of 550 µm ID requires about 1.5 PSI for water flow of 1 mL/min.

Ultra-sensitive Microliter UV/Vis Spectroscopy with 2 to 500 cm optical pathlength!



Mote Marine Laboratories in Sarasota, Florida has partnered with WPI, using the Company's waveguide technology. WPI customized the LWCC installed in the instrument package inside these AUV's (Autonomous Underwater Vehicles). The LWCCs are used to monitor the color of seawater in the Gulf of Mexico. One of the primary interest areas of this research is red tide algae blooms.





systems can be assembled using a TIDAS E spectrometer and a UV/VIS light source such as **D4H** and **FO-6000** (see pages 190-191).

Applications

LWCCs have been used in a variety of applications, such as liquid chromatography, stopped-flow and colormetric detection, drinking water analysis, as well as environmental and oceanographic monitoring systems (see References on www.wpiinc.com). WPI's Liquid Waveguide Capillary Cells are made of fused silica tubing with an outer coating of a low refractive index polymer. This results in high signal stability and easy removal of air bubbles trapped in the sensor cell due to the hydrophilic character of the cell wall.

LWCC Injection System

For flow analysis, including simple fluid injection analysis (FIA) setups, add WPI's LWCC injection system (WPI #89372). A selection valve provides baseline or cleaning solutions to the sample stream. The injection valve injects a sample into the stream, avoiding the introduction of air bubbles or changes of flow

Related Patents

Micro Chemical Analysis Employing Flow Through Detectors, 1995, U.S. Patent No. 5,444,807. Agueous Fluid Core Waveguide, 1996, U.S. Patent No. 5,507,447.

Long Capillary Waveguide Raman Cell, 1997, U.S. Patent No. 5,604,587

Chemical Sensing Techniques Employing Liquid-Core Optical Fibers, U.S. Patent No. 6,016,372

Waveguide Cleaning Kit (#501609), above, includes the most commonly needed cleaning solutions for the LWCC waveguides. The **LWCC** Start-up Kit (#KITLWCC), at right, includes two fiber optic cables (#FO-600-SMA1M), Sample *Injector Assembly (#58006), MiniStar™ Peristaltic* Pump, and WaveGuide Cleaning Kit (#501609).



LWCC-3050	Liquid Waveguide Capillary Cell, pathlength = 50 cm
LWCC-3100	Liquid Waveguide Capillary Cell, pathlength = 100 cm
LWCC-3250	Liquid Waveguide Capillary Cell, pathlength = 250 cm
LWCC-3500	Liquid Waveguide Capillary Cell, pathlength = 500 cm
LWCC-4010	Liquid Waveguide Capillary Cell, pathlength = 10 cm
LWCC-4050	Liquid Waveguide Capillary Cell, pathlength = 50cm, 2mm ID
LWCC-4100	Liquid Waveguide Capillary Cell, pathlength = 100cm, 2mm ID
· ·	

ACCESSORIES

A sample injector assembly can be used to conveniently fill an LWCC with sample solution using a peristaltic pump. Please note that the LWCC requires two optical fibers to connect to spectrophotometer system. Choose between anti-solarized 400 micron core or UV-enhanced cables (may be ordered in 1 or 3 meter lengths).

89372	LWCC Injection System
58006	Sample Injector Attachment
PERIPRO-4LS	Peri-Star™ Pro Peristaltic Pump (see page 162)
MINISTAR	Miniature Peristaltic Pump, 1-channel (see page 164)
FO-600-SMA1M	Fiber Optic cable, 1m, SMA, 600 µm core, UV-enhanced
501609	Waveguide Cleaning Kit (available only in USA)
KITLWCC	LWCC Start-up Kit*
58450	Kit, Adapter Syringe, LWCC

*includes FO-600-SMA1M (two), 58006, MINISTAR, 501609

Optical Glass and Quartz Cuvettes for

- High Quality Glass Cuvettes
- Standard Quartz Cuvettes
- Now offering an expanded line of standard, self-masking, and flow cuvettes!

WPI's glass and synthetic quartz cuvettes are ideal for UV/VIS/NIR absorbance or fluorescence experiments.

Synthetic quartz can be used in deep UV applications and is recommended for fluorescent applications, as it does not exhibit

background fluorescence. Quartz cuvettes (absorbance, fluorescence and flow) are shipped individually

packaged, glass cuvettes are shipped in packages of 10 cuvettes. These economic quartz and glass cuvettes are ideal for precision measurements because of their high quality materials used and their low

manufacturing tolerances. Typical transmission curves of glass and synthetic quartz cuvettes are shown in Fig. 1 (cuvettes were empty, thickness 1.25 mm x 2, including surface reflections, measured with a TIDAS II against air as reference).

TECHNICAL CHARACTERISTICS

Cuvette Material

Spectral Range (>80%)

Transmission Difference Between Different Cuvettes

Optical Glass 350 - 2500 nm Synthetic Quartz 200 - 2500 nm Less than 1% Less than 1% A complete transmission spectrum from 190 nm to 4 mm is shown in Fig. 2 (cuvettes were empty, thickness 1.25 mm x 2, including surface reflections).




Fig. 1—Transmission curves of Glass and Synthetic Quartz Cuvettes



Style A



Style B



Style C




Fig. 2—Complete transmission curves of Glass and Synthetic Quartz Cuvettes



Style D




Style E



Style F

Spectrophotometry and Fluorometry

WPI PN	Style	material	Polished windows	path [mm]	Dimen- sions [mm]	volume [mL]	Beam width [mm]	Price per cell
standard red	tangular	cuvettes						
CUV2101-1*	В	Quartz	2	1	3.5x12.5x45	0.35	10	
CUV2102-1*	В	Quartz	2	2	4.5x12.5x45	0.7	10	
CUV2011-1*	В	Quartz	2	5	7.5x12.5x45	1.7	10	
CUV1022-10	С	Optical Glass	2	10	12.5x12.5x45	3.5	10	pack of 10
CUV2012-1	C	Quartz	2	10	12.5x12.5x45	3.5	10	
CUV2105-1	C	Quartz	2	20	22.5x12.5x45	7	10	
CUV2106-1	С	Quartz	2	30	32.5x12.5x45	10.5	10	
CUV2107-1	C	Quartz	2	40	42.5x12.5x45	14	10	
CUV2108-1	C	Quartz	2	50	52.5x12.5x45	17.5	10	
*89337 Cu	vette space				2023-1, CUV2063-1)			
CUV2023-1*	D	Quartz	2	5	7.5x12.5x45	0.7	4	
CUV2031-1	D	Quartz	2	10	12.5x12.5x45	1.4	4	
CUV2025-1	D	Quartz	2	20	22.5x12.5x45	2.8	4	
CUV2028-1	D	Quartz	2	50	52.5x12.5x45	7	4	
CUV2032-1	D	Quartz	2	10	12.5x12.5x45	1	3	
CUV2032-1			2	10	12.5x12.5x45	0.7	2	
	D	Quartz						
CUV2033-1	D D	Quartz Quartz	2	10	12.5x12.5x45	0.35	1	
CUV2033-1 CUV2034-1	D	•	2	10	12.5x12.5x45	0.35	1	
CUV2033-1 CUV2034-1 Self masking	D	Quartz	2	10 5	12.5x12.5x45 7.4x12.5x45	0.35	1 Ø3	
CUV2033-1 CUV2034-1 Self masking CUV2063-1*	D g continuo	Quartz ous flowthrou	2 gh cell					
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1	D g continuc E	Quartz Dus flowthrou Quartz	2 igh cell 2	5	7.4x12.5x45	0.035	Ø3	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1 CUV2065-1	D g continuo E E	Quartz Pus flowthrou Quartz Quartz	2 gh cell 2 2	5 10	7.4x12.5x45 12.5x12.5x45	0.035 0.07	Ø3 Ø3	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1 CUV2065-1 CUV2066-1	D g continuo E E E	Quartz Pus flowthrou Quartz Quartz Quartz Quartz	2 2 2 2 2	5 10 20	7.4x12.5x45 12.5x12.5x45 22.6x12.5x45	0.035 0.07 0.14	Ø3 Ø3 Ø3	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1 CUV2065-1 CUV2066-1 CUV2062-1	D g continuo E E E E F	Quartz	2 2 2 2 2 2 2	5 10 20 30 10	7.4x12.5x45 12.5x12.5x45 22.6x12.5x45 32.6x12.4x45	0.035 0.07 0.14 0.21 0.48	Ø3 Ø3 Ø3 Ø3	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1 CUV2065-1 CUV2066-1 CUV2062-1 Self masking	D g continuo E E E E F	Quartz	2 2 2 2 2 2 2	5 10 20 30 10	7.4x12.5x45 12.5x12.5x45 22.6x12.5x45 32.6x12.4x45 12.5x12.5x45	0.035 0.07 0.14 0.21 0.48	Ø3 Ø3 Ø3 Ø3	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1 CUV2065-1 CUV2066-1 CUV2062-1 Self masking CUV2614-1	Continuo E E E E E F F g continuo	Quartz	gh cell 2 2 2 2 2 2 2 ugh cell, small	5 10 20 30 10	7.4x12.5x45 12.5x12.5x45 22.6x12.5x45 32.6x12.4x45 12.5x12.5x45 32.6x12.4x45	0.035 0.07 0.14 0.21 0.48	Ø 3 Ø 3 Ø 3 Ø 3 4x12	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1 CUV2065-1 CUV2062-1 Self masking CUV2614-1 Micro Cell w	Continuo E E E E E F F g continuo	Quartz	gh cell 2 2 2 2 2 2 2 ugh cell, small	5 10 20 30 10	7.4x12.5x45 12.5x12.5x45 22.6x12.5x45 32.6x12.4x45 12.5x12.5x45 32.6x12.4x45	0.035 0.07 0.14 0.21 0.48	Ø 3 Ø 3 Ø 3 Ø 3 4x12	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1 CUV2065-1 CUV2066-1 CUV2062-1 Self masking	E E E F s continuo H ith black	Quartz	gh cell 2 2 2 2 2 2 2 ugh cell, small 2	5 10 20 30 10 input, larg	7.4x12.5x45 12.5x12.5x45 22.6x12.5x45 32.6x12.4x45 12.5x12.5x45 12.5x12.5x45 12.4x12.4x35.6	0.035 0.07 0.14 0.21 0.48 n	Ø 3 Ø 3 Ø 3 Ø 3 4x12	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1 CUV2065-1 CUV2062-1 Self masking CUV2614-1 Micro Cell w CUV2674-1	E E E F s continuo H ith black	Quartz	gh cell 2 2 2 2 2 2 2 ugh cell, small 2	5 10 20 30 10 input, larg	7.4x12.5x45 12.5x12.5x45 22.6x12.5x45 32.6x12.4x45 12.5x12.5x45 12.5x12.5x45 12.4x12.4x35.6	0.035 0.07 0.14 0.21 0.48 n	Ø 3 Ø 3 Ø 3 Ø 3 4x12	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2061-1 CUV2065-1 CUV2062-1 Self masking CUV2614-1 Micro Cell w CUV2674-1 Fluorescence	E E E F S continuo H	Quartz Ous flow through the property of	2 2 2 2 2 2 2 2 3 agh cell, small 2	5 10 20 30 10 input, larg 10	7.4x12.5x45 12.5x12.5x45 22.6x12.5x45 32.6x12.4x45 12.5x12.5x45 ge output Z=8.5mn 12.4x12.4x35.6	0.035 0.07 0.14 0.21 0.48 n 0.03	Ø 3 Ø 3 Ø 3 Ø 3 4x12	
CUV2033-1 CUV2034-1 Self masking CUV2063-1* CUV2065-1 CUV2066-1 CUV2062-1 Self masking CUV2614-1 Micro Cell w CUV2674-1 Fluorescence CUV2051-1	E E E F G Continuo H ith black v J e A A	Quartz Ous flow throughout Quartz Quartz Quartz Quartz Quartz Quartz Quartz	2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 4	5 10 20 30 10 input, larg 10	7.4x12.5x45 12.5x12.5x45 22.6x12.5x45 32.6x12.4x45 12.5x12.5x45 ge output Z=8.5mn 12.4x12.4x35.6 12.5x12.5x45	0.035 0.07 0.14 0.21 0.48 n 0.03	Ø 3 Ø 3 Ø 3 Ø 3 4x12 Ø 2	







Style G

Style H

Style J

Accessories

Z-Dimensions Are Not Created Equal

Cuvettes come in a variety of shapes and sizes, but one of the most important specifications of a cuvette is its Z-dimension. The Z-dimension of an instrument (cuvette holder or spectrometer) is the distance from the bottom of the cuvette chamber floor to the center of its light beam (see image). A cuvette's Z-dimension must match the Z-dimension of the instrument with which it will be used.



Each manufacturer designs its instruments with a specific Z-dimension. Common Z-dimensions include 8.5 and 15mm, and sometimes 20mm. When purchasing small volume cuvettes, the correct Z-dimension becomes critical. Matching the Z-dimension of the cuvette to the Z-dimension of the instrument ensures that the light beam passes through the center of small samples. The table below shows the standard Z-dimension of the spectrometer sample compartments for many manufacturers.

Manufacturer	Z-Dimension
Agilent®	15 mm
Avantes®	15 mm
Beckman®	8.5 mm
Bio-Rad®	8.5 mm
Cecil®	15 mm
Eppendorf®	8.5 mm
Hewlett - Packard®	15 mm
Hitachi®	8.5 mm
Jasco®	11 mm
J & M®	8.5 mm
Ocean Optics®	15 mm
Perkin – Elmer®	15 mm
Pharmacia®	15 mm
Shimadzu®	15 mm
Spectronics®	8.5 mm
Stellarnet®	15 mm
Turner®	8.5 mm
Varian ®	20 mm
WPI	15mm

To determine the Z-dimension of a cuvette holder:

- Use strips of heavy paper that will fit neatly into a cuvette (for example, 12mm x 50mm) and not allow light to pass through the cuvette.
- Poke a tiny hole in each paper "sample." For example, one paper sample could have a hole at 8.5mm, one at 15mm, one at 20mm.
- One at a time, insert the paper samples into the cuvette and place the cuvette into the cuvette holder. The paper sample with the pin hole at the instrument's Z-dimension will allow light to pass. The other paper samples will not allow light to pass.

If you have an instrument that is not on the list and need to know its Z-dimension, please contact WPI at 941-371-1003 or technicalsupport@wpiinc.com.

Fiber Optic Collimator



WPI's Fiber Optic Collimator can be used for both collimating a light beam emitted by an optical fiber or coupling light from a collimated light beam into an optical fiber. The numerical aperture of the collimator is optimized for maximum coupling efficiency into typical fused silica fibers. The collimator can, for example, be used to guide a parallel light beam through a sample cuvette or an optical filter with virtually no optical losses. In this application, one collimator collimates the light into a parallel beam 5 mm in diameter, enabling it to pass a long distance without losing the energy. After the light passes the sample media, a second collimator can be used to collect the beam into the receiving fiber. A unique design feature of this collimator is that the distance between the lens and the optical fiber can be easily adjusted. This permits it to be used as a focusing device or for fine-tuning the color balance when coupling light from a light source into multimode fibers.

COLLIMATOR Specifications

LENS DIAMETER 5 mm LENS FOCAL DISTANCE

LENS MATERIAL Ultraviolet grade synthetic fused silica (KU-

WAVELENGTH RANGE 220 nm-2 μm MOUNTING THREADS 3/8-24 UNF

DIVERGENCE < 0.1 rad for 1 mm core fiber

FIBER CONNECTOR INTERFACE SMA or ST

300051	Fiber Optic Collimator (SMA)
300052	Fiber Optic Collimator (ST)

OPTIONAL ACCESSORIES

13395	SMA Bulkhead Feedthru connector/coupler, D-hole
13370	SMA half-length Bulkhead coupler/connector
CC-3-UV	Cosine Corrector







CC-3-UV

Filter Holder for Glass Fiber Filters

OFT1



Simple measurements for particulate absorption

WPI's filter holder for particulate absorption measurements is specially designed for field use. It is rugged and portable. It performs as well as a laboratory based spectrophotometer. It can be directly connected to WPI's line of fiber optic spectrometers and light sources. Instead of collecting your samples, transporting them to a laboratory, and accepting the loss of spectral information associated with it (Sosik, 1999), particulate absorption can now be measured on site.

How does it work?

Particulate absorption of fresh and seawater can be determined by filtering a known amount of sample through a Glass Fiber Filter (GF/F) and measuring the particulate absorption coefficient $a_{\text{P}}(\lambda)$ concentrated on the filter. This technique is called quantitative filter technique (QFT) and corrects for the pathlength amplification, an effect of scattering. The correction of the pathlength amplification and the correction of the non-linear relationship between the optical density of samples on a Whatman GF/F filter and in suspension are discussed in Mitchell (1990).

Detector and light source requirements

The optical throughput of QFT1 equipped with a classical GF/F filter is very low and requires a matched light source / spectrometer system. WPl's **TIDAS E** in combination with WPl's **FO-6000** tungsten light source or **D4H** deuterium/halogen light source can be used in the 380–730 nm and 280–730 nm wavelength range, respectively. The QFT1 can also be interfaced to any

other CCD, PDA or scanning type spectrometer with fiber optic capabilities.

Performance

A significant advantage of the filter holder is its large beam diameter of 5 mm, resulting in "averaging out" of larger non-organic particles frequently found on the filter pad when using natural samples. The removable filter fixture allows simple filter alternation and cleaning.

Specifications

GF/F Filter Diameter 25 mm

Wavelength Range 280-730 nm *

Fiber Optic Connection Ø 600 µm / SMA

Material in contact with filter pad Delrin

Weight0.5 kg (1 lb)

* Using a TIDAS E spectrometer and D4H UV/VIS light source.

References

Mitchell, B. G., "Algorithms for Determining the Absorption Coefficient of Aquatic Particles Using the Quantitative Filter Technique (QFT)", SPIE Vol. 1302 *Ocean Optics X* (1990), 137-148.

Sosik, H. M., "Storage of marine particulate samples for light-absorption measurements", *Limnol. Oceanogr.*, 44(4), 1999, 1139-1141

M. Belz, K. Larsen, K.-F. Klein, "Fiber optic sample cells for polychromatic detection of dissolved and particulate matter in natural waters", *Proc. SPIE*, Vol. 6377, Oct 2006, 63770X

89575 QFT1, Fiber Optic Holder for Glass Fiber Filters

In-Line Fiber Optic Filter Holder



This In-Line Fiber Optic Filter Holder allows the insertion of optical filters within a fiber optic pathway. The connectors of the filter holder assembly are compatible with WPI's range of fiber optic jumper cables and can be coupled using SMA or ST connectors.

Filters with outer diameters from 8 to 25.4 mm and thicknesses from 2 to 10 mm can be accomodated. The design limits lateral and axial movement of the filter when secured in the holder.

Two fiber optic collimators are internally mounted in the holder to pass collimated light through the filter and then refocus the filtered light into the aperture of the output fiber. Spectral range will be largely limited by the bandpass of the optical fibers (from UV to near IR using WPI UV-enhanced cables).

56200 In-Line Fiber Optic Filter Holder (SMA)

World Precision Instruments

www.wpiinc.com

V-Vette

Sample holder for spectroscopic analysis of microliter volume samples

- 2 microliter sample
- 1 mm pathlength
- No moving parts
- Baseline repeatability < 2 mAU
- Patent pending

V-Vette is a fiber optic sample cell with a pathlength of 1 mm for spectroscopic analysis of microliter volume samples. Light is coupled into and out of the sample cell via optical fibers. A 2 µL sample droplet can be conveniently placed into the v-shaped sample compartment from a pipetter. Absorbance of the sample is measured between the input and output fiber after a cover is placed on the sample compartment to minimize stray light. The sample can be picked up and reused or removed by blowing it off with dry air or wiping it off.

V-Vette Microliter Sample Holder



V-VETTE SPECIFICATIONS

FUNCTIONALITY	
PATHLENGTH	.1 mm
WAVELENGTH RANGE	.200 – 1000 nm
FIBER CONNECTION	.600 μm (SMA)
SAMPLE VOLUME	.2-5 µL
BASELINE REPEATABILITY	.< 2 mAU peak to peak

MicroLWCC

Low volume flow cell for FIA, HPLC and Process Analysis

MicroLWCC is a new fiber optic low volume flow cell for UV/VIS/NIR absorbance analysis. Based on WPI's established liquid core waveguide technology, the analyte solution functions as the core of a fluid filled light waveguide. Wetted parts in the sample cell light path are PEEK, fused silica and PTFE. Optical fibers are used to transport light to and from the sample cell. The cell can be used in biochemistry for DNA, RNA & protein quantification, colorimetric nutrient and trace metal analysis, drug discovery and dissolution testing, process control, and HPLC analysis.

LWCC-M SPECIFICATIONS

	LWCC-M-10	LWCC-M-50	LWCC-M-100
OPTICAL PATHLENGTH	10 mm	50mm	100mm
INTERNAL VOLUME	2.4 µL	12 μL	24 µL
REFRACTIVE INDEX @ 280 nm**	< 7 mAU	< 15 mAU	< 30 mAU
WAVELENGTH RANGE		200	0 – 1000 nm
FIBER CONNECTION [µm]			500 (SMA)
TRANSMISSION @ 254 nm	*		>40%
MAXIMUM PRESSURE			> 1000 psi
WETTED MATERIALS		PEEK, Fused	d Silica. PTFE

* Reference: 2 * 600 µm Fiber, butt-coupled

** Measured using ASTM E 685 - 93

WPI U.S. Patents: 5,444,807; 5,570,447; 5,604,587; 6,603,556; 6,385,380.



LWCC-M-10	Low Volume Flow Cell, 10 mm pathlength
LWCC-M-50	Low Volume Flow Cell, 50 mm pathlength
LWCC-M-100	Low Volume Flow Cell, 100 mm pathlength

References

M. Belz, "Simple and sensitive protein detection system using UV LEDs and liquid core waveguides", Advanced Environmental, Chemical, and Biological Sensing Technologies V, Optics East, Oct 2007, Proc. SPIE, Vol. 6755, 675505

M. Belz, F. A. Klein, H. S. Eckhardt, K. Klein, D. Dinges, K. T. V. Grattan, "Optical Detection Techniques and Light Delivery with UV LEDs and Optical Fibres", Third International Conference on Optical and Laser Diagnostics, Proc. IOP, City University, London, UK, May 2007.

M. Belz, P. Dress, A. Sukhitskiy, S. Liu, "Linearity and effective optical pathlength of liquid waveguide capillary cells", Part of the SPIE Conference on Internal Standardization and Calibration; Architectures for Chemical Sensors, Boston Massachusetts, September 1999, SPIE Vol. 3856, 271-281.

Optical Fibers

Bifurcated Fiber Optic Assemblies

Use to combine similar intensity light from differing sources or to split a light source into two fibers. For use with a dual spectrometer as a reference.



BIF22	Split or combine similar intensities (200/200)
BIF44	Split or combine similar intensities (400/400)
BIF41	Combine UV (400) + VIS (100)
BIF62	Combine UV (600) + VIS (200)
BIF66	Split or Combine Similar Intensities (600/600)

WPI can build custom fiber optic assemblies for many **UV/VIS/NIR** applications. Call for more information.

Plastic Fiber Optic Cables

More flexible than glass fibers, these inexpensive PMMA plastic fibers can be used for illumination and scientific applications. They are excellent for light transfer between 350 nm and 1000 nm. Their maximum temperature should be kept below 80° C.

PLASTIC FIBER OPTIC CABLES (NON UV), 400 TO 1000 NM


FOP1-SMA	Plastic Fiber Optic Cable, SMA connectors, 1 mm x 2 m
FOP1-SMA/ST	Plastic Fiber Optic Cable, ST/SMA connectors, 1 mm x 2 m
FOP1-ST	Plastic Fiber Optic Cable, ST connectors, 1 mm x 2 m

UV-enhanced fiber optic cables



generation of color centers in the fiber material. The lifetime of such a fiber, defined by the 1/e reduction in transmission at 240 nm, is normally less than 200 minutes. This effect renders them unsuitable for UV spectroscopy below 240 nm.

Anti-solarization fibers suitable for deep UV spectroscopy solve this problem. These fibers stabilize to less than 5% degradation over a period of 20,000 hours after an initial transmission "burn-in" loss of less than 25%. Additionally, this anti-solarization characteristic will not degrade over time.



Features

 Broad UV/Vis spectral range
 Laser damage resistant • High core to clad ratios • Broad temperature range • Bio-compatible materials Radiation resistant
 Sterilizable by ETO and gamma radiation • Higher transmission than PCS between 180-nm and 300 nm

Properties

- Multimode Pure silica core Numerical aperture: 0.22 ± 0.02 (standard)
- Standard prooftest: 70 kpsi Minimum bend radius: 100x clad radius (momentary), 600x clad radius (long term)

Anti-Solarization

The transmission of conventional UV-enhanced silica/silica fiber decreases rapidly at wavelengths below 240 nm when exposed to high intensities of a deuterium lamp. This effect is called "UV-solarization" and results from the

UV-ENHANCED FIBER OPTIC CABLES, 230 - 1000 NM

FO-50-SMA1M	Fiber Optic Cable, 1 m, SMA, 50 µm Core, UV-Enhanced
FO-50-SMA	Fiber Optic Cable, 3 m, SMA, 50 µm Core, UV-Enhanced
FO-100-SMA1M	Fiber Optic Cable, 1 m, SMA, 100 µm Core, UV-Enhanced
FO-100-SMA	Fiber Optic Cable, 3 m, SMA, 100 µm Core, UV-Enhanced
FO-200-SMA1M	Fiber Optic Cable, 1 m, SMA, 200 µm Core, UV-Enhanced
FO-200-SMA	Fiber Optic Cable, 3 m, SMA, 200 µm Core, UV-Enhanced
FO-400-SMA1M	Fiber Optic Cable, 1 m, SMA, 400 µm Core, UV-Enhanced
FO-400-SMA	Fiber Optic Cable, 3 m, SMA, 400 µm Core, UV-Enhanced
FO-400SMA/ST	Fiber Optic cable, 1 m, SMA/ST connector, 400 µm core, UV-Enhanced
FO-600-SMA1M	Fiber Optic Cable, 1 m, SMA, 600 µm Core, UV-Enhanced
FO-600-SMA	Fiber Optic Cable, 3 m, SMA, 600 µm Core, UV-Enhanced
FO-1000-SMA1M	Fiber Optic Cable, 1 m, SMA, 1000 µm Core, UV-Enhanced
FO-1000-SMA	Fiber Optic Cable, 3 m, SMA, 1000 µm Core, UV-Enhanced

ANTI SOLARIZATION FIBER OPTIC CABLES, 190 - 1000 NM

FO-200AS-SMA	Fiber Optic Cable, 1 m, SMA, 200 µm Core, Anti-Solarization
FO-400AS-SMA	Fiber Optic Cable, 1 m, SMA, 400 µm Core, Anti-Solarization
FO-600AS-SMA	Fiber Optic Cable, 1 m, SMA, 600 µm Core, Anti-Solarization

Product Index

Item	Page	Item	Page	Item	Page	Item	Page	Item	Page
Sym	ıbols		101		143, 147		52, 53, 54, 55		25
1B100-3	112		120		61, 101		171, 177		G 25
1B100-4			120	13370			128, 129, 132, 134, 135		
1B100-6			101, 106		101				
1B100F-3	112		85		198				
1B100F-4			85		101		18		
1B100F-6			101 99	13555 13620					17 20
1B120			85		101				
1B120-4			85		30				
1B120-6			81		101				
1B120F-3			177		99		22		101
1B120F-4		4070	177		122		22		
1B120F-6				13854			22		
1B150-3			94 94, 97				22		
1B150-4			138		20		121		
1B150-6			85		20		121		
1B150F-3			85 85						
1B150F-4					98	159/0	101		
1B150F-6			85				171		
1B200-4	112		85 85		99 99		121 121		169 169
1B200-6	112				99 99				
1B200F-4				14036-15			121		
1B200F-6			101	14036-10	99 99		121		
2B150F-4			101	14039-10	99		121		G 22 22
2B150F-6	114		101	14040-50			121		
3B120F-4	114		101		99		121		G 22
3B120F-6	114	53//	51, 54, 61		99		121		24
5B120F-4			51, 54, 61		99		121		G 24
5B120F-6	114				99		121		24
7B100F-4			52, 61		99		121		G 24
7B100F-6	114		172		99		199		
7B120F-4	114		172		99		195		
260	77		172		99		195		
705	69		172	14055-2			117, 119		G 27
705PF	69		172	14057-10			117, 119		50024827
711P		5435	51, 52		99	75090	117		27
712P	70, 71		51, 52		99		172		
715P	70, 71		115		99		172		101
773			109		102		172		101
900A	79				17		172		
900AP	79	5464	129, 132, 135		17		113		101
900APP			77, 104, 120		17				20
1300M	120		104		127				G 20
1350M	120		104		127		185		
1358	101		147		23		185		141
1571	94, 97		138		23		185		141
1965			104		27		185		141
1966	114		104		27		185		
1967	114		74		20		185		117, 119
2026-10	101	6820	94, 97		20		185		177
2478	120		94, 97		20		185		117
2479	120		51, 52		20		197		17
2505			51, 54, 61		20		197		17
2541	69		94, 97		20		197		
2547	68, 71		94, 96		23		187, 199		27
2851	101		94, 96		23		53		27
2932	103, 172		52		22		81		
2933	103, 172		99		22		81		50035627
2935	103		52		22		191		27
3161	101		103		22		151		27
3259			103		18		151		27
3260			117, 119, 177		18		51, 54		27
3294			99		101		51		27
3316			99		30		177		27
3414					22		101		27
3417-10			99		22		193, 198		27
3484			99		22		193, 198		
3485			99		22		75, 104		
3491			99		29		151	500778	177
3492			99		29	300305	58		145
3508			94, 97	15590	151	300443	38	500871	136
5500	IVI	13324	101	15623	101, 167	300521	177	500890	98

2005	Item	Page	Item	Page	Item	Page	Item	Page	Item	Page
2018	500895	98	501715	18	502122	63	503099	153	504280	50
2011										
50178										
501916-6. 23 501934 23 50195 77 50197 157 50588 146							503120	157. 159		
50720										
SHIZEA 2									504285	146
SICPIC 73 SOUTH 75 SOUTH 75 SOUTH 75 SOUTH 76 SO	501220-G	23								
Serice 2 Soff 50 50 50 50 50 50 50					502191	102				
Septiment Sept										
591222 20 501756- 23 502216- 47 503901 771 50444 145 501724-C. 20 501757- 23 502717- 44 503901 25 50444- 145 501724-C. 20 501757- 23 502717- 44 503901 25 504445- 145 501724-C. 30 501757- 73 502717- 45 503901 25 504445- 145 501725- 31 501757- 72 50227- 45 503901 25 504445- 145 501725- 31 50175- 72 50227- 45 503901 15 504401 177 501725- 31 50177- 72 50227- 45 503901 175 50476- 77 501275- 31 50177- 72 50227- 45 503901 144 50455- 77 501275- 31 50177- 31 50177- 20 50227- 45 503901 144 50455- 22 501278- 31 501778- 20 502228- 45 503901 144 50455- 22 501278- 31 50189- 20 502220- 45 503912- 144 50455- 27 501278- 31 50189- 20 502220- 45 503912- 144 50455- 27 501278- 31 50189- 20 50220- 45 503912- 144 50455- 27 501278- 31 50189- 20 50220- 45 503912- 144 50455- 27 501278- 31 50189- 20 50220- 45 503912- 144 50455- 27 501278- 31 50189- 20 50220- 45 503912- 144 50455- 27 501278- 31 50189- 20 50220- 45 503912- 144 50455- 27 501278- 31 50189- 20 50220- 45 503912- 144 50455- 27 501278- 31 50189- 20 50220- 45 503912- 144 50455- 27 501278- 31 50189- 20 50220- 45 503912- 144 50455- 27 50128- 31 50189- 20 50220- 47 50593- 103 50455- 27 50128- 31 50189- 20 50220- 47 50593- 103 50455- 27 50128- 31 50189- 20 50220- 47 50593- 103 50455- 27 50128- 31 50189- 20 50220- 47 50593- 103 50455- 27 50128- 31 50189- 20 50220- 47 50593- 103 50455- 27 50128- 31 50189- 20 50220- 47 50593- 103 50455- 27 50128- 31 50189- 20 50220- 47 50593- 103 50455- 27 50128- 31 50189- 20 50220- 47 50593- 103 50455- 27 50128- 31 50189- 20 50220- 47 50459- 47 50189- 31 50189- 20 50220- 47 50459- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 20 50220- 47 50189- 31 50189- 31 50189- 31 50189- 31 50189- 31 50189- 31 50189- 31 50189- 31 50189- 31 50189- 31 50189- 31 50189-							503260	22	504417	145
50 234			501756	23	502204	46	503301	171	504443	145
90141-	501234	20								
501241	501234-G	20	501757	23	502213	46, 47			504445	145
50/1246 38 50/7584 22 50/2255 47 50/356 155 50/459 27 50/7274 31 50/7594 22 50/2278 45 50/359 155 50/459 27 50/7274 31 50/7594 22 50/2278 45 50/359 155 50/459 27 50/7275 31 50/7594 22 50/2278 45 50/359 155 50/459 27 50/7277 31 50/7784 30 50/229 45 50/359 154 50/552 27 50/7277 31 50/7784 30 50/229 45 50/359 154 50/552 27 50/7277 31 50/7784 27 50/779 46 50/779 17 50/7279 31 50/7784 27 50/779 46 50/779 17 50/7270 31 50/7784 27 50/779 46 50/779 10 50/779 10 50/7280 31 50/7784 27 50/779 46 50/779 10 50/779 10 50/7280 31 50/7784 27 50/779 46 50/779 10 50/779 10 50/7280 31 50/7784 27 50/779 46 50/779 10 50/779 10 50/7280 31 50/7784 27 50/779 46 50/779 10 50/779 10 50/7280 31 50/7784 27 50/779 46 50/779 10 50/779 10 50/7280 31 50/7784 27 50/779 10 50/779 10 10 50/7280 31 50/7784 27 50/779 10 50/779 10 10 50/7280 31 50/7784 27 50/779 10 50/779 10 10 50/7280 31 50/7784 27 50/779 10 50/779 10 10 50/7280 31 50/7784 27 50/779 10 50/779 10 10 50/7280 31 50/7784 27 50/779 10 50/779 10 10 50/7280 31 50/7784 27 50/7784 27 50/779 10 50/779 10 10 50/7280 31 50/7784 27 50/7784 27 50/779 10 50/779 10 10 50/7280 31 50/7784 27 50/7784 27 50/779 27 50/7280 31 50/7784 27 50/7784 27 50/7784 27 50/7280 31 50/7784 27 50/7784 27 50/7784 27 50/7280 31 50/7784 27 50/7784 27 50/7784 27 50/7280 31 50/7784 27 50/7784 27 50/7784 27 50/7280 31 50/7784 27 50/7784 27 50/7784 27 50/7280 31 50/7784 27 50/7784 27 50/7784 27 50/7280 31 50/7784 27 50/7784 27 50/7784 27 50/7280 31 50/7784 27 50/7784 27	501241	18			502224	47			504459	28
901725-01326. 31 501736 22 502276. 46 508507 155 500500 27 7501777. 31 501736 22 502278. 45 508500 144 509550 27 7501777. 31 5017378 20 502258. 46 508510 144 509550 27 7501778. 31 5017378 20 502258. 46 508510 144 509552 27 7501788. 31 501838. 20 502251. 46 508510 144 509552 27 7501788. 31 501839. 20 502251. 46 508510 144 509552 27 7501789. 31 501820. 32 502251. 46 508510 149 509552 27 7501789. 31 501820. 32 502251. 46 508510 149 509552 27 7501789. 31 501820. 32 502251. 46 508510 149 509552 27 7501789. 31 501820. 32 502251. 46 508510 149 509552 27 7501789. 31 501820. 32 502251. 46 508510 149 509550. 32 7501820. 31 501820. 32 502251. 46 508510 149 509550. 32 7501820. 31 501820. 32 502251. 46 508510 149 509550. 32 7501820. 31 501820. 32 502251. 46 508510 149 509550. 32 7501820. 31 501820. 32 502251. 47 508512 149 509570. 122 7501820. 31 501821. 32 502251. 47 508512 149 509570. 152 7501820. 31 501821. 32 502251. 47 508512 149 509570. 152 7501820. 31 501821. 32 502251. 47 508512 149 509570. 152 7501820. 31 501821. 32 502251. 47 508512 149 509570. 164, 165 7501820. 31 501821. 32 502251. 47 508512 149 509570. 164, 165 7501820. 31 501821. 32 502251. 47 508517. 31 509577. 164, 165 7501820. 31 501821. 32 502251. 47 508587. 47 509589. 169 7501820. 31 501820. 32 502251. 47 508587. 47 509589. 169 7501820. 31 501820. 32 502251. 47 508587. 31 509579. 164, 165 7501820. 31 501820. 32 502251. 47 508599. 189 7501820. 31 501820. 32 502251. 47 508599. 189 7501820. 31 501820. 32 502251. 43 508599. 189 7501820. 31 501820. 32 502251. 43 508599. 189 7501820. 31 501820. 32 502251. 44 508599. 32 47 509589. 169 7501830. 31 501820. 32 502251. 43 508599. 32 47 509689. 169 7501830. 31 501840. 32 502251. 43 508599. 32 47 509689. 169 7501830. 31 501840. 32 502251. 34 508599. 32 47 509689. 149 7501830. 31 501850. 32 502251. 34 508599. 32 509779. 109 7501830. 31 501850. 32 502251. 34 508599. 32 509779. 109 7501830. 31 501850. 32 502251. 34 508599. 32 509779. 109 7501830. 31 501850. 32 502251. 34 508599. 32 509779. 109 7501830. 31 501850. 32 5022			501758	22	502225	47			504539	27
591274. 31 501796 22 501227 45 508590 155 504590 277 501277 31 501778 20 502229 45 508511 144 504552 77 501277 31 501778 20 502229 45 508511 144 504552 77 501277 31 501778 20 502229 45 508511 144 504552 77 501277 31 50178 20 502229 45 508511 144 504552 77 501279 31 50182 28 502231 46 508513 13, 15 504554 77 501279 31 50182 28 502231 46 508513 13, 15 504554 77 501282 31 50185 28 502231 46 508513 13, 15 504554 77 501282 31 50185 28 502231 46 508513 13, 15 504554 77 501282 31 50185 28 502231 46 508513 13, 15 504554 12 504554 77 501282 31 50185 28 502231 46 508513 13, 15 504554 12 504554 77 501285 31 501852 28 502236 47 508513 13, 15 504554 12 504554 78 501285 31 501855 28 502236 47 508513 13, 15 504556 15, 15 501285 13 501855 28 502236 47 508513 13, 15 504556 16, 16 501285 13 501855 28 502236 47 508513 13, 15 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 15 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 15 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 15 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 15 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 16 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 16 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 16 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 16 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 16 504556 16, 16 501285 31 501855 28 502236 47 508513 13, 16 504556 16, 16 501856 28 502236 47 508513 13, 16 504556 16, 16 501856 28 502236 47 508513 13, 16 504556 16, 16 501856 28 502236 48 502236 48 508513 13, 16 504556 16, 16 501856 28 502236 48 502236 48 50856 16, 16 501856 28 502236 48 502236	501273-501286.	31			502226	46	503507	155	504540	27
501275. 31 501734. 20 502279. 45 503511. 144 50455. 27 501278. 31 501639. 20 502290. 46 503511. 144 50455. 27 501278. 31 501639. 20 502290. 46 503511. 144 50455. 27 501278. 31 501634. 28 502211. 46 503512. 146 50455. 27 501279. 31 501642. 28 502211. 46 503512. 146 50455. 27 501279. 31 501642. 28 502212. 46 503512. 146 50455. 27 501279. 31 501642. 28 502213. 46 503512. 10 50455. 27 501280. 31 50165.01858. 28 5022134. 46 503512. 10 50455. 27 501281. 31 50165.01858. 28 5022134. 46 503532. 10 50 50455. 27 501281. 31 50165. 32 502215. 47 503532. 10 50 50455. 32 502215. 47 503532. 10 50 50455. 15 501281. 31 501651. 32 502215. 47 503532. 10 50 50455. 15 501281. 31 501651. 32 502215. 47 503532. 10 50 50455. 15 501281. 31 501651. 32 502215. 47 503532. 10 50 50455. 15 501281. 31 501651. 32 502217. 28 47 503532. 10 50 50455. 15 501281. 31 501655. 28 50224. 47 503547. 15 501281. 31 501655. 28 50224. 47 503567. 44 503532. 31 501552. 32 50245. 31 501651. 32 502217. 32 502217. 32 50 50231. 31 501651. 32 502217. 32 502217. 32 502217. 32 50 50231. 31 501651. 32 500455. 31 501651. 32 500454. 47 503567. 43 500689. 15 501281. 31 501651. 32 500244. 47 503567. 43 500689. 15 501281. 31 501651. 32 500244. 47 503567. 43 500689. 15 501281. 31 501652. 32 502242. 47 503567. 32 500491. 10 501651. 32 500491					502227	45	503509	155	504550	27
501278. 31 501839. 20 501229	501275	31			502228	45	503510	144	504551	27
501289			501778	20	502229	45	503511	144		
501280. 31 501846. 27 50222	501278	31	501839	20	502230	45	503512	144	504553	27
501781	501279	31	501842	28	502231	45	503513	119, 145	504554	27
501280 31 501851 28 502734 45 5028521 103 504877 152 501284 31 501852 28 502285 47 503532 103 504874 152 501284 31 501853 28 502285 47 503532 103 504874 161, 65 501286 31 501855 28 502285 47 503538 101, 152 504878 164, 165 501286 31 501855 28 50228 47 503538 63 504878 161, 161, 161, 161, 161, 161, 161, 161	501280	31	501846	27	502232	45	503529	103	504555	27
\$\frac{91738}{31}\$ \$\frac{31}{31}\$ \$\frac{91083}{31}\$ \$\frac{31}{31}\$ \$\frac{91083}{31}\$ \$\frac{31}{31}\$ \$\frac{91083}{31}\$ \$\frac{31}{31}\$ \$\frac{91085}{31}\$ \$\frac{31}{31}\$ \$\frac{91085}{316}\$ \$\frac{32}{31}\$ \$\frac{91085}{31}\$ \$\frac{31}{31}\$ \$\frac{91085}{31}\$ \$\frac{91085}{31}\$ \$\frac{91085}{31}\$ \$\frac{91085}{31}\$ \$\frac{91085}{31}\$ \$\frac{91085}{31}\$ \$\frac	501281	31					503530	103	504556	27
901284 31 501653 28 50226 47 503356 101, 152 50478, 161, 165 501286 31 501656 28 50274 46 503538 63 504578, 161, 165 501286 31 501655 28 50274 46 503538 63 504578, 161, 165, 161, 165 501286 31 501655 28 50274 47 503568, 36 504578, 161, 165, 161, 165 501288 31 501656 28 50224 47 503568 138 504698 169 501288 31 501657 28 50224 47 503568 138 504698 169 501281 31 501685 28 502744 47 503568 138 504698 169 501291 31 501685 28 502744 47 503569 138 504699 169 501291 31 501685 28 502744 47 503569 138 504699 169 501291 31 501685 28 502745 47 503570 138 504699 168 501312 24 501685 28 50225 43 503571 38 504596 141 501221 97 501683 28 50225 43 503571 38 504596 141 501221 97 501684 28 50225 43 503573 38 504597 141 501221 97 501684 29 50225 43 503573 38 504697 141 501221 31 501913 29 50226 43 503566 22 504713 101 501369 143 501913 29 50226 43 503666 22 504713 101 501369 143 501913 29 50226 43 503666 22 504713 101 501369 143 501913 29 50226 44 503666 22 504713 101 501370 143 501923 25 50338 44 503666 22 504713 101 501370 143 501923 25 50338 44 503666 22 504713 101 501371 443 501923 25 50338 44 503660 22 504713 101 501371 443 501923 24 503638 44 503660 22 504715 101 501371 443 501928 34 501933 44 503666 22 504713 101 501371 443 501928 34 501933 44 503666 22 504715 101 501371 443 501928 34 501928 34 503698 34 503			501851	28	502234	45	503531	103	504570	152
S01285	501283	31			502235	47			504574	152
501788. 31 501855. 28 50274 46 50553 3 3 504578. 162, 143, 144, 165 501288. 18 501856. 28 50274 47 50556 138 504599. 169 501288. 6 18 501857. 28 50224 47 50556 138 504599. 169 501291. 18 501888. 28 50224 47 50556 138 504599. 169 501291. 6 18 501888. 28 50224 47 50556 138 504599. 169 501291. 6 18 501886. 28 50224 47 50556 138 504599. 169 501291. 7 18 501862. 28 50224 47 50556 138 504591. 168 501371. 24 501862. 28 50225 47 50557. 138 504591. 168 501372. 97 501863. 28 50225 43 50557. 138 504591. 168 501372. 97 501863. 28 50225 43 50557. 38 504597. 141 501371. 97 501863. 28 50225 43 505599. 28,47 504606. 145 501352. 143 501913. 29 50226 43 505599. 28,47 504606. 145 501353. 143 501914. 29 50226. 43 505599. 28,47 504608. 47 501369. 143 501915. 29 50220. 44 503666. 22 50471 4. 101 501369. 143 501915. 29 50220. 44 503666. 22 50471 4. 101 501379. 143 501922. 25 50229 44 503666. 22 50471 4. 101 501371. 143 501986. 17 50239. 44 503666. 22 50471 6. 101 501372. 143 501986. 17 50239. 44 503666. 22 50471 6. 101 501373. 143 501986. 17 50239. 44 503660. 22 50471 6. 101 501373. 143 501986. 17 50239. 44 503660. 22 50471 6. 101 501373. 143 501986. 17 50239. 44 503670. 22 50471 6. 101 501373. 143 501986. 17 50239. 44 503670. 22 50471 6. 101 501373. 143 501980. 141 502610. 42 503671. 22 55500517. 19 501373. 143 501980. 141 502610. 42 503671. 22 55500517. 19 501373. 143 502000. 141 502610. 42 503671. 22 55500517. 19 501375. 143 502000. 141 502610. 42 503673. 20 55500771. 19 501376. 143 502001. 141 502610. 42 503673. 20 55500771. 19 501378. 143 502001. 141 502610. 42 503673. 20 55500771. 19 501379. 143 502001. 141 502610. 42 503673. 20 55500771. 19 501379. 143 502001. 141 502610. 42 503673. 20 55500771. 19 501379. 143 502001. 141 502610. 42 503673. 20 55500771. 19 501379. 143 502001. 141 502610. 42 503673. 20 55500771. 19 501379. 143 502001. 141 502610. 42 503673. 20 55500771. 19 501370. 143 502001. 141 502610. 42 503673. 20 55500771. 19 501378. 143 502001. 141 502610. 42 503673. 20 55500771. 19 501379. 143 502001. 141 502610. 141 502					502236	47	503536	101, 152	504576	
501288			501854	28	502237	28, 47	503537	103	504577	
501289 18 501857 28 502244 47 503569 138 504599 169										
501291	501288	18	501856	28	502242	47	503567	47	504588	169
501291-6 18	501288-G	18					503568	138		
501312	501291	18	501858	28	502244	47	503569	138	504590	169
501321	501291-G	18			502245	47	503570	138	504591	168
501356 76					502257	43				
501352							503573	38	504597	141
501553	501336	26	501864	28			503598	28, 47	504606	145
50186			501913	29			503599	28, 47		
501370					502262	43			504713	101
501377	501369	143	501915	29					504714	101
501372	501370	143			502303	44	503668	22	504715	101
501373							503669	22		
501375	501372	143	501986	96			503670	22	555001FT	19
501376			501989	24	502600	42			555005FT	19
									555007FT	19
Solition					502610	43	503738	30	555009FT	19
Solition	501377	143	502004	141, 143	502612	43		30	555041FT	19
	30.370	143	502005	141, 143	302011111111	43	503740	24	555047FT	19
					502616	43	503741	24	555402NT	25
501607 129 502010 141 503022 159 503844 87 555407NT 25 501608 48 502011 141 503023 159 503847 192 555408NT 25 501609 187,195 502012 141 503041 102 503848 192 555500S 21 501622 126,128,129 502013 141 503042 102 504045 31 555500S 21 501624 126 502016 141 503050 159 504046 31 555500S 21 501624 126 502017 141 503050 159 504046 31 555540S 21 501635 97 502017 141 503050 159 504046 31 555560S 21 501636 146 502018 141 503070 102 504049 31 555560S 21 501635 146 502019	501381	143	502007	141, 143	502650	42	503763	94	555403NT	25
Solidon Sol										
501609. 187, 195 502012. 141 503041 102 503848. 192 555500S. 21 501622. 126, 128, 129 502013. 141 503049. 159 504046. 31 55554S. 21 501624. 126 502016. 141, 147 503049. 159 504047. 31 55554S. 21 501624. 126 502016. 141 503050. 159 504047. 31 55554OS. 21 501635. 97 502017. 141 503051. 143 504048. 31 55556OS. 21 501637. 146 502018. 141 503070. 102 504122. 87 555580S. 21 501641. 62 502040. 155 503071. 102 504122. 87 555583S. 21 501642. 62 502040. 155 503071. 102 504128. 143 6000019. 95 501642. <t< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></t<>										
501622 126, 128, 129 502013 141 503042 102 504045 31 555524S 21 501623 126, 128, 129 502015 141, 147 503049 159 504046 31 5555305 21 501624 126 502016 141 503050 159 504047 31 555540S 21 501635 97 502017 141 503051 143 504048 31 555560S 21 501636 146 502018 141 503067 40 504049 31 555560S 21 501637 146 502019 141 503070 102 504122 87 555833 21 501641 62 502040 155 503071 102 504122 87 555833 21 501642 62 502041 155 503072 102 504128 143 600019 95 501643 62 502053			502011	141	503023	159				
501623 126, 128, 129 502015 141, 147 503049 159 504046 31 555530S 21			502012	141			503848	192	555500S	21
501624 126 502016 141 503050 159 504047 31 555540S 21 501635 97 502017 141 503051 143 504048 31 555560S 21 501636 146 502018 141 503067 40 504049 31 555560S 21 501637 146 502019 141 503070 102 504122 87 555583S 21 501641 62 502040 155 503071 102 504123 141, 143 555584S 21 501642 62 502041 155 503072 102 504128 143 600009 95 501643 62 502053 47 503073 102 504128 143 600011 51,54 501644 62,63 502054 47 503075 102 504134 147 600011 51,54 501650 138 502056 47<										
501635 97 502017 141 503051 143 504048 31 555560S 21 501636 146 502018 141 503067 40 504049 31 555562S 21 501637 146 502019 141 503070 102 504122 87 555583S 21 501641 62 502040 155 503071 102 504123 141, 143 555584S 21 501642 62 502041 155 503072 102 504128 143 600009 95 501643 62 502053 47 503075 102 504128 143 600011 51,54 501644 62,63 502054 47 503075 102 504134 147 600012 51,54 501651 138 502055 47 503076 102 504136 153 600015 51 501652 138 502056 47<	501623	126, 128, 129	502015	141, 147	503049	159				
501636 146 502018 141 503067 40 504049 31 555562S 21 501637 146 502019 141 503070 102 504122 87 555583S 21 501641 62 502040 155 503071 102 504123 141 43 555584S 21 501642 62 502041 155 503072 102 504128 143 600009 95 501643 62 502053 47 503073 102 504128 143 6000019 95 501644 62,63 502054 47 503075 102 504134 147 600012 51,54 501651 138 502055 47 503076 102 504136 153 600015 51 501652 138 502056 47 503077 102 504137 153 600016 51 501653 138 502066					503050	159				
501637 146 502019 141 503070 102 504122 87 5555835 21 501641 62 502040 155 503071 102 504123 141,143 5555845 21 501642 62 502041 155 503072 102 504128 143 600009 95 501643 62 502053 47 503073 102 504129 143 600011 51,54 501644 62,63 502054 47 503075 102 504134 147 600012 51,54 501651 138 502055 47 503076 102 504136 153 600015 51 501652 138 502056 47 503077 102 504137 153 600016 51 501653 138 502056 47 503077 102 504137 153 600016 51 501656 62 502062 46 <td></td> <td></td> <td>502017</td> <td> 141</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>			502017	141						
501641 62 502040 155 503071 102 504123 141, 143 555584S 21 501642 62 502041 155 503072 102 504128 143 600009 95 501643 62 502053 47 503073 102 504129 143 600011 51,54 501644 62,63 502054 47 503075 102 504134 147 600012 51,54 501651 138 502055 47 503076 102 504136 153 600015 51 501652 138 502056 47 503077 102 504137 153 600016 51 501653 138 502060 46 503078-4 102 504138 153 800100 67 501656- 62 502062 46 503089-4 102 504141 87 800292 119 501656-F 62 502063										
501642. 62 502041 155 503072. 102 504128. 143 600009. 95 501643. 62 502053. 47 503073. 102 504129. 143 600011. 51,54 501644. 62,63 502054. 47 503075. 102 504134. 147 600012. 51,54 501651. 138 502055. 47 503076. 102 504136. 153 600015. 51 501652. 138 502056. 47 503077. 102 504137. 153 600016. 51 501653. 138 502060. 46 503078-4. 102 504138. 153 800100. 67 501656. 62 502062. 46 503079-4. 102 504141. 87 800120. 193 501656. 62 502063. 46 503080-4. 102 504142. 87 800292. 119 501657. 62										
501643 62 502053 47 503073 102 504129 143 600011 51, 54 501644 62, 63 502054 47 503075 102 504134 147 600012 51, 54 501651 138 502055 47 503076 102 504136 153 600015 51 501652 138 502056 47 503077 102 504137 153 600016 51 501653 138 502060 46 503079-4 102 504138 153 800100 67 501656 62 502062 46 503079-4 102 504141 87 800120 193 501656-C 62 502063 46 503080-4 102 504141 87 800120 193 501656-F 62 502067 47 503081-4 102 504142 87 800292 119 501658-F 62 502068 <										
501644. 62,63 502054. 47 503075. 102 504134. 147 600012. 51,54 501651. 138 502055. 47 503076. 102 504136. 153 600015. 51 501652. 138 502056. 47 503077. 102 504137. 153 600016. 51 501653. 138 502060. 46 503078.4 102 504138. 153 800100. 67 501656. 62 502062. 46 503079.4 102 504141. 87 800120. 193 501656. 62 502063. 46 503089.4 102 504142. 87 800292. 119 501656.F. 62 502067. 47 503081.4 102 504164. 155 801566. 173 501657. 62 502068. 47 503082.4 102 504165. 155 801963. 173 501658.F. 62										
501651 138 502055 47 503076 102 504136 153 600015 51 501652 138 502056 47 503077 102 504137 153 600016 51 501653 138 502060 46 503078-4 102 504138 153 800100 67 501656 62 502062 46 503079-4 102 504141 87 800120 193 501656-C 62 502063 46 503084-4 102 504142 87 800292 119 501656-F 62 502067 47 503081-4 102 504164 155 801566 173 501657- 62 502068 47 503082-4 102 504165 155 801963 173 501658-F 62 502070 47 503083 102 504166 155 μTip 116 501658-F 62 502102 126										,
501652 138 502056 47 503077 102 504137 153 600016 51 501653 138 502060 46 503078-4 102 504138 153 800100 67 501656 62 502062 46 503079-4 102 504141 87 800120 193 501656-C 62 502063 46 503080-4 102 504142 87 800292 119 501656-F 62 502067 47 503081-4 102 504164 155 801566 173 501657 62 502068 47 503082-4 102 504165 155 801963 173 501658 62 502070 47 503083 102 504166 155 μTip 116 501658.F 62 502102 126 503084 102 504187 29 501705 18 502105 132 503085 102 504188										
501653 138 502060 46 503078-4 102 504138 153 800100 67 501656 62 502062 46 503079-4 102 504141 87 800120 193 501656-C 62 502063 46 503080-4 102 504142 87 800292 119 501656-F 62 502067 47 503081-4 102 504164 155 801566 173 501657 62 502068 47 503082-4 102 504165 155 801963 173 501658 62 502070 47 503083 102 504166 155 μTip 116 501658-F 62 502102 126 503084 102 504187 29 501705 18 502105 132 503085 102 504188 29 A 501705-G 18 502109 160 503086 102 504221 145										
501656 62 502062 46 503079-4 102 504141 87 800120 193 501656-C 62 502063 46 503080-4 102 504142 87 800292 119 501656-F 62 502067 47 503081-4 102 504164 155 801566 173 501657 62 502068 47 503082-4 102 504165 155 801963 173 501658 62 502070 47 503083 102 504166 155 μTip 116 501658-F 62 502102 126 503084 102 504187 29 501705 18 502105 132 503085 102 504188 29 A 501705-G 18 502109 160 503086 102 504221 145 A300 Pulsemaster 89 501714 18 502110 160 503097 153 504278 50										
501656-C 62 502063 46 503080-4 102 504142 87 800292 119 501656-F 62 502067 47 503081-4 102 504164 155 801566 173 501657 62 502068 47 503082-4 102 504165 155 801963 173 501658 62 502070 47 503083 102 504166 155 μTip 116 501658-F 62 502102 126 503084 102 504187 29 29 501705 18 502105 132 503085 102 504188 29 A 501705-G 18 502109 160 503086 102 504221 145 A300 Pulsemaster 89 501714 18 502110 160 503097 153 504278 50 A310 Accupulser 88										
501656-F. 62 502067. 47 503081-4. 102 504164. 155 801566. 173 501657. 62 502068. 47 503082-4. 102 504165. 155 801963. 173 501658. 62 502070. 47 503083. 102 504166. 155 μTip. 116 501658-F. 62 502102. 126 503084. 102 504187. 29 501705. 18 502105. 132 503085. 102 504188. 29 A 501705-G. 18 502109. 160 503086. 102 504221. 145 A300 Pulsemaster. 89 501714. 18 502110. 160 503097. 153 504278. 50 A310 Accupulser. 88										
501657 62 502068 47 503082-4 102 504165 155 801963 173 501658 62 502070 47 503083 102 504166 155 μTip 116 501658-F 62 502102 126 503084 102 504187 29 501705 18 502105 132 503085 102 504188 29 A 501705-G 18 502109 160 503086 102 504221 145 A300 Pulsemaster 89 501714 18 502110 160 503097 153 504278 50 A310 Accupulser 88										
501658 62 502070 47 503083 102 504166 155 μTip 116 501658-F 62 502102 126 503084 102 504187 29 501705 18 502105 132 503085 102 504188 29 A 501705-G 18 502109 160 503086 102 504221 145 A300 Pulsemaster 89 501714 18 502110 160 503097 153 504278 50 A310 Accupulser 88										
501658-F. 62 502102. 126 503084. 102 504187. 29 501705. 18 502105. 132 503085. 102 504188. 29 A 501705-G. 18 502109. 160 503086. 102 504221. 145 A300 Pulsemaster. 89 501714. 18 502110. 160 503097. 153 504278. 50 A310 Accupulser. 88										
501705 18 502105 132 503085 102 504188 29 A 501705-G 18 502109 160 503086 102 504221 145 A300 Pulsemaster 89 501714 18 502110 160 503097 153 504278 50 A310 Accupulser 88									μTip	116
501705-G 18 502109 160 503086 102 504221 145 A300 Pulsemaster 89 501714 18 502110 160 503097 153 504278 50 A310 Accupulser 88										_
501714	501705	18								A
501714									A300 Pulser	naster 89
501714-G										
	501714-G	18	502120	63	503098	153	504279	50	A320D	90

Item

Page

Item

Page

Item

Page

Item

Page

Page

		9	G	G
F	GN-NET7161	ILS025LT 171	K	M9 136
FD35-100 154	GN-NET25 161	ILS gas-tight syringes 171		M10 136
1833 10011111111111111111111111111111111	GN-PC7 161	implantable probe 62	Kelly Hemostatic Forceps	M10L136
FD35PDL-100	GN-PC25	IMPLANTABLE PROBES 38, 49	KG Transducers	M11 136
FD223A	GN-TTL	INC22 50	KITE-L 132	M325 135
FD223AP	GO1-100	INC24 50	KITE-R	M3301 132
FD3510-100 154	G02-100	INC26 50	KITE-TB-L 132	M3301EH . 128, 129, 132, 134, 135
FD5040-100 154			KITE-TB-R 132	
fiber optic cables 201	GO3-100	Incapacitance Meter	KITLWCC 195	M3301L
Fiber Optic Cables 201	G04-100	Incremental Hot Cold Plate Analgesia	KIT-PHYSIO-I26	M3301-M3-L 132
Fiber Optic Collimator	Gold-plated pins, sockets 104	Meter	KIT-UPUV-STARTUP 187	M3301-M3-R
Fiber Optic Filter Holder 199	Gold wire 107	Induction Chambers	KIT-UPVIS-STARTUP 187	M3301R 132
Fiber Optic Illuminator , 120	Gold Wire 107	Infusion Pump 166, 167	KWIK-2 66	Macro sensors 51
Filling Solution 51	GR100-4114	Infusion/Withdrawal Pump 167	KWIKCAL-2 66	Magnetic Ball Joint 136
Filter Holder	GR100-6114	injector 170, 174	Kwik-Cast	Magnetic Holding Devices 136, 138
Filter Holder, Inline 193	Grip Strength Meter 35	inline filter holder 193	KWIK-CAST	Manometer 181
Fire-Polished 112	GSNO 67	installation 151	Kwik-Cast™95	Manual Micromanipulator 132
Five-barrel glass 114	GSNO-5067	Internal Cannula 50	KWIKGARD	Manual Microsyringe Pump 181
	GSNO-10067	Intracellular Amplifier 69, 70	Kwik CardM 04	Mask, Gas Anesthesia
Flexible NO Sensor 53	g-SPIN 103	intracellular injection 170	Kwik-Gard™94	MAX3820107
FLEXREF65	G-SPIN6	INV-101144	KWIKGLUE	MAX4020
flow analysis 185	G-SPIN6-220	Inverted Trinocular Microscope 144	KWIKGUN94	Mayo-Hegar Needle Holder 23, 25
Flow Cell, Low Volume 200	Guide Cannula 50	in vitro fertilization	KWIKH-2 66	MB2 137
Flow Cells 195			KWIKMIX 94	
FLOX-PATCH60	Guillotine 28	IO-KIT	KWIKPOT-2 66	MBMPH5-7 108, 109
FLOX-PROBE60	н	ion selective electrodes 66	Kwik-Sil 94	MBS 121
Fluorescence-based optical sensor	п	iontophoresis	KWIK-SIL 95	MCL3-HS6 129
60	Half-Cells 108	Iridium Microelectrodes 104	Kwik-Sil™95	MCL3-SM325 128
FluoroDish	Halogen Lamp 192	Iris Forceps 18	Kwik-Tip™	McPherson Forceps 19
FO-50-SMA	Heater Controller 151	Iris Scissors 21, 22	KWIKTPP-2	McPherson-Vannas Scissors 20
FO-50-SMA1M	Heating Filament 117, 119	IRM23E01KT 105	KZ1101 40, 116	Mcpherson-Westcott Scissors 21
FO-100-SMA	Heating Plate	IRM23E10 105	KZ1106 40	Mcpherson-Westcott Stitch Scissors
	high throughput screening 80	IRM23E15 105	KZ1100 40	21
FO-100-SMA1M	Horizontal Tissue Bath 5	IRM23E20KT 105	L	MD4L 135
FO-200AS-SMA		IRM23E25 105	_	MD4-M3-L
FO-200-SMA	Hot Plate Analgesia Meter 37	IRM23E25KT	Lancets, Animal 27	MD4-M3-R
FO-200-SMA1M	hot spots	IRM23E30	Lapping Film, Alumina; 121	MD4R 135
FO-400AS-SMA201	HS6	IRM23E30KT	Lapping Film, Diamond; 121	
FO-400-SMA 201	HS6-3		laser 10	mechanical allodynia
FO-400-SMA1M 201	HS6-3M 129	ISO2	LED light source 193	MechanoCulture
FO-400SMA/ST	HT-1 49	ISO2 Filling Solution	LED-Lite 193	megohm meter
FO-600AS-SMA 201	HT-2 49	ISO-80	LED-LITE 193	MEH1F45 108
FO-600-SMA 201	Humidity Probe 151	ISO-80P 76	LED Ringlight	MEH1F45xx 109
FO-600-SMA1M 195, 201	Hydrogen Electrode 66	ISO-H2S-2 51	LEDspec	MEH1R 108
FO-1000-SMA 201	Hydrogen Peroxide Sensors 51, 54	ISO-HPO-2 51, 54	LEDSpec-2	MEH1RF 110
FO-1000-SMA1M	Hydrogen Sulfide Sensor 51	ISO-HPO-10054		MEH1RFxx109
FO-6000 193	Hydrostatic pressure 181	ISO-HPO-100-L53	LEDSpec-4	MEH1Rxx 109
FO-6000-FILT 193	.,, -, -, -, -, -, -, -, -, -, -, -, -, -	Isolated Bioamplifier 76	Light guides	MEH1S 110
	I	Isolated Differential Amplifier 76	Light Sources	MEH1SF 110
FOIMPH		ISO-NOP 51, 52	Linear Motor Control	MEH1SFxx
FOIMPH-LF	IE010	ISO-NOP00753	Linear Stimulus Isolator 93	MEH1Sxx
Footswitch (OmniDrill35) 28	IE190	ISO-NOP30L 53	Liquid Ion Exchangers 66	MEH2R
FOP1-SMA 201	IE200	ISO-NOP70L 53	Liquid Waveguide Capillary Cell 194	MEH2RF 110
FOP1-SMA/ST	IGS100 55		Live cell imaging 151	
FOP1-ST 201	II-39 37	ISO-NOP3005 53	L-shaped sensors 53	MEH2RFW
Forceps 18	II-336T 37	ISO-NOP302053	Luer-to-Tubing Coupler 98	MEH2RFWxx
force transducers 78	II-336TG 37	ISO-NOPF100 53	Luer Valve Kit	MEH2RFxx
Force Transducers 13	II-390 37	ISO-NOPF200 53	LWCC 194	MEH2RW 110
FORT10g 78	II-390G 37	ISO-NOPF200-L10 53	LWCC-3050 194. 195	MEH2RWxx 109
FORT25 78	II-520MR35	ISO-NOPF500 53	LWCC-3100 194, 195	MEH2Rxx109
FORT100	II-600MR36	ISO-NOPNM 52	LWCC-3250	MEH2S 110
FORT250	II-2200 35	ISO-NOP Rejuvenator 67	LWCC-3500	MEH2SF 110
FORT1000	II-239034	ISO-OXY-2 51, 54		MEH2SFW 110
FORT5000	II-239134	Isostim 90	LWCC-4010	MEH2SFWxx109
	II-2392 34	ISO-TEMP-254	LWCC-4050 195	MEH2SFxx 109
FORT transducers		IT-1E	LWCC-4100 195	MEH2SW110
Frame Clamp 102	II-2393 34	IT-14	LWCC-M-10	MEH2SWxx
FrameWorks 102	II-2394 34	IT-18	LWCC-M-50 200	MEH2Sxx
Free Radical Analyzer 56	II-2395 34	IT-18EXLONG	LWCC-M-100 200	MEH3F45 110
_	II-2396			
G	II-2397 34	IT-21	M	MEH3F45xx 109
Gaskets for PicoNozzle 172	II-2400 34	IT-23 38, 49	M1 137	MEH3FW45110
GC2250	II-2500 35	IV-ECU-H5	M1L	MEH3FW45xx
GC24 50	II-2888 34	IV-ECU-HC 150	M2 133, 135	MEH3R
GC26	II-2393134	IV-ECU-HCP 150		MEH3RF110
Glass Cuvettes	II-MRBP-M	IV-ECU-HOC	M-3 127, 132, 135	MEH3RFW 110
	II-MRBP-R	_	M4C 129, 132, 133, 134, 135	MEH3RFWxx 109
Glass Handling Forceps	II-PE34	J	M5 133, 135, 137	MEH3RFxx109
glass micropipettes	ILS005LT	Joystick control 128, 134	M6 133, 135	MEH3RW 110
glass rod	ILS010LT	JUV 67	M7	MEH3RWxx 109
Glue Gun , 97	.2301021	•	M8 137	

Page

Item

Page

Item

Item

Page

Item

MEH3Rxx 109 Micro pH Electrodes 64 Multi-barrel micropipette coupling. Patch Clamp 113 PSMT5. MEH3SS 110 Micropipette Beveler 120 115 Paw Pressure Meter 35 PT0110 MEH3SBW.x 109 Micropipette Holders 108 Multi-Channel Stimulator 89 PB150F-4 114 PT0203 MEH3SFW.x 109 Micropipette Storage 114 Multi-Channel Stimulator 89 PB150F-6 114 PT0402 MEH3SFW.x 101 Micropositioner 130 muscle hyperalgesia 35 perfusion chambers 84 PTM23CC01NS MEH3SFW.x 109 Microprobe Thermometers 48 PERIPRO-2HS 159 PTP401 MEH3SWx.x 109 Micro-Scale Compression 15 N PERIPRO-2HS 159 PTP406 MEH3SWx.x 109 microscope stage 155 NANOFIL 178 PERIPRO-8LS 159 PTPxxx MEH3SW45 111 microscope stage 133 N	107 107 107 105 105 107 107 107 107 107 107 107 107 122 123 88, 89 Stimu 89 Stimu 156 167 173
MEH3S. 110 Micropipette Beveler. 120 115 Paw Pressure Meter. 35 PT0110 MEH3SBW. 111 Micropipette Calibration 118 Multi-Channel Stimulator. 89 PB150F-4 114 PT0402 MEH3SFW. 109 Micropipette Storage. 114 Multi-Channel Stimulator. 89 PB150F-4 114 PT0402 MEH3SFW. 111 Micropipette Storage. 114 Multipipette Puller. 123 PB150F-6 114 PTM3CCO2INS MEH3SFW. 110 Microprobet Thermometers. 48 Muscle Research System. 4 Perfusion System. 160 PTP201 MEH3SFW.x. 109 Micro-Reference Electrode. 65 N PERIPRO-2HS. 159 PTP401 MEH3SW.x. 109 Micro-Scale Compression. 15 NANOFIL. 178 PERIPRO-4HS. 159 PTP201 MEH3SWx. 109 microscope stage. 133 NanoFil100 178 PERIPRO-4LS. 159 PTD406 MEH3SWx.	107 107 107 105 105 107 107 107 107 107 107 107 107 122 123 88, 89 Stimu 89 Stimu 156 167 173
MEH3SBW 111 Micropipette Calibration 118 Multichannel Perfusion System 160 Paw Volume Meter 35 PT0203 MEH3SBWxx 109 Micropipette Holders 108 Multi-Channel Stimulator 89 PB150F-4 114 PT0402 MEH3SFW 111 Micropositioner 130 Multipipette Puller 123 PB150F-6 114 PTM3CC02INS MEH3SFWX 109 Microprosesure System 79 Muscle Research System 4 Perfusion chambers 84 PTM23CC001NON MEH3SFWX 109 Microprobe Thermometers 48 PERIPRO-2HS 159 PTP401 MEH3SWX 109 Micro-Reference Electrode 65 N PERIPRO-2HS 159 PTP406 MEH3SWX 109 Micro-Scale Compression 15 NANOFIL 178 PERIPRO-4LS 159 PTPxxx MEH3SWX 109 microscope Stage 133 NANOFIL-100 178 PERIPRO-4LS 159 PTT0110 MEH3W45xx 109 Micr	107 107 105 105 107 107 107 107 107 107 107 107 122 123 88, 89 Stimu 156 156 167 173
MEH3SBWxx 109 Micropipette Holders 108 Multi-Channel Stimulator 89 PB150F-4 114 PT0402 MEH3SF 111 Micropipette Storage 114 Multipipette Puller 123 PB150F-6 114 PTM3CC02INS MEH3SFWx 109 Microprossure System 79 Muscle Research System 4 Perfusion Chambers 84 PTM23CC02INS MEH3SFWx 109 Microprobe Thermometers 48 PERIPRO-2HS 159 PTP401 MEH3SW. 111 Micro-Reference Electrode 65 N PERIPRO-4HS 159 PTP406 MEH3SWx 109 Micro-Scale Compression 15 NANOFIL 178 PERIPRO-4LS 159 PTP406 MEH3SWx 109 microscope stage 133 NANOFIL 178 PERIPRO-4LS 159, 185, 195 PTPxxxx MEH3W45xx 109 Microscope Stage Adapter 133 NanoFil Application Kits 180 Peristatic pump 158 PTTxxxx MEH6RFW 111	107 105 105 107 107 107 107 107 107 107 107 122 123 88, 89 Stimu 156 156 167 173
MEH3SF. 111 Micropipette Storage. 114 Multipipette Puller. 123 PB150F-6. 114 PTM3CC02INS MEH3SFW. 111 Micropositioner	105 107 107 107 107 107 107 107 107 107 122 123 88,89 Stimu 89 156 167 173
MEH3SFW 111 Micropositioner 130 muscle hyperalgesia 35 perfusion chambers 84 PTM23CC001NON MEH3SFWxx 109 Microprobe Thermometers 48 PERIPRO-2HS 159 PTP401 MEH3SFxx 109 Micro-Reference Electrode 65 N PERIPRO-2HS 159 PTP401 MEH3SWxx 109 Micro-Scale Compression 15 NANOFIL 178 PERIPRO-4HS 159 PTPxxx MEH3SWx 109 microscope slides 155 NANOFIL-100 178 PERIPRO-8LS 159, 185, 195 PTPxxx MEH3W45 111 microscope stage 133 NanoFil Application Kits 180 Peristaltic pump 158 PTTxxxx MEH6RF 111 Microscope Stage Adapter 133 Nanoliter 2010 174, 177 Peri-Star Pro 158 PTTxxxx MEH6RFWx 111 MicroScquisher 15 NanOliter Injector 177 PE Tubing, Flexible 50 PTxxxxx MEH6RFWx 109 Micr	105 107 107 107 107 107 107 107 107 107 122 123 88, 89 Stimu 89 156 167 173 172
MEH3SFWxx. 109 Micropressure System	107 107 107 107 107 107 107 107 107 122 123 88, 89 Stimu 89 156 167 173 172
MEH3SFxx. 109 Microprobe Thermometers. 48 N PERIPRO-2HS. 159 PTP401 MEH3SW. 111 Micro-Scale Compression. 15 N PERIPRO-4HS. 159 PTP406 MEH3SWxx. 109 Micro-Scale Compression. 15 NANOFIL. 178 PERIPRO-4LS. 159, 185, 195 PTPxxx. MEH3SWx. 109 microscope slage. 155 NANOFIL. 178 PERIPRO-4LS. 159, 185, 195 PTPxxx. MEH3W45x. 109 Microscope stage Adapter. 133 NanoFil Application Kits. 180 Peristaltic pump. 158 PTTxxxx. MEH6RF. 111 Microscope Stage Adapter. 133 NanoIlter 2010. 174, 177 Peri-Star Pro. 158 PTTxxxx. MEH6RFW. 111 MicroSquisher. 15 NANOLITER2010. 177 PE Tubing, Flexible. 50 PTxxxx. MEH6RFWx. 109 MicroVolume Syringes. 15 Nanoliter Injector. 177 PG10150-4 113 PUL-1000. MEH6S	107 107 107 107 107 107 122 123 88, 89 Stimu 89 156 167 173
MEH3SW. 111 Micro-Reference Electrode 65 N PERIPRO-4HS 159 PTP406 MEH3SWxx. 109 Micro-Scale Compression. 15 NANOFIL 178 PERIPRO-4LS 159, 185, 195 PTPxxx MEH3Sxx. 109 microscope slades 155 NANOFIL-100 178 PERIPRO-8LS 159 PTT0110 MEH3W45xx. 109 Microscope Stage Adapter 133 NanoFil Application Kits 180 Peri-Star Pro 158 PTT0203 MEH6RF 111 Micro Sensors 51 NANOLITER2010 177 PE Tubing, Flexible 50 PTxxxx MEH6RFWx 111 MicroSquisher 15 Nanoliter Injector 177 PG Tubing, Flexible 50 PTxxxx MEH6RFWx 109 MicroValume Syringes 171 Nanoliter Injector 177 PG10150-4 113 PUL-1000 MEH6RFX 109 Microvolume Syringes 171 Needle Holders 23 PG52151-4 113 Puller, Programmable MEH	107 107 107 107 107 122 123 88, 89 Stimu 89 156 167 173
MEH3Sxx. 109 microscope slides 155 NANOFIL-100 178 PERIPRO-8LS 159 PTT0110 MEH3W45. 111 microscope stage. 133 NanoFil Application Kits. 180 Peristaltic pump. 158 PTT0203 MEH3W45xx. 109 Microscope Stage Adapter. 133 Nanoliter 2010 174, 177 Peri-Star Pro. 158 PTTxxxx. MEH6RF. 111 Micro Sensors. 51 NANOLITER2010 177 PE Tubing, Flexible 50 PTxxxx. MEH6RFWx. 111 MicroSquisher. 15 Nanoliter Injector 177 PG10150-4 113 PUL-1000 MEH6RFWxx. 109 MicroVolume Syringes. 171 Needle Holders. 23 PG52151-4 113 Puller, Programmable. MEH6SFW. 111 Millivolt & megohm meter 77 NEEDLE MICROPROBES. 38, 49 PG52165-4 113 Pulsemaster™Multi-Channe MEH6SFW. 111 Miniature Micropositioners. 131 Neurotransmitter Detection. 58 pH Meter. <td> 107 107 107 107 122 123 88, 89 Stimu 89 156 167 173</td>	107 107 107 107 122 123 88, 89 Stimu 89 156 167 173
MEH3Sxx. 109 microscope slides. 155 NANOFIL-100. 178 PERIPRO-8LS. 159 PTT0110. MEH3W45. 111 microscope stage. 133 NanoFil Application Kits. 180 Peristaltic pump. 158 PTT0203. MEH3W45xx. 109 Microscope Stage Adapter. 133 Nanoliter 2010. 174, 177 Peri-Star Pro. 158 PTTxxxx. MEH6RF. 111 Micro Sensors. 51 NANOLITER2010. 177 PE Tubing, Flexible. 50 PTxxxx. MEH6RFWx. 111 MicroSquisher. 15 NanoSensor. 52 PG10150-4 113 PUL-1000. MEH6RFWx. 109 Microvolume Syringes. 171 Needle Holders. 23 PG52151-4 113 Puller, Programmable. MEH6SFW. 111 Miniature Micropositioners. 131 Needle Holders. 23 PG52165-4 113 Pulsemaster™Multi-Channo. MEH6SFW. 111 Miniature Micropositioners. 131 Neurotransmitter Detection. 58 pH Meter.	107 107 107 122 123 88, 89 Stimu- 89 156 167 173
MEH3W45 111 microscope stage 133 NanoFil Application Kits 180 Peristaltic pump 158 PTT0203 MEH3W45xx 109 Microscope Stage Adapter 133 Nanoliter 2010 174, 177 Peri-Star Pro 158 PTTxxxx MEH6RF 111 Micro Sensors 51 NANOLITER2010 177 PE Tubing, Flexible 50 PTxxxx MEH6RFW 111 MicroSquisher 15 Nanoliter Injector 177 PG10150-4 113 PUL-1000 MEH6RFW.x 109 Microvolume Syringes 171 Needle Holders 23 PG52151-4 113 Puller, Programmable MEH6SF 111 Millivolt & megohm meter 77 NEEDLE MICROPROBES 38, 49 PG52165-4 113 Pulsemaster™Multi-Channe MEH6SFW 111 Miniature Micropositioners 131 Neurotransmitter Detection 58 pH Meter 63 lator MEH6SFW.x 109 Miniature Peristaltic Pump	107 107 122 123 88, 89 Stimu 89 156 167 173
MEH6RF. 111 Micro Sensors 51 NANOLITER2010 177 PE Tubing, Flexible 50 PTxxxx MEH6RFW. 111 MicroSquisher. 15 Nanoliter Injector 177 PG10150-4 113 PUL-1000 MEH6RFW.xx. 109 MicroVolume Syringes 171 Needle Holders 23 PG52151-4 113 Puller, Programmable MEH6SF. 111 Millivolt & megohm meter 77 NEEDLE MICROPROBES 38, 49 PG52165-4 113 Pulsemaster™Multi-Channe MEH6SFW. 111 Miniature Micropositioners 131 Neurotransmitter Detection 58 pH Meter 63 lator MEH6SFWxx 109 Miniature Peristaltic Pump 157 NF26BV-2 178 pH Meter, Fiber Optic 63 Pumps & Fluid Handling MEH6SFxx 109 Mini DipTip 191 NF33-36BL 179 pH mini sensors 63 Push-Pull Pump	107 122 123 88, 89 Stimu 89 156 167 173 172
MEH6RF	122 123 88, 89 Stimu- 89 156 167 173 172
MEH6RFWxx. 109 MicroTip. 62 NanoSensor 52 PG10165-4 113 Puller, Programmable MEH6RFxx. 109 Microvolume Syringes 171 Needle Holders 23 PG52151-4 113 pulse generator MEH6SF. 111 Millivolt & megohm meter 77 NEEDLE MICROPROBES 38, 49 PG52165-4 113 Pulsemaster™Multi-Channe MEH6SFW. 111 Miniature Micropositioners 131 Neurotransmitter Detection 58 pH Meter 63 lator MEH6SFWxx 109 Miniature Peristaltic Pump 157 NF26BV-2 178 pH Meter, Fiber Optic 63 Pumps & Fluid Handling MEH6SFxx 109 Mini DipTip 191 NF33-36BL 179 pH mini sensors 63 Push-Pull Pump	123 88, 89 Stimu- 89 156 167 173
MEH6RFxx	88, 89 Stimu- 89 156 167 173 172
MEH6SF	Stimu- 89 156 167 173 172
MEH6SFW	89 156 167 173 172
MEH6SFWxx	156 167 173 172
MEH6SFxx	167 173 172
	173 172
MEH/ 111 MiniFlow 67 MEGG 26DV 170 PH-OVITCA-MICVO 62 DV/970	172
111 33 3001	
MEH7W	147
MEH7Wxx	
MEH7xx	
MEH8	
MEH8xx	
MEH145	
MEH145xx	
MEH345	
MEH345xx	
177 00	
ALFILOROG ALA ALA ALA ALA ALA ALA ALA ALA ALA AL	
MES	
MF28G-5	
MF28G67-5	
MF34G-5	140
MF200-1	
MF200-2	199
MF200-H2	
MF200-H3	
MF200-H4	5.
MF200-M1	
MF200-M2	<i>A</i> 1 1 <i>A</i> 7
MEZUU MICROTORGE	
MICRO4	Meter
Micro4 Controller	WICCCI.
MicroBeveler	26
MicroC	
Micro Cannula	
Micro-Cautery	87
Microcentrifuge	87
Microcentrifuge tube,	
MICROC-P	
Microdialysis Pumps	87
Microelectrode Beveler	
Microelectrodes	87
MicroFil	
MicroFlow	
Microforge	
MicroImplant	
Microiontophoresis	
Microliter Sample Holder	
MicroLWCC	
Micromanipulator 128, 129, 132 MT-26/4	
Missansky AF NT 20/4 20 40 Paradalasa C7	
ON/ACII 2013013	106
0/1-WICKO-AO1	
ACC TOTAL TO THE TAX A STATE OF	
NECCUTIFICATION AT D. 20 40	
Nulti Parrel Class 114	
24 Multi-barrel Glass 114 Parallel Rail	

Item	Page	Item	Page	Item	Page	Item	Page	Item	Page
Reflex Clip Applier	27	SI-KG2B	12	SSM33A70	105	TAXIC-350	11	TM33A10KT	105
REG2		SI-KG4		SSM33A120		TAXIC-600		TM33A20	
REG5K		SI-KGxx	13	SST30407-50		TAXIC-603		TM33B01	
REG10		SILFLEX-2		SSTxxxxx		TAXIC-650		TM33B01KT	
REG10K		Silicone RTV		Stage Micrometer		TAXIC-653		TM33B05	
REG20		Silicone Tubing		Stainless Steel Rods		TB-1		TM33B05KT	
REG50		Silver Epoxy		Stainless Steel Wire		TBM4M		TM33B10	
REG100		Silver Wire		Steel Base Plate		TBR1025		TM33B10KT	
REG200		SI-MKB		Stereotaxic	. 42, 43	TBR4100		TM33B20	
REG1000		SI-MKBM-OXY		Stereotaxic Accessories	46	TBR4100-416		TM33BxxKT	
REGAL Pipetters		SI-MKB-OXY		Stereotaxic Frame		TBR4100-424T		TM33C05	
REGS5		SI-MOTDB	12	Stereo Zoom Microscope	140, 142	TBS	129, 133	TM33C10	105
REGS6		SI-MOT-MKB	12	Sterilizer	29	TEER measurement		TM33CCINS	105
REGS7	168	SI-MOT-MT	12	Sterilizing Trays	29	Temperature Controller.	38	TM33CCNON	105
REMS-24	80	Single Barrel	112	stimulator	. 88, 89	Temperature Control Uni		TM53CCINS	105
REMS-96	80	SI-PF100		Stimulator/Isolator	90	temperature probe		TPP Electrode	66
REMS AutoSampler	80	SI-SARCAM		Stimulus Isolator 9	1, 92, 93	Temperature Probes		Transbridge	
Replacement Electrode Holde		SI-TCM2		STM3		Temperature Sensor		Transducer Amplifier	
Replacement Needles		SI-TM1		Stopcock		TGW0515		Transmission probes	
Replacement Sleeves		SI-TM2		Straight Forceps		TGW1510		Trio	
Replacement Tubing		SI-TM3		STX2		TGWxxxx		Troutman Needle Holder	25
RET-2		SI-TM4		STX3		Thin Wall		TST33A001KT	
RET-3		SI-TM5		STX100		Three-barrel glass		TST33A05KT	
Retractors		SI-TM6		STX100C		Tidas	100	TST33A10KT	
Ring Light		SI-TM7		STX100C96		TIDAS-D2		TST33A20KT	
Ding Light Adaptor	1.47								
Ring Light Adapter		SI-TM8		STX100F		TidasDAQ		TST33AxKT	
Ring Light Guide		SI-TM9	13	STX100M		TIDAS-H		TST33C05KT	
RN001		SI-TM10	13	Sub-microliter injection		TIDAS-H2		TST53A10KT	
RN0005	1/1	SI-TM11		Super Adhesive		Tilt Base	,	TST150-6	
RN010	1/1	SM325		SuperCut		Tilt Base Assembly		TTL Control Module	
RN025		SM325-M		SUPER-Dri-Ref		TIP1TW1		Tubing Cartridge	
Rodent Brain Matrices		SNAP		Sure-Seal Chambers		TIP01TW1F		tungsten light source	
Roger Wirecutting Scissors		SNAP25		SurgioScope	146	TIP01TW1F-L		Tungsten Microelectrodes	
RPE-KIT	180	SNAP50 52, 5		SYLG184	94	TIP1TW1-L		Tungsten Wire	
RTD sensor		SNAP100	67	Sylgard	96	TIP2TW1	116	TW100-3	112
RTV Coating	97	S-nitrosoglutathione	67	Syringe Pump	162, 164	TIP02TW1F		TW100-4	
RTV Prime Coat,	97	S-Nitroso-N-acetyl-D,L-peni	icillamine	Syringe Pumps	161, 166	TIP02TW1F-L	116	TW100-6	
RTV Sealant	97	67		Syringes with Luers		TIP2TW1-L	116	TW100F-3	112
		SN-PCZ-50L	125	SYS-260	77	TIP03TW1F	116	TW100F-4	
S		SN-PCZ-50R	125	SYS-705	69	TIP03TW1F-L	116	TW100F-6	112
sample cell	186	SN-PZ-50	125	SYS-773	71	TIP04TW1F	116	TW120-3	
Sample holder, Microliter		SN-PZ-50L	125	SYS-900A	79	TIP5TW1		TW120-4	112
Sarcomere Spacing Module		SN-PZ-50R	125	SYS-A300	89	TIP05TW1F	116	TW120-6	112
Scalpel Blades		SP100i		SYS-A310	88	TIP05TW1F-L	116	TW120F-3	112
Scalpel Handle		SP100iZ	167	SYS-A320D	90	TIP5TW1-L	116	TW120F-4	112
Scissors	20		166, 167	SYS-A320R	90	TIP5TW1I S01	116	TW120F-6	112
Scotch-Weld		SP101iZ	167	SYS-A3629	0 91 93	TIP5TW1LS02	116	TW150-3	112
Screwdriver Set		SP120p		SYS-A365D		TIP10TW1		TW150-4	
SDR2		SP120pZ		SYS-A365R		TIP10TW1-L		TW150-6	
Sensor Guide		SP200i		SYS-A382		TIP10TW1LS01		TW150F-3	
		SP200iZ		SYS-A385R		TIP10TW1LS02		TW150F-4	
Septum Theta		SP210c		SYS-A385RC		TIP10XV119		TW150F-6	
Seven-barrel glass		SP210cZ		SYS-BP1		TIP30TW1		Two-barrel glass	
SGE001RN		SP210iw		SYS-DAM50		TIP30TW1-L		TWU-Dattel glass	114
SGE0005RN		SP210iwZ		SYS-DAM80		TIP30TW1LS01		U	
SGE010RNS						TIP30TW1LS02		•	20
SGE025RN		SP220i		SYS-EVC4000				UBATH	
SGE050RN		SP220iZ		SYS-ISO2		TIPCA		UltraMicroPump III	
SGE050TLL		SP230iw		SYS-MICRO4		TIPH		UltraPath	
SGE100RN		SP230iwZ		SYS-MICROC		TIPK		Ultrasonic Cleaner	
SGE100TLL	171	SP250i		SYS-MS314		TIPMIX01-05		UMP3	
SGE250TLL	171	SP250iZ		SYS-OMEGAZ		TIPMIX01-05-L		UMP3-2	
SGE syringes		SP260p		SYS-PM01D		TIPMIX05-10		UMP3-3	171
SI-BAM21-LC	9	SP260pZ		SYS-PM01R		TIPMIX05-10-L		UMP3-4	171
SI-BAM21-LCB	8,9	spectrometer	188	SYS-PM015D		TIPTPP		UMS	
SI-BF-100		SPLG100		SYS-PM015R		Tissue Mounting Hooks.		Uniaxial Tissue Tester	14
SI-BMFA		SPLG101	163	SYS-PM100D	181	Titanium Forceps		Universal Focus Mount	143
SI-COLUB		SPLG110	163	SYS-PM100R		Titanium wire	107	Universal Manipulator Stand	
SI-CTS200		SPLG200	165	SYS-PV820		TM31/33Axx		Universal Pipette Tips	
Signal Conditioning Amplifiers		SPLG210		SYS-PV830		TM31/33AxxKT		UPUV	
Signal Generator		SPLG212		SYS-REMS		TM31A10		UPVIS	
SI-HTB2		SPLG270		SYS-TBM4M		TM31A10KT		USBCAM33	
SI-HTB2M		SPLG272		-		TM31A20	105	USBCAM50	
SI-HTB4		SS316xx		T		TM31C05		USBCAM133	
SI-HTB4M		SS31614		Tail Flick Test	27	TM33A05		USBCAM152	
SI-KG2		SSM33A20KT		TAXIC-300		TM33A10		USBCAM202	
JI*NUZ	13	J314133/120111	103	11 N/IC-2000	44		105	ODDCVINITATE	132

Item Page	Item Page	e Item	Page Item	Page	Item Page
USS1L	USS7L	V-Clamp Ventilator Connecti Vetbond	ion Kit	200 V	X xxxxD 85
USS2S	USS8S	VFP Vibration-Free Tab Vibration Isolation Video Camera	les	119, 145	Z ZBEECAL
USS4S 85 USS5L 85 USS5S 85 USS6L 85 USS6S 85	V V503382	Voltage Clamp Voltage Electrode von Frey Anesthesi			Z-LITE-186, 120, 147 Z-LITE-Z, 120 Z-LITE-Z186

Terms & Conditions

WPI values its customers and will make every effort to maintain the highest professional and ethical standards. However, WPI regards as invalid any special terms and conditions set forth by any other party, including purchasing, quotations, warranties or other statements. WPI's response to a quotation, purchase order or inquiry shall not be construed as acceptance of any conditions stated therein.

Credit Terms

Net 30 days from date of invoice to customers with satisfactory credit references.

Standard Equipment & Accessories

To insure against misunderstandings or delays, always include model number, name, and significant specifications of systems and accessories when ordering. Also specify any options desired. WPI retains the right to change specifications or discontinue products without prior notice and without responsibility for incorporating such changes into equipment previously sold. Please state the desired operating line voltage and frequency.

Warranty

WPI strives to maintain the highest quality standards in all products. We warranty these products against defects in workmanship or materials, and all WPI products retain a warranty that is stated in the product's manual. (Typically, one year for parts and labor except for consumable items such as electrodes, glass, etc., for which warranty varies from 30 days to six months.) However, WPI hereby specifies that no use or application of any product may be suitable or successful for a specific application and therefore waives all liability in that regard. It is the sole responsibility of the buyer/user to apply the product in a manner consistent with its intended use.

WPI products are not approved for human use unless specifically stated in writing with accompanying FDA (or applicable) documents. WPI assumes no liability or legal, moral, ethical, or fiduciary responsibility for a use of any product in the WPI catalog in treating or studying humans.

Prices

WPI maintains competitive world wide pricing. Prices listed for standard items as described are net and subject to change without notice. Prices (in U.S. dollars) do not include shipping charges. Formal quotations are valid for a period of 30 days unless otherwise specified. Purchases may be charged to your VISA, MasterCard, or American Express account.

Taxes

We reserve the right to add any taxes required by local law or ordinance. Institutions operating under tax-free conditions must state

applicable license or contract number and be prepared to furnish necessary exemption certificate.

Packaging and Shipping

Equipment is shipped in adequate commercial packing. All packages should be checked thoroughly before being discarded. Shipping and handling expenses, including adequate insurance coverage, are paid by the customer.

Claims and Returns

Inspect all shipments upon receipt. Missing cartons or obvious damage to cartons should be noted on the delivery receipt before signing. Concealed loss or damage should be reported at once to the carrier and an inspection requested. All claims for shortage or damage must be made within ten (10) days after receipt of shipment. Claims for lost shipments must be made within thirty (30) days of receipt of invoice or other notification of shipment. Please save damaged or pilfered cartons until claim is settled. In some instances, photographic documentation may be required. Some items are time-sensitive; WPI assumes no extended warranty or any liability for use beyond the date specified on the container

Do not return any goods to us without obtaining prior approval and instructions from our Returns Department. Goods returned unauthorized or by collect freight may be refused. Goods accepted for restocking will be exchanged or credited to your WPI account and are subject to a restock fee of up to 25%. Goods returned which were ordered by customers in error are subject to a 25% restock fee. Equipment which was built as a special order cannot be returned.

Repairs

Contact our Returns Department for assistance in the repair of apparatus. Do not return goods until instructions have been received. Returned items must be securely packed to prevent further damage in transit. The Customer is responsible for paying shipping expenses, including adequate insurance on all items returned for repairs. Identification of the item(s) by model number, name, serial number and proof of purchase (packing slip number or invoice number) as well as complete description of the difficulties experienced should be written on the RMA request form and a copy of the form must be included with the item.

Ordering Information

Orders from **North America** should be directed to WPI's main office:

World Precision Instruments, Inc. 175 Sarasota Center Boulevard Sarasota, Florida 34240-9258

Telephone: 941-371-1003 (collect calls not accepted). Facsimile copies of purchase orders, inquiries or other correspondence may be transmitted to WPI's Fax number: 941-377-5428. When sending written purchase orders to confirm a telephone order, please clearly mark "CONFIRMING" on the order to avoid it being shipped a second time. Our Sales Department's e-mail address is: sales@wpiinc.com

North American customers may also order many products on-line through WPI's Web site — www.wpiinc.com. All customers may use the on-line Quote Request Form and on-line Purchase Order Form.

Orders from Belgium, Denmark, Eire, Finland, Luxembourg, Norway, Portugal, Spain, Sweden, Switzerland (FR) and the United Kingdom should be directed to

WPI UK: 1 Hunting Gate

Hitchin, Hertfordshire SG4 0TJ England

Tel: +44 (0)1462 424700 Fax: +44 (0)1462 424701

E-mail: wpiuk@wpi-europe.com

Customers in **France** may call +33 (0)970 44 9000 or e-mail wpifr@wpi-europe.com.

Orders from Austria, Germany, Greece, Italy, Malta, Netherlands, Russian States, Switzerland (DE) and Turkey should be directed to WPI Germany: Zossener Strasse 55-58 D-10961 Berlin, Germany

Tel: +49 (0)30-6188845 Fax: +49 (0)30-6188670 E-mail: wpide@wpi-europe.com

Orders from China, Hong Kong and Taiwan should be directed to

WPI China: WPI Shanghai Trading Co., Ltd., Rm 20a No8 Dong Fang Rd., Lu Jia Zui Financial District Shanghai PRC

Tel: + 86 688 85517

Email: chinasales@china.wpiinc.com

Orders from Brazil should be directed to

WPI Brazil:

Av: Conselheiro Nebias, 756 sala 2611 Santos - Sao Paulo - Brasil

CEP: 11045-002

Fone: +011 55 13 40629703 Email: Info@brazil.wpiinc.com

Orders from the countries listed below should be directed to the appropriate WPI distributors. All other international orders should be directed to WPI's main office and factory in Sarasota, Florida.

To expedite shipment of orders outside North America (other than to areas served by the international offices above), payment must be made in advance via wire transfer, check payable in U.S. dollars, or charged to a credit card. Proforma Invoices will gladly be furnished upon request. Please specify proper electrical current for your order.

WPI Distributors

Argentina

Summit Research S.A

Rodrigo de Ibarrola 3140, 1D C1419CGB Ciudad de Buenos Aires, Argentina

Telefono: (5411)-4566-8246 Fax: (5411) 4638-5049 E-mail: rb@summit-research.com

Australia, New Zealand Coherent Life Sciences Pty. Ltd.

116 Sir Donald Bradman Dr. Hilton, SA 5033

Tel: (61) 8-8150-5200 Fax: (61) 8-8352-2020 E-mail: sales@coherent.com.au

Rrasi

SanBio Científica Ltda.

Av. Mauro Krepski nr. 107 – CEP. 13140-542 - Jardim Ype Paulinia - SP

> Fone: 19 3884-1293 E-mail: vendas@sanbiocientifica.com.br Website: Sanbiocientifica.com.br

Chile

Equilab Ltd

DíazYCompañía Ltd. El Quillay 627 sitio 96 Empresarial Valle Grande, Santiago Chile

Tel: (56-2)437-0217 Fax: (56-2)465-0066 Email: secvtas3@equilab.cl Web: www.equilab.cl Contact: Leonardo Martinez

Hungary, Poland MDE Heidelberg

Altrottstrasse 31 Walldorf 69190 Germany Tel: +36-87-703-001 Fax: +36-87-703-002

Website: www.experimetria.com E-mail: info@experimetria.com This Distributor also covers:

Albania, Belarus, Bosnia, Bulgaria, Croatia, Czech Republic, Estonia, Herzegovina, Latvia, Lithuania, Republic of Makedonija, Romania, Serbia and Montenegro, Slovakia, Slovenia

India

Labindia Instruments Pvt. Ltd.

Head Office: 201, Nand Chambers, L B S Marg, Near Vandana Cinema, Thane 400 602

Tel.: 022 – 2598 6062

Email: borgaonkarhv@labindia.com Website: www.labindia.com

Other Labindia offices: Chandigarh, Delhi, Gurgaon, Lucknow, Guwahati, Baroda, Bhopal, Kolkata, Pune, Hyderabad, Bangalore, Chennai, and Thiruvananthapuram

Israel

NBT NewBiotechnology Ltd.

Uri Schechter, Managing Director 3 Mekor Haim St. P.O. Box 8662 Jerusalem 91086 Tel: 972-2-6732001 Fax: 972-2-6731611

E-mail: nbtsales@nbtltd.com

Website: www.nbtltd.com

Japan LMS Co., Ltd.

Tanaka Bldg. 3-6-7 Hongo Bunkyo-Ku, Tokyo 113-0033 Tel: 03-5842-4171 Fax: 03-3818-7095 E-mail: intldpt@lms.co.jp Website: www.lms.co.jp

Physio-Tech Co., Ltd.

1-6-3 lwamoto-cho, Chiyoda-ku Tokyo 100-0032 Tel: 81-3-38642781 Fax: 81-3-38642787 E-mail: sales@physio-tech.co.jp

Korea

Scitech Korea Inc,

RM 801, 74 Gangbuk Electronics Industrial Complex, Deongneung-ro 40-gil, Gangbuk-gu, Seoul 142-705, Korea

Tel: +82 2 986 4413~5 Fax: +82 2 986 4429 Email: scitech00@korea.com

México

Alta Tecnologia en Laboratorios S.A. de C.

Sucursales:

GUADALAJARA, JAL.

Av. Inglaterra #2588
Col. Arcos Vallarta
Guadalajara, Jalisco C.P. 44150
Telefono: 52 (33) 36 16 52 04
52 (33) 36 16 24 69
Fax: 52 (33) 36 16 97 30
E-mail: altatecenlabs@prodigy.net.mx
Website: www.atlsadecv.com.mx

MEXICO D.F.

Comoporis #43, Col el Caracol Coyoacan, Mexico D.F. C.P. 04739 Telefono: 52 555 606 24 69 Y 52 04

Fax: 52 555 665 79 68 E-mail: atl@prodigy.net.mx Website: www.atlsadecv.com.mx

Singapore, Malaysia, Indonesia, Thailand & Pakistan

Prime Bioscience Pte Ltd

Block 196 Pandan Loop #07-02 Pantech Business Hub Singapore 128384

Tel: + 65 6364 0391 Fax: + 65 6269 0745 E-mail: info@primebioscience.com

Website: www.primebioscience.com

South Africa

Axiology Labs

Johannesburg, South Africa Tel.: +27 82778 5889 E-mail: info@axiologylabs.com Web: axiologylabs.com Axiology also includes: **Kenya**, **Rwanda**, **Uganda**, **Tanzania**, **Zambia**, **Mozambique**, **Zimbabwe**, **Botswana**, **Namibia**, **Nigeria**, **Ghana**

UAE & GCC

Laboratory Automation Services Middle East FZ-LLC

Dubai Biotechnology and Research Park (DuBiotech) PBU No. D025, P.O.Box 345806 Dubai,

United Arab Emirates
Tel: 971-(0)4-8180 700
Fax: 971-(0)4-8180 790
Email: info@las-me.com
Website: www.las-me.com
GCC also includes: Qatar Om

GCC also includes: **Qatar, Oman, Bahrain, Kuwait, Saudi Arabia**

Labgulf FZC

P.O.Box 9668 Saif Zone

Sharjah, United Arab Emirates Tel.: +9716 5579893 Fax: +9716 5579894 E-mail: Jahsales@eim ae

E-mail: labsales@eim.ae Website: www.labgulf.com



Instrumenting scientific ideas
International Edition